



EN_ MATERIALS COMPOSING ADLER ORTHO® IMPLANTS
IT_ MATERIALI CHE COMPONGONO GLI IMPIANTI ADLER
FR_ MATÉRIAUX COMPOSANT LES IMPLANTS ADLER ORTHO®
DE_ MATERIALIEN AUS ADLER ORTHO® IMPLANTATEN
ES_ MATERIALES QUE COMPONEN LOS IMPLANTES ADLER ORTHO®
PT_ MATERIAIS QUE COMPÕEM OS IMPLANTES ADLER ORTHO®
NL_ MATERIALEN WAARUIT ADLER ORTHO® IMPLANTATEN ZIJN OPGEBOUWD
EL_ ΥΛΙΚΑ ΠΟΥ ΣΥΝΘΕΤΟΥΝ ΤΑ ΕΜΦΥΤΕΥΜΑΤΑ ADLER ORTHO®

Material	Specifications
Stainless Steel M30NW	ISO 5832-9
Inox Steel AISI 316L	ISO 5932-1
Alumina Ceramics	ISO 6474-1 and ISO 6474-2
PMMA (Polymethyl Methacrylate)	ISO 5833
CoCrMo Alloy	ISO 5832-4 and ISO 5832-12 ASTM F75
UHMWPE (Ultra High Molecular Weight Polyethylene)	ISO 5834-1 and ISO 5834-2
Ti6Al4V Alloy	ISO 5832-3 ASTM F2924
Unalloyed Titanium	ISO 5832-2
Hydroxyapatite coating	ISO 13779-1
Titanium Niobium Nitride Coatings (TiNbN)	ISO 5832-2

EN_To find out the analytical composition of the constituent materials of the implants refer to the following:

IT_Per conoscere la composizione analitica dei materiali costitutivi impianti si rimanda a quanto segue:

FR_Pour connaître la composition analytique des matières constitutives des implants, se référer à ce qui suit:

DE_Die analytische Zusammensetzung der Inhaltsstoffe der implantate kann wie folgt ermittelt werden:

ES_Para conocer la composición analítica de los materiales constitutivos de las implantes, consúltese lo siguiente:

PT_Para descobrir a composição analítica dos materiais constituintes das implantes, consulte o seguinte:

NL_Om de analytische samenstelling van de samenstellende materialen van de implanten te weten te komen, raadpleeg je het volgende:

EL_Για την αναλυτική σύνθεση των υλικών κατασκευής των εμφυτευμάτων, παρακαλούμε ανατρέξτε στα ακόλουθα:

Material

Stainless Steel M30NW

Specifications

ISO 5832-9

Table 1 — Chemical composition

Element	Mass fraction %
Carbon	0,08 maximum
Silicon	0,75 maximum
Manganese	2,00 to 4,25
Nickel	9,0 to 11,0
Chromium	19,5 to 22,0
Molybdenum	2,0 to 3,0
Niobium	0,25 to 0,80
Sulfur	0,01 maximum
Phosphorus	0,025 maximum
Copper	0,25 maximum
Nitrogen	0,25 to 0,50
Iron	Balance
Residuals	—
Each	0,1 maximum
Total	0,4 maximum

Material

Inox Steel AISI 316L

Specifications

ISO 5932-1

Table 1 — Chemical Composition

Element	Mass fraction %
Carbon	0,030 max.
Silicon	1,0 max.
Manganese	2,0 max.
Phosphorus	0,025 max.
Sulfur	0,010 max.
Nitrogen	0,10 max.
Chromium	17,0 to 19,0 max.
Molybdenum	2,25 to 3,0
Nickel	13,0 to 15,0
Copper	0,50 max.
Iron	Balance

Material
Alumina Ceramics

Specifications
ISO 6474-1 and ISO 6474-2

Table 1 — Limits for material properties

Property	Unit	Property category	Requirement		Subclause	References
			Type A	Type B		
Average bulk density	kg/m ³ × 10 ³	1	≥ 3,94	≥ 3,90	6.1	ISO 18754 EN 623-2
Chemical composition:						
Basic material, Al ₂ O ₃	% mass fraction	1	≥ 99,7	≥ 99,5		ISO 12677 EN 725-1
Sintering additive, MgO	% mass fraction	1	≤ 0,2	≤ 0,2	6.2	
Limits of impurities, total amount of SiO ₂ + CaO + Na ₂ O	% mass fraction	1	≤ 0,1	≤ 0,3		

Table 1 — Limits for material property category 1

Property	Unit	Property category	Requirement		Subclause	References
			Type X	Type S		
Average relative bulk density	%	1	≥99	≥99	6.1	ISO 18754 EN 623-2
Chemical composition:						
Alumina, Al ₂ O ₃	% mass fraction	1	60 to 90	60 to 90	6.2	ISO 12677
Zirconia, ZrO ₂ + HfO ₂	% mass fraction	1	10 to 30	10 to 30		
Amount of HfO ₂ in ZrO ₂	% mass fraction	1	≤5	≤5		
Intended additives	% mass fraction	1	≤10	≤10		
Total amount of impurities	% mass fraction	1	≤0,2	≤0,2		

Table 2 — Limits for material property category 2

Property	Unit	Property category	Requirement		Subclause	References
			Type X	Type S		
Radioactivity (measured on raw materials)						
Zirconia Other intended additives	Bq/kg	2 See 6.5	≤200	≤200	6.5	ISO 13356

Material

PMMA (Polymethyl Methacrylate)

Specifications

ISO 5833

Table 1 — Requirements and test methods for setting properties of liquid-powder mixtures

Mixture	Doughing time			Setting time		Maximum temperature		
	Average min	Maximum deviation from average min	Test method	Average min	Test method	Average °C	Maximum deviation from average °C	Test method
Syringe usage (see clause 5)	—	—	—	6,5 to 15	Annex C	90	± 5	Annex C
Dough state usage (see 6.1)	≤ 5	1,5	Annex B	3 to 15	Annex C	90	± 5	Annex C

Table 2 — Requirements and test methods for set and polymerized cement

Average compressive strength		Bending modulus		Bending strength	
MPa	Test method	MPa	Test method	MPa	Test method
≥ 70	Annex E	≥ 1 800	Annex F	≥ 50	Annex F

Material
CoCrMo Alloy

Specifications
**ISO 5832-4 (equivalent to ASTM F75)
and ISO 5832-12**

COLED

Table 1 — Chemical composition

Element	Compositional limits % (m/m)
Chromium	26,5 to 30,0
Molybdenum	4,5 to 7,0
Nickel	1,0 max.
Iron	1,0 max.
Carbon	0,35 max.
Manganese	1,0 max.
Silicon	1,0 max.
Cobalt	Balance

FORGED

Chemical composition

Element	Mass fraction %	
	Alloy 1 – low carbon	Alloy 2 – high carbon
Chromium	26,3 – 30,0	26,3 – 30,0
Molybdenum	5,0 – 7,0	5,0 – 7,0
iron	<0,75	<0,75
Manganese	<1,0	<1,0
Silicon	<1,0	<1,0
Carbon	<0,14	0,15 to 0,35
Nickel	<1,0	<1,0
Nitrogen	<0,25	<0,25
Cobalt	Balance	Balance

Material
CoCrMo Alloy

*Cobalt-28 Chromium-6 Molybdenum Alloy Castings
and Casting Alloy for Surgical Implants*

Specifications
ASTM F75

(equivalent to ISO 5832-4)

TABLE 1 Chemical Composition

Element	Composition, % (Mass/Mass)	
	min	max
Chromium	27.00	30.00
Molybdenum	5.00	7.00
Nickel	...	0.50
Iron	...	0.75
Carbon	...	0.35
Silicon	...	1.00
Manganese	...	1.00
Tungsten	...	0.20
Phosphorous	...	0.020
Sulfur	...	0.010
Nitrogen	...	0.25
Aluminum	...	0.10
Titanium	...	0.10
Boron	...	0.010
Cobalt ^A	balance	balance

^A Approximately equal to the difference of 100 % and the sum percentage of the other specified elements. The percentage of the cobalt difference is not required to be reported.

TABLE 2 Product Analysis Tolerances^{A,B}

Element	Tolerance Under the Minimum or Over the Maximum Limit, % (Mass/Mass) ^C
Chromium	0.30
Molybdenum	0.15
Nickel	0.05
Iron	0.03
Carbon	0.02
Silicon	0.05
Manganese	0.03
Tungsten	0.04
Phosphorous	0.005
Sulfur	0.003
Nitrogen	0.02 ^D
Aluminum	0.02
Titanium	0.02
Boron	0.002

^A See Test Methods **E354**.

^B Refer to AMS Standard 2269 for chemical check analysis limits (except nitrogen).

^C For elements in which only a maximum percentage is indicated, the "under minimum limit" is not applicable.

^D Refer to AMS 2248 for chemical check analysis limits.

Material

UHMWPE (Ultra High Molecular Weight Polyethylene)

Specifications

ISO 5834-1 and ISO 5834-2

Table 2 — Maximum ash and trace element content

Element	Maximum quantity permitted mg/kg			Test method according to subclause
	Type 1	Type 2	Type 3 ^a	
Ash	125	125	300	8.3
Titanium	40	40	150	8.4
Calcium	5	5	50	8.4
Chlorine	30	30	90	8.4
Aluminium	20	20	100	8.4
^a Type 3 polymer is no longer manufactured. However, in order to cover existing supplies held in stockpile, this Type 3 material is retained in this document until the next revision.				

Material
Ti6Al4V Alloy

Specifications
ISO 5832-3

Table 1 — Chemical composition

Element	Compositional limits mass fraction %
Aluminium	5,5 to 6,75
Vanadium	3,5 to 4,5
Iron	0,3 max.
Oxygen	0,2 max.
Carbon	0,08 max.
Nitrogen	0,05 max.
Hydrogen	0,015 max. ^a
Titanium	Balance
^a Except for billets, for which the maximum hydrogen content shall be of a mass fraction of 0,010 %.	

Material
Ti6Al4V Alloy

Specifications
ASTM F2924

*Additive Manufacturing Titanium-6 Aluminum-4 Vanadium with Powder Bed
Fusion such as electron beam melting and laser melting*

TABLE 1 Composition

Element	min	max
Aluminum	5.50	6.75
Vanadium	3.50	4.50
Iron	—	0.30
Oxygen	—	0.20
Carbon	—	0.08
Nitrogen	—	0.05
Hydrogen	—	0.015
Yttrium	—	0.005
Other elements, each	—	0.10
Other elements, total	—	0.40
Titanium	remainder	

TABLE 2 Check Analysis Tolerances

Element	Permissible Variation in Check Analysis
Aluminum	± 0.40
Vanadium	± 0.15
Iron	± 0.10
Oxygen	± 0.02
Carbon	± 0.02
Nitrogen	± 0.02
Hydrogen	± 0.002
Yttrium	± 0.0006
Other Elements, each	± 0.02

Material
Unalloyed Titanium

Specifications
ISO 5832-2

Table 1 — Chemical composition

Element	Maximum compositional limits				
	percent mass fraction				
	Grade 1 ELI	Grade 1	Grade 2	Grade 3	Grades 4A and 4B
Nitrogen	0,012	0,03	0,03	0,05	0,05
Carbon	0,03	0,08	0,08	0,08	0,08
Hydrogen	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a
Iron	0,10	0,20	0,30	0,30	0,50
Oxygen	0,10	0,18	0,25	0,35	0,40
Titanium	Balance	Balance	Balance	Balance	Balance

^a Except for billets, for which the maximum hydrogen content shall be 0,010 0 % (mass fraction) and for flat products for which the maximum hydrogen content shall be 0,015 % (mass fraction).

Material

Hydroxyapatite coating

Specifications

ISO 13779-1

ASTM F 1185

Table 1-Limits of specific trace elements

trace element	maximum limit (mg/kg)
Arsenic	3 ppm
Cadmium	5 ppm
Mercury	5 ppm
Lead	30 ppm
Total heavy metal	50 ppm

Material

Titanium Niobium Nitride Coatings (TiNbN)

Specifications

ISO 5832-2

Table 1 — Chemical composition

Element	Maximum compositional limits				
	percent mass fraction				
	Grade 1 ELI	Grade 1	Grade 2	Grade 3	Grades 4A and 4B
Nitrogen	0,012	0,03	0,03	0,05	0,05
Carbon	0,03	0,08	0,08	0,08	0,08
Hydrogen	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a	0,012 5 ^a
Iron	0,10	0,20	0,30	0,30	0,50
Oxygen	0,10	0,18	0,25	0,35	0,40
Titanium	Balance	Balance	Balance	Balance	Balance

^a Except for billets, for which the maximum hydrogen content shall be 0,010 0 % (mass fraction) and for flat products for which the maximum hydrogen content shall be 0,015 % (mass fraction).

Devices reference codes and related material.

Codici di riferimento dei dispositivi e relativo materiale.

Codes de référence des dispositifs et du matériel associé.

Referenzcodes der Geräte und des zugehörigen Materials.

Códigos de referencia de los dispositivos y material relacionado.

Códigos de referência de dispositivos e material relacionado.

Referentiecodes en aanverwant materiaal voor apparaten.

Κωδικοί αναφοράς συσκευών και σχετικό υλικό.

Reference Code	Material
0101101	Ti6Al4V Alloy + Hydroxyapatite
0101102	Ti6Al4V Alloy + Hydroxyapatite
0101103	Ti6Al4V Alloy + Hydroxyapatite
0101104	Ti6Al4V Alloy + Hydroxyapatite
0101105	Ti6Al4V Alloy + Hydroxyapatite
0101106	Ti6Al4V Alloy + Hydroxyapatite
0101107	Ti6Al4V Alloy + Hydroxyapatite
0101108	Ti6Al4V Alloy + Hydroxyapatite
0101111	Ti6Al4V Alloy + Hydroxyapatite
0101112	Ti6Al4V Alloy + Hydroxyapatite
0101113	Ti6Al4V Alloy + Hydroxyapatite
0101114	Ti6Al4V Alloy + Hydroxyapatite
0101115	Ti6Al4V Alloy + Hydroxyapatite
0101116	Ti6Al4V Alloy + Hydroxyapatite
0101117	Ti6Al4V Alloy + Hydroxyapatite
0101118	Ti6Al4V Alloy + Hydroxyapatite
0101201	Ti6Al4V Alloy + Hydroxyapatite
0101202	Ti6Al4V Alloy + Hydroxyapatite
0101203	Ti6Al4V Alloy + Hydroxyapatite
0101204	Ti6Al4V Alloy + Hydroxyapatite
0101205	Ti6Al4V Alloy + Hydroxyapatite
0101206	Ti6Al4V Alloy + Hydroxyapatite
0101207	Ti6Al4V Alloy + Hydroxyapatite
0101208	Ti6Al4V Alloy + Hydroxyapatite
0101211	Ti6Al4V Alloy + Hydroxyapatite
0101212	Ti6Al4V Alloy + Hydroxyapatite
0101213	Ti6Al4V Alloy + Hydroxyapatite
0101214	Ti6Al4V Alloy + Hydroxyapatite
0101215	Ti6Al4V Alloy + Hydroxyapatite
0101216	Ti6Al4V Alloy + Hydroxyapatite
0101217	Ti6Al4V Alloy + Hydroxyapatite

0101218	Ti6Al4V Alloy + Hydroxyapatite
0101401	Ti6Al4V Alloy + Hydroxyapatite
0101402	Ti6Al4V Alloy + Hydroxyapatite
0101403	Ti6Al4V Alloy + Hydroxyapatite
0101404	Ti6Al4V Alloy + Hydroxyapatite
0101405	Ti6Al4V Alloy + Hydroxyapatite
0101406	Ti6Al4V Alloy + Hydroxyapatite
0101407	Ti6Al4V Alloy + Hydroxyapatite
0101408	Ti6Al4V Alloy + Hydroxyapatite
0101411	Ti6Al4V Alloy + Hydroxyapatite
0101412	Ti6Al4V Alloy + Hydroxyapatite
0101413	Ti6Al4V Alloy + Hydroxyapatite
0101414	Ti6Al4V Alloy + Hydroxyapatite
0101415	Ti6Al4V Alloy + Hydroxyapatite
0101416	Ti6Al4V Alloy + Hydroxyapatite
0101417	Ti6Al4V Alloy + Hydroxyapatite
0101418	Ti6Al4V Alloy + Hydroxyapatite
0101501	Ti6Al4V Alloy + Hydroxyapatite
0101502	Ti6Al4V Alloy + Hydroxyapatite
0101503	Ti6Al4V Alloy + Hydroxyapatite
0101504	Ti6Al4V Alloy + Hydroxyapatite
0101505	Ti6Al4V Alloy + Hydroxyapatite
0101506	Ti6Al4V Alloy + Hydroxyapatite
0101507	Ti6Al4V Alloy + Hydroxyapatite
0101508	Ti6Al4V Alloy + Hydroxyapatite
0101511	Ti6Al4V Alloy + Hydroxyapatite
0101512	Ti6Al4V Alloy + Hydroxyapatite
0101513	Ti6Al4V Alloy + Hydroxyapatite
0101514	Ti6Al4V Alloy + Hydroxyapatite
0101515	Ti6Al4V Alloy + Hydroxyapatite
0101516	Ti6Al4V Alloy + Hydroxyapatite
0101517	Ti6Al4V Alloy + Hydroxyapatite
0101518	Ti6Al4V Alloy + Hydroxyapatite
0101601	Ti6Al4V Alloy
0101602	Ti6Al4V Alloy
0101603	Ti6Al4V Alloy
0101604	Ti6Al4V Alloy
0101605	Ti6Al4V Alloy
0101606	Ti6Al4V Alloy
0101607	Ti6Al4V Alloy
0101608	Ti6Al4V Alloy
0101611	Ti6Al4V Alloy
0101612	Ti6Al4V Alloy

0101613	Ti6Al4V Alloy
0101614	Ti6Al4V Alloy
0101615	Ti6Al4V Alloy
0101616	Ti6Al4V Alloy
0101617	Ti6Al4V Alloy
0101618	Ti6Al4V Alloy
0101701	Ti6Al4V Alloy + Hydroxyapatite
0101702	Ti6Al4V Alloy + Hydroxyapatite
0101703	Ti6Al4V Alloy + Hydroxyapatite
0101704	Ti6Al4V Alloy + Hydroxyapatite
0101705	Ti6Al4V Alloy + Hydroxyapatite
0101706	Ti6Al4V Alloy + Hydroxyapatite
0101707	Ti6Al4V Alloy + Hydroxyapatite
0101708	Ti6Al4V Alloy + Hydroxyapatite
0101711	Ti6Al4V Alloy + Hydroxyapatite
0101712	Ti6Al4V Alloy + Hydroxyapatite
0101713	Ti6Al4V Alloy + Hydroxyapatite
0101714	Ti6Al4V Alloy + Hydroxyapatite
0101715	Ti6Al4V Alloy + Hydroxyapatite
0101716	Ti6Al4V Alloy + Hydroxyapatite
0101717	Ti6Al4V Alloy + Hydroxyapatite
0101718	Ti6Al4V Alloy + Hydroxyapatite
0101801	Ti6Al4V Alloy + Hydroxyapatite
0101802	Ti6Al4V Alloy + Hydroxyapatite
0101803	Ti6Al4V Alloy + Hydroxyapatite
0101804	Ti6Al4V Alloy + Hydroxyapatite
0101805	Ti6Al4V Alloy + Hydroxyapatite
0101806	Ti6Al4V Alloy + Hydroxyapatite
0101807	Ti6Al4V Alloy + Hydroxyapatite
0101808	Ti6Al4V Alloy + Hydroxyapatite
0101811	Ti6Al4V Alloy + Hydroxyapatite
0101812	Ti6Al4V Alloy + Hydroxyapatite
0101813	Ti6Al4V Alloy + Hydroxyapatite
0101814	Ti6Al4V Alloy + Hydroxyapatite
0101815	Ti6Al4V Alloy + Hydroxyapatite
0101816	Ti6Al4V Alloy + Hydroxyapatite
0101817	Ti6Al4V Alloy + Hydroxyapatite
0101818	Ti6Al4V Alloy + Hydroxyapatite
0101901	Ti6Al4V Alloy + Hydroxyapatite
0101902	Ti6Al4V Alloy + Hydroxyapatite
0101903	Ti6Al4V Alloy + Hydroxyapatite
0101904	Ti6Al4V Alloy + Hydroxyapatite
0101905	Ti6Al4V Alloy + Hydroxyapatite

0101906	Ti6Al4V Alloy + Hydroxyapatite
0101907	Ti6Al4V Alloy + Hydroxyapatite
0101908	Ti6Al4V Alloy + Hydroxyapatite
0101911	Ti6Al4V Alloy + Hydroxyapatite
0101912	Ti6Al4V Alloy + Hydroxyapatite
0101913	Ti6Al4V Alloy + Hydroxyapatite
0101914	Ti6Al4V Alloy + Hydroxyapatite
0101915	Ti6Al4V Alloy + Hydroxyapatite
0101916	Ti6Al4V Alloy + Hydroxyapatite
0101917	Ti6Al4V Alloy + Hydroxyapatite
0101918	Ti6Al4V Alloy + Hydroxyapatite
0102101	Ti6Al4V Alloy
0102102	Ti6Al4V Alloy
0102103	Ti6Al4V Alloy
0102104	Ti6Al4V Alloy
0102105	Ti6Al4V Alloy
0102106	Ti6Al4V Alloy
0102107	Ti6Al4V Alloy
0102108	Ti6Al4V Alloy
0102109	Ti6Al4V Alloy
0102300	Ti6Al4V Alloy
0102301	Ti6Al4V Alloy
0102302	Ti6Al4V Alloy
0102303	Ti6Al4V Alloy
0102304	Ti6Al4V Alloy
0102305	Ti6Al4V Alloy
0102306	Ti6Al4V Alloy
0102307	Ti6Al4V Alloy
0102308	Ti6Al4V Alloy
0102309	Ti6Al4V Alloy
0102391	Ti6Al4V Alloy
0102401	Ti6Al4V Alloy
0102402	Ti6Al4V Alloy
0102403	Ti6Al4V Alloy
0102404	Ti6Al4V Alloy
0102405	Ti6Al4V Alloy
0102406	Ti6Al4V Alloy
0102407	Ti6Al4V Alloy
0102408	Ti6Al4V Alloy
0102409	Ti6Al4V Alloy
0105113	Ti6Al4V Alloy
0105114	Ti6Al4V Alloy
0105115	Ti6Al4V Alloy

0105116	Ti6Al4V Alloy
0105117	Ti6Al4V Alloy
0105118	Ti6Al4V Alloy
0105119	Ti6Al4V Alloy
0105120	Ti6Al4V Alloy
0105121	Ti6Al4V Alloy
0105122	Ti6Al4V Alloy
0105124	Ti6Al4V Alloy
0105126	Ti6Al4V Alloy
0106008	Ti6Al4V Alloy
0106009	Ti6Al4V Alloy
0106010	Ti6Al4V Alloy
0106011	Ti6Al4V Alloy
0106012	Ti6Al4V Alloy
0106013	Ti6Al4V Alloy
0106014	Ti6Al4V Alloy
0106015	Ti6Al4V Alloy
0106016	Ti6Al4V Alloy
0106017	Ti6Al4V Alloy
0106018	Ti6Al4V Alloy
0106019	Ti6Al4V Alloy
0106020	Ti6Al4V Alloy
0106201	Ti6Al4V Alloy
0106202	Ti6Al4V Alloy
0106203	Ti6Al4V Alloy
0106204	Ti6Al4V Alloy
0106205	Ti6Al4V Alloy
0106206	Ti6Al4V Alloy
0106207	Ti6Al4V Alloy
0106208	Ti6Al4V Alloy
0106209	Ti6Al4V Alloy
0106210	Ti6Al4V Alloy
0106211	Ti6Al4V Alloy
0106401	Ti6Al4V Alloy
0106402	Ti6Al4V Alloy
0106403	Ti6Al4V Alloy
0106404	Ti6Al4V Alloy
0106405	Ti6Al4V Alloy
0106406	Ti6Al4V Alloy
0106407	Ti6Al4V Alloy
0106408	Ti6Al4V Alloy
0106409	Ti6Al4V Alloy
0106410	Ti6Al4V Alloy

0106411	Ti6Al4V Alloy
0106606	Ti6Al4V Alloy
0106607	Ti6Al4V Alloy
0106706	Ti6Al4V Alloy
0106707	Ti6Al4V Alloy
0106808	Ti6Al4V Alloy + Hydroxyapatite
0106809	Ti6Al4V Alloy + Hydroxyapatite
0106810	Ti6Al4V Alloy + Hydroxyapatite
0106811	Ti6Al4V Alloy + Hydroxyapatite
0106812	Ti6Al4V Alloy + Hydroxyapatite
0106813	Ti6Al4V Alloy + Hydroxyapatite
0106814	Ti6Al4V Alloy + Hydroxyapatite
0106815	Ti6Al4V Alloy + Hydroxyapatite
0106816	Ti6Al4V Alloy + Hydroxyapatite
0106817	Ti6Al4V Alloy + Hydroxyapatite
0106818	Ti6Al4V Alloy + Hydroxyapatite
0106908	Ti6Al4V Alloy + Hydroxyapatite
0106909	Ti6Al4V Alloy + Hydroxyapatite
0106910	Ti6Al4V Alloy + Hydroxyapatite
0106911	Ti6Al4V Alloy + Hydroxyapatite
0106912	Ti6Al4V Alloy + Hydroxyapatite
0106913	Ti6Al4V Alloy + Hydroxyapatite
0106914	Ti6Al4V Alloy + Hydroxyapatite
0106915	Ti6Al4V Alloy + Hydroxyapatite
0106916	Ti6Al4V Alloy + Hydroxyapatite
0106917	Ti6Al4V Alloy + Hydroxyapatite
0106918	Ti6Al4V Alloy + Hydroxyapatite
0107100	Ti6Al4V Alloy + Hydroxyapatite
0107101	Ti6Al4V Alloy + Hydroxyapatite
0107102	Ti6Al4V Alloy + Hydroxyapatite
0107103	Ti6Al4V Alloy + Hydroxyapatite
0107104	Ti6Al4V Alloy + Hydroxyapatite
0107105	Ti6Al4V Alloy + Hydroxyapatite
0107106	Ti6Al4V Alloy + Hydroxyapatite
0107107	Ti6Al4V Alloy + Hydroxyapatite
0107108	Ti6Al4V Alloy + Hydroxyapatite
0107109	Ti6Al4V Alloy + Hydroxyapatite
0107110	Ti6Al4V Alloy + Hydroxyapatite
0107700	Ti6Al4V Alloy + Hydroxyapatite
0107701	Ti6Al4V Alloy + Hydroxyapatite
0107702	Ti6Al4V Alloy + Hydroxyapatite
0107703	Ti6Al4V Alloy + Hydroxyapatite
0107704	Ti6Al4V Alloy + Hydroxyapatite

0107705	Ti6Al4V Alloy + Hydroxyapatite
0107706	Ti6Al4V Alloy + Hydroxyapatite
0107707	Ti6Al4V Alloy + Hydroxyapatite
0107708	Ti6Al4V Alloy + Hydroxyapatite
0107709	Ti6Al4V Alloy + Hydroxyapatite
0107710	Ti6Al4V Alloy + Hydroxyapatite
0107756	Ti6Al4V Alloy + Hydroxyapatite
0107757	Ti6Al4V Alloy + Hydroxyapatite
0107800	Ti6Al4V Alloy + Hydroxyapatite
0107801	Ti6Al4V Alloy + Hydroxyapatite
0107802	Ti6Al4V Alloy + Hydroxyapatite
0107803	Ti6Al4V Alloy + Hydroxyapatite
0107804	Ti6Al4V Alloy + Hydroxyapatite
0107805	Ti6Al4V Alloy + Hydroxyapatite
0107806	Ti6Al4V Alloy + Hydroxyapatite
0107807	Ti6Al4V Alloy + Hydroxyapatite
0107808	Ti6Al4V Alloy + Hydroxyapatite
0107809	Ti6Al4V Alloy + Hydroxyapatite
0107810	Ti6Al4V Alloy + Hydroxyapatite
0107856	Ti6Al4V Alloy + Hydroxyapatite
0107857	Ti6Al4V Alloy + Hydroxyapatite
0107900	Ti6Al4V Alloy + Hydroxyapatite
0107901	Ti6Al4V Alloy + Hydroxyapatite
0107902	Ti6Al4V Alloy + Hydroxyapatite
0107903	Ti6Al4V Alloy + Hydroxyapatite
0107904	Ti6Al4V Alloy + Hydroxyapatite
0107905	Ti6Al4V Alloy + Hydroxyapatite
0107906	Ti6Al4V Alloy + Hydroxyapatite
0107907	Ti6Al4V Alloy + Hydroxyapatite
0107908	Ti6Al4V Alloy + Hydroxyapatite
0107909	Ti6Al4V Alloy + Hydroxyapatite
0107910	Ti6Al4V Alloy + Hydroxyapatite
0108000	Ti6Al4V Alloy + Hydroxyapatite
0108001	Ti6Al4V Alloy + Hydroxyapatite
0108002	Ti6Al4V Alloy + Hydroxyapatite
0108003	Ti6Al4V Alloy + Hydroxyapatite
0108004	Ti6Al4V Alloy + Hydroxyapatite
0108005	Ti6Al4V Alloy + Hydroxyapatite
0108006	Ti6Al4V Alloy + Hydroxyapatite
0108007	Ti6Al4V Alloy + Hydroxyapatite
0108008	Ti6Al4V Alloy + Hydroxyapatite
0108009	Ti6Al4V Alloy + Hydroxyapatite
0108010	Ti6Al4V Alloy + Hydroxyapatite

0108209	Ti6Al4V Alloy
0108210	Ti6Al4V Alloy
0108211	Ti6Al4V Alloy
0108212	Ti6Al4V Alloy
0108213	Ti6Al4V Alloy
0108214	Ti6Al4V Alloy
0108215	Ti6Al4V Alloy
0108216	Ti6Al4V Alloy
0108217	Ti6Al4V Alloy
0108218	Ti6Al4V Alloy
0108219	Ti6Al4V Alloy
0108220	Ti6Al4V Alloy
0108308	Ti6Al4V Alloy
0108309	Ti6Al4V Alloy
0108310	Ti6Al4V Alloy
0108311	Ti6Al4V Alloy
0108312	Ti6Al4V Alloy
0108313	Ti6Al4V Alloy
0108314	Ti6Al4V Alloy
0108315	Ti6Al4V Alloy
0108316	Ti6Al4V Alloy
0108317	Ti6Al4V Alloy
0108318	Ti6Al4V Alloy
0108319	Ti6Al4V Alloy
0108320	Ti6Al4V Alloy
0108408	Ti6Al4V Alloy
0108409	Ti6Al4V Alloy
0108410	Ti6Al4V Alloy
0108411	Ti6Al4V Alloy
0108412	Ti6Al4V Alloy
0108413	Ti6Al4V Alloy
0108414	Ti6Al4V Alloy
0108415	Ti6Al4V Alloy
0108416	Ti6Al4V Alloy
0108417	Ti6Al4V Alloy
0108418	Ti6Al4V Alloy
0108419	Ti6Al4V Alloy
0108420	Ti6Al4V Alloy
0108500	Ti6Al4V Alloy + Hydroxyapatite
0108501	Ti6Al4V Alloy + Hydroxyapatite
0108502	Ti6Al4V Alloy + Hydroxyapatite
0108503	Ti6Al4V Alloy + Hydroxyapatite
0108504	Ti6Al4V Alloy + Hydroxyapatite

0108505	Ti6Al4V Alloy + Hydroxyapatite
0108506	Ti6Al4V Alloy + Hydroxyapatite
0108507	Ti6Al4V Alloy + Hydroxyapatite
0108508	Ti6Al4V Alloy + Hydroxyapatite
0108509	Ti6Al4V Alloy + Hydroxyapatite
0108510	Ti6Al4V Alloy + Hydroxyapatite
0108511	Ti6Al4V Alloy + Hydroxyapatite
0108512	Ti6Al4V Alloy + Hydroxyapatite
0108602	Ti6Al4V Alloy + Hydroxyapatite
0108603	Ti6Al4V Alloy + Hydroxyapatite
0108604	Ti6Al4V Alloy + Hydroxyapatite
0108605	Ti6Al4V Alloy + Hydroxyapatite
0108606	Ti6Al4V Alloy + Hydroxyapatite
0108607	Ti6Al4V Alloy + Hydroxyapatite
0108608	Ti6Al4V Alloy + Hydroxyapatite
0108609	Ti6Al4V Alloy + Hydroxyapatite
0108610	Ti6Al4V Alloy + Hydroxyapatite
0108611	Ti6Al4V Alloy + Hydroxyapatite
0108612	Ti6Al4V Alloy + Hydroxyapatite
0108700	Ti6Al4V Alloy + Hydroxyapatite
0108701	Ti6Al4V Alloy + Hydroxyapatite
0108702	Ti6Al4V Alloy + Hydroxyapatite
0108703	Ti6Al4V Alloy + Hydroxyapatite
0108704	Ti6Al4V Alloy + Hydroxyapatite
0108705	Ti6Al4V Alloy + Hydroxyapatite
0108706	Ti6Al4V Alloy + Hydroxyapatite
0108707	Ti6Al4V Alloy + Hydroxyapatite
0108708	Ti6Al4V Alloy + Hydroxyapatite
0108709	Ti6Al4V Alloy + Hydroxyapatite
0108710	Ti6Al4V Alloy + Hydroxyapatite
0108711	Ti6Al4V Alloy + Hydroxyapatite
0108712	Ti6Al4V Alloy + Hydroxyapatite
0108802	Ti6Al4V Alloy + Hydroxyapatite
0108803	Ti6Al4V Alloy + Hydroxyapatite
0108804	Ti6Al4V Alloy + Hydroxyapatite
0108805	Ti6Al4V Alloy + Hydroxyapatite
0108806	Ti6Al4V Alloy + Hydroxyapatite
0108807	Ti6Al4V Alloy + Hydroxyapatite
0108808	Ti6Al4V Alloy + Hydroxyapatite
0108809	Ti6Al4V Alloy + Hydroxyapatite
0108810	Ti6Al4V Alloy + Hydroxyapatite
0108811	Ti6Al4V Alloy + Hydroxyapatite
0108812	Ti6Al4V Alloy + Hydroxyapatite

0108901	Ti6Al4V Alloy
0108902	Ti6Al4V Alloy
0108903	Ti6Al4V Alloy
0108904	Ti6Al4V Alloy
0108905	Ti6Al4V Alloy
0108906	Ti6Al4V Alloy
0108907	Ti6Al4V Alloy
0108908	Ti6Al4V Alloy
0108909	Ti6Al4V Alloy
0108910	Ti6Al4V Alloy
0108911	Ti6Al4V Alloy
0109001	Ti6Al4V Alloy
0109002	Ti6Al4V Alloy
0109003	Ti6Al4V Alloy
0109004	Ti6Al4V Alloy
0109005	Ti6Al4V Alloy
0109006	Ti6Al4V Alloy
0109007	Ti6Al4V Alloy
0109008	Ti6Al4V Alloy
0109009	Ti6Al4V Alloy
0109010	Ti6Al4V Alloy
0109011	Ti6Al4V Alloy
0109101	Ti6Al4V Alloy
0109102	Ti6Al4V Alloy
0109103	Ti6Al4V Alloy
0109104	Ti6Al4V Alloy
0109105	Ti6Al4V Alloy
0109106	Ti6Al4V Alloy
0109107	Ti6Al4V Alloy
0109108	Ti6Al4V Alloy
0109109	Ti6Al4V Alloy
0109110	Ti6Al4V Alloy
0109400	Ti6Al4V Alloy
0109401	Ti6Al4V Alloy
0109402	Ti6Al4V Alloy
0109403	Ti6Al4V Alloy
0109404	Ti6Al4V Alloy
0109405	Ti6Al4V Alloy
0109406	Ti6Al4V Alloy
0109407	Ti6Al4V Alloy
0109408	Ti6Al4V Alloy
0109409	Ti6Al4V Alloy
0109410	Ti6Al4V Alloy

0109490	Ti6Al4V Alloy
0109500	Ti6Al4V Alloy
0109501	Ti6Al4V Alloy
0109502	Ti6Al4V Alloy
0109503	Ti6Al4V Alloy
0109504	Ti6Al4V Alloy
0109505	Ti6Al4V Alloy
0109506	Ti6Al4V Alloy
0109507	Ti6Al4V Alloy
0109508	Ti6Al4V Alloy
0109509	Ti6Al4V Alloy
0109510	Ti6Al4V Alloy
0109590	Ti6Al4V Alloy
0109601	Ti6Al4V Alloy
0109602	Ti6Al4V Alloy
0109603	Ti6Al4V Alloy
0109604	Ti6Al4V Alloy
0109605	Ti6Al4V Alloy
0109606	Ti6Al4V Alloy
0109607	Ti6Al4V Alloy
0109608	Ti6Al4V Alloy
0109609	Ti6Al4V Alloy
0109610	Ti6Al4V Alloy
0109611	Ti6Al4V Alloy
0109701	Ti6Al4V Alloy
0109702	Ti6Al4V Alloy
0109703	Ti6Al4V Alloy
0109704	Ti6Al4V Alloy
0109705	Ti6Al4V Alloy
0109706	Ti6Al4V Alloy
0109707	Ti6Al4V Alloy
0109708	Ti6Al4V Alloy
0109709	Ti6Al4V Alloy
0109710	Ti6Al4V Alloy
0109711	Ti6Al4V Alloy
0110001	Ti6Al4V Alloy
0110002	Ti6Al4V Alloy
0110003	Ti6Al4V Alloy
0110004	Ti6Al4V Alloy
0110005	Ti6Al4V Alloy
0110006	Ti6Al4V Alloy
0110007	Ti6Al4V Alloy
0110008	Ti6Al4V Alloy

0110009	Ti6Al4V Alloy
0110010	Ti6Al4V Alloy
0110011	Ti6Al4V Alloy
0110101	Ti6Al4V Alloy
0110102	Ti6Al4V Alloy
0110103	Ti6Al4V Alloy
0110104	Ti6Al4V Alloy
0110105	Ti6Al4V Alloy
0110106	Ti6Al4V Alloy
0110107	Ti6Al4V Alloy
0110108	Ti6Al4V Alloy
0110109	Ti6Al4V Alloy
0110110	Ti6Al4V Alloy
0110111	Ti6Al4V Alloy
0110201	Ti6Al4V Alloy
0110202	Ti6Al4V Alloy
0110203	Ti6Al4V Alloy
0110204	Ti6Al4V Alloy
0110205	Ti6Al4V Alloy
0110206	Ti6Al4V Alloy
0110207	Ti6Al4V Alloy
0110208	Ti6Al4V Alloy
0110209	Ti6Al4V Alloy
0110210	Ti6Al4V Alloy
0110211	Ti6Al4V Alloy
0110301	Ti6Al4V Alloy
0110302	Ti6Al4V Alloy
0110303	Ti6Al4V Alloy
0110304	Ti6Al4V Alloy
0110305	Ti6Al4V Alloy
0110306	Ti6Al4V Alloy
0110307	Ti6Al4V Alloy
0110308	Ti6Al4V Alloy
0110309	Ti6Al4V Alloy
0110310	Ti6Al4V Alloy
0110311	Ti6Al4V Alloy
0110401	Ti6Al4V Alloy
0110402	Ti6Al4V Alloy
0110403	Ti6Al4V Alloy
0110404	Ti6Al4V Alloy
0110405	Ti6Al4V Alloy
0110406	Ti6Al4V Alloy
0110407	Ti6Al4V Alloy

0110408	Ti6Al4V Alloy
0110409	Ti6Al4V Alloy
0110410	Ti6Al4V Alloy
0110411	Ti6Al4V Alloy
0110501	Ti6Al4V Alloy
0110502	Ti6Al4V Alloy
0110503	Ti6Al4V Alloy
0110504	Ti6Al4V Alloy
0110505	Ti6Al4V Alloy
0110506	Ti6Al4V Alloy
0110507	Ti6Al4V Alloy
0110508	Ti6Al4V Alloy
0110509	Ti6Al4V Alloy
0110510	Ti6Al4V Alloy
0110511	Ti6Al4V Alloy
0181201	Ti6Al4V Alloy
0181202	Ti6Al4V Alloy
0181203	Ti6Al4V Alloy
0181204	Ti6Al4V Alloy
0181205	Ti6Al4V Alloy
0181206	Ti6Al4V Alloy
0181207	Ti6Al4V Alloy
0181208	Ti6Al4V Alloy
0181211	Ti6Al4V Alloy
0181212	Ti6Al4V Alloy
0181213	Ti6Al4V Alloy
0181214	Ti6Al4V Alloy
0181215	Ti6Al4V Alloy
0181216	Ti6Al4V Alloy
0181217	Ti6Al4V Alloy
0181218	Ti6Al4V Alloy
0181401	Ti6Al4V Alloy
0181402	Ti6Al4V Alloy
0181403	Ti6Al4V Alloy
0181404	Ti6Al4V Alloy
0181405	Ti6Al4V Alloy
0181406	Ti6Al4V Alloy
0181407	Ti6Al4V Alloy
0181408	Ti6Al4V Alloy
0181411	Ti6Al4V Alloy
0181412	Ti6Al4V Alloy
0181413	Ti6Al4V Alloy
0181414	Ti6Al4V Alloy

0181415	Ti6Al4V Alloy
0181416	Ti6Al4V Alloy
0181417	Ti6Al4V Alloy
0181418	Ti6Al4V Alloy
0181501	Ti6Al4V Alloy
0181502	Ti6Al4V Alloy
0181503	Ti6Al4V Alloy
0181504	Ti6Al4V Alloy
0181505	Ti6Al4V Alloy
0181506	Ti6Al4V Alloy
0181507	Ti6Al4V Alloy
0181508	Ti6Al4V Alloy
0181511	Ti6Al4V Alloy
0181512	Ti6Al4V Alloy
0181513	Ti6Al4V Alloy
0181514	Ti6Al4V Alloy
0181515	Ti6Al4V Alloy
0181516	Ti6Al4V Alloy
0181517	Ti6Al4V Alloy
0181518	Ti6Al4V Alloy
0181601	STAINLESS STEEL
0181602	STAINLESS STEEL
0181603	STAINLESS STEEL
0181604	STAINLESS STEEL
0181605	STAINLESS STEEL
0181606	STAINLESS STEEL
0181607	STAINLESS STEEL
0181608	STAINLESS STEEL
0181611	STAINLESS STEEL
0181612	STAINLESS STEEL
0181613	STAINLESS STEEL
0181614	STAINLESS STEEL
0181615	STAINLESS STEEL
0181616	STAINLESS STEEL
0181617	STAINLESS STEEL
0181618	STAINLESS STEEL
0181701	STAINLESS STEEL
0181702	STAINLESS STEEL
0181703	STAINLESS STEEL
0181704	STAINLESS STEEL
0181705	STAINLESS STEEL
0181706	STAINLESS STEEL
0181707	STAINLESS STEEL

0181708	STAINLESS STEEL
0181711	STAINLESS STEEL
0181712	STAINLESS STEEL
0181713	STAINLESS STEEL
0181714	STAINLESS STEEL
0181715	STAINLESS STEEL
0181716	STAINLESS STEEL
0181717	STAINLESS STEEL
0181718	STAINLESS STEEL
0184501	STAINLESS STEEL
0184502	STAINLESS STEEL
0184503	STAINLESS STEEL
0184504	STAINLESS STEEL
0185001	STAINLESS STEEL
0185002	STAINLESS STEEL
0185003	STAINLESS STEEL
0185004	STAINLESS STEEL
0185005	STAINLESS STEEL
0185011	STAINLESS STEEL
0185012	STAINLESS STEEL
0185013	STAINLESS STEEL
0185014	STAINLESS STEEL
0185015	STAINLESS STEEL
0185501	STAINLESS STEEL
0185502	STAINLESS STEEL
0185503	STAINLESS STEEL
0185504	STAINLESS STEEL
0185505	STAINLESS STEEL
0185901	PMMA
0185902	PMMA
0185903	PMMA
0185904	PMMA
0185905	PMMA
0187300	Ti6Al4V Alloy
0187301	Ti6Al4V Alloy
0187302	Ti6Al4V Alloy
0187303	Ti6Al4V Alloy
0187304	Ti6Al4V Alloy
0187305	Ti6Al4V Alloy
0187306	Ti6Al4V Alloy
0187307	Ti6Al4V Alloy
0187308	Ti6Al4V Alloy
0187309	Ti6Al4V Alloy

0187310	Ti6Al4V Alloy
0187700	STAINLESS STEEL
0187701	STAINLESS STEEL
0187702	STAINLESS STEEL
0187703	STAINLESS STEEL
0187704	STAINLESS STEEL
0187705	STAINLESS STEEL
0187706	STAINLESS STEEL
0187707	STAINLESS STEEL
0187708	STAINLESS STEEL
0187709	STAINLESS STEEL
0187710	STAINLESS STEEL
0187800	STAINLESS STEEL
0187801	STAINLESS STEEL
0187802	STAINLESS STEEL
0187803	STAINLESS STEEL
0187804	STAINLESS STEEL
0187805	STAINLESS STEEL
0187806	STAINLESS STEEL
0187807	STAINLESS STEEL
0187808	STAINLESS STEEL
0187809	STAINLESS STEEL
0187810	STAINLESS STEEL
0188501	STAINLESS STEEL
0188502	STAINLESS STEEL
0188503	STAINLESS STEEL
0188504	STAINLESS STEEL
0188505	STAINLESS STEEL
0188506	STAINLESS STEEL
0188507	STAINLESS STEEL
0188508	STAINLESS STEEL
0188509	STAINLESS STEEL
0188510	STAINLESS STEEL
0188511	STAINLESS STEEL
0188512	STAINLESS STEEL
0188602	STAINLESS STEEL
0188603	STAINLESS STEEL
0188604	STAINLESS STEEL
0188605	STAINLESS STEEL
0188606	STAINLESS STEEL
0188607	STAINLESS STEEL
0188608	STAINLESS STEEL
0188609	STAINLESS STEEL

0188610	STAINLESS STEEL
0188611	STAINLESS STEEL
0188612	STAINLESS STEEL
0188701	STAINLESS STEEL
0188702	STAINLESS STEEL
0188703	STAINLESS STEEL
0188704	STAINLESS STEEL
0188705	STAINLESS STEEL
0188706	STAINLESS STEEL
0188707	STAINLESS STEEL
0188708	STAINLESS STEEL
0188709	STAINLESS STEEL
0188710	STAINLESS STEEL
0188711	STAINLESS STEEL
0188712	STAINLESS STEEL
0188802	STAINLESS STEEL
0188803	STAINLESS STEEL
0188804	STAINLESS STEEL
0188805	STAINLESS STEEL
0188806	STAINLESS STEEL
0188807	STAINLESS STEEL
0188808	STAINLESS STEEL
0188809	STAINLESS STEEL
0188810	STAINLESS STEEL
0188811	STAINLESS STEEL
0188812	STAINLESS STEEL
0188901	STAINLESS STEEL
0188902	STAINLESS STEEL
0188903	STAINLESS STEEL
0188904	STAINLESS STEEL
0188905	STAINLESS STEEL
0188906	STAINLESS STEEL
0188907	STAINLESS STEEL
0188908	STAINLESS STEEL
0188909	STAINLESS STEEL
0188910	STAINLESS STEEL
0188911	STAINLESS STEEL
0189503	STAINLESS STEEL
0189504	STAINLESS STEEL
0189505	STAINLESS STEEL
0189506	STAINLESS STEEL
0189507	STAINLESS STEEL
0189508	STAINLESS STEEL

0189509	STAINLESS STEEL
0189510	STAINLESS STEEL
0189511	STAINLESS STEEL
0189803	STAINLESS STEEL
0189804	STAINLESS STEEL
0189805	STAINLESS STEEL
0189806	STAINLESS STEEL
0189807	STAINLESS STEEL
0189808	STAINLESS STEEL
0189809	STAINLESS STEEL
0189810	STAINLESS STEEL
0189811	STAINLESS STEEL
0189903	STAINLESS STEEL
0189904	STAINLESS STEEL
0189905	STAINLESS STEEL
0189906	STAINLESS STEEL
0189907	STAINLESS STEEL
0189908	STAINLESS STEEL
0189909	STAINLESS STEEL
0189910	STAINLESS STEEL
0189911	STAINLESS STEEL
0193113	Ti6Al4V Alloy
0193114	Ti6Al4V Alloy
0193115	Ti6Al4V Alloy
0193116	Ti6Al4V Alloy
0193117	Ti6Al4V Alloy
0193118	Ti6Al4V Alloy
0193119	Ti6Al4V Alloy
0193120	Ti6Al4V Alloy
0193121	Ti6Al4V Alloy
0193122	Ti6Al4V Alloy
0193213	Ti6Al4V Alloy
0193214	Ti6Al4V Alloy
0193215	Ti6Al4V Alloy
0193216	Ti6Al4V Alloy
0193217	Ti6Al4V Alloy
0193218	Ti6Al4V Alloy
0193219	Ti6Al4V Alloy
0193220	Ti6Al4V Alloy
0193221	Ti6Al4V Alloy
0193222	Ti6Al4V Alloy
0200010	Ti6Al4V Alloy
0200015	Ti6Al4V Alloy

0200020	Ti6Al4V Alloy
0200025	Ti6Al4V Alloy
0200030	Ti6Al4V Alloy
0200035	Ti6Al4V Alloy
0200090	Ti6Al4V Alloy
0200095	Ti6Al4V Alloy
0200110	Ti6Al4V Alloy
0200115	Ti6Al4V Alloy
0200120	Ti6Al4V Alloy
0200125	Ti6Al4V Alloy
0200130	Ti6Al4V Alloy
0200135	Ti6Al4V Alloy
0200210	Ti6Al4V Alloy
0200215	Ti6Al4V Alloy
0200220	Ti6Al4V Alloy
0200225	Ti6Al4V Alloy
0200230	Ti6Al4V Alloy
0200235	Ti6Al4V Alloy
0200310	Ti6Al4V Alloy
0200315	Ti6Al4V Alloy
0200320	Ti6Al4V Alloy
0200325	Ti6Al4V Alloy
0200330	Ti6Al4V Alloy
0200335	Ti6Al4V Alloy
0200410	Ti6Al4V Alloy
0200415	Ti6Al4V Alloy
0200420	Ti6Al4V Alloy
0200425	Ti6Al4V Alloy
0200430	Ti6Al4V Alloy
0200435	Ti6Al4V Alloy
0200511	Ti6Al4V Alloy
0200512	Ti6Al4V Alloy
0200513	Ti6Al4V Alloy
0200514	Ti6Al4V Alloy
0200515	Ti6Al4V Alloy
0200516	Ti6Al4V Alloy
0200517	Ti6Al4V Alloy
0200518	Ti6Al4V Alloy
0200519	Ti6Al4V Alloy
0200520	Ti6Al4V Alloy
0200522	Ti6Al4V Alloy
0200524	Ti6Al4V Alloy
0200605	Ti6Al4V Alloy

0200630	Titanium
0200631	Titanium + Hydroxyapatite
0200632	Ti6Al4V Alloy
0200650	Titanium
0200651	Titanium + Hydroxyapatite
0200652	Ti6Al4V Alloy
0200670	Titanium
0200671	Titanium + Hydroxyapatite
0200672	Ti6Al4V Alloy
0200700	Ti6Al4V Alloy
0200705	Ti6Al4V Alloy
0200710	Ti6Al4V Alloy
0200715	Ti6Al4V Alloy
0200720	Ti6Al4V Alloy
0200740	Ti6Al4V Alloy
0200745	Ti6Al4V Alloy + TiNbN
0200750	Ti6Al4V Alloy
0200755	Ti6Al4V Alloy + TiNbN
0200760	Ti6Al4V Alloy
0200765	Ti6Al4V Alloy + TiNbN
0200770	Ti6Al4V Alloy
0200775	Ti6Al4V Alloy + TiNbN
0200811	Ti6Al4V Alloy + Hydroxyapatite
0200812	Ti6Al4V Alloy + Hydroxyapatite
0200813	Ti6Al4V Alloy + Hydroxyapatite
0200814	Ti6Al4V Alloy + Hydroxyapatite
0200815	Ti6Al4V Alloy + Hydroxyapatite
0200816	Ti6Al4V Alloy + Hydroxyapatite
0200817	Ti6Al4V Alloy + Hydroxyapatite
0200818	Ti6Al4V Alloy + Hydroxyapatite
0200819	Ti6Al4V Alloy + Hydroxyapatite
0200820	Ti6Al4V Alloy + Hydroxyapatite
0200822	Ti6Al4V Alloy + Hydroxyapatite
0200824	Ti6Al4V Alloy + Hydroxyapatite
0200910	Ti6Al4V Alloy
0200920	Ti6Al4V Alloy
0200930	Ti6Al4V Alloy
0200940	Ti6Al4V Alloy
0200950	Ti6Al4V Alloy
0200960	Ti6Al4V Alloy
0200965	Ti6Al4V Alloy
0200970	Ti6Al4V Alloy
0200980	Ti6Al4V Alloy

0201000	Ti6Al4V Alloy
0201010	Ti6Al4V Alloy
0201020	Ti6Al4V Alloy
0201030	Ti6Al4V Alloy
0201040	Ti6Al4V Alloy
0201050	Ti6Al4V Alloy
0201060	Ti6Al4V Alloy
0201070	Ti6Al4V Alloy
0201080	Ti6Al4V Alloy
0201085	Ti6Al4V Alloy
0201090	Ti6Al4V Alloy
0201110	Ti6Al4V Alloy
0201120	Ti6Al4V Alloy
0201130	Ti6Al4V Alloy
0201140	Ti6Al4V Alloy
0201150	Ti6Al4V Alloy
0201160	Ti6Al4V Alloy
0201170	Ti6Al4V Alloy
0201180	Ti6Al4V Alloy
0201185	Ti6Al4V Alloy
0201190	Ti6Al4V Alloy
0201193	Ti6Al4V Alloy
0201196	Ti6Al4V Alloy
0201210	Ti6Al4V Alloy
0201220	Ti6Al4V Alloy
0201230	Ti6Al4V Alloy
0201240	Ti6Al4V Alloy
0201250	Ti6Al4V Alloy
0201260	Ti6Al4V Alloy
0201270	Ti6Al4V Alloy
0201280	Ti6Al4V Alloy
0201285	Ti6Al4V Alloy
0201290	Ti6Al4V Alloy
0201293	Ti6Al4V Alloy
0201296	Ti6Al4V Alloy
0201511	Ti6Al4V Alloy
0201512	Ti6Al4V Alloy
0201513	Ti6Al4V Alloy
0201514	Ti6Al4V Alloy
0201515	Ti6Al4V Alloy
0201516	Ti6Al4V Alloy
0201517	Ti6Al4V Alloy
0201518	Ti6Al4V Alloy

0201519	Ti6Al4V Alloy
0201520	Ti6Al4V Alloy
0201522	Ti6Al4V Alloy
0201524	Ti6Al4V Alloy
0201811	Ti6Al4V Alloy + Hydroxyapatite
0201812	Ti6Al4V Alloy + Hydroxyapatite
0201813	Ti6Al4V Alloy + Hydroxyapatite
0201814	Ti6Al4V Alloy + Hydroxyapatite
0201815	Ti6Al4V Alloy + Hydroxyapatite
0201816	Ti6Al4V Alloy + Hydroxyapatite
0201817	Ti6Al4V Alloy + Hydroxyapatite
0201818	Ti6Al4V Alloy + Hydroxyapatite
0201819	Ti6Al4V Alloy + Hydroxyapatite
0201820	Ti6Al4V Alloy + Hydroxyapatite
0201822	Ti6Al4V Alloy + Hydroxyapatite
0201824	Ti6Al4V Alloy + Hydroxyapatite
0202000	Ti6Al4V Alloy + Hydroxyapatite
0202010	Ti6Al4V Alloy + Hydroxyapatite
0202020	Ti6Al4V Alloy + Hydroxyapatite
0202030	Ti6Al4V Alloy + Hydroxyapatite
0202040	Ti6Al4V Alloy + Hydroxyapatite
0202050	Ti6Al4V Alloy + Hydroxyapatite
0202060	Ti6Al4V Alloy + Hydroxyapatite
0202070	Ti6Al4V Alloy + Hydroxyapatite
0202080	Ti6Al4V Alloy + Hydroxyapatite
0202085	Ti6Al4V Alloy + Hydroxyapatite
0202090	Ti6Al4V Alloy + Hydroxyapatite
0202110	Ti6Al4V Alloy + Hydroxyapatite
0202120	Ti6Al4V Alloy + Hydroxyapatite
0202130	Ti6Al4V Alloy + Hydroxyapatite
0202140	Ti6Al4V Alloy + Hydroxyapatite
0202150	Ti6Al4V Alloy + Hydroxyapatite
0202160	Ti6Al4V Alloy + Hydroxyapatite
0202170	Ti6Al4V Alloy + Hydroxyapatite
0202180	Ti6Al4V Alloy + Hydroxyapatite
0202185	Ti6Al4V Alloy + Hydroxyapatite
0202190	Ti6Al4V Alloy + Hydroxyapatite
0202193	Ti6Al4V Alloy + Hydroxyapatite
0202196	Ti6Al4V Alloy + Hydroxyapatite
0202210	Ti6Al4V Alloy + Hydroxyapatite
0202220	Ti6Al4V Alloy + Hydroxyapatite
0202230	Ti6Al4V Alloy + Hydroxyapatite
0202240	Ti6Al4V Alloy + Hydroxyapatite

0202250	Ti6Al4V Alloy + Hydroxyapatite
0202260	Ti6Al4V Alloy + Hydroxyapatite
0202270	Ti6Al4V Alloy + Hydroxyapatite
0202280	Ti6Al4V Alloy + Hydroxyapatite
0202285	Ti6Al4V Alloy + Hydroxyapatite
0202290	Ti6Al4V Alloy + Hydroxyapatite
0202293	Ti6Al4V Alloy + Hydroxyapatite
0202296	Ti6Al4V Alloy + Hydroxyapatite
0203000	Ti6Al4V Alloy
0203010	Ti6Al4V Alloy
0203020	Ti6Al4V Alloy
0203030	Ti6Al4V Alloy
0203040	Ti6Al4V Alloy
0203050	Ti6Al4V Alloy
0203060	Ti6Al4V Alloy
0203070	Ti6Al4V Alloy
0203080	Ti6Al4V Alloy
0203085	Ti6Al4V Alloy
0203090	Ti6Al4V Alloy
0203110	Ti6Al4V Alloy
0203120	Ti6Al4V Alloy
0203130	Ti6Al4V Alloy
0203140	Ti6Al4V Alloy
0203150	Ti6Al4V Alloy
0203160	Ti6Al4V Alloy
0203170	Ti6Al4V Alloy
0203180	Ti6Al4V Alloy
0203185	Ti6Al4V Alloy
0203190	Ti6Al4V Alloy
0203193	Ti6Al4V Alloy
0203196	Ti6Al4V Alloy
0203210	Ti6Al4V Alloy
0203220	Ti6Al4V Alloy
0203230	Ti6Al4V Alloy
0203240	Ti6Al4V Alloy
0203250	Ti6Al4V Alloy
0203260	Ti6Al4V Alloy
0203270	Ti6Al4V Alloy
0203280	Ti6Al4V Alloy
0203285	Ti6Al4V Alloy
0203290	Ti6Al4V Alloy
0203293	Ti6Al4V Alloy
0203296	Ti6Al4V Alloy

0204000	Ti6Al4V Alloy + Hydroxyapatite
0204010	Ti6Al4V Alloy + Hydroxyapatite
0204020	Ti6Al4V Alloy + Hydroxyapatite
0204030	Ti6Al4V Alloy + Hydroxyapatite
0204040	Ti6Al4V Alloy + Hydroxyapatite
0204050	Ti6Al4V Alloy + Hydroxyapatite
0204060	Ti6Al4V Alloy + Hydroxyapatite
0204070	Ti6Al4V Alloy + Hydroxyapatite
0204080	Ti6Al4V Alloy + Hydroxyapatite
0204085	Ti6Al4V Alloy + Hydroxyapatite
0204090	Ti6Al4V Alloy + Hydroxyapatite
0204110	Ti6Al4V Alloy + Hydroxyapatite
0204120	Ti6Al4V Alloy + Hydroxyapatite
0204130	Ti6Al4V Alloy + Hydroxyapatite
0204140	Ti6Al4V Alloy + Hydroxyapatite
0204150	Ti6Al4V Alloy + Hydroxyapatite
0204160	Ti6Al4V Alloy + Hydroxyapatite
0204170	Ti6Al4V Alloy + Hydroxyapatite
0204180	Ti6Al4V Alloy + Hydroxyapatite
0204185	Ti6Al4V Alloy + Hydroxyapatite
0204190	Ti6Al4V Alloy + Hydroxyapatite
0204193	Ti6Al4V Alloy + Hydroxyapatite
0204196	Ti6Al4V Alloy + Hydroxyapatite
0204210	Ti6Al4V Alloy + Hydroxyapatite
0204220	Ti6Al4V Alloy + Hydroxyapatite
0204230	Ti6Al4V Alloy + Hydroxyapatite
0204240	Ti6Al4V Alloy + Hydroxyapatite
0204250	Ti6Al4V Alloy + Hydroxyapatite
0204260	Ti6Al4V Alloy + Hydroxyapatite
0204270	Ti6Al4V Alloy + Hydroxyapatite
0204280	Ti6Al4V Alloy + Hydroxyapatite
0204285	Ti6Al4V Alloy + Hydroxyapatite
0204290	Ti6Al4V Alloy + Hydroxyapatite
0204293	Ti6Al4V Alloy + Hydroxyapatite
0204296	Ti6Al4V Alloy + Hydroxyapatite
0205000	Ti6Al4V Alloy
0205010	Ti6Al4V Alloy
0205020	Ti6Al4V Alloy
0205100	Ti6Al4V Alloy
0205110	Ti6Al4V Alloy
0205120	Ti6Al4V Alloy
0205130	Ti6Al4V Alloy
0205200	Ti6Al4V Alloy

0205210	Ti6Al4V Alloy
0205220	Ti6Al4V Alloy
0205230	Ti6Al4V Alloy
0206000	Ti6Al4V Alloy
0206010	Ti6Al4V Alloy + Hydroxyapatite
0206020	Ti6Al4V Alloy + Hydroxyapatite
0206100	Ti6Al4V Alloy + Hydroxyapatite
0206110	Ti6Al4V Alloy + Hydroxyapatite
0206120	Ti6Al4V Alloy + Hydroxyapatite
0206130	Ti6Al4V Alloy + Hydroxyapatite
0206200	Ti6Al4V Alloy + Hydroxyapatite
0206210	Ti6Al4V Alloy + Hydroxyapatite
0206220	Ti6Al4V Alloy + Hydroxyapatite
0206230	Ti6Al4V Alloy + Hydroxyapatite
0207000	Ti6Al4V Alloy
0207010	Ti6Al4V Alloy
0207020	Ti6Al4V Alloy
0207100	Ti6Al4V Alloy
0207110	Ti6Al4V Alloy
0207120	Ti6Al4V Alloy
0207130	Ti6Al4V Alloy
0207200	Ti6Al4V Alloy
0207210	Ti6Al4V Alloy
0207220	Ti6Al4V Alloy
0207230	Ti6Al4V Alloy
0208000	Ti6Al4V Alloy
0208010	Ti6Al4V Alloy + Hydroxyapatite
0208020	Ti6Al4V Alloy + Hydroxyapatite
0208100	Ti6Al4V Alloy + Hydroxyapatite
0208110	Ti6Al4V Alloy + Hydroxyapatite
0208120	Ti6Al4V Alloy + Hydroxyapatite
0208130	Ti6Al4V Alloy + Hydroxyapatite
0208200	Ti6Al4V Alloy + Hydroxyapatite
0208210	Ti6Al4V Alloy + Hydroxyapatite
0208220	Ti6Al4V Alloy + Hydroxyapatite
0208230	Ti6Al4V Alloy + Hydroxyapatite
0280510	Ti6Al4V Alloy
0280511	Ti6Al4V Alloy
0280512	Ti6Al4V Alloy
0280513	Ti6Al4V Alloy
0280514	Ti6Al4V Alloy
0280515	Ti6Al4V Alloy
0280516	Ti6Al4V Alloy

0280517	Ti6Al4V Alloy
0280518	Ti6Al4V Alloy
0280519	Ti6Al4V Alloy
0280520	Ti6Al4V Alloy
0280522	Ti6Al4V Alloy
0460110	Ti6Al4V Alloy
0460210	Ti6Al4V Alloy
0460220	Ti6Al4V Alloy
0460310	Ti6Al4V Alloy
0460320	Ti6Al4V Alloy
0460330	Ti6Al4V Alloy
0469110	Ti6Al4V Alloy
0469120	Ti6Al4V Alloy
0469130	Ti6Al4V Alloy
0469210	Ti6Al4V Alloy
0469220	Ti6Al4V Alloy
0469230	Ti6Al4V Alloy
0469310	Ti6Al4V Alloy
0469320	Ti6Al4V Alloy
0469330	Ti6Al4V Alloy
0470110	Ti6Al4V Alloy
0470210	Ti6Al4V Alloy
0470220	Ti6Al4V Alloy
0470310	Ti6Al4V Alloy
0470320	Ti6Al4V Alloy
0470330	Ti6Al4V Alloy
0479110	Ti6Al4V Alloy
0479120	Ti6Al4V Alloy
0479130	Ti6Al4V Alloy
0479210	Ti6Al4V Alloy
0479220	Ti6Al4V Alloy
0479230	Ti6Al4V Alloy
0479310	Ti6Al4V Alloy
0479320	Ti6Al4V Alloy
0479330	Ti6Al4V Alloy
0510281	Alumina Ceramics
0510282	Alumina Ceramics
0510283	Alumina Ceramics
0510321	Alumina Ceramics
0510322	Alumina Ceramics
0510323	Alumina Ceramics
0510361	Alumina Ceramics
0510362	Alumina Ceramics

0510363	Alumina Ceramics
0511321	Alumina Ceramics
0511322	Alumina Ceramics
0511323	Alumina Ceramics
0511361	Alumina Ceramics
0511362	Alumina Ceramics
0511363	Alumina Ceramics
0512221	Alumina Ceramics
0512222	Alumina Ceramics
0512223	Alumina Ceramics
0512281	Alumina Ceramics
0512282	Alumina Ceramics
0512283	Alumina Ceramics
0512321	Alumina Ceramics
0512322	Alumina Ceramics
0512323	Alumina Ceramics
0512361	Alumina Ceramics
0512362	Alumina Ceramics
0512363	Alumina Ceramics
0513001	Ti6Al4V Alloy
0513002	Ti6Al4V Alloy
0513003	Ti6Al4V Alloy
0513004	Ti6Al4V Alloy
0513280	Alumina Ceramics
0513281	Alumina Ceramics
0513282	Alumina Ceramics
0513283	Alumina Ceramics
0513284	Alumina Ceramics
0513320	Alumina Ceramics
0513321	Alumina Ceramics
0513322	Alumina Ceramics
0513323	Alumina Ceramics
0513324	Alumina Ceramics
0513360	Alumina Ceramics
0513361	Alumina Ceramics
0513362	Alumina Ceramics
0513363	Alumina Ceramics
0513364	Alumina Ceramics
0513400	Alumina Ceramics
0513401	Alumina Ceramics
0513402	Alumina Ceramics
0513403	Alumina Ceramics
0513404	Alumina Ceramics

0513440	Alumina Ceramics
0513441	Alumina Ceramics
0513442	Alumina Ceramics
0513443	Alumina Ceramics
0513444	Alumina Ceramics
0513480	Alumina Ceramics
0513481	Alumina Ceramics
0513482	Alumina Ceramics
0513483	Alumina Ceramics
0513484	Alumina Ceramics
0514181	Alumina Ceramics
0514182	Alumina Ceramics
0514183	Alumina Ceramics
0514281	Alumina Ceramics
0514282	Alumina Ceramics
0514283	Alumina Ceramics
0514321	Alumina Ceramics
0514322	Alumina Ceramics
0514323	Alumina Ceramics
0514324	Alumina Ceramics
0514361	Alumina Ceramics
0514362	Alumina Ceramics
0514363	Alumina Ceramics
0514364	Alumina Ceramics
0514401	Alumina Ceramics
0514402	Alumina Ceramics
0514403	Alumina Ceramics
0514404	Alumina Ceramics
0514441	Alumina Ceramics
0514442	Alumina Ceramics
0514443	Alumina Ceramics
0514444	Alumina Ceramics
0520281	Co-Cr-Mo Alloy
0520282	Co-Cr-Mo Alloy
0520283	Co-Cr-Mo Alloy
0522221	Co-Cr-Mo Alloy
0522222	Co-Cr-Mo Alloy
0522223	Co-Cr-Mo Alloy
0524221	Co-Cr-Mo Alloy
0524222	Co-Cr-Mo Alloy
0524223	Co-Cr-Mo Alloy
0524224	Co-Cr-Mo Alloy
0524284	Co-Cr-Mo Alloy

0524285	Co-Cr-Mo Alloy
0524321	Co-Cr-Mo Alloy
0524322	Co-Cr-Mo Alloy
0524323	Co-Cr-Mo Alloy
0524324	Co-Cr-Mo Alloy
0524325	Co-Cr-Mo Alloy
0524361	Co-Cr-Mo Alloy
0524362	Co-Cr-Mo Alloy
0524363	Co-Cr-Mo Alloy
0524364	Co-Cr-Mo Alloy
0524365	Co-Cr-Mo Alloy
0524401	Co-Cr-Mo Alloy
0524402	Co-Cr-Mo Alloy
0524403	Co-Cr-Mo Alloy
0524404	Co-Cr-Mo Alloy
0524405	Co-Cr-Mo Alloy
0550221	Co-Cr-Mo Alloy + TiNbN
0550222	Co-Cr-Mo Alloy + TiNbN
0550223	Co-Cr-Mo Alloy + TiNbN
0550224	Co-Cr-Mo Alloy + TiNbN
0550281	Co-Cr-Mo Alloy + TiNbN
0550282	Co-Cr-Mo Alloy + TiNbN
0550283	Co-Cr-Mo Alloy + TiNbN
0550284	Co-Cr-Mo Alloy + TiNbN
0550285	Co-Cr-Mo Alloy + TiNbN
0550321	Co-Cr-Mo Alloy + TiNbN
0550322	Co-Cr-Mo Alloy + TiNbN
0550323	Co-Cr-Mo Alloy + TiNbN
0550324	Co-Cr-Mo Alloy + TiNbN
0550325	Co-Cr-Mo Alloy + TiNbN
0550361	Co-Cr-Mo Alloy + TiNbN
0550362	Co-Cr-Mo Alloy + TiNbN
0550363	Co-Cr-Mo Alloy + TiNbN
0550364	Co-Cr-Mo Alloy + TiNbN
0550365	Co-Cr-Mo Alloy + TiNbN
0550401	Co-Cr-Mo Alloy + TiNbN
0550402	Co-Cr-Mo Alloy + TiNbN
0550403	Co-Cr-Mo Alloy + TiNbN
0550404	Co-Cr-Mo Alloy + TiNbN
0550405	Co-Cr-Mo Alloy + TiNbN
0601015	Ti6Al4V Alloy
0601020	Ti6Al4V Alloy
0601025	Ti6Al4V Alloy

0601030	Ti6Al4V Alloy
0601035	Ti6Al4V Alloy
0601040	Ti6Al4V Alloy
0601045	Ti6Al4V Alloy
0601050	Ti6Al4V Alloy
0602020	Ti6Al4V Alloy
0602025	Ti6Al4V Alloy
0602030	Ti6Al4V Alloy
0602035	Ti6Al4V Alloy
0602040	Ti6Al4V Alloy
0602045	Ti6Al4V Alloy
0602050	Ti6Al4V Alloy
0700000	Ti6Al4V Alloy
0700001	Ti6Al4V Alloy
0700003	Ti6Al4V Alloy
0700004	Ti6Al4V Alloy
0700005	Ti6Al4V Alloy
0700006	Ti6Al4V Alloy
0700007	Ti6Al4V Alloy
0700010	Ti6Al4V Alloy
0702020	Co-Cr-Mo Alloy
0702025	Co-Cr-Mo Alloy
0702030	Co-Cr-Mo Alloy
0702035	Co-Cr-Mo Alloy
0702040	Co-Cr-Mo Alloy
0702045	Co-Cr-Mo Alloy
0702050	Co-Cr-Mo Alloy
0710042	Ti6Al4V Alloy + Hydroxyapatite
0710044	Ti6Al4V Alloy + Hydroxyapatite
0710146	Ti6Al4V Alloy + Hydroxyapatite
0710148	Ti6Al4V Alloy + Hydroxyapatite
0710250	Ti6Al4V Alloy + Hydroxyapatite
0710252	Ti6Al4V Alloy + Hydroxyapatite
0710354	Ti6Al4V Alloy + Hydroxyapatite
0710356	Ti6Al4V Alloy + Hydroxyapatite
0710458	Ti6Al4V Alloy + Hydroxyapatite
0710460	Ti6Al4V Alloy + Hydroxyapatite
0710462	Ti6Al4V Alloy + Hydroxyapatite
0710564	Ti6Al4V Alloy + Hydroxyapatite
0710566	Ti6Al4V Alloy + Hydroxyapatite
0733252	Ti6Al4V Alloy
0733254	Ti6Al4V Alloy
0733256	Ti6Al4V Alloy

0733258	Ti6Al4V Alloy
0733260	Ti6Al4V Alloy
0733262	Ti6Al4V Alloy
0733264	Ti6Al4V Alloy
0733266	Ti6Al4V Alloy
0734144	Ti6Al4V Alloy
0734246	Ti6Al4V Alloy
0734248	Ti6Al4V Alloy
0734250	Ti6Al4V Alloy
0734252	Ti6Al4V Alloy
0734254	Ti6Al4V Alloy
0734256	Ti6Al4V Alloy
0734258	Ti6Al4V Alloy
0734260	Ti6Al4V Alloy
0734262	Ti6Al4V Alloy
0734264	Ti6Al4V Alloy
0734266	Ti6Al4V Alloy
0734650	Ti6Al4V Alloy
0734652	Ti6Al4V Alloy
0738252	Ti6Al4V Alloy
0738254	Ti6Al4V Alloy
0738256	Ti6Al4V Alloy
0738258	Ti6Al4V Alloy
0738260	Ti6Al4V Alloy
0738262	Ti6Al4V Alloy
0738264	Ti6Al4V Alloy
0738266	Ti6Al4V Alloy
0739354	Ti6Al4V Alloy
0739356	Ti6Al4V Alloy
0739358	Ti6Al4V Alloy
0739360	Ti6Al4V Alloy
0739362	Ti6Al4V Alloy
0739364	Ti6Al4V Alloy
0739366	Ti6Al4V Alloy
0740352	Ti6Al4V Alloy
0740354	Ti6Al4V Alloy
0740356	Ti6Al4V Alloy
0740358	Ti6Al4V Alloy
0740360	Ti6Al4V Alloy
0740362	Ti6Al4V Alloy
0740364	Ti6Al4V Alloy
0740366	Ti6Al4V Alloy
0742356	Ti6Al4V Alloy
0742358	Ti6Al4V Alloy
0742360	Ti6Al4V Alloy

0742362	Ti6Al4V Alloy
0742364	Ti6Al4V Alloy
0742366	Ti6Al4V Alloy
0745142	Ti6Al4V Alloy + Ti
0745144	Ti6Al4V Alloy + Ti
0745146	Ti6Al4V Alloy + Ti
0745148	Ti6Al4V Alloy + Ti
0745350	Ti6Al4V Alloy + Ti
0745352	Ti6Al4V Alloy + Ti
0745454	Ti6Al4V Alloy + Ti
0745456	Ti6Al4V Alloy + Ti
0745558	Ti6Al4V Alloy + Ti
0745560	Ti6Al4V Alloy + Ti
0745562	Ti6Al4V Alloy + Ti
0745564	Ti6Al4V Alloy + Ti
0745566	Ti6Al4V Alloy + Ti
0745568	Ti6Al4V Alloy + Ti
0746142	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746144	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746146	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746148	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746350	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746352	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746454	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746456	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746558	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746560	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746562	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746564	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746566	Ti6Al4V Alloy + Ti + Hydroxyapatite
0746568	Ti6Al4V Alloy + Ti + Hydroxyapatite
0753142	Ti6Al4V Alloy
0753144	Ti6Al4V Alloy
0753146	Ti6Al4V Alloy
0753148	Ti6Al4V Alloy
0753350	Ti6Al4V Alloy
0753352	Ti6Al4V Alloy
0753454	Ti6Al4V Alloy
0753456	Ti6Al4V Alloy
0753558	Ti6Al4V Alloy
0753560	Ti6Al4V Alloy
0753562	Ti6Al4V Alloy
0753564	Ti6Al4V Alloy

0753566	Ti6Al4V Alloy
0753568	Ti6Al4V Alloy
0753570	Ti6Al4V Alloy
0753572	Ti6Al4V Alloy
0753574	Ti6Al4V Alloy
0755142	Ti6Al4V Alloy
0755144	Ti6Al4V Alloy
0755246	Ti6Al4V Alloy
0755248	Ti6Al4V Alloy
0755350	Ti6Al4V Alloy
0755352	Ti6Al4V Alloy
0755454	Ti6Al4V Alloy
0755456	Ti6Al4V Alloy
0755558	Ti6Al4V Alloy
0755560	Ti6Al4V Alloy
0755562	Ti6Al4V Alloy
0755564	Ti6Al4V Alloy
0755566	Ti6Al4V Alloy
0756142	Ti6Al4V Alloy + Alumina Ceramics
0756144	Ti6Al4V Alloy + Alumina Ceramics
0756246	Ti6Al4V Alloy + Alumina Ceramics
0756248	Ti6Al4V Alloy + Alumina Ceramics
0756350	Ti6Al4V Alloy + Alumina Ceramics
0756352	Ti6Al4V Alloy + Alumina Ceramics
0756454	Ti6Al4V Alloy + Alumina Ceramics
0756456	Ti6Al4V Alloy + Alumina Ceramics
0756558	Ti6Al4V Alloy + Alumina Ceramics
0756560	Ti6Al4V Alloy + Alumina Ceramics
0756562	Ti6Al4V Alloy + Alumina Ceramics
0756564	Ti6Al4V Alloy + Alumina Ceramics
0756566	Ti6Al4V Alloy + Alumina Ceramics
0757142	Ti6Al4V Alloy
0757144	Ti6Al4V Alloy
0757246	Ti6Al4V Alloy
0757248	Ti6Al4V Alloy
0757350	Ti6Al4V Alloy
0757352	Ti6Al4V Alloy
0757454	Ti6Al4V Alloy
0757456	Ti6Al4V Alloy
0757558	Ti6Al4V Alloy
0757560	Ti6Al4V Alloy
0757562	Ti6Al4V Alloy
0757564	Ti6Al4V Alloy

0757566	Ti6Al4V Alloy
0758142	Ti6Al4V Alloy + Alumina Ceramics
0758144	Ti6Al4V Alloy + Alumina Ceramics
0758246	Ti6Al4V Alloy + Alumina Ceramics
0758248	Ti6Al4V Alloy + Alumina Ceramics
0758350	Ti6Al4V Alloy + Alumina Ceramics
0758352	Ti6Al4V Alloy + Alumina Ceramics
0758454	Ti6Al4V Alloy + Alumina Ceramics
0758456	Ti6Al4V Alloy + Alumina Ceramics
0758558	Ti6Al4V Alloy + Alumina Ceramics
0758560	Ti6Al4V Alloy + Alumina Ceramics
0758562	Ti6Al4V Alloy + Alumina Ceramics
0758564	Ti6Al4V Alloy + Alumina Ceramics
0758566	Ti6Al4V Alloy + Alumina Ceramics
0760000	Ti6Al4V Alloy
0760005	Ti6Al4V Alloy
0760006	Ti6Al4V Alloy
0760007	Ti6Al4V Alloy
0760008	Ti6Al4V Alloy
0760015	Ti6Al4V Alloy + Hydroxyapatite
0760025	Ti6Al4V Alloy + Hydroxyapatite
0760035	Ti6Al4V Alloy + Hydroxyapatite
0760045	Ti6Al4V Alloy + Hydroxyapatite
0760130	Ti6Al4V Alloy + Hydroxyapatite
0760150	Ti6Al4V Alloy + Hydroxyapatite
0760170	Ti6Al4V Alloy + Hydroxyapatite
0761000	Ti6Al4V Alloy
0761015	Ti6Al4V Alloy + Hydroxyapatite
0761020	Ti6Al4V Alloy + Hydroxyapatite
0761025	Ti6Al4V Alloy + Hydroxyapatite
0761030	Ti6Al4V Alloy + Hydroxyapatite
0761035	Ti6Al4V Alloy + Hydroxyapatite
0761040	Ti6Al4V Alloy + Hydroxyapatite
0761045	Ti6Al4V Alloy + Hydroxyapatite
0761050	Ti6Al4V Alloy + Hydroxyapatite
0761055	Ti6Al4V Alloy + Hydroxyapatite
0761060	Ti6Al4V Alloy + Hydroxyapatite
0761065	Ti6Al4V Alloy + Hydroxyapatite
0761115	Ti6Al4V Alloy + Hydroxyapatite
0761120	Ti6Al4V Alloy + Hydroxyapatite
0761125	Ti6Al4V Alloy + Hydroxyapatite
0761130	Ti6Al4V Alloy + Hydroxyapatite
0761135	Ti6Al4V Alloy + Hydroxyapatite

0761140	Ti6Al4V Alloy + Hydroxyapatite
0761145	Ti6Al4V Alloy + Hydroxyapatite
0761150	Ti6Al4V Alloy + Hydroxyapatite
0761155	Ti6Al4V Alloy + Hydroxyapatite
0761160	Ti6Al4V Alloy + Hydroxyapatite
0761165	Ti6Al4V Alloy + Hydroxyapatite
0762150	Ti6Al4V Alloy
0762152	Ti6Al4V Alloy
0762154	Ti6Al4V Alloy
0762356	Ti6Al4V Alloy
0762358	Ti6Al4V Alloy
0762460	Ti6Al4V Alloy
0762462	Ti6Al4V Alloy
0762564	Ti6Al4V Alloy
0762566	Ti6Al4V Alloy
0762852	Ti6Al4V Alloy + Hydroxyapatite
0762854	Ti6Al4V Alloy + Hydroxyapatite
0762856	Ti6Al4V Alloy + Hydroxyapatite
0762858	Ti6Al4V Alloy + Hydroxyapatite
0762860	Ti6Al4V Alloy + Hydroxyapatite
0762864	Ti6Al4V Alloy + Hydroxyapatite
0762868	Ti6Al4V Alloy + Hydroxyapatite
0763000	Ti6Al4V Alloy
0763001	Ti6Al4V Alloy
0763002	Ti6Al4V Alloy
0763003	Ti6Al4V Alloy
0763004	Ti6Al4V Alloy
0763101	Ti6Al4V Alloy
0763102	Ti6Al4V Alloy
0763103	Ti6Al4V Alloy
0763111	Ti6Al4V Alloy
0763112	Ti6Al4V Alloy
0763113	Ti6Al4V Alloy
0764000	Ti6Al4V Alloy
0764001	Ti6Al4V Alloy
0764111	Ti6Al4V Alloy
0764113	Ti6Al4V Alloy
0764115	Ti6Al4V Alloy
0764117	Ti6Al4V Alloy
0764118	Ti6Al4V Alloy
0764141	Ti6Al4V Alloy
0764143	Ti6Al4V Alloy
0764145	Ti6Al4V Alloy

0764147	Ti6Al4V Alloy
0764148	Ti6Al4V Alloy
0764171	Ti6Al4V Alloy
0764173	Ti6Al4V Alloy
0764175	Ti6Al4V Alloy
0764177	Ti6Al4V Alloy
0764178	Ti6Al4V Alloy
0765152	Ti6Al4V Alloy
0765154	Ti6Al4V Alloy
0765156	Ti6Al4V Alloy
0765358	Ti6Al4V Alloy
0765360	Ti6Al4V Alloy
0765462	Ti6Al4V Alloy
0765464	Ti6Al4V Alloy
0765566	Ti6Al4V Alloy
0765568	Ti6Al4V Alloy
0767142	Ti6Al4V Alloy + UHMWPE
0767144	Ti6Al4V Alloy + UHMWPE
0767146	Ti6Al4V Alloy + UHMWPE
0767148	Ti6Al4V Alloy + UHMWPE
0767150	Ti6Al4V Alloy + UHMWPE
0767152	Ti6Al4V Alloy + UHMWPE
0767154	Ti6Al4V Alloy + UHMWPE
0767156	Ti6Al4V Alloy + UHMWPE
0767158	Ti6Al4V Alloy + UHMWPE
0767160	Ti6Al4V Alloy + UHMWPE
0767162	Ti6Al4V Alloy + UHMWPE
0767164	Ti6Al4V Alloy + UHMWPE
0767166	Ti6Al4V Alloy + UHMWPE
0767168	Ti6Al4V Alloy + UHMWPE
0767242	Ti6Al4V Alloy + UHMWPE
0767244	Ti6Al4V Alloy + UHMWPE
0767246	Ti6Al4V Alloy + UHMWPE
0767248	Ti6Al4V Alloy + UHMWPE
0767250	Ti6Al4V Alloy + UHMWPE
0767252	Ti6Al4V Alloy + UHMWPE
0767254	Ti6Al4V Alloy + UHMWPE
0767256	Ti6Al4V Alloy + UHMWPE
0767258	Ti6Al4V Alloy + UHMWPE
0767260	Ti6Al4V Alloy + UHMWPE
0767262	Ti6Al4V Alloy + UHMWPE
0767264	Ti6Al4V Alloy + UHMWPE
0767266	Ti6Al4V Alloy + UHMWPE

0767268	Ti6Al4V Alloy + UHMWPE
0767346	Ti6Al4V Alloy + UHMWPE
0767348	Ti6Al4V Alloy + UHMWPE
0767350	Ti6Al4V Alloy + UHMWPE
0767352	Ti6Al4V Alloy + UHMWPE
0767354	Ti6Al4V Alloy + UHMWPE
0767356	Ti6Al4V Alloy + UHMWPE
0767358	Ti6Al4V Alloy + UHMWPE
0767360	Ti6Al4V Alloy + UHMWPE
0767362	Ti6Al4V Alloy + UHMWPE
0767364	Ti6Al4V Alloy + UHMWPE
0767366	Ti6Al4V Alloy + UHMWPE
0767368	Ti6Al4V Alloy + UHMWPE
0767450	Ti6Al4V Alloy + UHMWPE
0767452	Ti6Al4V Alloy + UHMWPE
0767454	Ti6Al4V Alloy + UHMWPE
0767456	Ti6Al4V Alloy + UHMWPE
0767458	Ti6Al4V Alloy + UHMWPE
0767460	Ti6Al4V Alloy + UHMWPE
0767462	Ti6Al4V Alloy + UHMWPE
0767464	Ti6Al4V Alloy + UHMWPE
0767466	Ti6Al4V Alloy + UHMWPE
0767468	Ti6Al4V Alloy + UHMWPE
0767554	Ti6Al4V Alloy + UHMWPE
0767556	Ti6Al4V Alloy + UHMWPE
0767558	Ti6Al4V Alloy + UHMWPE
0767560	Ti6Al4V Alloy + UHMWPE
0767562	Ti6Al4V Alloy + UHMWPE
0767564	Ti6Al4V Alloy + UHMWPE
0767566	Ti6Al4V Alloy + UHMWPE
0767568	Ti6Al4V Alloy + UHMWPE
0767658	Ti6Al4V Alloy + UHMWPE
0767660	Ti6Al4V Alloy + UHMWPE
0767662	Ti6Al4V Alloy + UHMWPE
0767664	Ti6Al4V Alloy + UHMWPE
0767666	Ti6Al4V Alloy + UHMWPE
0767668	Ti6Al4V Alloy + UHMWPE
0767762	Ti6Al4V Alloy + UHMWPE
0767764	Ti6Al4V Alloy + UHMWPE
0767766	Ti6Al4V Alloy + UHMWPE
0767768	Ti6Al4V Alloy + UHMWPE
0768142	Ti6Al4V Alloy + UHMWPE
0768144	Ti6Al4V Alloy + UHMWPE

0768146	Ti6Al4V Alloy + UHMWPE
0768148	Ti6Al4V Alloy + UHMWPE
0768150	Ti6Al4V Alloy + UHMWPE
0768152	Ti6Al4V Alloy + UHMWPE
0768154	Ti6Al4V Alloy + UHMWPE
0768156	Ti6Al4V Alloy + UHMWPE
0768158	Ti6Al4V Alloy + UHMWPE
0768160	Ti6Al4V Alloy + UHMWPE
0768162	Ti6Al4V Alloy + UHMWPE
0768164	Ti6Al4V Alloy + UHMWPE
0768166	Ti6Al4V Alloy + UHMWPE
0768168	Ti6Al4V Alloy + UHMWPE
0768242	Ti6Al4V Alloy + UHMWPE
0768244	Ti6Al4V Alloy + UHMWPE
0768246	Ti6Al4V Alloy + UHMWPE
0768248	Ti6Al4V Alloy + UHMWPE
0768250	Ti6Al4V Alloy + UHMWPE
0768252	Ti6Al4V Alloy + UHMWPE
0768254	Ti6Al4V Alloy + UHMWPE
0768256	Ti6Al4V Alloy + UHMWPE
0768258	Ti6Al4V Alloy + UHMWPE
0768260	Ti6Al4V Alloy + UHMWPE
0768262	Ti6Al4V Alloy + UHMWPE
0768264	Ti6Al4V Alloy + UHMWPE
0768266	Ti6Al4V Alloy + UHMWPE
0768268	Ti6Al4V Alloy + UHMWPE
0768346	Ti6Al4V Alloy + UHMWPE
0768348	Ti6Al4V Alloy + UHMWPE
0768350	Ti6Al4V Alloy + UHMWPE
0768352	Ti6Al4V Alloy + UHMWPE
0768354	Ti6Al4V Alloy + UHMWPE
0768356	Ti6Al4V Alloy + UHMWPE
0768358	Ti6Al4V Alloy + UHMWPE
0768360	Ti6Al4V Alloy + UHMWPE
0768362	Ti6Al4V Alloy + UHMWPE
0768364	Ti6Al4V Alloy + UHMWPE
0768366	Ti6Al4V Alloy + UHMWPE
0768368	Ti6Al4V Alloy + UHMWPE
0768450	Ti6Al4V Alloy + UHMWPE
0768452	Ti6Al4V Alloy + UHMWPE
0768454	Ti6Al4V Alloy + UHMWPE
0768456	Ti6Al4V Alloy + UHMWPE
0768458	Ti6Al4V Alloy + UHMWPE

0768460	Ti6Al4V Alloy + UHMWPE
0768462	Ti6Al4V Alloy + UHMWPE
0768464	Ti6Al4V Alloy + UHMWPE
0768466	Ti6Al4V Alloy + UHMWPE
0768468	Ti6Al4V Alloy + UHMWPE
0768554	Ti6Al4V Alloy + UHMWPE
0768556	Ti6Al4V Alloy + UHMWPE
0768558	Ti6Al4V Alloy + UHMWPE
0768560	Ti6Al4V Alloy + UHMWPE
0768562	Ti6Al4V Alloy + UHMWPE
0768564	Ti6Al4V Alloy + UHMWPE
0768566	Ti6Al4V Alloy + UHMWPE
0768568	Ti6Al4V Alloy + UHMWPE
0768658	Ti6Al4V Alloy + UHMWPE
0768660	Ti6Al4V Alloy + UHMWPE
0768662	Ti6Al4V Alloy + UHMWPE
0768664	Ti6Al4V Alloy + UHMWPE
0768666	Ti6Al4V Alloy + UHMWPE
0768668	Ti6Al4V Alloy + UHMWPE
0768762	Ti6Al4V Alloy + UHMWPE
0768764	Ti6Al4V Alloy + UHMWPE
0768766	Ti6Al4V Alloy + UHMWPE
0768768	Ti6Al4V Alloy + UHMWPE
0769142	Ti6Al4V Alloy + UHMWPE
0769144	Ti6Al4V Alloy + UHMWPE
0769146	Ti6Al4V Alloy + UHMWPE
0769148	Ti6Al4V Alloy + UHMWPE
0769150	Ti6Al4V Alloy + UHMWPE
0769152	Ti6Al4V Alloy + UHMWPE
0769154	Ti6Al4V Alloy + UHMWPE
0769156	Ti6Al4V Alloy + UHMWPE
0769158	Ti6Al4V Alloy + UHMWPE
0769160	Ti6Al4V Alloy + UHMWPE
0769162	Ti6Al4V Alloy + UHMWPE
0769164	Ti6Al4V Alloy + UHMWPE
0769166	Ti6Al4V Alloy + UHMWPE
0769168	Ti6Al4V Alloy + UHMWPE
0769242	Ti6Al4V Alloy + UHMWPE
0769244	Ti6Al4V Alloy + UHMWPE
0769246	Ti6Al4V Alloy + UHMWPE
0769248	Ti6Al4V Alloy + UHMWPE
0769250	Ti6Al4V Alloy + UHMWPE
0769252	Ti6Al4V Alloy + UHMWPE

0769254	Ti6Al4V Alloy + UHMWPE
0769256	Ti6Al4V Alloy + UHMWPE
0769258	Ti6Al4V Alloy + UHMWPE
0769260	Ti6Al4V Alloy + UHMWPE
0769262	Ti6Al4V Alloy + UHMWPE
0769264	Ti6Al4V Alloy + UHMWPE
0769266	Ti6Al4V Alloy + UHMWPE
0769268	Ti6Al4V Alloy + UHMWPE
0769346	Ti6Al4V Alloy + UHMWPE
0769348	Ti6Al4V Alloy + UHMWPE
0769350	Ti6Al4V Alloy + UHMWPE
0769352	Ti6Al4V Alloy + UHMWPE
0769354	Ti6Al4V Alloy + UHMWPE
0769356	Ti6Al4V Alloy + UHMWPE
0769358	Ti6Al4V Alloy + UHMWPE
0769360	Ti6Al4V Alloy + UHMWPE
0769362	Ti6Al4V Alloy + UHMWPE
0769364	Ti6Al4V Alloy + UHMWPE
0769366	Ti6Al4V Alloy + UHMWPE
0769368	Ti6Al4V Alloy + UHMWPE
0769450	Ti6Al4V Alloy + UHMWPE
0769452	Ti6Al4V Alloy + UHMWPE
0769454	Ti6Al4V Alloy + UHMWPE
0769456	Ti6Al4V Alloy + UHMWPE
0769458	Ti6Al4V Alloy + UHMWPE
0769460	Ti6Al4V Alloy + UHMWPE
0769462	Ti6Al4V Alloy + UHMWPE
0769464	Ti6Al4V Alloy + UHMWPE
0769466	Ti6Al4V Alloy + UHMWPE
0769468	Ti6Al4V Alloy + UHMWPE
0769554	Ti6Al4V Alloy + UHMWPE
0769556	Ti6Al4V Alloy + UHMWPE
0769558	Ti6Al4V Alloy + UHMWPE
0769560	Ti6Al4V Alloy + UHMWPE
0769562	Ti6Al4V Alloy + UHMWPE
0769564	Ti6Al4V Alloy + UHMWPE
0769566	Ti6Al4V Alloy + UHMWPE
0769568	Ti6Al4V Alloy + UHMWPE
0769658	Ti6Al4V Alloy + UHMWPE
0769660	Ti6Al4V Alloy + UHMWPE
0769662	Ti6Al4V Alloy + UHMWPE
0769664	Ti6Al4V Alloy + UHMWPE
0769666	Ti6Al4V Alloy + UHMWPE

0769668	Ti6Al4V Alloy + UHMWPE
0769762	Ti6Al4V Alloy + UHMWPE
0769764	Ti6Al4V Alloy + UHMWPE
0769766	Ti6Al4V Alloy + UHMWPE
0769768	Ti6Al4V Alloy + UHMWPE
0770142	Ti6Al4V Alloy + UHMWPE
0770144	Ti6Al4V Alloy + UHMWPE
0770146	Ti6Al4V Alloy + UHMWPE
0770148	Ti6Al4V Alloy + UHMWPE
0770150	Ti6Al4V Alloy + UHMWPE
0770152	Ti6Al4V Alloy + UHMWPE
0770154	Ti6Al4V Alloy + UHMWPE
0770156	Ti6Al4V Alloy + UHMWPE
0770158	Ti6Al4V Alloy + UHMWPE
0770160	Ti6Al4V Alloy + UHMWPE
0770162	Ti6Al4V Alloy + UHMWPE
0770164	Ti6Al4V Alloy + UHMWPE
0770166	Ti6Al4V Alloy + UHMWPE
0770168	Ti6Al4V Alloy + UHMWPE
0770242	Ti6Al4V Alloy + UHMWPE
0770244	Ti6Al4V Alloy + UHMWPE
0770246	Ti6Al4V Alloy + UHMWPE
0770248	Ti6Al4V Alloy + UHMWPE
0770250	Ti6Al4V Alloy + UHMWPE
0770252	Ti6Al4V Alloy + UHMWPE
0770254	Ti6Al4V Alloy + UHMWPE
0770256	Ti6Al4V Alloy + UHMWPE
0770258	Ti6Al4V Alloy + UHMWPE
0770260	Ti6Al4V Alloy + UHMWPE
0770262	Ti6Al4V Alloy + UHMWPE
0770264	Ti6Al4V Alloy + UHMWPE
0770266	Ti6Al4V Alloy + UHMWPE
0770268	Ti6Al4V Alloy + UHMWPE
0770346	Ti6Al4V Alloy + UHMWPE
0770348	Ti6Al4V Alloy + UHMWPE
0770350	Ti6Al4V Alloy + UHMWPE
0770352	Ti6Al4V Alloy + UHMWPE
0770354	Ti6Al4V Alloy + UHMWPE
0770356	Ti6Al4V Alloy + UHMWPE
0770358	Ti6Al4V Alloy + UHMWPE
0770360	Ti6Al4V Alloy + UHMWPE
0770362	Ti6Al4V Alloy + UHMWPE
0770364	Ti6Al4V Alloy + UHMWPE

0770366	Ti6Al4V Alloy + UHMWPE
0770368	Ti6Al4V Alloy + UHMWPE
0770450	Ti6Al4V Alloy + UHMWPE
0770452	Ti6Al4V Alloy + UHMWPE
0770454	Ti6Al4V Alloy + UHMWPE
0770456	Ti6Al4V Alloy + UHMWPE
0770458	Ti6Al4V Alloy + UHMWPE
0770460	Ti6Al4V Alloy + UHMWPE
0770462	Ti6Al4V Alloy + UHMWPE
0770464	Ti6Al4V Alloy + UHMWPE
0770466	Ti6Al4V Alloy + UHMWPE
0770468	Ti6Al4V Alloy + UHMWPE
0770554	Ti6Al4V Alloy + UHMWPE
0770556	Ti6Al4V Alloy + UHMWPE
0770558	Ti6Al4V Alloy + UHMWPE
0770560	Ti6Al4V Alloy + UHMWPE
0770562	Ti6Al4V Alloy + UHMWPE
0770564	Ti6Al4V Alloy + UHMWPE
0770566	Ti6Al4V Alloy + UHMWPE
0770568	Ti6Al4V Alloy + UHMWPE
0770658	Ti6Al4V Alloy + UHMWPE
0770660	Ti6Al4V Alloy + UHMWPE
0770662	Ti6Al4V Alloy + UHMWPE
0770664	Ti6Al4V Alloy + UHMWPE
0770666	Ti6Al4V Alloy + UHMWPE
0770668	Ti6Al4V Alloy + UHMWPE
0770762	Ti6Al4V Alloy + UHMWPE
0770764	Ti6Al4V Alloy + UHMWPE
0770766	Ti6Al4V Alloy + UHMWPE
0770768	Ti6Al4V Alloy + UHMWPE
0771142	Co-Cr-Mo Alloy + Hydroxyapatite
0771144	Co-Cr-Mo Alloy + Hydroxyapatite
0771146	Co-Cr-Mo Alloy + Hydroxyapatite
0771248	Co-Cr-Mo Alloy + Hydroxyapatite
0771250	Co-Cr-Mo Alloy + Hydroxyapatite
0771252	Co-Cr-Mo Alloy + Hydroxyapatite
0771254	Co-Cr-Mo Alloy + Hydroxyapatite
0771256	Co-Cr-Mo Alloy + Hydroxyapatite
0771258	Co-Cr-Mo Alloy + Hydroxyapatite
0771260	Co-Cr-Mo Alloy + Hydroxyapatite
0771262	Co-Cr-Mo Alloy + Hydroxyapatite
0771264	Co-Cr-Mo Alloy + Hydroxyapatite
0771266	Co-Cr-Mo Alloy + Hydroxyapatite

0771268	Co-Cr-Mo Alloy + Hydroxyapatite
0771600	Ti6Al4V Alloy + Hydroxyapatite
0771800	Co-Cr-Mo Alloy + Hydroxyapatite
0771842	Co-Cr-Mo Alloy + Hydroxyapatite
0771844	Co-Cr-Mo Alloy + Hydroxyapatite
0771846	Co-Cr-Mo Alloy + Hydroxyapatite
0771948	Co-Cr-Mo Alloy + Hydroxyapatite
0771950	Co-Cr-Mo Alloy + Hydroxyapatite
0771952	Co-Cr-Mo Alloy + Hydroxyapatite
0771954	Co-Cr-Mo Alloy + Hydroxyapatite
0771956	Co-Cr-Mo Alloy + Hydroxyapatite
0771958	Co-Cr-Mo Alloy + Hydroxyapatite
0771960	Co-Cr-Mo Alloy + Hydroxyapatite
0771962	Co-Cr-Mo Alloy + Hydroxyapatite
0771964	Co-Cr-Mo Alloy + Hydroxyapatite
0771966	Co-Cr-Mo Alloy + Hydroxyapatite
0771968	Co-Cr-Mo Alloy + Hydroxyapatite
0772142	STAINLESS STEEL
0772144	STAINLESS STEEL
0772146	STAINLESS STEEL
0772248	STAINLESS STEEL
0772250	STAINLESS STEEL
0772252	STAINLESS STEEL
0772254	STAINLESS STEEL
0772256	STAINLESS STEEL
0772258	STAINLESS STEEL
0772260	STAINLESS STEEL
0772262	STAINLESS STEEL
0772264	STAINLESS STEEL
0772266	STAINLESS STEEL
0772268	STAINLESS STEEL
0775044	Ti6Al4V Alloy
0775046	Ti6Al4V Alloy
0775048	Ti6Al4V Alloy
0775050	Ti6Al4V Alloy
0775052	Ti6Al4V Alloy
0775054	Ti6Al4V Alloy
0775056	Ti6Al4V Alloy
0775058	Ti6Al4V Alloy
0775060	Ti6Al4V Alloy
0775062	Ti6Al4V Alloy
0775064	Ti6Al4V Alloy
0775066	Ti6Al4V Alloy

0775068	Ti6Al4V Alloy
0775070	Ti6Al4V Alloy
0775072	Ti6Al4V Alloy
0775074	Ti6Al4V Alloy
0775076	Ti6Al4V Alloy
0775078	Ti6Al4V Alloy
0775080	Ti6Al4V Alloy
0775082	Ti6Al4V Alloy
0775172	Ti6Al4V Alloy
0775174	Ti6Al4V Alloy
0775176	Ti6Al4V Alloy
0775178	Ti6Al4V Alloy
0775180	Ti6Al4V Alloy
0775182	Ti6Al4V Alloy
0776110	Ti6Al4V Alloy
0776115	Ti6Al4V Alloy
0776120	Ti6Al4V Alloy
0776150	Ti6Al4V Alloy
0776155	Ti6Al4V Alloy
0776160	Ti6Al4V Alloy
0776210	Ti6Al4V Alloy
0776215	Ti6Al4V Alloy
0776220	Ti6Al4V Alloy
0776250	Ti6Al4V Alloy
0776255	Ti6Al4V Alloy
0776260	Ti6Al4V Alloy
0776310	Ti6Al4V Alloy
0776315	Ti6Al4V Alloy
0776320	Ti6Al4V Alloy
0776350	Ti6Al4V Alloy
0776355	Ti6Al4V Alloy
0776360	Ti6Al4V Alloy
0776410	Ti6Al4V Alloy
0776415	Ti6Al4V Alloy
0776420	Ti6Al4V Alloy
0776450	Ti6Al4V Alloy
0776455	Ti6Al4V Alloy
0776460	Ti6Al4V Alloy
0776510	Ti6Al4V Alloy
0776515	Ti6Al4V Alloy
0776520	Ti6Al4V Alloy
0776550	Ti6Al4V Alloy
0776555	Ti6Al4V Alloy

0776560	Ti6Al4V Alloy
0777010	Titanium
0777020	Titanium
0777030	Titanium
0778015	Ti6Al4V Alloy
0778020	Ti6Al4V Alloy
0778025	Ti6Al4V Alloy
0778030	Ti6Al4V Alloy
0778035	Ti6Al4V Alloy
0778040	Ti6Al4V Alloy
0778045	Ti6Al4V Alloy
0778050	Ti6Al4V Alloy
0778115	Ti6Al4V Alloy
0778120	Ti6Al4V Alloy
0778125	Ti6Al4V Alloy
0778130	Ti6Al4V Alloy
0778135	Ti6Al4V Alloy
0778140	Ti6Al4V Alloy
0778145	Ti6Al4V Alloy
0778150	Ti6Al4V Alloy
0779015	Ti6Al4V Alloy
0779020	Ti6Al4V Alloy
0779025	Ti6Al4V Alloy
0779030	Ti6Al4V Alloy
0779035	Ti6Al4V Alloy
0779040	Ti6Al4V Alloy
0779045	Ti6Al4V Alloy
0779050	Ti6Al4V Alloy
0780010	Ti6Al4V Alloy
0780020	Ti6Al4V Alloy
0780030	Ti6Al4V Alloy
0780090	Ti6Al4V Alloy
0781044	Ti6Al4V Alloy
0781046	Ti6Al4V Alloy
0781048	Ti6Al4V Alloy
0781050	Ti6Al4V Alloy
0781152	Ti6Al4V Alloy
0781154	Ti6Al4V Alloy
0781256	Ti6Al4V Alloy
0781258	Ti6Al4V Alloy
0781360	Ti6Al4V Alloy
0781362	Ti6Al4V Alloy
0782044	Ti6Al4V Alloy

0782046	Ti6Al4V Alloy
0782048	Ti6Al4V Alloy
0782050	Ti6Al4V Alloy
0782152	Ti6Al4V Alloy
0782154	Ti6Al4V Alloy
0782256	Ti6Al4V Alloy
0782258	Ti6Al4V Alloy
0782360	Ti6Al4V Alloy
0782362	Ti6Al4V Alloy
0783044	Ti6Al4V Alloy
0783046	Ti6Al4V Alloy
0783048	Ti6Al4V Alloy
0783050	Ti6Al4V Alloy
0783152	Ti6Al4V Alloy
0783154	Ti6Al4V Alloy
0783256	Ti6Al4V Alloy
0783258	Ti6Al4V Alloy
0783360	Ti6Al4V Alloy
0783362	Ti6Al4V Alloy
0810000	Alumina Ceramics
0812800	Alumina Ceramics
0813201	Alumina Ceramics
0813202	Alumina Ceramics
0813203	Alumina Ceramics
0813204	Alumina Ceramics
0813205	Alumina Ceramics
0813603	Alumina Ceramics
0813604	Alumina Ceramics
0813605	Alumina Ceramics
0814004	Alumina Ceramics
0814005	Alumina Ceramics
0823201	Alumina Ceramics
0823602	Alumina Ceramics
0824003	Alumina Ceramics
0824404	Alumina Ceramics
0824805	Alumina Ceramics
0832200	UHMWPE
0832201	UHMWPE
0832202	UHMWPE
0832203	UHMWPE
0832204	UHMWPE
0832205	UHMWPE
0832230	UHMWPE

0832231	UHMWPE
0832232	UHMWPE
0832233	UHMWPE
0832234	UHMWPE
0832235	UHMWPE
0832240	UHMWPE
0832241	UHMWPE
0832242	UHMWPE
0832243	UHMWPE
0832244	UHMWPE
0832245	UHMWPE
0832260	UHMWPE
0832261	UHMWPE
0832262	UHMWPE
0832263	UHMWPE
0832264	UHMWPE
0832265	UHMWPE
0832280	UHMWPE
0832281	UHMWPE
0832282	UHMWPE
0832283	UHMWPE
0832284	UHMWPE
0832285	UHMWPE
0832290	UHMWPE
0832291	UHMWPE
0832292	UHMWPE
0832293	UHMWPE
0832294	UHMWPE
0832295	UHMWPE
0832800	UHMWPE
0832801	UHMWPE
0832802	UHMWPE
0832803	UHMWPE
0832804	UHMWPE
0832805	UHMWPE
0832811	UHMWPE
0832813	UHMWPE
0832814	UHMWPE
0832815	UHMWPE
0832821	UHMWPE
0832823	UHMWPE
0832824	UHMWPE
0832825	UHMWPE

0832830	UHMWPE
0832831	UHMWPE
0832832	UHMWPE
0832833	UHMWPE
0832834	UHMWPE
0832835	UHMWPE
0832840	UHMWPE
0832841	UHMWPE
0832842	UHMWPE
0832843	UHMWPE
0832844	UHMWPE
0832845	UHMWPE
0832851	UHMWPE
0832853	UHMWPE
0832854	UHMWPE
0832855	UHMWPE
0832860	UHMWPE
0832861	UHMWPE
0832862	UHMWPE
0832863	UHMWPE
0832864	UHMWPE
0832865	UHMWPE
0832871	UHMWPE
0832873	UHMWPE
0832874	UHMWPE
0832875	UHMWPE
0832880	UHMWPE
0832881	UHMWPE
0832882	UHMWPE
0832883	UHMWPE
0832884	UHMWPE
0832885	UHMWPE
0832886	UHMWPE
0832887	UHMWPE
0832888	UHMWPE
0832889	UHMWPE
0832890	UHMWPE
0832891	UHMWPE
0832892	UHMWPE
0832893	UHMWPE
0832894	UHMWPE
0832895	UHMWPE
0832896	UHMWPE

0832897	UHMWPE
0832898	UHMWPE
0832899	UHMWPE
0833221	UHMWPE
0833223	UHMWPE
0833224	UHMWPE
0833225	UHMWPE
0833230	UHMWPE
0833231	UHMWPE
0833232	UHMWPE
0833233	UHMWPE
0833234	UHMWPE
0833235	UHMWPE
0833260	UHMWPE
0833261	UHMWPE
0833262	UHMWPE
0833263	UHMWPE
0833264	UHMWPE
0833265	UHMWPE
0833271	UHMWPE
0833273	UHMWPE
0833274	UHMWPE
0833275	UHMWPE
0833280	UHMWPE
0833281	UHMWPE
0833282	UHMWPE
0833283	UHMWPE
0833284	UHMWPE
0833285	UHMWPE
0833286	UHMWPE
0833287	UHMWPE
0833288	UHMWPE
0833289	UHMWPE
0833290	UHMWPE
0833291	UHMWPE
0833292	UHMWPE
0833293	UHMWPE
0833294	UHMWPE
0833295	UHMWPE
0833296	UHMWPE
0833297	UHMWPE
0833298	UHMWPE
0833299	UHMWPE

0833623	UHMWPE
0833624	UHMWPE
0833625	UHMWPE
0833633	UHMWPE
0833634	UHMWPE
0833635	UHMWPE
0833663	UHMWPE
0833664	UHMWPE
0833665	UHMWPE
0833673	UHMWPE
0833674	UHMWPE
0833675	UHMWPE
0833683	UHMWPE
0833684	UHMWPE
0833685	UHMWPE
0833687	UHMWPE
0833688	UHMWPE
0833689	UHMWPE
0833693	UHMWPE
0833694	UHMWPE
0833695	UHMWPE
0833697	UHMWPE
0833698	UHMWPE
0833699	UHMWPE
0834024	UHMWPE
0834025	UHMWPE
0834034	UHMWPE
0834035	UHMWPE
0834064	UHMWPE
0834065	UHMWPE
0834074	UHMWPE
0834075	UHMWPE
0834084	UHMWPE
0834085	UHMWPE
0834088	UHMWPE
0834089	UHMWPE
0834094	UHMWPE
0834095	UHMWPE
0834098	UHMWPE
0834099	UHMWPE
0842200	UHMWPE
0842201	UHMWPE
0842202	UHMWPE

0842203	UHMWPE
0842204	UHMWPE
0842205	UHMWPE
0842230	UHMWPE
0842231	UHMWPE
0842232	UHMWPE
0842233	UHMWPE
0842234	UHMWPE
0842235	UHMWPE
0842240	UHMWPE
0842241	UHMWPE
0842242	UHMWPE
0842243	UHMWPE
0842244	UHMWPE
0842245	UHMWPE
0842260	UHMWPE
0842261	UHMWPE
0842262	UHMWPE
0842263	UHMWPE
0842264	UHMWPE
0842265	UHMWPE
0842280	UHMWPE
0842281	UHMWPE
0842282	UHMWPE
0842283	UHMWPE
0842284	UHMWPE
0842285	UHMWPE
0842290	UHMWPE
0842291	UHMWPE
0842292	UHMWPE
0842293	UHMWPE
0842294	UHMWPE
0842295	UHMWPE
0842800	UHMWPE
0842801	UHMWPE
0842802	UHMWPE
0842803	UHMWPE
0842804	UHMWPE
0842805	UHMWPE
0842811	UHMWPE
0842813	UHMWPE
0842814	UHMWPE
0842815	UHMWPE

0842821	UHMWPE
0842823	UHMWPE
0842824	UHMWPE
0842825	UHMWPE
0842830	UHMWPE
0842831	UHMWPE
0842832	UHMWPE
0842833	UHMWPE
0842834	UHMWPE
0842835	UHMWPE
0842840	UHMWPE
0842841	UHMWPE
0842842	UHMWPE
0842843	UHMWPE
0842844	UHMWPE
0842845	UHMWPE
0842851	UHMWPE
0842852	UHMWPE
0842853	UHMWPE
0842854	UHMWPE
0842860	UHMWPE
0842861	UHMWPE
0842862	UHMWPE
0842863	UHMWPE
0842864	UHMWPE
0842865	UHMWPE
0842871	UHMWPE
0842873	UHMWPE
0842874	UHMWPE
0842875	UHMWPE
0842880	UHMWPE
0842881	UHMWPE
0842882	UHMWPE
0842883	UHMWPE
0842884	UHMWPE
0842885	UHMWPE
0842886	UHMWPE
0842887	UHMWPE
0842888	UHMWPE
0842889	UHMWPE
0842890	UHMWPE
0842891	UHMWPE
0842892	UHMWPE

0842893	UHMWPE
0842894	UHMWPE
0842895	UHMWPE
0842896	UHMWPE
0842897	UHMWPE
0842898	UHMWPE
0842899	UHMWPE
0843221	UHMWPE
0843223	UHMWPE
0843224	UHMWPE
0843225	UHMWPE
0843230	UHMWPE
0843231	UHMWPE
0843232	UHMWPE
0843233	UHMWPE
0843234	UHMWPE
0843235	UHMWPE
0843260	UHMWPE
0843261	UHMWPE
0843262	UHMWPE
0843263	UHMWPE
0843264	UHMWPE
0843265	UHMWPE
0843271	UHMWPE
0843273	UHMWPE
0843274	UHMWPE
0843275	UHMWPE
0843280	UHMWPE
0843281	UHMWPE
0843282	UHMWPE
0843283	UHMWPE
0843284	UHMWPE
0843285	UHMWPE
0843286	UHMWPE
0843287	UHMWPE
0843288	UHMWPE
0843289	UHMWPE
0843290	UHMWPE
0843291	UHMWPE
0843292	UHMWPE
0843293	UHMWPE
0843294	UHMWPE
0843295	UHMWPE

0843296	UHMWPE
0843297	UHMWPE
0843298	UHMWPE
0843299	UHMWPE
0843623	UHMWPE
0843624	UHMWPE
0843625	UHMWPE
0843633	UHMWPE
0843634	UHMWPE
0843635	UHMWPE
0843663	UHMWPE
0843664	UHMWPE
0843665	UHMWPE
0843673	UHMWPE
0843674	UHMWPE
0843675	UHMWPE
0843683	UHMWPE
0843684	UHMWPE
0843685	UHMWPE
0843687	UHMWPE
0843688	UHMWPE
0843689	UHMWPE
0843693	UHMWPE
0843694	UHMWPE
0843695	UHMWPE
0843697	UHMWPE
0843698	UHMWPE
0843699	UHMWPE
0844024	UHMWPE
0844025	UHMWPE
0844034	UHMWPE
0844035	UHMWPE
0844064	UHMWPE
0844065	UHMWPE
0844074	UHMWPE
0844075	UHMWPE
0844084	UHMWPE
0844085	UHMWPE
0844088	UHMWPE
0844089	UHMWPE
0844094	UHMWPE
0844095	UHMWPE
0844098	UHMWPE

0844099	UHMWPE
0852800	UHMWPE
0852801	UHMWPE
0852802	UHMWPE
0852803	UHMWPE
0852804	UHMWPE
0852805	UHMWPE
0852811	UHMWPE
0852813	UHMWPE
0852814	UHMWPE
0852815	UHMWPE
0852821	UHMWPE
0852823	UHMWPE
0852824	UHMWPE
0852825	UHMWPE
0852830	UHMWPE
0852831	UHMWPE
0852832	UHMWPE
0852833	UHMWPE
0852834	UHMWPE
0852835	UHMWPE
0852840	UHMWPE
0852841	UHMWPE
0852842	UHMWPE
0852843	UHMWPE
0852844	UHMWPE
0852845	UHMWPE
0852851	UHMWPE
0852852	UHMWPE
0852853	UHMWPE
0852854	UHMWPE
0852860	UHMWPE
0852861	UHMWPE
0852862	UHMWPE
0852863	UHMWPE
0852864	UHMWPE
0852865	UHMWPE
0852871	UHMWPE
0852873	UHMWPE
0852874	UHMWPE
0852875	UHMWPE
0852880	UHMWPE
0852881	UHMWPE

0852882	UHMWPE
0852883	UHMWPE
0852884	UHMWPE
0852885	UHMWPE
0852886	UHMWPE
0852887	UHMWPE
0852888	UHMWPE
0852889	UHMWPE
0852890	UHMWPE
0852891	UHMWPE
0852892	UHMWPE
0852893	UHMWPE
0852894	UHMWPE
0852895	UHMWPE
0852896	UHMWPE
0852897	UHMWPE
0852898	UHMWPE
0852899	UHMWPE
0853221	UHMWPE
0853223	UHMWPE
0853224	UHMWPE
0853225	UHMWPE
0853230	UHMWPE
0853231	UHMWPE
0853232	UHMWPE
0853233	UHMWPE
0853234	UHMWPE
0853235	UHMWPE
0853260	UHMWPE
0853261	UHMWPE
0853262	UHMWPE
0853263	UHMWPE
0853264	UHMWPE
0853265	UHMWPE
0853271	UHMWPE
0853273	UHMWPE
0853274	UHMWPE
0853275	UHMWPE
0853280	UHMWPE
0853281	UHMWPE
0853282	UHMWPE
0853283	UHMWPE
0853284	UHMWPE

0853285	UHMWPE
0853286	UHMWPE
0853287	UHMWPE
0853288	UHMWPE
0853289	UHMWPE
0853290	UHMWPE
0853291	UHMWPE
0853292	UHMWPE
0853293	UHMWPE
0853294	UHMWPE
0853295	UHMWPE
0853296	UHMWPE
0853297	UHMWPE
0853298	UHMWPE
0853299	UHMWPE
0853623	UHMWPE
0853624	UHMWPE
0853625	UHMWPE
0853633	UHMWPE
0853634	UHMWPE
0853635	UHMWPE
0853663	UHMWPE
0853664	UHMWPE
0853665	UHMWPE
0853673	UHMWPE
0853674	UHMWPE
0853675	UHMWPE
0853683	UHMWPE
0853684	UHMWPE
0853685	UHMWPE
0853687	UHMWPE
0853688	UHMWPE
0853689	UHMWPE
0853693	UHMWPE
0853694	UHMWPE
0853695	UHMWPE
0853697	UHMWPE
0853698	UHMWPE
0853699	UHMWPE
0854024	UHMWPE
0854025	UHMWPE
0854034	UHMWPE
0854035	UHMWPE

0854064	UHMWPE
0854065	UHMWPE
0854074	UHMWPE
0854075	UHMWPE
0854084	UHMWPE
0854085	UHMWPE
0854088	UHMWPE
0854089	UHMWPE
0854094	UHMWPE
0854095	UHMWPE
0854098	UHMWPE
0854099	UHMWPE
0871001	UHMWPE
0871002	UHMWPE
0871003	UHMWPE
0871004	UHMWPE
0871005	UHMWPE
0871006	UHMWPE
0871007	UHMWPE
0871008	UHMWPE
0871009	UHMWPE
0871010	UHMWPE
0871011	UHMWPE
0871012	UHMWPE
0871013	UHMWPE
0871014	UHMWPE
0871015	UHMWPE
0871016	UHMWPE
0871017	UHMWPE
0871018	UHMWPE
0871019	UHMWPE
0871020	UHMWPE
0871101	UHMWPE
0871102	UHMWPE
0871103	UHMWPE
0871104	UHMWPE
0871105	UHMWPE
0871106	UHMWPE
0871107	UHMWPE
0871108	UHMWPE
0871109	UHMWPE
0871110	UHMWPE
0871111	UHMWPE

0871112	UHMWPE
0871113	UHMWPE
0871114	UHMWPE
0871115	UHMWPE
0871116	UHMWPE
0871117	UHMWPE
0871118	UHMWPE
0871119	UHMWPE
0871120	UHMWPE
0871201	UHMWPE
0871202	UHMWPE
0871203	UHMWPE
0871204	UHMWPE
0871205	UHMWPE
0871206	UHMWPE
0871207	UHMWPE
0871208	UHMWPE
0871209	UHMWPE
0871210	UHMWPE
0871211	UHMWPE
0871212	UHMWPE
0871213	UHMWPE
0871214	UHMWPE
0871215	UHMWPE
0871216	UHMWPE
0871217	UHMWPE
0871218	UHMWPE
0871219	UHMWPE
0871220	UHMWPE
0871301	Co-Cr-Mo Alloy + TiNbN
0871302	Co-Cr-Mo Alloy + TiNbN
0871303	Co-Cr-Mo Alloy + TiNbN
0871304	Co-Cr-Mo Alloy + TiNbN
0871305	Co-Cr-Mo Alloy + TiNbN
0871311	Co-Cr-Mo Alloy + TiNbN
0871312	Co-Cr-Mo Alloy + TiNbN
0871313	Co-Cr-Mo Alloy + TiNbN
0871314	Co-Cr-Mo Alloy + TiNbN
0871315	Co-Cr-Mo Alloy + TiNbN
0871321	Co-Cr-Mo Alloy + TiNbN
0871322	Co-Cr-Mo Alloy + TiNbN
0871323	Co-Cr-Mo Alloy + TiNbN
0871324	Co-Cr-Mo Alloy + TiNbN

0871325	Co-Cr-Mo Alloy + TiNbN
0871331	Co-Cr-Mo Alloy + TiNbN
0871332	Co-Cr-Mo Alloy + TiNbN
0871333	Co-Cr-Mo Alloy + TiNbN
0871334	Co-Cr-Mo Alloy + TiNbN
0871335	Co-Cr-Mo Alloy + TiNbN
0871341	Co-Cr-Mo Alloy + TiNbN
0871342	Co-Cr-Mo Alloy + TiNbN
0871343	Co-Cr-Mo Alloy + TiNbN
0871344	Co-Cr-Mo Alloy + TiNbN
0871345	Co-Cr-Mo Alloy + TiNbN
0871401	Co-Cr-Mo Alloy + TiNbN
0871402	Co-Cr-Mo Alloy + TiNbN
0871403	Co-Cr-Mo Alloy + TiNbN
0871404	Co-Cr-Mo Alloy+ TiNbN
0871405	Co-Cr-Mo Alloy + TiNbN
0871411	Co-Cr-Mo Alloy + TiNbN
0871412	Co-Cr-Mo Alloy + TiNbN
0871413	Co-Cr-Mo Alloy + TiNbN
0871414	Co-Cr-Mo Alloy + TiNbN
0871415	Co-Cr-Mo Alloy + TiNbN
0871421	Co-Cr-Mo Alloy + TiNbN
0871422	Co-Cr-Mo Alloy + TiNbN
0871423	Co-Cr-Mo Alloy + TiNbN
0871424	Co-Cr-Mo Alloy + TiNbN
0871425	Co-Cr-Mo Alloy + TiNbN
0871431	Co-Cr-Mo Alloy + TiNbN
0871432	Co-Cr-Mo Alloy + TiNbN
0871433	Co-Cr-Mo Alloy + TiNbN
0871434	Co-Cr-Mo Alloy + TiNbN
0871435	Co-Cr-Mo Alloy + TiNbN
0871441	Co-Cr-Mo Alloy + TiNbN
0871442	Co-Cr-Mo Alloy + TiNbN
0871443	Co-Cr-Mo Alloy + TiNbN
0871444	Co-Cr-Mo Alloy + TiNbN
0871445	Co-Cr-Mo Alloy + TiNbN
0872142	UHMWPE
0872144	UHMWPE
0872146	UHMWPE
0872148	UHMWPE
0872150	UHMWPE
0872152	UHMWPE
0872154	UHMWPE

0872156	UHMWPE
0872158	UHMWPE
0872160	UHMWPE
0872162	UHMWPE
0872164	UHMWPE
0872166	UHMWPE
0872168	UHMWPE
0872246	UHMWPE
0872248	UHMWPE
0872250	UHMWPE
0872252	UHMWPE
0872254	UHMWPE
0872256	UHMWPE
0872258	UHMWPE
0872260	UHMWPE
0872262	UHMWPE
0872264	UHMWPE
0872266	UHMWPE
0872268	UHMWPE
0873142	UHMWPE
0873144	UHMWPE
0873146	UHMWPE
0873148	UHMWPE
0873150	UHMWPE
0873152	UHMWPE
0873154	UHMWPE
0873156	UHMWPE
0873158	UHMWPE
0873160	UHMWPE
0873162	UHMWPE
0873164	UHMWPE
0873166	UHMWPE
0873168	UHMWPE
0873246	UHMWPE
0873248	UHMWPE
0873250	UHMWPE
0873252	UHMWPE
0873254	UHMWPE
0873256	UHMWPE
0873258	UHMWPE
0873260	UHMWPE
0873262	UHMWPE
0873264	UHMWPE

0873266	UHMWPE
0873268	UHMWPE
0874142	UHMWPE
0874144	UHMWPE
0874146	UHMWPE
0874148	UHMWPE
0874150	UHMWPE
0874152	UHMWPE
0874154	UHMWPE
0874156	UHMWPE
0874158	UHMWPE
0874160	UHMWPE
0874162	UHMWPE
0874164	UHMWPE
0874166	UHMWPE
0874168	UHMWPE
0874246	UHMWPE
0874248	UHMWPE
0874250	UHMWPE
0874252	UHMWPE
0874254	UHMWPE
0874256	UHMWPE
0874258	UHMWPE
0874260	UHMWPE
0874262	UHMWPE
0874264	UHMWPE
0874266	UHMWPE
0874268	UHMWPE
0875142	UHMWPE
0875144	UHMWPE
0875146	UHMWPE
0875148	UHMWPE
0875150	UHMWPE
0875152	UHMWPE
0875154	UHMWPE
0875156	UHMWPE
0875158	UHMWPE
0875160	UHMWPE
0875162	UHMWPE
0875164	UHMWPE
0875166	UHMWPE
0875168	UHMWPE
0875246	UHMWPE

0875248	UHMWPE
0875250	UHMWPE
0875252	UHMWPE
0875254	UHMWPE
0875256	UHMWPE
0875258	UHMWPE
0875260	UHMWPE
0875262	UHMWPE
0875264	UHMWPE
0875266	UHMWPE
0875268	UHMWPE
0876142	UHMWPE
0876144	UHMWPE
0876146	UHMWPE
0876148	UHMWPE
0876150	UHMWPE
0876152	UHMWPE
0876154	UHMWPE
0876156	UHMWPE
0876158	UHMWPE
0876160	UHMWPE
0876162	UHMWPE
0876164	UHMWPE
0876166	UHMWPE
0876168	UHMWPE
0876246	UHMWPE
0876248	UHMWPE
0876250	UHMWPE
0876252	UHMWPE
0876254	UHMWPE
0876256	UHMWPE
0876258	UHMWPE
0876260	UHMWPE
0876262	UHMWPE
0876264	UHMWPE
0876266	UHMWPE
0876268	UHMWPE
0877142	UHMWPE
0877144	UHMWPE
0877146	UHMWPE
0877148	UHMWPE
0877150	UHMWPE
0877152	UHMWPE

0877154	UHMWPE
0877156	UHMWPE
0877158	UHMWPE
0877160	UHMWPE
0877162	UHMWPE
0877164	UHMWPE
0877166	UHMWPE
0877168	UHMWPE
0877246	UHMWPE
0877248	UHMWPE
0877250	UHMWPE
0877252	UHMWPE
0877254	UHMWPE
0877256	UHMWPE
0877258	UHMWPE
0877260	UHMWPE
0877262	UHMWPE
0877264	UHMWPE
0877266	UHMWPE
0877268	UHMWPE
0878010	UHMWPE
0878020	UHMWPE
0878030	UHMWPE
0878040	UHMWPE
0879010	UHMWPE
0879020	UHMWPE
0879030	UHMWPE
0879040	UHMWPE
0880010	UHMWPE
0880020	UHMWPE
0880030	UHMWPE
0880040	UHMWPE
0881239	Alumina Ceramics
0881244	Alumina Ceramics
0881344	Alumina Ceramics
0881348	Alumina Ceramics
0881352	Alumina Ceramics
0881448	Alumina Ceramics
0881452	Alumina Ceramics
0895001	Co-Cr-Mo Alloy
0895002	Co-Cr-Mo Alloy
0895003	Co-Cr-Mo Alloy
0895004	Co-Cr-Mo Alloy

0895005	Co-Cr-Mo Alloy
0895011	Co-Cr-Mo Alloy
0895012	Co-Cr-Mo Alloy
0895013	Co-Cr-Mo Alloy
0895014	Co-Cr-Mo Alloy
0895015	Co-Cr-Mo Alloy
0895021	Co-Cr-Mo Alloy
0895022	Co-Cr-Mo Alloy
0895023	Co-Cr-Mo Alloy
0895024	Co-Cr-Mo Alloy
0895025	Co-Cr-Mo Alloy
0895031	Co-Cr-Mo Alloy
0895032	Co-Cr-Mo Alloy
0895033	Co-Cr-Mo Alloy
0895034	Co-Cr-Mo Alloy
0895035	Co-Cr-Mo Alloy
0895041	Co-Cr-Mo Alloy
0895042	Co-Cr-Mo Alloy
0895043	Co-Cr-Mo Alloy
0895044	Co-Cr-Mo Alloy
0895045	Co-Cr-Mo Alloy
0895101	Co-Cr-Mo Alloy
0895102	Co-Cr-Mo Alloy
0895103	Co-Cr-Mo Alloy
0895104	Co-Cr-Mo Alloy
0895105	Co-Cr-Mo Alloy
0895111	Co-Cr-Mo Alloy
0895112	Co-Cr-Mo Alloy
0895113	Co-Cr-Mo Alloy
0895114	Co-Cr-Mo Alloy
0895115	Co-Cr-Mo Alloy
0895121	Co-Cr-Mo Alloy
0895122	Co-Cr-Mo Alloy
0895123	Co-Cr-Mo Alloy
0895124	Co-Cr-Mo Alloy
0895125	Co-Cr-Mo Alloy
0895131	Co-Cr-Mo Alloy
0895132	Co-Cr-Mo Alloy
0895133	Co-Cr-Mo Alloy
0895134	Co-Cr-Mo Alloy
0895135	Co-Cr-Mo Alloy
0895141	Co-Cr-Mo Alloy
0895142	Co-Cr-Mo Alloy

0895143	Co-Cr-Mo Alloy
0895144	Co-Cr-Mo Alloy
0895145	Co-Cr-Mo Alloy
0896001	Co-Cr-Mo Alloy
0896002	Co-Cr-Mo Alloy
0896003	Co-Cr-Mo Alloy
0896004	Co-Cr-Mo Alloy
0896005	Co-Cr-Mo Alloy
0896010	Co-Cr-Mo Alloy
0896011	Co-Cr-Mo Alloy
0896012	Co-Cr-Mo Alloy
0896013	Co-Cr-Mo Alloy
0896014	Co-Cr-Mo Alloy
0896015	Co-Cr-Mo Alloy
0896020	Co-Cr-Mo Alloy
0896021	Co-Cr-Mo Alloy
0896022	Co-Cr-Mo Alloy
0896023	Co-Cr-Mo Alloy
0896024	Co-Cr-Mo Alloy
0896025	Co-Cr-Mo Alloy
0896030	Co-Cr-Mo Alloy
0896031	Co-Cr-Mo Alloy
0896032	Co-Cr-Mo Alloy
0896033	Co-Cr-Mo Alloy
0896034	Co-Cr-Mo Alloy
0896035	Co-Cr-Mo Alloy
0896040	Co-Cr-Mo Alloy
0896041	Co-Cr-Mo Alloy
0896042	Co-Cr-Mo Alloy
0896043	Co-Cr-Mo Alloy
0896044	Co-Cr-Mo Alloy
0896045	Co-Cr-Mo Alloy
0896101	Co-Cr-Mo Alloy
0896102	Co-Cr-Mo Alloy
0896103	Co-Cr-Mo Alloy
0896104	Co-Cr-Mo Alloy
0896105	Co-Cr-Mo Alloy
0896111	Co-Cr-Mo Alloy
0896112	Co-Cr-Mo Alloy
0896113	Co-Cr-Mo Alloy
0896114	Co-Cr-Mo Alloy
0896115	Co-Cr-Mo Alloy
0896121	Co-Cr-Mo Alloy

0896122	Co-Cr-Mo Alloy
0896123	Co-Cr-Mo Alloy
0896124	Co-Cr-Mo Alloy
0896125	Co-Cr-Mo Alloy
0896131	Co-Cr-Mo Alloy
0896132	Co-Cr-Mo Alloy
0896133	Co-Cr-Mo Alloy
0896134	Co-Cr-Mo Alloy
0896135	Co-Cr-Mo Alloy
0896141	Co-Cr-Mo Alloy
0896142	Co-Cr-Mo Alloy
0896143	Co-Cr-Mo Alloy
0896144	Co-Cr-Mo Alloy
0896145	Co-Cr-Mo Alloy
0896501	Co-Cr-Mo Alloy + TiNbN
0896502	Co-Cr-Mo Alloy + TiNbN
0896503	Co-Cr-Mo Alloy + TiNbN
0896504	Co-Cr-Mo Alloy + TiNbN
0896505	Co-Cr-Mo Alloy + TiNbN
0896511	Co-Cr-Mo Alloy + TiNbN
0896512	Co-Cr-Mo Alloy + TiNbN
0896513	Co-Cr-Mo Alloy + TiNbN
0896514	Co-Cr-Mo Alloy + TiNbN
0896515	Co-Cr-Mo Alloy + TiNbN
0896521	Co-Cr-Mo Alloy + TiNbN
0896522	Co-Cr-Mo Alloy + TiNbN
0896523	Co-Cr-Mo Alloy + TiNbN
0896524	Co-Cr-Mo Alloy + TiNbN
0896525	Co-Cr-Mo Alloy + TiNbN
0896531	Co-Cr-Mo Alloy + TiNbN
0896532	Co-Cr-Mo Alloy + TiNbN
0896533	Co-Cr-Mo Alloy + TiNbN
0896534	Co-Cr-Mo Alloy + TiNbN
0896535	Co-Cr-Mo Alloy + TiNbN
0896541	Co-Cr-Mo Alloy + TiNbN
0896542	Co-Cr-Mo Alloy + TiNbN
0896543	Co-Cr-Mo Alloy + TiNbN
0896544	Co-Cr-Mo Alloy + TiNbN
0896545	Co-Cr-Mo Alloy + TiNbN
0896601	Co-Cr-Mo Alloy + TiNbN
0896602	Co-Cr-Mo Alloy + TiNbN
0896603	Co-Cr-Mo Alloy + TiNbN
0896604	Co-Cr-Mo Alloy + TiNbN

0896605	Co-Cr-Mo Alloy + TiNbN
0896611	Co-Cr-Mo Alloy + TiNbN
0896612	Co-Cr-Mo Alloy + TiNbN
0896613	Co-Cr-Mo Alloy + TiNbN
0896614	Co-Cr-Mo Alloy + TiNbN
0896615	Co-Cr-Mo Alloy + TiNbN
0896621	Co-Cr-Mo Alloy + TiNbN
0896622	Co-Cr-Mo Alloy + TiNbN
0896623	Co-Cr-Mo Alloy + TiNbN
0896624	Co-Cr-Mo Alloy + TiNbN
0896625	Co-Cr-Mo Alloy + TiNbN
0896631	Co-Cr-Mo Alloy + TiNbN
0896632	Co-Cr-Mo Alloy + TiNbN
0896633	Co-Cr-Mo Alloy + TiNbN
0896634	Co-Cr-Mo Alloy + TiNbN
0896635	Co-Cr-Mo Alloy + TiNbN
0896641	Co-Cr-Mo Alloy + TiNbN
0896642	Co-Cr-Mo Alloy + TiNbN
0896643	Co-Cr-Mo Alloy + TiNbN
0896644	Co-Cr-Mo Alloy + TiNbN
0896645	Co-Cr-Mo Alloy + TiNbN
0902800	UHMWPE
0902801	UHMWPE
0902802	UHMWPE
0902803	UHMWPE
0902804	UHMWPE
0902810	UHMWPE
0902811	UHMWPE
0902812	UHMWPE
0902813	UHMWPE
0902814	UHMWPE
0902820	UHMWPE
0902821	UHMWPE
0902822	UHMWPE
0902823	UHMWPE
0902824	UHMWPE
0902830	UHMWPE
0902831	UHMWPE
0902832	UHMWPE
0902833	UHMWPE
0902834	UHMWPE
0902840	UHMWPE
0902841	UHMWPE

0902842	UHMWPE
0902843	UHMWPE
0902844	UHMWPE
0902850	UHMWPE
0902851	UHMWPE
0902852	UHMWPE
0902853	UHMWPE
0902854	UHMWPE
0903210	UHMWPE
0903211	UHMWPE
0903212	UHMWPE
0903213	UHMWPE
0903214	UHMWPE
0903230	UHMWPE
0903231	UHMWPE
0903232	UHMWPE
0903233	UHMWPE
0903234	UHMWPE
0903240	UHMWPE
0903241	UHMWPE
0903242	UHMWPE
0903243	UHMWPE
0903244	UHMWPE
0903250	UHMWPE
0903251	UHMWPE
0903252	UHMWPE
0903253	UHMWPE
0903254	UHMWPE
0903610	UHMWPE
0903611	UHMWPE
0903612	UHMWPE
0903613	UHMWPE
0903614	UHMWPE
0903630	UHMWPE
0903631	UHMWPE
0903632	UHMWPE
0903633	UHMWPE
0903634	UHMWPE
0903640	UHMWPE
0903641	UHMWPE
0903642	UHMWPE
0903643	UHMWPE
0903644	UHMWPE

0903650	UHMWPE
0903651	UHMWPE
0903652	UHMWPE
0903653	UHMWPE
0903654	UHMWPE
0904010	UHMWPE
0904011	UHMWPE
0904012	UHMWPE
0904013	UHMWPE
0904014	UHMWPE
0904030	UHMWPE
0904031	UHMWPE
0904032	UHMWPE
0904033	UHMWPE
0904034	UHMWPE
0904040	UHMWPE
0904041	UHMWPE
0904042	UHMWPE
0904043	UHMWPE
0904044	UHMWPE
0904050	UHMWPE
0904051	UHMWPE
0904052	UHMWPE
0904053	UHMWPE
0904054	UHMWPE
0912800	UHMWPE
0912801	UHMWPE
0912802	UHMWPE
0912803	UHMWPE
0912804	UHMWPE
0912810	UHMWPE
0912811	UHMWPE
0912812	UHMWPE
0912813	UHMWPE
0912814	UHMWPE
0912820	UHMWPE
0912821	UHMWPE
0912822	UHMWPE
0912823	UHMWPE
0912824	UHMWPE
0912830	UHMWPE
0912831	UHMWPE
0912832	UHMWPE

0912833	UHMWPE
0912834	UHMWPE
0912840	UHMWPE
0912841	UHMWPE
0912842	UHMWPE
0912843	UHMWPE
0912844	UHMWPE
0912850	UHMWPE
0912851	UHMWPE
0912852	UHMWPE
0912853	UHMWPE
0912854	UHMWPE
0913210	UHMWPE
0913211	UHMWPE
0913212	UHMWPE
0913213	UHMWPE
0913214	UHMWPE
0913230	UHMWPE
0913231	UHMWPE
0913232	UHMWPE
0913233	UHMWPE
0913234	UHMWPE
0913240	UHMWPE
0913241	UHMWPE
0913242	UHMWPE
0913243	UHMWPE
0913244	UHMWPE
0913250	UHMWPE
0913251	UHMWPE
0913252	UHMWPE
0913253	UHMWPE
0913254	UHMWPE
0913610	UHMWPE
0913611	UHMWPE
0913612	UHMWPE
0913613	UHMWPE
0913614	UHMWPE
0913630	UHMWPE
0913631	UHMWPE
0913632	UHMWPE
0913633	UHMWPE
0913634	UHMWPE
0913640	UHMWPE

0913641	UHMWPE
0913642	UHMWPE
0913643	UHMWPE
0913644	UHMWPE
0913650	UHMWPE
0913651	UHMWPE
0913652	UHMWPE
0913653	UHMWPE
0913654	UHMWPE
0914010	UHMWPE
0914011	UHMWPE
0914012	UHMWPE
0914013	UHMWPE
0914014	UHMWPE
0914030	UHMWPE
0914031	UHMWPE
0914032	UHMWPE
0914033	UHMWPE
0914034	UHMWPE
0914040	UHMWPE
0914041	UHMWPE
0914042	UHMWPE
0914043	UHMWPE
0914044	UHMWPE
0914050	UHMWPE
0914051	UHMWPE
0914052	UHMWPE
0914053	UHMWPE
0914054	UHMWPE
0922800	UHMWPE
0922801	UHMWPE
0922802	UHMWPE
0922803	UHMWPE
0922804	UHMWPE
0922810	UHMWPE
0922811	UHMWPE
0922812	UHMWPE
0922813	UHMWPE
0922814	UHMWPE
0922820	UHMWPE
0922821	UHMWPE
0922822	UHMWPE
0922823	UHMWPE

0922824	UHMWPE
0922830	UHMWPE
0922831	UHMWPE
0922832	UHMWPE
0922833	UHMWPE
0922834	UHMWPE
0922840	UHMWPE
0922841	UHMWPE
0922842	UHMWPE
0922843	UHMWPE
0922844	UHMWPE
0922850	UHMWPE
0922851	UHMWPE
0922852	UHMWPE
0922853	UHMWPE
0922854	UHMWPE
0923210	UHMWPE
0923211	UHMWPE
0923212	UHMWPE
0923213	UHMWPE
0923214	UHMWPE
0923230	UHMWPE
0923231	UHMWPE
0923232	UHMWPE
0923233	UHMWPE
0923234	UHMWPE
0923240	UHMWPE
0923241	UHMWPE
0923242	UHMWPE
0923243	UHMWPE
0923244	UHMWPE
0923250	UHMWPE
0923251	UHMWPE
0923252	UHMWPE
0923253	UHMWPE
0923254	UHMWPE
0923610	UHMWPE
0923611	UHMWPE
0923612	UHMWPE
0923613	UHMWPE
0923614	UHMWPE
0923630	UHMWPE
0923631	UHMWPE

0923632	UHMWPE
0923633	UHMWPE
0923634	UHMWPE
0923640	UHMWPE
0923641	UHMWPE
0923642	UHMWPE
0923643	UHMWPE
0923644	UHMWPE
0923650	UHMWPE
0923651	UHMWPE
0923652	UHMWPE
0923653	UHMWPE
0923654	UHMWPE
0924010	UHMWPE
0924011	UHMWPE
0924012	UHMWPE
0924013	UHMWPE
0924014	UHMWPE
0924030	UHMWPE
0924031	UHMWPE
0924032	UHMWPE
0924033	UHMWPE
0924034	UHMWPE
0924040	UHMWPE
0924041	UHMWPE
0924042	UHMWPE
0924043	UHMWPE
0924044	UHMWPE
0924050	UHMWPE
0924051	UHMWPE
0924052	UHMWPE
0924053	UHMWPE
0924054	UHMWPE
0942200	UHMWPE
0942201	UHMWPE
0942202	UHMWPE
0942203	UHMWPE
0942204	UHMWPE
0942205	UHMWPE
0942230	UHMWPE
0942231	UHMWPE
0942232	UHMWPE
0942233	UHMWPE

0942234	UHMWPE
0942235	UHMWPE
0942800	UHMWPE
0942801	UHMWPE
0942802	UHMWPE
0942803	UHMWPE
0942804	UHMWPE
0942805	UHMWPE
0942811	UHMWPE
0942813	UHMWPE
0942814	UHMWPE
0942815	UHMWPE
0942821	UHMWPE
0942823	UHMWPE
0942824	UHMWPE
0942825	UHMWPE
0942830	UHMWPE
0942831	UHMWPE
0942832	UHMWPE
0942833	UHMWPE
0942834	UHMWPE
0942835	UHMWPE
0943223	UHMWPE
0943224	UHMWPE
0943225	UHMWPE
0943233	UHMWPE
0943234	UHMWPE
0943235	UHMWPE
0943624	UHMWPE
0943625	UHMWPE
0943634	UHMWPE
0943635	UHMWPE
1101001	Co-Cr-Mo Alloy + Hydroxyapatite
1101002	Co-Cr-Mo Alloy + Hydroxyapatite
1101003	Co-Cr-Mo Alloy + Hydroxyapatite
1101004	Co-Cr-Mo Alloy + Hydroxyapatite
1101005	Co-Cr-Mo Alloy + Hydroxyapatite
1101006	Co-Cr-Mo Alloy + Hydroxyapatite
1101007	Co-Cr-Mo Alloy + Hydroxyapatite
1101011	Co-Cr-Mo Alloy + Hydroxyapatite
1101012	Co-Cr-Mo Alloy + Hydroxyapatite
1101013	Co-Cr-Mo Alloy + Hydroxyapatite
1101014	Co-Cr-Mo Alloy + Hydroxyapatite

1101015	Co-Cr-Mo Alloy + Hydroxyapatite
1101016	Co-Cr-Mo Alloy + Hydroxyapatite
1101017	Co-Cr-Mo Alloy + Hydroxyapatite
1101022	Co-Cr-Mo Alloy + Hydroxyapatite
1101023	Co-Cr-Mo Alloy + Hydroxyapatite
1101024	Co-Cr-Mo Alloy + Hydroxyapatite
1101025	Co-Cr-Mo Alloy + Hydroxyapatite
1101026	Co-Cr-Mo Alloy + Hydroxyapatite
1101027	Co-Cr-Mo Alloy + Hydroxyapatite
1101032	Co-Cr-Mo Alloy + Hydroxyapatite
1101033	Co-Cr-Mo Alloy + Hydroxyapatite
1101034	Co-Cr-Mo Alloy + Hydroxyapatite
1101035	Co-Cr-Mo Alloy + Hydroxyapatite
1101036	Co-Cr-Mo Alloy + Hydroxyapatite
1101037	Co-Cr-Mo Alloy + Hydroxyapatite
1102001	Co-Cr-Mo Alloy
1102002	Co-Cr-Mo Alloy
1102003	Co-Cr-Mo Alloy
1102004	Co-Cr-Mo Alloy
1102005	Co-Cr-Mo Alloy
1102006	Co-Cr-Mo Alloy
1102007	Co-Cr-Mo Alloy
1102011	Co-Cr-Mo Alloy
1102012	Co-Cr-Mo Alloy
1102013	Co-Cr-Mo Alloy
1102014	Co-Cr-Mo Alloy
1102015	Co-Cr-Mo Alloy
1102016	Co-Cr-Mo Alloy
1102017	Co-Cr-Mo Alloy
1102022	Co-Cr-Mo Alloy
1102023	Co-Cr-Mo Alloy
1102024	Co-Cr-Mo Alloy
1102025	Co-Cr-Mo Alloy
1102026	Co-Cr-Mo Alloy
1102027	Co-Cr-Mo Alloy
1102032	Co-Cr-Mo Alloy
1102033	Co-Cr-Mo Alloy
1102034	Co-Cr-Mo Alloy
1102035	Co-Cr-Mo Alloy
1102036	Co-Cr-Mo Alloy
1102037	Co-Cr-Mo Alloy
1102101	Co-Cr-Mo Alloy
1102102	Co-Cr-Mo Alloy

1102103	Co-Cr-Mo Alloy
1102104	Co-Cr-Mo Alloy
1102105	Co-Cr-Mo Alloy
1102106	Co-Cr-Mo Alloy
1102107	Co-Cr-Mo Alloy
1102111	Co-Cr-Mo Alloy
1102112	Co-Cr-Mo Alloy
1102113	Co-Cr-Mo Alloy
1102114	Co-Cr-Mo Alloy
1102115	Co-Cr-Mo Alloy
1102116	Co-Cr-Mo Alloy
1102117	Co-Cr-Mo Alloy
1103001	Co-Cr-Mo Alloy
1103002	Co-Cr-Mo Alloy
1103003	Co-Cr-Mo Alloy
1103004	Co-Cr-Mo Alloy
1103005	Co-Cr-Mo Alloy
1103006	Co-Cr-Mo Alloy
1103007	Co-Cr-Mo Alloy
1103011	Co-Cr-Mo Alloy
1103012	Co-Cr-Mo Alloy
1103013	Co-Cr-Mo Alloy
1103014	Co-Cr-Mo Alloy
1103015	Co-Cr-Mo Alloy
1103016	Co-Cr-Mo Alloy
1103017	Co-Cr-Mo Alloy
1103022	Co-Cr-Mo Alloy
1103023	Co-Cr-Mo Alloy
1103024	Co-Cr-Mo Alloy
1103025	Co-Cr-Mo Alloy
1103026	Co-Cr-Mo Alloy
1103027	Co-Cr-Mo Alloy
1103032	Co-Cr-Mo Alloy
1103033	Co-Cr-Mo Alloy
1103034	Co-Cr-Mo Alloy
1103035	Co-Cr-Mo Alloy
1103036	Co-Cr-Mo Alloy
1103037	Co-Cr-Mo Alloy
1104001	Co-Cr-Mo Alloy + Hydroxyapatite
1104002	Co-Cr-Mo Alloy + Hydroxyapatite
1104003	Co-Cr-Mo Alloy + Hydroxyapatite
1104004	Co-Cr-Mo Alloy + Hydroxyapatite
1104005	Co-Cr-Mo Alloy + Hydroxyapatite

1104006	Co-Cr-Mo Alloy + Hydroxyapatite
1104007	Co-Cr-Mo Alloy + Hydroxyapatite
1104011	Co-Cr-Mo Alloy + Hydroxyapatite
1104012	Co-Cr-Mo Alloy + Hydroxyapatite
1104013	Co-Cr-Mo Alloy + Hydroxyapatite
1104014	Co-Cr-Mo Alloy + Hydroxyapatite
1104015	Co-Cr-Mo Alloy + Hydroxyapatite
1104016	Co-Cr-Mo Alloy + Hydroxyapatite
1104017	Co-Cr-Mo Alloy + Hydroxyapatite
1104022	Co-Cr-Mo Alloy + Hydroxyapatite
1104023	Co-Cr-Mo Alloy + Hydroxyapatite
1104024	Co-Cr-Mo Alloy + Hydroxyapatite
1104025	Co-Cr-Mo Alloy + Hydroxyapatite
1104026	Co-Cr-Mo Alloy + Hydroxyapatite
1104027	Co-Cr-Mo Alloy + Hydroxyapatite
1104032	Co-Cr-Mo Alloy + Hydroxyapatite
1104033	Co-Cr-Mo Alloy + Hydroxyapatite
1104034	Co-Cr-Mo Alloy + Hydroxyapatite
1104035	Co-Cr-Mo Alloy + Hydroxyapatite
1104036	Co-Cr-Mo Alloy + Hydroxyapatite
1104037	Co-Cr-Mo Alloy + Hydroxyapatite
1104062	Co-Cr-Mo Alloy + Hydroxyapatite
1104063	Co-Cr-Mo Alloy + Hydroxyapatite
1104064	Co-Cr-Mo Alloy + Hydroxyapatite
1104065	Co-Cr-Mo Alloy + Hydroxyapatite
1104066	Co-Cr-Mo Alloy + Hydroxyapatite
1104071	Co-Cr-Mo Alloy + Hydroxyapatite
1104072	Co-Cr-Mo Alloy + Hydroxyapatite
1104073	Co-Cr-Mo Alloy + Hydroxyapatite
1104074	Co-Cr-Mo Alloy + Hydroxyapatite
1104075	Co-Cr-Mo Alloy + Hydroxyapatite
1104076	Co-Cr-Mo Alloy + Hydroxyapatite
1104081	Co-Cr-Mo Alloy + Hydroxyapatite
1104082	Co-Cr-Mo Alloy + Hydroxyapatite
1104083	Co-Cr-Mo Alloy + Hydroxyapatite
1104084	Co-Cr-Mo Alloy + Hydroxyapatite
1104085	Co-Cr-Mo Alloy + Hydroxyapatite
1104086	Co-Cr-Mo Alloy + Hydroxyapatite
1104092	Co-Cr-Mo Alloy + Hydroxyapatite
1104093	Co-Cr-Mo Alloy + Hydroxyapatite
1104094	Co-Cr-Mo Alloy + Hydroxyapatite
1104095	Co-Cr-Mo Alloy + Hydroxyapatite
1104096	Co-Cr-Mo Alloy + Hydroxyapatite

1104262	Co-Cr-Mo Alloy + Hydroxyapatite
1104263	Co-Cr-Mo Alloy + Hydroxyapatite
1104264	Co-Cr-Mo Alloy + Hydroxyapatite
1104265	Co-Cr-Mo Alloy + Hydroxyapatite
1104266	Co-Cr-Mo Alloy + Hydroxyapatite
1104267	Co-Cr-Mo Alloy + Hydroxyapatite
1104272	Co-Cr-Mo Alloy + Hydroxyapatite
1104273	Co-Cr-Mo Alloy + Hydroxyapatite
1104274	Co-Cr-Mo Alloy + Hydroxyapatite
1104275	Co-Cr-Mo Alloy + Hydroxyapatite
1104276	Co-Cr-Mo Alloy + Hydroxyapatite
1104277	Co-Cr-Mo Alloy + Hydroxyapatite
1105001	Co-Cr-Mo Alloy + Hydroxyapatite
1105002	Co-Cr-Mo Alloy + Hydroxyapatite
1105003	Co-Cr-Mo Alloy + Hydroxyapatite
1105004	Co-Cr-Mo Alloy + Hydroxyapatite
1105005	Co-Cr-Mo Alloy + Hydroxyapatite
1105006	Co-Cr-Mo Alloy + Hydroxyapatite
1105007	Co-Cr-Mo Alloy + Hydroxyapatite
1105011	Co-Cr-Mo Alloy + Hydroxyapatite
1105012	Co-Cr-Mo Alloy + Hydroxyapatite
1105013	Co-Cr-Mo Alloy + Hydroxyapatite
1105014	Co-Cr-Mo Alloy + Hydroxyapatite
1105015	Co-Cr-Mo Alloy + Hydroxyapatite
1105016	Co-Cr-Mo Alloy + Hydroxyapatite
1105017	Co-Cr-Mo Alloy + Hydroxyapatite
1105022	Co-Cr-Mo Alloy + Hydroxyapatite
1105023	Co-Cr-Mo Alloy + Hydroxyapatite
1105024	Co-Cr-Mo Alloy + Hydroxyapatite
1105025	Co-Cr-Mo Alloy + Hydroxyapatite
1105026	Co-Cr-Mo Alloy + Hydroxyapatite
1105027	Co-Cr-Mo Alloy + Hydroxyapatite
1105032	Co-Cr-Mo Alloy + Hydroxyapatite
1105033	Co-Cr-Mo Alloy + Hydroxyapatite
1105034	Co-Cr-Mo Alloy + Hydroxyapatite
1105035	Co-Cr-Mo Alloy + Hydroxyapatite
1105036	Co-Cr-Mo Alloy + Hydroxyapatite
1105037	Co-Cr-Mo Alloy + Hydroxyapatite
1105062	Co-Cr-Mo Alloy + Hydroxyapatite
1105063	Co-Cr-Mo Alloy + Hydroxyapatite
1105064	Co-Cr-Mo Alloy + Hydroxyapatite
1105065	Co-Cr-Mo Alloy + Hydroxyapatite
1105066	Co-Cr-Mo Alloy + Hydroxyapatite

1105071	Co-Cr-Mo Alloy + Hydroxyapatite
1105072	Co-Cr-Mo Alloy + Hydroxyapatite
1105073	Co-Cr-Mo Alloy + Hydroxyapatite
1105074	Co-Cr-Mo Alloy + Hydroxyapatite
1105075	Co-Cr-Mo Alloy + Hydroxyapatite
1105076	Co-Cr-Mo Alloy + Hydroxyapatite
1105081	Co-Cr-Mo Alloy + Hydroxyapatite
1105082	Co-Cr-Mo Alloy + Hydroxyapatite
1105083	Co-Cr-Mo Alloy + Hydroxyapatite
1105084	Co-Cr-Mo Alloy + Hydroxyapatite
1105085	Co-Cr-Mo Alloy + Hydroxyapatite
1105086	Co-Cr-Mo Alloy + Hydroxyapatite
1105092	Co-Cr-Mo Alloy + Hydroxyapatite
1105093	Co-Cr-Mo Alloy + Hydroxyapatite
1105094	Co-Cr-Mo Alloy + Hydroxyapatite
1105095	Co-Cr-Mo Alloy + Hydroxyapatite
1105096	Co-Cr-Mo Alloy + Hydroxyapatite
1105101	Co-Cr-Mo Alloy + Hydroxyapatite
1105102	Co-Cr-Mo Alloy + Hydroxyapatite
1105103	Co-Cr-Mo Alloy + Hydroxyapatite
1105104	Co-Cr-Mo Alloy + Hydroxyapatite
1105105	Co-Cr-Mo Alloy + Hydroxyapatite
1105106	Co-Cr-Mo Alloy + Hydroxyapatite
1105107	Co-Cr-Mo Alloy + Hydroxyapatite
1105111	Co-Cr-Mo Alloy + Hydroxyapatite
1105112	Co-Cr-Mo Alloy + Hydroxyapatite
1105113	Co-Cr-Mo Alloy + Hydroxyapatite
1105114	Co-Cr-Mo Alloy + Hydroxyapatite
1105115	Co-Cr-Mo Alloy + Hydroxyapatite
1105116	Co-Cr-Mo Alloy + Hydroxyapatite
1105117	Co-Cr-Mo Alloy + Hydroxyapatite
1105262	Co-Cr-Mo Alloy + Hydroxyapatite
1105263	Co-Cr-Mo Alloy + Hydroxyapatite
1105264	Co-Cr-Mo Alloy + Hydroxyapatite
1105265	Co-Cr-Mo Alloy + Hydroxyapatite
1105266	Co-Cr-Mo Alloy + Hydroxyapatite
1105267	Co-Cr-Mo Alloy + Hydroxyapatite
1105272	Co-Cr-Mo Alloy + Hydroxyapatite
1105273	Co-Cr-Mo Alloy + Hydroxyapatite
1105274	Co-Cr-Mo Alloy + Hydroxyapatite
1105275	Co-Cr-Mo Alloy + Hydroxyapatite
1105276	Co-Cr-Mo Alloy + Hydroxyapatite
1105277	Co-Cr-Mo Alloy + Hydroxyapatite

1107001	Ti6Al4V Alloy + TiNbN
1107002	Ti6Al4V Alloy + TiNbN
1107003	Ti6Al4V Alloy + TiNbN
1107004	Ti6Al4V Alloy + TiNbN
1107005	Ti6Al4V Alloy + TiNbN
1107006	Ti6Al4V Alloy + TiNbN
1107007	Ti6Al4V Alloy + TiNbN
1107011	Ti6Al4V Alloy + TiNbN
1107012	Ti6Al4V Alloy + TiNbN
1107013	Ti6Al4V Alloy + TiNbN
1107014	Ti6Al4V Alloy + TiNbN
1107015	Ti6Al4V Alloy + TiNbN
1107016	Ti6Al4V Alloy + TiNbN
1107017	Ti6Al4V Alloy + TiNbN
1107051	Ti6Al4V Alloy + TiNbN
1107052	Ti6Al4V Alloy + TiNbN
1107053	Ti6Al4V Alloy + TiNbN
1107054	Ti6Al4V Alloy + TiNbN
1107055	Ti6Al4V Alloy + TiNbN
1107056	Ti6Al4V Alloy + TiNbN
1107061	Ti6Al4V Alloy + TiNbN
1107062	Ti6Al4V Alloy + TiNbN
1107063	Ti6Al4V Alloy + TiNbN
1107064	Ti6Al4V Alloy + TiNbN
1107065	Ti6Al4V Alloy + TiNbN
1107066	Ti6Al4V Alloy + TiNbN
1107101	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107102	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107103	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107104	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107105	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107106	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107107	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107111	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107112	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107113	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107114	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107115	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107116	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107117	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107151	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107152	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107153	Ti6Al4V Alloy + TiNbN + Hydroxyapatite

1107154	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107155	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107156	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107161	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107162	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107163	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107164	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107165	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107166	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107202	Ti6Al4V Alloy + TiNbN
1107203	Ti6Al4V Alloy + TiNbN
1107204	Ti6Al4V Alloy + TiNbN
1107205	Ti6Al4V Alloy + TiNbN
1107206	Ti6Al4V Alloy + TiNbN
1107207	Ti6Al4V Alloy + TiNbN
1107212	Ti6Al4V Alloy + TiNbN
1107213	Ti6Al4V Alloy + TiNbN
1107214	Ti6Al4V Alloy + TiNbN
1107215	Ti6Al4V Alloy + TiNbN
1107216	Ti6Al4V Alloy + TiNbN
1107217	Ti6Al4V Alloy + TiNbN
1107252	Ti6Al4V Alloy + TiNbN
1107253	Ti6Al4V Alloy + TiNbN
1107254	Ti6Al4V Alloy + TiNbN
1107255	Ti6Al4V Alloy + TiNbN
1107256	Ti6Al4V Alloy + TiNbN
1107262	Ti6Al4V Alloy + TiNbN
1107263	Ti6Al4V Alloy + TiNbN
1107264	Ti6Al4V Alloy + TiNbN
1107265	Ti6Al4V Alloy + TiNbN
1107266	Ti6Al4V Alloy + TiNbN
1107302	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107303	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107304	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107305	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107306	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107307	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107312	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107313	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107314	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107315	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107316	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107317	Ti6Al4V Alloy + TiNbN + Hydroxyapatite

1107352	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107353	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107354	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107355	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107356	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107362	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107363	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107364	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107365	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107366	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107402	Ti6Al4V Alloy + TiNbN
1107403	Ti6Al4V Alloy + TiNbN
1107404	Ti6Al4V Alloy + TiNbN
1107405	Ti6Al4V Alloy + TiNbN
1107406	Ti6Al4V Alloy + TiNbN
1107407	Ti6Al4V Alloy + TiNbN
1107412	Ti6Al4V Alloy + TiNbN
1107413	Ti6Al4V Alloy + TiNbN
1107414	Ti6Al4V Alloy + TiNbN
1107415	Ti6Al4V Alloy + TiNbN
1107416	Ti6Al4V Alloy + TiNbN
1107417	Ti6Al4V Alloy + TiNbN
1107502	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107503	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107504	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107505	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107506	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107507	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107512	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107513	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107514	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107515	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107516	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1107517	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1111001	Co-Cr-Mo Alloy + Ti
1111002	Co-Cr-Mo Alloy + Ti
1111003	Co-Cr-Mo Alloy + Ti
1111004	Co-Cr-Mo Alloy + Ti
1111005	Co-Cr-Mo Alloy + Ti
1111006	Co-Cr-Mo Alloy + Ti
1111011	Co-Cr-Mo Alloy + Ti
1111012	Co-Cr-Mo Alloy + Ti
1111013	Co-Cr-Mo Alloy + Ti

1111014	Co-Cr-Mo Alloy + Ti
1111015	Co-Cr-Mo Alloy + Ti
1111016	Co-Cr-Mo Alloy + Ti
1112001	Co-Cr-Mo Alloy
1112002	Co-Cr-Mo Alloy
1112003	Co-Cr-Mo Alloy
1112004	Co-Cr-Mo Alloy
1112005	Co-Cr-Mo Alloy
1112006	Co-Cr-Mo Alloy
1112011	Co-Cr-Mo Alloy
1112012	Co-Cr-Mo Alloy
1112013	Co-Cr-Mo Alloy
1112014	Co-Cr-Mo Alloy
1112015	Co-Cr-Mo Alloy
1112016	Co-Cr-Mo Alloy
1113001	Co-Cr-Mo Alloy
1113002	Co-Cr-Mo Alloy
1113003	Co-Cr-Mo Alloy
1113004	Co-Cr-Mo Alloy
1113005	Co-Cr-Mo Alloy
1113006	Co-Cr-Mo Alloy
1113011	Co-Cr-Mo Alloy
1113012	Co-Cr-Mo Alloy
1113013	Co-Cr-Mo Alloy
1113014	Co-Cr-Mo Alloy
1113015	Co-Cr-Mo Alloy
1113016	Co-Cr-Mo Alloy
1114001	Co-Cr-Mo Alloy+ Hydroxyapatite
1114002	Co-Cr-Mo Alloy+ Hydroxyapatite
1114003	Co-Cr-Mo Alloy+ Hydroxyapatite
1114004	Co-Cr-Mo Alloy+ Hydroxyapatite
1114005	Co-Cr-Mo Alloy+ Hydroxyapatite
1114006	Co-Cr-Mo Alloy+ Hydroxyapatite
1114011	Co-Cr-Mo Alloy+ Hydroxyapatite
1114012	Co-Cr-Mo Alloy+ Hydroxyapatite
1114013	Co-Cr-Mo Alloy+ Hydroxyapatite
1114014	Co-Cr-Mo Alloy+ Hydroxyapatite
1114015	Co-Cr-Mo Alloy+ Hydroxyapatite
1114016	Co-Cr-Mo Alloy+ Hydroxyapatite
1115001	Co-Cr-Mo Alloy + Hydroxyapatite
1115002	Co-Cr-Mo Alloy + Hydroxyapatite
1115003	Co-Cr-Mo Alloy + Hydroxyapatite
1115004	Co-Cr-Mo Alloy + Hydroxyapatite

1115005	Co-Cr-Mo Alloy + Hydroxyapatite
1115006	Co-Cr-Mo Alloy + Hydroxyapatite
1115011	Co-Cr-Mo Alloy + Hydroxyapatite
1115012	Co-Cr-Mo Alloy + Hydroxyapatite
1115013	Co-Cr-Mo Alloy + Hydroxyapatite
1115014	Co-Cr-Mo Alloy + Hydroxyapatite
1115015	Co-Cr-Mo Alloy + Hydroxyapatite
1115016	Co-Cr-Mo Alloy + Hydroxyapatite
1116001	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116002	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116003	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116004	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116005	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116006	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116011	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116012	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116013	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116014	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116015	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1116016	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1120001	Co-Cr-Mo Alloy
1120002	Co-Cr-Mo Alloy
1120003	Co-Cr-Mo Alloy
1120004	Co-Cr-Mo Alloy
1120005	Co-Cr-Mo Alloy
1120006	Co-Cr-Mo Alloy
1120011	Co-Cr-Mo Alloy
1120012	Co-Cr-Mo Alloy
1120013	Co-Cr-Mo Alloy
1120014	Co-Cr-Mo Alloy
1120015	Co-Cr-Mo Alloy
1120016	Co-Cr-Mo Alloy
1120051	Co-Cr-Mo Alloy
1120052	Co-Cr-Mo Alloy
1120053	Co-Cr-Mo Alloy
1120054	Co-Cr-Mo Alloy
1120055	Co-Cr-Mo Alloy
1120061	Co-Cr-Mo Alloy
1120062	Co-Cr-Mo Alloy
1120063	Co-Cr-Mo Alloy
1120064	Co-Cr-Mo Alloy
1120065	Co-Cr-Mo Alloy
1121001	Co-Cr-Mo Alloy + Hydroxyapatite

1121002	Co-Cr-Mo Alloy + Hydroxyapatite
1121003	Co-Cr-Mo Alloy + Hydroxyapatite
1121004	Co-Cr-Mo Alloy + Hydroxyapatite
1121005	Co-Cr-Mo Alloy + Hydroxyapatite
1121006	Co-Cr-Mo Alloy + Hydroxyapatite
1121011	Co-Cr-Mo Alloy + Hydroxyapatite
1121012	Co-Cr-Mo Alloy + Hydroxyapatite
1121013	Co-Cr-Mo Alloy + Hydroxyapatite
1121014	Co-Cr-Mo Alloy + Hydroxyapatite
1121015	Co-Cr-Mo Alloy + Hydroxyapatite
1121016	Co-Cr-Mo Alloy + Hydroxyapatite
1121051	Co-Cr-Mo Alloy + Hydroxyapatite
1121052	Co-Cr-Mo Alloy + Hydroxyapatite
1121053	Co-Cr-Mo Alloy + Hydroxyapatite
1121054	Co-Cr-Mo Alloy + Hydroxyapatite
1121055	Co-Cr-Mo Alloy + Hydroxyapatite
1121061	Co-Cr-Mo Alloy + Hydroxyapatite
1121062	Co-Cr-Mo Alloy + Hydroxyapatite
1121063	Co-Cr-Mo Alloy + Hydroxyapatite
1121064	Co-Cr-Mo Alloy + Hydroxyapatite
1121065	Co-Cr-Mo Alloy + Hydroxyapatite
1122001	Co-Cr-Mo Alloy
1122002	Co-Cr-Mo Alloy
1122003	Co-Cr-Mo Alloy
1122004	Co-Cr-Mo Alloy
1122005	Co-Cr-Mo Alloy
1122006	Co-Cr-Mo Alloy
1122011	Co-Cr-Mo Alloy
1122012	Co-Cr-Mo Alloy
1122013	Co-Cr-Mo Alloy
1122014	Co-Cr-Mo Alloy
1122015	Co-Cr-Mo Alloy
1122016	Co-Cr-Mo Alloy
1122051	Co-Cr-Mo Alloy
1122052	Co-Cr-Mo Alloy
1122053	Co-Cr-Mo Alloy
1122054	Co-Cr-Mo Alloy
1122055	Co-Cr-Mo Alloy
1122061	Co-Cr-Mo Alloy
1122062	Co-Cr-Mo Alloy
1122063	Co-Cr-Mo Alloy
1122064	Co-Cr-Mo Alloy
1122065	Co-Cr-Mo Alloy

1123001	Co-Cr-Mo Alloy + Hydroxyapatite
1123002	Co-Cr-Mo Alloy + Hydroxyapatite
1123003	Co-Cr-Mo Alloy + Hydroxyapatite
1123004	Co-Cr-Mo Alloy + Hydroxyapatite
1123005	Co-Cr-Mo Alloy + Hydroxyapatite
1123006	Co-Cr-Mo Alloy + Hydroxyapatite
1123011	Co-Cr-Mo Alloy + Hydroxyapatite
1123012	Co-Cr-Mo Alloy + Hydroxyapatite
1123013	Co-Cr-Mo Alloy + Hydroxyapatite
1123014	Co-Cr-Mo Alloy + Hydroxyapatite
1123015	Co-Cr-Mo Alloy + Hydroxyapatite
1123016	Co-Cr-Mo Alloy + Hydroxyapatite
1123051	Co-Cr-Mo Alloy + Hydroxyapatite
1123052	Co-Cr-Mo Alloy + Hydroxyapatite
1123053	Co-Cr-Mo Alloy + Hydroxyapatite
1123054	Co-Cr-Mo Alloy + Hydroxyapatite
1123055	Co-Cr-Mo Alloy + Hydroxyapatite
1123061	Co-Cr-Mo Alloy + Hydroxyapatite
1123062	Co-Cr-Mo Alloy + Hydroxyapatite
1123063	Co-Cr-Mo Alloy + Hydroxyapatite
1123064	Co-Cr-Mo Alloy + Hydroxyapatite
1123065	Co-Cr-Mo Alloy + Hydroxyapatite
1124001	Ti6Al4V Alloy
1124002	Ti6Al4V Alloy
1124003	Ti6Al4V Alloy
1124004	Ti6Al4V Alloy
1124005	Ti6Al4V Alloy
1124006	Ti6Al4V Alloy
1124011	Ti6Al4V Alloy
1124012	Ti6Al4V Alloy
1124013	Ti6Al4V Alloy
1124014	Ti6Al4V Alloy
1124015	Ti6Al4V Alloy
1124016	Ti6Al4V Alloy
1124051	Ti6Al4V Alloy
1124053	Ti6Al4V Alloy
1124055	Ti6Al4V Alloy
1124061	Ti6Al4V Alloy
1124063	Ti6Al4V Alloy
1124065	Ti6Al4V Alloy
1124101	Ti6Al4V Alloy
1124102	Ti6Al4V Alloy
1124103	Ti6Al4V Alloy

1124104	Ti6Al4V Alloy
1124105	Ti6Al4V Alloy
1124106	Ti6Al4V Alloy
1124111	Ti6Al4V Alloy
1124112	Ti6Al4V Alloy
1124113	Ti6Al4V Alloy
1124114	Ti6Al4V Alloy
1124115	Ti6Al4V Alloy
1124116	Ti6Al4V Alloy
1124151	Ti6Al4V Alloy
1124153	Ti6Al4V Alloy
1124155	Ti6Al4V Alloy
1124161	Ti6Al4V Alloy
1124163	Ti6Al4V Alloy
1124165	Ti6Al4V Alloy
1125001	Ti6Al4V Alloy
1125002	Ti6Al4V Alloy
1125003	Ti6Al4V Alloy
1125004	Ti6Al4V Alloy
1125005	Ti6Al4V Alloy
1125006	Ti6Al4V Alloy
1125011	Ti6Al4V Alloy
1125012	Ti6Al4V Alloy
1125013	Ti6Al4V Alloy
1125014	Ti6Al4V Alloy
1125015	Ti6Al4V Alloy
1125016	Ti6Al4V Alloy
1125051	Ti6Al4V Alloy
1125053	Ti6Al4V Alloy
1125055	Ti6Al4V Alloy
1125061	Ti6Al4V Alloy
1125063	Ti6Al4V Alloy
1125065	Ti6Al4V Alloy
1125101	Ti6Al4V Alloy
1125102	Ti6Al4V Alloy
1125103	Ti6Al4V Alloy
1125104	Ti6Al4V Alloy
1125105	Ti6Al4V Alloy
1125106	Ti6Al4V Alloy
1125111	Ti6Al4V Alloy
1125112	Ti6Al4V Alloy
1125113	Ti6Al4V Alloy
1125114	Ti6Al4V Alloy

1125115	Ti6Al4V Alloy
1125116	Ti6Al4V Alloy
1125151	Ti6Al4V Alloy
1125153	Ti6Al4V Alloy
1125155	Ti6Al4V Alloy
1125161	Ti6Al4V Alloy
1125163	Ti6Al4V Alloy
1125165	Ti6Al4V Alloy
1125201	Ti6Al4V Alloy
1125202	Ti6Al4V Alloy
1125203	Ti6Al4V Alloy
1125204	Ti6Al4V Alloy
1125205	Ti6Al4V Alloy
1125206	Ti6Al4V Alloy
1125211	Ti6Al4V Alloy
1125212	Ti6Al4V Alloy
1125213	Ti6Al4V Alloy
1125214	Ti6Al4V Alloy
1125215	Ti6Al4V Alloy
1125216	Ti6Al4V Alloy
1125251	Ti6Al4V Alloy
1125253	Ti6Al4V Alloy
1125255	Ti6Al4V Alloy
1125261	Ti6Al4V Alloy
1125263	Ti6Al4V Alloy
1125265	Ti6Al4V Alloy
1125301	Ti6Al4V Alloy
1125302	Ti6Al4V Alloy
1125303	Ti6Al4V Alloy
1125304	Ti6Al4V Alloy
1125305	Ti6Al4V Alloy
1125306	Ti6Al4V Alloy
1125311	Ti6Al4V Alloy
1125312	Ti6Al4V Alloy
1125313	Ti6Al4V Alloy
1125314	Ti6Al4V Alloy
1125315	Ti6Al4V Alloy
1125316	Ti6Al4V Alloy
1125351	Ti6Al4V Alloy
1125353	Ti6Al4V Alloy
1125355	Ti6Al4V Alloy
1125361	Ti6Al4V Alloy
1125363	Ti6Al4V Alloy

1125365	Ti6Al4V Alloy
1125401	Ti6Al4V Alloy
1125402	Ti6Al4V Alloy
1125403	Ti6Al4V Alloy
1125404	Ti6Al4V Alloy
1125405	Ti6Al4V Alloy
1125406	Ti6Al4V Alloy
1125411	Ti6Al4V Alloy
1125412	Ti6Al4V Alloy
1125413	Ti6Al4V Alloy
1125414	Ti6Al4V Alloy
1125415	Ti6Al4V Alloy
1125416	Ti6Al4V Alloy
1125451	Ti6Al4V Alloy
1125453	Ti6Al4V Alloy
1125455	Ti6Al4V Alloy
1125461	Ti6Al4V Alloy
1125463	Ti6Al4V Alloy
1125465	Ti6Al4V Alloy
1125501	Ti6Al4V Alloy
1125502	Ti6Al4V Alloy
1125503	Ti6Al4V Alloy
1125504	Ti6Al4V Alloy
1125505	Ti6Al4V Alloy
1125506	Ti6Al4V Alloy
1125511	Ti6Al4V Alloy
1125512	Ti6Al4V Alloy
1125513	Ti6Al4V Alloy
1125514	Ti6Al4V Alloy
1125515	Ti6Al4V Alloy
1125516	Ti6Al4V Alloy
1125551	Ti6Al4V Alloy
1125553	Ti6Al4V Alloy
1125555	Ti6Al4V Alloy
1125561	Ti6Al4V Alloy
1125563	Ti6Al4V Alloy
1125565	Ti6Al4V Alloy
1125905	Ti6Al4V Alloy
1125910	Ti6Al4V Alloy
1126110	UHMWPE
1126112	UHMWPE
1126114	UHMWPE
1126116	UHMWPE

1126117	UHMWPE
1126118	UHMWPE
1126120	UHMWPE
1126122	UHMWPE
1126124	UHMWPE
1126126	UHMWPE
1126128	UHMWPE
1126130	UHMWPE
1126140	UHMWPE
1126142	UHMWPE
1126144	UHMWPE
1126146	UHMWPE
1126147	UHMWPE
1126148	UHMWPE
1126150	UHMWPE
1126152	UHMWPE
1126154	UHMWPE
1126156	UHMWPE
1126158	UHMWPE
1126160	UHMWPE
1126170	UHMWPE
1126172	UHMWPE
1126174	UHMWPE
1126176	UHMWPE
1126177	UHMWPE
1126178	UHMWPE
1126180	UHMWPE
1126182	UHMWPE
1126184	UHMWPE
1126186	UHMWPE
1126188	UHMWPE
1126190	UHMWPE
1126210	UHMWPE
1126212	UHMWPE
1126214	UHMWPE
1126216	UHMWPE
1126217	UHMWPE
1126218	UHMWPE
1126220	UHMWPE
1126222	UHMWPE
1126224	UHMWPE
1126226	UHMWPE
1126228	UHMWPE

1126230	UHMWPE
1126240	UHMWPE
1126242	UHMWPE
1126244	UHMWPE
1126246	UHMWPE
1126247	UHMWPE
1126248	UHMWPE
1126250	UHMWPE
1126252	UHMWPE
1126254	UHMWPE
1126256	UHMWPE
1126258	UHMWPE
1126260	UHMWPE
1126270	UHMWPE
1126272	UHMWPE
1126274	UHMWPE
1126276	UHMWPE
1126277	UHMWPE
1126278	UHMWPE
1126280	UHMWPE
1126282	UHMWPE
1126284	UHMWPE
1126286	UHMWPE
1126288	UHMWPE
1126290	UHMWPE
1126310	UHMWPE
1126312	UHMWPE
1126314	UHMWPE
1126316	UHMWPE
1126317	UHMWPE
1126318	UHMWPE
1126320	UHMWPE
1126322	UHMWPE
1126324	UHMWPE
1126326	UHMWPE
1126328	UHMWPE
1126330	UHMWPE
1126340	UHMWPE
1126342	UHMWPE
1126344	UHMWPE
1126346	UHMWPE
1126347	UHMWPE
1126348	UHMWPE

1126350	UHMWPE
1126352	UHMWPE
1126354	UHMWPE
1126356	UHMWPE
1126358	UHMWPE
1126360	UHMWPE
1126370	UHMWPE
1126372	UHMWPE
1126374	UHMWPE
1126376	UHMWPE
1126377	UHMWPE
1126378	UHMWPE
1126380	UHMWPE
1126382	UHMWPE
1126384	UHMWPE
1126386	UHMWPE
1126388	UHMWPE
1126390	UHMWPE
1126410	UHMWPE
1126412	UHMWPE
1126414	UHMWPE
1126416	UHMWPE
1126417	UHMWPE
1126418	UHMWPE
1126420	UHMWPE
1126422	UHMWPE
1126424	UHMWPE
1126426	UHMWPE
1126428	UHMWPE
1126430	UHMWPE
1126440	UHMWPE
1126442	UHMWPE
1126444	UHMWPE
1126446	UHMWPE
1126447	UHMWPE
1126448	UHMWPE
1126450	UHMWPE
1126452	UHMWPE
1126454	UHMWPE
1126456	UHMWPE
1126458	UHMWPE
1126460	UHMWPE
1126470	UHMWPE

1126472	UHMWPE
1126474	UHMWPE
1126476	UHMWPE
1126477	UHMWPE
1126478	UHMWPE
1126480	UHMWPE
1126482	UHMWPE
1126484	UHMWPE
1126486	UHMWPE
1126488	UHMWPE
1126490	UHMWPE
1126510	UHMWPE
1126512	UHMWPE
1126514	UHMWPE
1126516	UHMWPE
1126517	UHMWPE
1126518	UHMWPE
1126520	UHMWPE
1126522	UHMWPE
1126524	UHMWPE
1126526	UHMWPE
1126528	UHMWPE
1126530	UHMWPE
1126540	UHMWPE
1126542	UHMWPE
1126544	UHMWPE
1126546	UHMWPE
1126547	UHMWPE
1126548	UHMWPE
1126550	UHMWPE
1126552	UHMWPE
1126554	UHMWPE
1126556	UHMWPE
1126558	UHMWPE
1126560	UHMWPE
1126570	UHMWPE
1126572	UHMWPE
1126574	UHMWPE
1126576	UHMWPE
1126577	UHMWPE
1126578	UHMWPE
1126580	UHMWPE
1126582	UHMWPE

1126584	UHMWPE
1126586	UHMWPE
1126588	UHMWPE
1126590	UHMWPE
1126610	UHMWPE
1126612	UHMWPE
1126614	UHMWPE
1126616	UHMWPE
1126617	UHMWPE
1126618	UHMWPE
1126620	UHMWPE
1126622	UHMWPE
1126624	UHMWPE
1126626	UHMWPE
1126628	UHMWPE
1126630	UHMWPE
1126640	UHMWPE
1126642	UHMWPE
1126644	UHMWPE
1126646	UHMWPE
1126647	UHMWPE
1126648	UHMWPE
1126650	UHMWPE
1126652	UHMWPE
1126654	UHMWPE
1126656	UHMWPE
1126658	UHMWPE
1126660	UHMWPE
1126670	UHMWPE
1126672	UHMWPE
1126674	UHMWPE
1126676	UHMWPE
1126677	UHMWPE
1126678	UHMWPE
1126680	UHMWPE
1126682	UHMWPE
1126684	UHMWPE
1126686	UHMWPE
1126688	UHMWPE
1126690	UHMWPE
1127110	UHMWPE
1127112	UHMWPE
1127114	UHMWPE

1127116	UHMWPE
1127117	UHMWPE
1127118	UHMWPE
1127120	UHMWPE
1127122	UHMWPE
1127124	UHMWPE
1127126	UHMWPE
1127128	UHMWPE
1127130	UHMWPE
1127140	UHMWPE
1127142	UHMWPE
1127144	UHMWPE
1127146	UHMWPE
1127147	UHMWPE
1127148	UHMWPE
1127150	UHMWPE
1127152	UHMWPE
1127154	UHMWPE
1127156	UHMWPE
1127158	UHMWPE
1127160	UHMWPE
1127170	UHMWPE
1127172	UHMWPE
1127174	UHMWPE
1127176	UHMWPE
1127177	UHMWPE
1127178	UHMWPE
1127180	UHMWPE
1127182	UHMWPE
1127184	UHMWPE
1127186	UHMWPE
1127188	UHMWPE
1127190	UHMWPE
1127210	UHMWPE
1127212	UHMWPE
1127214	UHMWPE
1127216	UHMWPE
1127217	UHMWPE
1127218	UHMWPE
1127220	UHMWPE
1127222	UHMWPE
1127224	UHMWPE
1127226	UHMWPE

1127228	UHMWPE
1127230	UHMWPE
1127240	UHMWPE
1127242	UHMWPE
1127244	UHMWPE
1127246	UHMWPE
1127247	UHMWPE
1127248	UHMWPE
1127250	UHMWPE
1127252	UHMWPE
1127254	UHMWPE
1127256	UHMWPE
1127258	UHMWPE
1127260	UHMWPE
1127270	UHMWPE
1127272	UHMWPE
1127274	UHMWPE
1127276	UHMWPE
1127277	UHMWPE
1127278	UHMWPE
1127280	UHMWPE
1127282	UHMWPE
1127284	UHMWPE
1127286	UHMWPE
1127288	UHMWPE
1127290	UHMWPE
1127310	UHMWPE
1127312	UHMWPE
1127314	UHMWPE
1127316	UHMWPE
1127317	UHMWPE
1127318	UHMWPE
1127320	UHMWPE
1127322	UHMWPE
1127324	UHMWPE
1127326	UHMWPE
1127328	UHMWPE
1127330	UHMWPE
1127340	UHMWPE
1127342	UHMWPE
1127344	UHMWPE
1127346	UHMWPE
1127347	UHMWPE

1127348	UHMWPE
1127350	UHMWPE
1127352	UHMWPE
1127354	UHMWPE
1127356	UHMWPE
1127358	UHMWPE
1127360	UHMWPE
1127370	UHMWPE
1127372	UHMWPE
1127374	UHMWPE
1127376	UHMWPE
1127377	UHMWPE
1127378	UHMWPE
1127380	UHMWPE
1127382	UHMWPE
1127384	UHMWPE
1127386	UHMWPE
1127388	UHMWPE
1127390	UHMWPE
1127410	UHMWPE
1127412	UHMWPE
1127414	UHMWPE
1127416	UHMWPE
1127417	UHMWPE
1127418	UHMWPE
1127420	UHMWPE
1127422	UHMWPE
1127424	UHMWPE
1127426	UHMWPE
1127428	UHMWPE
1127430	UHMWPE
1127440	UHMWPE
1127442	UHMWPE
1127444	UHMWPE
1127446	UHMWPE
1127447	UHMWPE
1127448	UHMWPE
1127450	UHMWPE
1127452	UHMWPE
1127454	UHMWPE
1127456	UHMWPE
1127458	UHMWPE
1127460	UHMWPE

1127470	UHMWPE
1127472	UHMWPE
1127474	UHMWPE
1127476	UHMWPE
1127477	UHMWPE
1127478	UHMWPE
1127480	UHMWPE
1127482	UHMWPE
1127484	UHMWPE
1127486	UHMWPE
1127488	UHMWPE
1127490	UHMWPE
1127510	UHMWPE
1127512	UHMWPE
1127514	UHMWPE
1127516	UHMWPE
1127517	UHMWPE
1127518	UHMWPE
1127520	UHMWPE
1127522	UHMWPE
1127524	UHMWPE
1127526	UHMWPE
1127528	UHMWPE
1127530	UHMWPE
1127540	UHMWPE
1127542	UHMWPE
1127544	UHMWPE
1127546	UHMWPE
1127547	UHMWPE
1127548	UHMWPE
1127550	UHMWPE
1127552	UHMWPE
1127554	UHMWPE
1127556	UHMWPE
1127558	UHMWPE
1127560	UHMWPE
1127570	UHMWPE
1127572	UHMWPE
1127574	UHMWPE
1127576	UHMWPE
1127577	UHMWPE
1127578	UHMWPE
1127580	UHMWPE

1127582	UHMWPE
1127584	UHMWPE
1127586	UHMWPE
1127588	UHMWPE
1127590	UHMWPE
1127610	UHMWPE
1127612	UHMWPE
1127614	UHMWPE
1127616	UHMWPE
1127617	UHMWPE
1127618	UHMWPE
1127620	UHMWPE
1127622	UHMWPE
1127624	UHMWPE
1127626	UHMWPE
1127628	UHMWPE
1127630	UHMWPE
1127640	UHMWPE
1127642	UHMWPE
1127644	UHMWPE
1127646	UHMWPE
1127647	UHMWPE
1127648	UHMWPE
1127650	UHMWPE
1127652	UHMWPE
1127654	UHMWPE
1127656	UHMWPE
1127658	UHMWPE
1127660	UHMWPE
1127670	UHMWPE
1127672	UHMWPE
1127674	UHMWPE
1127676	UHMWPE
1127677	UHMWPE
1127678	UHMWPE
1127680	UHMWPE
1127682	UHMWPE
1127684	UHMWPE
1127686	UHMWPE
1127688	UHMWPE
1127690	UHMWPE
1128001	Co-Cr-Mo Alloy
1128002	Co-Cr-Mo Alloy

1128003	Co-Cr-Mo Alloy
1128004	Co-Cr-Mo Alloy
1128005	Co-Cr-Mo Alloy
1128006	Co-Cr-Mo Alloy
1128011	Co-Cr-Mo Alloy
1128012	Co-Cr-Mo Alloy
1128013	Co-Cr-Mo Alloy
1128014	Co-Cr-Mo Alloy
1128015	Co-Cr-Mo Alloy
1128016	Co-Cr-Mo Alloy
1128051	Co-Cr-Mo Alloy
1128052	Co-Cr-Mo Alloy
1128053	Co-Cr-Mo Alloy
1128054	Co-Cr-Mo Alloy
1128055	Co-Cr-Mo Alloy
1128061	Co-Cr-Mo Alloy
1128062	Co-Cr-Mo Alloy
1128063	Co-Cr-Mo Alloy
1128064	Co-Cr-Mo Alloy
1128065	Co-Cr-Mo Alloy
1128100	UHMWPE
1128200	Ti6Al4V Alloy
1128901	Co-Cr-Mo Alloy + TiNbN
1128902	Co-Cr-Mo Alloy + TiNbN
1128903	Co-Cr-Mo Alloy + TiNbN
1128904	Co-Cr-Mo Alloy + TiNbN
1128905	Co-Cr-Mo Alloy + TiNbN
1128906	Co-Cr-Mo Alloy + TiNbN
1128911	Co-Cr-Mo Alloy + TiNbN
1128912	Co-Cr-Mo Alloy + TiNbN
1128913	Co-Cr-Mo Alloy + TiNbN
1128914	Co-Cr-Mo Alloy + TiNbN
1128915	Co-Cr-Mo Alloy + TiNbN
1128916	Co-Cr-Mo Alloy + TiNbN
1128951	Co-Cr-Mo Alloy + TiNbN
1128952	Co-Cr-Mo Alloy + TiNbN
1128953	Co-Cr-Mo Alloy + TiNbN
1128954	Co-Cr-Mo Alloy + TiNbN
1128955	Co-Cr-Mo Alloy + TiNbN
1128961	Co-Cr-Mo Alloy + TiNbN
1128962	Co-Cr-Mo Alloy + TiNbN
1128963	Co-Cr-Mo Alloy + TiNbN
1128964	Co-Cr-Mo Alloy + TiNbN

1128965	Co-Cr-Mo Alloy + TiNbN
1129001	Ti6Al4V Alloy
1129002	Ti6Al4V Alloy
1129003	Ti6Al4V Alloy
1129004	Ti6Al4V Alloy
1129005	Ti6Al4V Alloy
1129006	Ti6Al4V Alloy
1129011	Ti6Al4V Alloy
1129012	Ti6Al4V Alloy
1129013	Ti6Al4V Alloy
1129014	Ti6Al4V Alloy
1129015	Ti6Al4V Alloy
1129016	Ti6Al4V Alloy
1129051	Ti6Al4V Alloy
1129053	Ti6Al4V Alloy
1129055	Ti6Al4V Alloy
1129061	Ti6Al4V Alloy
1129063	Ti6Al4V Alloy
1129065	Ti6Al4V Alloy
1129101	Ti6Al4V Alloy
1129102	Ti6Al4V Alloy
1129103	Ti6Al4V Alloy
1129104	Ti6Al4V Alloy
1129105	Ti6Al4V Alloy
1129106	Ti6Al4V Alloy
1129111	Ti6Al4V Alloy
1129112	Ti6Al4V Alloy
1129113	Ti6Al4V Alloy
1129114	Ti6Al4V Alloy
1129115	Ti6Al4V Alloy
1129116	Ti6Al4V Alloy
1129151	Ti6Al4V Alloy
1129153	Ti6Al4V Alloy
1129155	Ti6Al4V Alloy
1129161	Ti6Al4V Alloy
1129163	Ti6Al4V Alloy
1129165	Ti6Al4V Alloy
1130001	Co-Cr-Mo Alloy + Hydroxyapatite
1130002	Co-Cr-Mo Alloy + Hydroxyapatite
1130003	Co-Cr-Mo Alloy + Hydroxyapatite
1130004	Co-Cr-Mo Alloy + Hydroxyapatite
1130005	Co-Cr-Mo Alloy + Hydroxyapatite
1130006	Co-Cr-Mo Alloy + Hydroxyapatite

1130007	Co-Cr-Mo Alloy + Hydroxyapatite
1130011	Co-Cr-Mo Alloy + Hydroxyapatite
1130012	Co-Cr-Mo Alloy + Hydroxyapatite
1130013	Co-Cr-Mo Alloy + Hydroxyapatite
1130014	Co-Cr-Mo Alloy + Hydroxyapatite
1130015	Co-Cr-Mo Alloy + Hydroxyapatite
1130016	Co-Cr-Mo Alloy + Hydroxyapatite
1130017	Co-Cr-Mo Alloy + Hydroxyapatite
1130053	Co-Cr-Mo Alloy + Hydroxyapatite
1130054	Co-Cr-Mo Alloy + Hydroxyapatite
1130063	Co-Cr-Mo Alloy + Hydroxyapatite
1130064	Co-Cr-Mo Alloy + Hydroxyapatite
1130101	Co-Cr-Mo Alloy + Hydroxyapatite
1130102	Co-Cr-Mo Alloy + Hydroxyapatite
1130103	Co-Cr-Mo Alloy + Hydroxyapatite
1130104	Co-Cr-Mo Alloy + Hydroxyapatite
1130105	Co-Cr-Mo Alloy + Hydroxyapatite
1130106	Co-Cr-Mo Alloy + Hydroxyapatite
1130107	Co-Cr-Mo Alloy + Hydroxyapatite
1130111	Co-Cr-Mo Alloy + Hydroxyapatite
1130112	Co-Cr-Mo Alloy + Hydroxyapatite
1130113	Co-Cr-Mo Alloy + Hydroxyapatite
1130114	Co-Cr-Mo Alloy + Hydroxyapatite
1130115	Co-Cr-Mo Alloy + Hydroxyapatite
1130116	Co-Cr-Mo Alloy + Hydroxyapatite
1130117	Co-Cr-Mo Alloy + Hydroxyapatite
1130153	Co-Cr-Mo Alloy + Hydroxyapatite
1130154	Co-Cr-Mo Alloy + Hydroxyapatite
1130163	Co-Cr-Mo Alloy + Hydroxyapatite
1130164	Co-Cr-Mo Alloy + Hydroxyapatite
1130201	Co-Cr-Mo Alloy + Hydroxyapatite
1130202	Co-Cr-Mo Alloy + Hydroxyapatite
1130203	Co-Cr-Mo Alloy + Hydroxyapatite
1130204	Co-Cr-Mo Alloy + Hydroxyapatite
1130205	Co-Cr-Mo Alloy + Hydroxyapatite
1130206	Co-Cr-Mo Alloy + Hydroxyapatite
1130207	Co-Cr-Mo Alloy + Hydroxyapatite
1130211	Co-Cr-Mo Alloy + Hydroxyapatite
1130212	Co-Cr-Mo Alloy + Hydroxyapatite
1130213	Co-Cr-Mo Alloy + Hydroxyapatite
1130214	Co-Cr-Mo Alloy + Hydroxyapatite
1130215	Co-Cr-Mo Alloy + Hydroxyapatite
1130216	Co-Cr-Mo Alloy + Hydroxyapatite

1130217	Co-Cr-Mo Alloy + Hydroxyapatite
1130222	Co-Cr-Mo Alloy + Hydroxyapatite
1130223	Co-Cr-Mo Alloy + Hydroxyapatite
1130224	Co-Cr-Mo Alloy + Hydroxyapatite
1130225	Co-Cr-Mo Alloy + Hydroxyapatite
1130226	Co-Cr-Mo Alloy + Hydroxyapatite
1130227	Co-Cr-Mo Alloy + Hydroxyapatite
1130232	Co-Cr-Mo Alloy + Hydroxyapatite
1130233	Co-Cr-Mo Alloy + Hydroxyapatite
1130234	Co-Cr-Mo Alloy + Hydroxyapatite
1130235	Co-Cr-Mo Alloy + Hydroxyapatite
1130236	Co-Cr-Mo Alloy + Hydroxyapatite
1130237	Co-Cr-Mo Alloy + Hydroxyapatite
1130242	Co-Cr-Mo Alloy + Hydroxyapatite
1130243	Co-Cr-Mo Alloy + Hydroxyapatite
1130244	Co-Cr-Mo Alloy + Hydroxyapatite
1130245	Co-Cr-Mo Alloy + Hydroxyapatite
1130246	Co-Cr-Mo Alloy + Hydroxyapatite
1130247	Co-Cr-Mo Alloy + Hydroxyapatite
1130251	Co-Cr-Mo Alloy + Hydroxyapatite
1130252	Co-Cr-Mo Alloy + Hydroxyapatite
1130253	Co-Cr-Mo Alloy + Hydroxyapatite
1130254	Co-Cr-Mo Alloy + Hydroxyapatite
1130255	Co-Cr-Mo Alloy + Hydroxyapatite
1130256	Co-Cr-Mo Alloy + Hydroxyapatite
1130261	Co-Cr-Mo Alloy + Hydroxyapatite
1130262	Co-Cr-Mo Alloy + Hydroxyapatite
1130263	Co-Cr-Mo Alloy + Hydroxyapatite
1130264	Co-Cr-Mo Alloy + Hydroxyapatite
1130265	Co-Cr-Mo Alloy + Hydroxyapatite
1130266	Co-Cr-Mo Alloy + Hydroxyapatite
1130272	Co-Cr-Mo Alloy + Hydroxyapatite
1130273	Co-Cr-Mo Alloy + Hydroxyapatite
1130274	Co-Cr-Mo Alloy + Hydroxyapatite
1130275	Co-Cr-Mo Alloy + Hydroxyapatite
1130276	Co-Cr-Mo Alloy + Hydroxyapatite
1130282	Co-Cr-Mo Alloy + Hydroxyapatite
1130283	Co-Cr-Mo Alloy + Hydroxyapatite
1130284	Co-Cr-Mo Alloy + Hydroxyapatite
1130285	Co-Cr-Mo Alloy + Hydroxyapatite
1130286	Co-Cr-Mo Alloy + Hydroxyapatite
1130292	Co-Cr-Mo Alloy + Hydroxyapatite
1130293	Co-Cr-Mo Alloy + Hydroxyapatite

1130294	Co-Cr-Mo Alloy + Hydroxyapatite
1130295	Co-Cr-Mo Alloy + Hydroxyapatite
1130296	Co-Cr-Mo Alloy + Hydroxyapatite
1130297	Co-Cr-Mo Alloy + Hydroxyapatite
1130301	Co-Cr-Mo Alloy + Hydroxyapatite
1130302	Co-Cr-Mo Alloy + Hydroxyapatite
1130303	Co-Cr-Mo Alloy + Hydroxyapatite
1130304	Co-Cr-Mo Alloy + Hydroxyapatite
1130305	Co-Cr-Mo Alloy + Hydroxyapatite
1130306	Co-Cr-Mo Alloy + Hydroxyapatite
1130307	Co-Cr-Mo Alloy + Hydroxyapatite
1130311	Co-Cr-Mo Alloy + Hydroxyapatite
1130312	Co-Cr-Mo Alloy + Hydroxyapatite
1130313	Co-Cr-Mo Alloy + Hydroxyapatite
1130314	Co-Cr-Mo Alloy + Hydroxyapatite
1130315	Co-Cr-Mo Alloy + Hydroxyapatite
1130316	Co-Cr-Mo Alloy + Hydroxyapatite
1130317	Co-Cr-Mo Alloy + Hydroxyapatite
1130353	Co-Cr-Mo Alloy + Hydroxyapatite
1130354	Co-Cr-Mo Alloy + Hydroxyapatite
1130363	Co-Cr-Mo Alloy + Hydroxyapatite
1130364	Co-Cr-Mo Alloy + Hydroxyapatite
1130401	Co-Cr-Mo Alloy + Hydroxyapatite
1130402	Co-Cr-Mo Alloy + Hydroxyapatite
1130403	Co-Cr-Mo Alloy + Hydroxyapatite
1130404	Co-Cr-Mo Alloy + Hydroxyapatite
1130405	Co-Cr-Mo Alloy + Hydroxyapatite
1130406	Co-Cr-Mo Alloy + Hydroxyapatite
1130407	Co-Cr-Mo Alloy + Hydroxyapatite
1130411	Co-Cr-Mo Alloy + Hydroxyapatite
1130412	Co-Cr-Mo Alloy + Hydroxyapatite
1130413	Co-Cr-Mo Alloy + Hydroxyapatite
1130414	Co-Cr-Mo Alloy + Hydroxyapatite
1130415	Co-Cr-Mo Alloy + Hydroxyapatite
1130416	Co-Cr-Mo Alloy + Hydroxyapatite
1130417	Co-Cr-Mo Alloy + Hydroxyapatite
1130453	Co-Cr-Mo Alloy + Hydroxyapatite
1130454	Co-Cr-Mo Alloy + Hydroxyapatite
1130463	Co-Cr-Mo Alloy + Hydroxyapatite
1130464	Co-Cr-Mo Alloy + Hydroxyapatite
1130501	Co-Cr-Mo Alloy + Hydroxyapatite
1130502	Co-Cr-Mo Alloy + Hydroxyapatite
1130503	Co-Cr-Mo Alloy + Hydroxyapatite

1130504	Co-Cr-Mo Alloy + Hydroxyapatite
1130505	Co-Cr-Mo Alloy + Hydroxyapatite
1130506	Co-Cr-Mo Alloy + Hydroxyapatite
1130507	Co-Cr-Mo Alloy + Hydroxyapatite
1130511	Co-Cr-Mo Alloy + Hydroxyapatite
1130512	Co-Cr-Mo Alloy + Hydroxyapatite
1130513	Co-Cr-Mo Alloy + Hydroxyapatite
1130514	Co-Cr-Mo Alloy + Hydroxyapatite
1130515	Co-Cr-Mo Alloy + Hydroxyapatite
1130516	Co-Cr-Mo Alloy + Hydroxyapatite
1130517	Co-Cr-Mo Alloy + Hydroxyapatite
1130553	Co-Cr-Mo Alloy + Hydroxyapatite
1130554	Co-Cr-Mo Alloy + Hydroxyapatite
1130563	Co-Cr-Mo Alloy + Hydroxyapatite
1130564	Co-Cr-Mo Alloy + Hydroxyapatite
1131001	Co-Cr-Mo Alloy
1131002	Co-Cr-Mo Alloy
1131003	Co-Cr-Mo Alloy
1131004	Co-Cr-Mo Alloy
1131005	Co-Cr-Mo Alloy
1131006	Co-Cr-Mo Alloy
1131007	Co-Cr-Mo Alloy
1131011	Co-Cr-Mo Alloy
1131012	Co-Cr-Mo Alloy
1131013	Co-Cr-Mo Alloy
1131014	Co-Cr-Mo Alloy
1131015	Co-Cr-Mo Alloy
1131016	Co-Cr-Mo Alloy
1131017	Co-Cr-Mo Alloy
1131053	Co-Cr-Mo Alloy
1131054	Co-Cr-Mo Alloy
1131063	Co-Cr-Mo Alloy
1131064	Co-Cr-Mo Alloy
1131101	Co-Cr-Mo Alloy
1131102	Co-Cr-Mo Alloy
1131103	Co-Cr-Mo Alloy
1131104	Co-Cr-Mo Alloy
1131105	Co-Cr-Mo Alloy
1131106	Co-Cr-Mo Alloy
1131107	Co-Cr-Mo Alloy
1131111	Co-Cr-Mo Alloy
1131112	Co-Cr-Mo Alloy
1131113	Co-Cr-Mo Alloy

1131114	Co-Cr-Mo Alloy
1131115	Co-Cr-Mo Alloy
1131116	Co-Cr-Mo Alloy
1131117	Co-Cr-Mo Alloy
1131153	Co-Cr-Mo Alloy
1131154	Co-Cr-Mo Alloy
1131163	Co-Cr-Mo Alloy
1131164	Co-Cr-Mo Alloy
1131201	Co-Cr-Mo Alloy
1131202	Co-Cr-Mo Alloy
1131203	Co-Cr-Mo Alloy
1131204	Co-Cr-Mo Alloy
1131205	Co-Cr-Mo Alloy
1131206	Co-Cr-Mo Alloy
1131207	Co-Cr-Mo Alloy
1131211	Co-Cr-Mo Alloy
1131212	Co-Cr-Mo Alloy
1131213	Co-Cr-Mo Alloy
1131214	Co-Cr-Mo Alloy
1131215	Co-Cr-Mo Alloy
1131216	Co-Cr-Mo Alloy
1131217	Co-Cr-Mo Alloy
1131222	Co-Cr-Mo Alloy
1131223	Co-Cr-Mo Alloy
1131224	Co-Cr-Mo Alloy
1131225	Co-Cr-Mo Alloy
1131226	Co-Cr-Mo Alloy
1131227	Co-Cr-Mo Alloy
1131232	Co-Cr-Mo Alloy
1131233	Co-Cr-Mo Alloy
1131234	Co-Cr-Mo Alloy
1131235	Co-Cr-Mo Alloy
1131236	Co-Cr-Mo Alloy
1131237	Co-Cr-Mo Alloy
1131253	Co-Cr-Mo Alloy
1131254	Co-Cr-Mo Alloy
1131255	Co-Cr-Mo Alloy
1131263	Co-Cr-Mo Alloy
1131264	Co-Cr-Mo Alloy
1131265	Co-Cr-Mo Alloy
1131273	Co-Cr-Mo Alloy
1131274	Co-Cr-Mo Alloy
1131283	Co-Cr-Mo Alloy

1131284	Co-Cr-Mo Alloy
1131301	Co-Cr-Mo Alloy
1131302	Co-Cr-Mo Alloy
1131303	Co-Cr-Mo Alloy
1131304	Co-Cr-Mo Alloy
1131305	Co-Cr-Mo Alloy
1131306	Co-Cr-Mo Alloy
1131307	Co-Cr-Mo Alloy
1131311	Co-Cr-Mo Alloy
1131312	Co-Cr-Mo Alloy
1131313	Co-Cr-Mo Alloy
1131314	Co-Cr-Mo Alloy
1131315	Co-Cr-Mo Alloy
1131316	Co-Cr-Mo Alloy
1131317	Co-Cr-Mo Alloy
1131353	Co-Cr-Mo Alloy
1131354	Co-Cr-Mo Alloy
1131363	Co-Cr-Mo Alloy
1131364	Co-Cr-Mo Alloy
1131401	Co-Cr-Mo Alloy
1131402	Co-Cr-Mo Alloy
1131403	Co-Cr-Mo Alloy
1131404	Co-Cr-Mo Alloy
1131405	Co-Cr-Mo Alloy
1131406	Co-Cr-Mo Alloy
1131407	Co-Cr-Mo Alloy
1131411	Co-Cr-Mo Alloy
1131412	Co-Cr-Mo Alloy
1131413	Co-Cr-Mo Alloy
1131414	Co-Cr-Mo Alloy
1131415	Co-Cr-Mo Alloy
1131416	Co-Cr-Mo Alloy
1131417	Co-Cr-Mo Alloy
1131453	Co-Cr-Mo Alloy
1131454	Co-Cr-Mo Alloy
1131463	Co-Cr-Mo Alloy
1131464	Co-Cr-Mo Alloy
1131501	Co-Cr-Mo Alloy
1131502	Co-Cr-Mo Alloy
1131503	Co-Cr-Mo Alloy
1131504	Co-Cr-Mo Alloy
1131505	Co-Cr-Mo Alloy
1131506	Co-Cr-Mo Alloy

1131507	Co-Cr-Mo Alloy
1131511	Co-Cr-Mo Alloy
1131512	Co-Cr-Mo Alloy
1131513	Co-Cr-Mo Alloy
1131514	Co-Cr-Mo Alloy
1131515	Co-Cr-Mo Alloy
1131516	Co-Cr-Mo Alloy
1131517	Co-Cr-Mo Alloy
1131553	Co-Cr-Mo Alloy
1131554	Co-Cr-Mo Alloy
1131563	Co-Cr-Mo Alloy
1131564	Co-Cr-Mo Alloy
1132001	Co-Cr-Mo Alloy + Hydroxyapatite
1132002	Co-Cr-Mo Alloy + Hydroxyapatite
1132003	Co-Cr-Mo Alloy + Hydroxyapatite
1132004	Co-Cr-Mo Alloy + Hydroxyapatite
1132005	Co-Cr-Mo Alloy + Hydroxyapatite
1132006	Co-Cr-Mo Alloy + Hydroxyapatite
1132007	Co-Cr-Mo Alloy + Hydroxyapatite
1132011	Co-Cr-Mo Alloy + Hydroxyapatite
1132012	Co-Cr-Mo Alloy + Hydroxyapatite
1132013	Co-Cr-Mo Alloy + Hydroxyapatite
1132014	Co-Cr-Mo Alloy + Hydroxyapatite
1132015	Co-Cr-Mo Alloy + Hydroxyapatite
1132016	Co-Cr-Mo Alloy + Hydroxyapatite
1132017	Co-Cr-Mo Alloy + Hydroxyapatite
1132053	Co-Cr-Mo Alloy + Hydroxyapatite
1132054	Co-Cr-Mo Alloy + Hydroxyapatite
1132063	Co-Cr-Mo Alloy + Hydroxyapatite
1132064	Co-Cr-Mo Alloy + Hydroxyapatite
1132101	Co-Cr-Mo Alloy + Hydroxyapatite
1132102	Co-Cr-Mo Alloy + Hydroxyapatite
1132103	Co-Cr-Mo Alloy + Hydroxyapatite
1132104	Co-Cr-Mo Alloy + Hydroxyapatite
1132105	Co-Cr-Mo Alloy + Hydroxyapatite
1132106	Co-Cr-Mo Alloy + Hydroxyapatite
1132107	Co-Cr-Mo Alloy + Hydroxyapatite
1132111	Co-Cr-Mo Alloy + Hydroxyapatite
1132112	Co-Cr-Mo Alloy + Hydroxyapatite
1132113	Co-Cr-Mo Alloy + Hydroxyapatite
1132114	Co-Cr-Mo Alloy + Hydroxyapatite
1132115	Co-Cr-Mo Alloy + Hydroxyapatite
1132116	Co-Cr-Mo Alloy + Hydroxyapatite

1132117	Co-Cr-Mo Alloy + Hydroxyapatite
1132153	Co-Cr-Mo Alloy + Hydroxyapatite
1132154	Co-Cr-Mo Alloy + Hydroxyapatite
1132163	Co-Cr-Mo Alloy + Hydroxyapatite
1132164	Co-Cr-Mo Alloy + Hydroxyapatite
1132201	Co-Cr-Mo Alloy + Hydroxyapatite
1132202	Co-Cr-Mo Alloy + Hydroxyapatite
1132203	Co-Cr-Mo Alloy + Hydroxyapatite
1132204	Co-Cr-Mo Alloy + Hydroxyapatite
1132205	Co-Cr-Mo Alloy + Hydroxyapatite
1132206	Co-Cr-Mo Alloy + Hydroxyapatite
1132207	Co-Cr-Mo Alloy + Hydroxyapatite
1132211	Co-Cr-Mo Alloy + Hydroxyapatite
1132212	Co-Cr-Mo Alloy + Hydroxyapatite
1132213	Co-Cr-Mo Alloy + Hydroxyapatite
1132214	Co-Cr-Mo Alloy + Hydroxyapatite
1132215	Co-Cr-Mo Alloy + Hydroxyapatite
1132216	Co-Cr-Mo Alloy + Hydroxyapatite
1132217	Co-Cr-Mo Alloy + Hydroxyapatite
1132222	Co-Cr-Mo Alloy + Hydroxyapatite
1132223	Co-Cr-Mo Alloy + Hydroxyapatite
1132224	Co-Cr-Mo Alloy + Hydroxyapatite
1132225	Co-Cr-Mo Alloy + Hydroxyapatite
1132226	Co-Cr-Mo Alloy + Hydroxyapatite
1132227	Co-Cr-Mo Alloy + Hydroxyapatite
1132232	Co-Cr-Mo Alloy + Hydroxyapatite
1132233	Co-Cr-Mo Alloy + Hydroxyapatite
1132234	Co-Cr-Mo Alloy + Hydroxyapatite
1132235	Co-Cr-Mo Alloy + Hydroxyapatite
1132236	Co-Cr-Mo Alloy + Hydroxyapatite
1132237	Co-Cr-Mo Alloy + Hydroxyapatite
1132242	Co-Cr-Mo Alloy + Hydroxyapatite
1132243	Co-Cr-Mo Alloy + Hydroxyapatite
1132244	Co-Cr-Mo Alloy + Hydroxyapatite
1132245	Co-Cr-Mo Alloy + Hydroxyapatite
1132246	Co-Cr-Mo Alloy + Hydroxyapatite
1132247	Co-Cr-Mo Alloy + Hydroxyapatite
1132251	Co-Cr-Mo Alloy + Hydroxyapatite
1132252	Co-Cr-Mo Alloy + Hydroxyapatite
1132253	Co-Cr-Mo Alloy + Hydroxyapatite
1132254	Co-Cr-Mo Alloy + Hydroxyapatite
1132255	Co-Cr-Mo Alloy + Hydroxyapatite
1132256	Co-Cr-Mo Alloy + Hydroxyapatite

1132261	Co-Cr-Mo Alloy + Hydroxyapatite
1132262	Co-Cr-Mo Alloy + Hydroxyapatite
1132263	Co-Cr-Mo Alloy + Hydroxyapatite
1132264	Co-Cr-Mo Alloy + Hydroxyapatite
1132265	Co-Cr-Mo Alloy + Hydroxyapatite
1132266	Co-Cr-Mo Alloy + Hydroxyapatite
1132272	Co-Cr-Mo Alloy + Hydroxyapatite
1132273	Co-Cr-Mo Alloy + Hydroxyapatite
1132274	Co-Cr-Mo Alloy + Hydroxyapatite
1132275	Co-Cr-Mo Alloy + Hydroxyapatite
1132276	Co-Cr-Mo Alloy + Hydroxyapatite
1132282	Co-Cr-Mo Alloy + Hydroxyapatite
1132283	Co-Cr-Mo Alloy + Hydroxyapatite
1132284	Co-Cr-Mo Alloy + Hydroxyapatite
1132285	Co-Cr-Mo Alloy + Hydroxyapatite
1132286	Co-Cr-Mo Alloy + Hydroxyapatite
1132292	Co-Cr-Mo Alloy + Hydroxyapatite
1132293	Co-Cr-Mo Alloy + Hydroxyapatite
1132294	Co-Cr-Mo Alloy + Hydroxyapatite
1132295	Co-Cr-Mo Alloy + Hydroxyapatite
1132296	Co-Cr-Mo Alloy + Hydroxyapatite
1132297	Co-Cr-Mo Alloy + Hydroxyapatite
1132401	Co-Cr-Mo Alloy + Hydroxyapatite
1132402	Co-Cr-Mo Alloy + Hydroxyapatite
1132403	Co-Cr-Mo Alloy + Hydroxyapatite
1132404	Co-Cr-Mo Alloy + Hydroxyapatite
1132405	Co-Cr-Mo Alloy + Hydroxyapatite
1132406	Co-Cr-Mo Alloy + Hydroxyapatite
1132407	Co-Cr-Mo Alloy + Hydroxyapatite
1132411	Co-Cr-Mo Alloy + Hydroxyapatite
1132412	Co-Cr-Mo Alloy + Hydroxyapatite
1132413	Co-Cr-Mo Alloy + Hydroxyapatite
1132414	Co-Cr-Mo Alloy + Hydroxyapatite
1132415	Co-Cr-Mo Alloy + Hydroxyapatite
1132416	Co-Cr-Mo Alloy + Hydroxyapatite
1132417	Co-Cr-Mo Alloy + Hydroxyapatite
1132453	Co-Cr-Mo Alloy + Hydroxyapatite
1132454	Co-Cr-Mo Alloy + Hydroxyapatite
1132463	Co-Cr-Mo Alloy + Hydroxyapatite
1132464	Co-Cr-Mo Alloy + Hydroxyapatite
1133001	Co-Cr-Mo Alloy
1133002	Co-Cr-Mo Alloy
1133003	Co-Cr-Mo Alloy

1133004	Co-Cr-Mo Alloy
1133005	Co-Cr-Mo Alloy
1133006	Co-Cr-Mo Alloy
1133007	Co-Cr-Mo Alloy
1133011	Co-Cr-Mo Alloy
1133012	Co-Cr-Mo Alloy
1133013	Co-Cr-Mo Alloy
1133014	Co-Cr-Mo Alloy
1133015	Co-Cr-Mo Alloy
1133016	Co-Cr-Mo Alloy
1133017	Co-Cr-Mo Alloy
1133053	Co-Cr-Mo Alloy
1133054	Co-Cr-Mo Alloy
1133063	Co-Cr-Mo Alloy
1133064	Co-Cr-Mo Alloy
1133101	Co-Cr-Mo Alloy
1133102	Co-Cr-Mo Alloy
1133103	Co-Cr-Mo Alloy
1133104	Co-Cr-Mo Alloy
1133105	Co-Cr-Mo Alloy
1133106	Co-Cr-Mo Alloy
1133107	Co-Cr-Mo Alloy
1133111	Co-Cr-Mo Alloy
1133112	Co-Cr-Mo Alloy
1133113	Co-Cr-Mo Alloy
1133114	Co-Cr-Mo Alloy
1133115	Co-Cr-Mo Alloy
1133116	Co-Cr-Mo Alloy
1133117	Co-Cr-Mo Alloy
1133153	Co-Cr-Mo Alloy
1133154	Co-Cr-Mo Alloy
1133163	Co-Cr-Mo Alloy
1133164	Co-Cr-Mo Alloy
1133201	Co-Cr-Mo Alloy
1133202	Co-Cr-Mo Alloy
1133203	Co-Cr-Mo Alloy
1133204	Co-Cr-Mo Alloy
1133205	Co-Cr-Mo Alloy
1133206	Co-Cr-Mo Alloy
1133207	Co-Cr-Mo Alloy
1133211	Co-Cr-Mo Alloy
1133212	Co-Cr-Mo Alloy
1133213	Co-Cr-Mo Alloy

1133214	Co-Cr-Mo Alloy
1133215	Co-Cr-Mo Alloy
1133216	Co-Cr-Mo Alloy
1133217	Co-Cr-Mo Alloy
1133222	Co-Cr-Mo Alloy
1133223	Co-Cr-Mo Alloy
1133224	Co-Cr-Mo Alloy
1133225	Co-Cr-Mo Alloy
1133226	Co-Cr-Mo Alloy
1133227	Co-Cr-Mo Alloy
1133232	Co-Cr-Mo Alloy
1133233	Co-Cr-Mo Alloy
1133234	Co-Cr-Mo Alloy
1133235	Co-Cr-Mo Alloy
1133236	Co-Cr-Mo Alloy
1133237	Co-Cr-Mo Alloy
1133253	Co-Cr-Mo Alloy
1133254	Co-Cr-Mo Alloy
1133255	Co-Cr-Mo Alloy
1133263	Co-Cr-Mo Alloy
1133264	Co-Cr-Mo Alloy
1133265	Co-Cr-Mo Alloy
1133273	Co-Cr-Mo Alloy
1133274	Co-Cr-Mo Alloy
1133283	Co-Cr-Mo Alloy
1133284	Co-Cr-Mo Alloy
1133401	Co-Cr-Mo Alloy
1133402	Co-Cr-Mo Alloy
1133403	Co-Cr-Mo Alloy
1133404	Co-Cr-Mo Alloy
1133405	Co-Cr-Mo Alloy
1133406	Co-Cr-Mo Alloy
1133407	Co-Cr-Mo Alloy
1133411	Co-Cr-Mo Alloy
1133412	Co-Cr-Mo Alloy
1133413	Co-Cr-Mo Alloy
1133414	Co-Cr-Mo Alloy
1133415	Co-Cr-Mo Alloy
1133416	Co-Cr-Mo Alloy
1133417	Co-Cr-Mo Alloy
1133453	Co-Cr-Mo Alloy
1133454	Co-Cr-Mo Alloy
1133463	Co-Cr-Mo Alloy

1133464	Co-Cr-Mo Alloy
1134001	Ti6Al4V Alloy + TiNbN
1134002	Ti6Al4V Alloy + TiNbN
1134003	Ti6Al4V Alloy + TiNbN
1134004	Ti6Al4V Alloy + TiNbN
1134005	Ti6Al4V Alloy + TiNbN
1134006	Ti6Al4V Alloy + TiNbN
1134007	Ti6Al4V Alloy + TiNbN
1134011	Ti6Al4V Alloy + TiNbN
1134012	Ti6Al4V Alloy + TiNbN
1134013	Ti6Al4V Alloy + TiNbN
1134014	Ti6Al4V Alloy + TiNbN
1134015	Ti6Al4V Alloy + TiNbN
1134016	Ti6Al4V Alloy + TiNbN
1134017	Ti6Al4V Alloy + TiNbN
1134022	Ti6Al4V Alloy + TiNbN
1134023	Ti6Al4V Alloy + TiNbN
1134024	Ti6Al4V Alloy + TiNbN
1134025	Ti6Al4V Alloy + TiNbN
1134026	Ti6Al4V Alloy + TiNbN
1134027	Ti6Al4V Alloy + TiNbN
1134032	Ti6Al4V Alloy + TiNbN
1134033	Ti6Al4V Alloy + TiNbN
1134034	Ti6Al4V Alloy + TiNbN
1134035	Ti6Al4V Alloy + TiNbN
1134036	Ti6Al4V Alloy + TiNbN
1134037	Ti6Al4V Alloy + TiNbN
1134051	Ti6Al4V Alloy + TiNbN
1134052	Ti6Al4V Alloy + TiNbN
1134053	Ti6Al4V Alloy + TiNbN
1134054	Ti6Al4V Alloy + TiNbN
1134055	Ti6Al4V Alloy + TiNbN
1134056	Ti6Al4V Alloy + TiNbN
1134061	Ti6Al4V Alloy + TiNbN
1134062	Ti6Al4V Alloy + TiNbN
1134063	Ti6Al4V Alloy + TiNbN
1134064	Ti6Al4V Alloy + TiNbN
1134065	Ti6Al4V Alloy + TiNbN
1134066	Ti6Al4V Alloy + TiNbN
1134072	Ti6Al4V Alloy + TiNbN
1134073	Ti6Al4V Alloy + TiNbN
1134074	Ti6Al4V Alloy + TiNbN
1134075	Ti6Al4V Alloy + TiNbN

1134076	Ti6Al4V Alloy + TiNbN
1134082	Ti6Al4V Alloy + TiNbN
1134083	Ti6Al4V Alloy + TiNbN
1134084	Ti6Al4V Alloy + TiNbN
1134085	Ti6Al4V Alloy + TiNbN
1134086	Ti6Al4V Alloy + TiNbN
1134101	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134102	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134103	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134104	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134105	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134106	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134107	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134111	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134112	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134113	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134114	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134115	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134116	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134117	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134122	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134123	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134124	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134125	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134126	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134127	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134132	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134133	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134134	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134135	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134136	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134137	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134151	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134152	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134153	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134154	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134155	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134156	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134161	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134162	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134163	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134164	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134165	Ti6Al4V Alloy + TiNbN + Hydroxyapatite

1134166	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134172	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134173	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134174	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134175	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134176	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134182	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134183	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134184	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134185	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1134186	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1138001	Co-Cr-Mo Alloy
1138002	Co-Cr-Mo Alloy
1138003	Co-Cr-Mo Alloy
1138004	Co-Cr-Mo Alloy
1138005	Co-Cr-Mo Alloy
1138006	Co-Cr-Mo Alloy
1138007	Co-Cr-Mo Alloy
1138011	Co-Cr-Mo Alloy
1138012	Co-Cr-Mo Alloy
1138013	Co-Cr-Mo Alloy
1138014	Co-Cr-Mo Alloy
1138015	Co-Cr-Mo Alloy
1138016	Co-Cr-Mo Alloy
1138017	Co-Cr-Mo Alloy
1138053	Co-Cr-Mo Alloy
1138054	Co-Cr-Mo Alloy
1138063	Co-Cr-Mo Alloy
1138064	Co-Cr-Mo Alloy
1138101	Co-Cr-Mo Alloy
1138102	Co-Cr-Mo Alloy
1138103	Co-Cr-Mo Alloy
1138104	Co-Cr-Mo Alloy
1138105	Co-Cr-Mo Alloy
1138106	Co-Cr-Mo Alloy
1138107	Co-Cr-Mo Alloy
1138111	Co-Cr-Mo Alloy
1138112	Co-Cr-Mo Alloy
1138113	Co-Cr-Mo Alloy
1138114	Co-Cr-Mo Alloy
1138115	Co-Cr-Mo Alloy
1138116	Co-Cr-Mo Alloy
1138117	Co-Cr-Mo Alloy

1138153	Co-Cr-Mo Alloy
1138154	Co-Cr-Mo Alloy
1138163	Co-Cr-Mo Alloy
1138164	Co-Cr-Mo Alloy
1138201	Co-Cr-Mo Alloy
1138202	Co-Cr-Mo Alloy
1138203	Co-Cr-Mo Alloy
1138204	Co-Cr-Mo Alloy
1138205	Co-Cr-Mo Alloy
1138206	Co-Cr-Mo Alloy
1138207	Co-Cr-Mo Alloy
1138211	Co-Cr-Mo Alloy
1138212	Co-Cr-Mo Alloy
1138213	Co-Cr-Mo Alloy
1138214	Co-Cr-Mo Alloy
1138215	Co-Cr-Mo Alloy
1138216	Co-Cr-Mo Alloy
1138217	Co-Cr-Mo Alloy
1138222	Co-Cr-Mo Alloy
1138223	Co-Cr-Mo Alloy
1138224	Co-Cr-Mo Alloy
1138225	Co-Cr-Mo Alloy
1138226	Co-Cr-Mo Alloy
1138227	Co-Cr-Mo Alloy
1138232	Co-Cr-Mo Alloy
1138233	Co-Cr-Mo Alloy
1138234	Co-Cr-Mo Alloy
1138235	Co-Cr-Mo Alloy
1138236	Co-Cr-Mo Alloy
1138237	Co-Cr-Mo Alloy
1138242	Co-Cr-Mo Alloy
1138243	Co-Cr-Mo Alloy
1138244	Co-Cr-Mo Alloy
1138245	Co-Cr-Mo Alloy
1138246	Co-Cr-Mo Alloy
1138247	Co-Cr-Mo Alloy
1138251	Co-Cr-Mo Alloy
1138252	Co-Cr-Mo Alloy
1138253	Co-Cr-Mo Alloy
1138254	Co-Cr-Mo Alloy
1138255	Co-Cr-Mo Alloy
1138256	Co-Cr-Mo Alloy
1138261	Co-Cr-Mo Alloy

1138262	Co-Cr-Mo Alloy
1138263	Co-Cr-Mo Alloy
1138264	Co-Cr-Mo Alloy
1138265	Co-Cr-Mo Alloy
1138266	Co-Cr-Mo Alloy
1138272	Co-Cr-Mo Alloy
1138273	Co-Cr-Mo Alloy
1138274	Co-Cr-Mo Alloy
1138275	Co-Cr-Mo Alloy
1138276	Co-Cr-Mo Alloy
1138282	Co-Cr-Mo Alloy
1138283	Co-Cr-Mo Alloy
1138284	Co-Cr-Mo Alloy
1138285	Co-Cr-Mo Alloy
1138286	Co-Cr-Mo Alloy
1138292	Co-Cr-Mo Alloy
1138293	Co-Cr-Mo Alloy
1138294	Co-Cr-Mo Alloy
1138295	Co-Cr-Mo Alloy
1138296	Co-Cr-Mo Alloy
1138297	Co-Cr-Mo Alloy
1138301	Co-Cr-Mo Alloy
1138302	Co-Cr-Mo Alloy
1138303	Co-Cr-Mo Alloy
1138304	Co-Cr-Mo Alloy
1138305	Co-Cr-Mo Alloy
1138306	Co-Cr-Mo Alloy
1138307	Co-Cr-Mo Alloy
1138311	Co-Cr-Mo Alloy
1138312	Co-Cr-Mo Alloy
1138313	Co-Cr-Mo Alloy
1138314	Co-Cr-Mo Alloy
1138315	Co-Cr-Mo Alloy
1138316	Co-Cr-Mo Alloy
1138317	Co-Cr-Mo Alloy
1138353	Co-Cr-Mo Alloy
1138354	Co-Cr-Mo Alloy
1138363	Co-Cr-Mo Alloy
1138364	Co-Cr-Mo Alloy
1138401	Co-Cr-Mo Alloy
1138402	Co-Cr-Mo Alloy
1138403	Co-Cr-Mo Alloy
1138404	Co-Cr-Mo Alloy

1138405	Co-Cr-Mo Alloy
1138406	Co-Cr-Mo Alloy
1138407	Co-Cr-Mo Alloy
1138411	Co-Cr-Mo Alloy
1138412	Co-Cr-Mo Alloy
1138413	Co-Cr-Mo Alloy
1138414	Co-Cr-Mo Alloy
1138415	Co-Cr-Mo Alloy
1138416	Co-Cr-Mo Alloy
1138417	Co-Cr-Mo Alloy
1138453	Co-Cr-Mo Alloy
1138454	Co-Cr-Mo Alloy
1138463	Co-Cr-Mo Alloy
1138464	Co-Cr-Mo Alloy
1138501	Co-Cr-Mo Alloy
1138502	Co-Cr-Mo Alloy
1138503	Co-Cr-Mo Alloy
1138504	Co-Cr-Mo Alloy
1138505	Co-Cr-Mo Alloy
1138506	Co-Cr-Mo Alloy
1138507	Co-Cr-Mo Alloy
1138511	Co-Cr-Mo Alloy
1138512	Co-Cr-Mo Alloy
1138513	Co-Cr-Mo Alloy
1138514	Co-Cr-Mo Alloy
1138515	Co-Cr-Mo Alloy
1138516	Co-Cr-Mo Alloy
1138517	Co-Cr-Mo Alloy
1138553	Co-Cr-Mo Alloy
1138554	Co-Cr-Mo Alloy
1138563	Co-Cr-Mo Alloy
1138564	Co-Cr-Mo Alloy
1138901	Co-Cr-Mo Alloy + TiNbN
1138902	Co-Cr-Mo Alloy + TiNbN
1138903	Co-Cr-Mo Alloy + TiNbN
1138904	Co-Cr-Mo Alloy + TiNbN
1138905	Co-Cr-Mo Alloy + TiNbN
1138906	Co-Cr-Mo Alloy + TiNbN
1138907	Co-Cr-Mo Alloy + TiNbN
1138911	Co-Cr-Mo Alloy + TiNbN
1138912	Co-Cr-Mo Alloy + TiNbN
1138913	Co-Cr-Mo Alloy + TiNbN
1138914	Co-Cr-Mo Alloy + TiNbN

1138915	Co-Cr-Mo Alloy + TiNbN
1138916	Co-Cr-Mo Alloy + TiNbN
1138917	Co-Cr-Mo Alloy + TiNbN
1138922	Co-Cr-Mo Alloy + TiNbN
1138923	Co-Cr-Mo Alloy + TiNbN
1138924	Co-Cr-Mo Alloy + TiNbN
1138925	Co-Cr-Mo Alloy + TiNbN
1138926	Co-Cr-Mo Alloy + TiNbN
1138927	Co-Cr-Mo Alloy + TiNbN
1138932	Co-Cr-Mo Alloy + TiNbN
1138933	Co-Cr-Mo Alloy + TiNbN
1138934	Co-Cr-Mo Alloy + TiNbN
1138935	Co-Cr-Mo Alloy + TiNbN
1138936	Co-Cr-Mo Alloy + TiNbN
1138937	Co-Cr-Mo Alloy + TiNbN
1138942	Co-Cr-Mo Alloy + TiNbN
1138943	Co-Cr-Mo Alloy + TiNbN
1138944	Co-Cr-Mo Alloy + TiNbN
1138945	Co-Cr-Mo Alloy + TiNbN
1138946	Co-Cr-Mo Alloy + TiNbN
1138947	Co-Cr-Mo Alloy + TiNbN
1138951	Co-Cr-Mo Alloy + TiNbN
1138952	Co-Cr-Mo Alloy + TiNbN
1138953	Co-Cr-Mo Alloy + TiNbN
1138954	Co-Cr-Mo Alloy + TiNbN
1138955	Co-Cr-Mo Alloy + TiNbN
1138956	Co-Cr-Mo Alloy + TiNbN
1138961	Co-Cr-Mo Alloy + TiNbN
1138962	Co-Cr-Mo Alloy + TiNbN
1138963	Co-Cr-Mo Alloy + TiNbN
1138964	Co-Cr-Mo Alloy + TiNbN
1138965	Co-Cr-Mo Alloy + TiNbN
1138966	Co-Cr-Mo Alloy + TiNbN
1138972	Co-Cr-Mo Alloy + TiNbN
1138973	Co-Cr-Mo Alloy + TiNbN
1138974	Co-Cr-Mo Alloy + TiNbN
1138975	Co-Cr-Mo Alloy + TiNbN
1138976	Co-Cr-Mo Alloy + TiNbN
1138982	Co-Cr-Mo Alloy + TiNbN
1138983	Co-Cr-Mo Alloy + TiNbN
1138984	Co-Cr-Mo Alloy + TiNbN
1138985	Co-Cr-Mo Alloy + TiNbN
1138986	Co-Cr-Mo Alloy + TiNbN

1138992	Co-Cr-Mo Alloy + TiNbN
1138993	Co-Cr-Mo Alloy + TiNbN
1138994	Co-Cr-Mo Alloy + TiNbN
1138995	Co-Cr-Mo Alloy + TiNbN
1138996	Co-Cr-Mo Alloy + TiNbN
1138997	Co-Cr-Mo Alloy + TiNbN
1139042	Ti6Al4V Alloy + TiNbN
1139043	Ti6Al4V Alloy + TiNbN
1139044	Ti6Al4V Alloy + TiNbN
1139045	Ti6Al4V Alloy + TiNbN
1139046	Ti6Al4V Alloy + TiNbN
1139047	Ti6Al4V Alloy + TiNbN
1139092	Ti6Al4V Alloy + TiNbN
1139093	Ti6Al4V Alloy + TiNbN
1139094	Ti6Al4V Alloy + TiNbN
1139095	Ti6Al4V Alloy + TiNbN
1139096	Ti6Al4V Alloy + TiNbN
1139097	Ti6Al4V Alloy + TiNbN
1139142	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139143	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139144	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139145	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139146	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139147	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139192	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139193	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139194	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139195	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139196	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1139197	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1181001	Co-Cr-Mo Alloy
1181002	Co-Cr-Mo Alloy
1181003	Co-Cr-Mo Alloy
1181004	Co-Cr-Mo Alloy
1181005	Co-Cr-Mo Alloy
1181006	Co-Cr-Mo Alloy
1181007	Co-Cr-Mo Alloy
1181011	Co-Cr-Mo Alloy
1181012	Co-Cr-Mo Alloy
1181013	Co-Cr-Mo Alloy
1181014	Co-Cr-Mo Alloy
1181015	Co-Cr-Mo Alloy
1181016	Co-Cr-Mo Alloy

1181017	Co-Cr-Mo Alloy
1181022	Co-Cr-Mo Alloy
1181023	Co-Cr-Mo Alloy
1181024	Co-Cr-Mo Alloy
1181025	Co-Cr-Mo Alloy
1181026	Co-Cr-Mo Alloy
1181027	Co-Cr-Mo Alloy
1181032	Co-Cr-Mo Alloy
1181033	Co-Cr-Mo Alloy
1181034	Co-Cr-Mo Alloy
1181035	Co-Cr-Mo Alloy
1181036	Co-Cr-Mo Alloy
1181037	Co-Cr-Mo Alloy
1181062	Co-Cr-Mo Alloy
1181063	Co-Cr-Mo Alloy
1181064	Co-Cr-Mo Alloy
1181065	Co-Cr-Mo Alloy
1181066	Co-Cr-Mo Alloy
1181071	Co-Cr-Mo Alloy
1181072	Co-Cr-Mo Alloy
1181073	Co-Cr-Mo Alloy
1181074	Co-Cr-Mo Alloy
1181075	Co-Cr-Mo Alloy
1181076	Co-Cr-Mo Alloy
1181081	Co-Cr-Mo Alloy
1181082	Co-Cr-Mo Alloy
1181083	Co-Cr-Mo Alloy
1181084	Co-Cr-Mo Alloy
1181085	Co-Cr-Mo Alloy
1181086	Co-Cr-Mo Alloy
1181092	Co-Cr-Mo Alloy
1181093	Co-Cr-Mo Alloy
1181094	Co-Cr-Mo Alloy
1181095	Co-Cr-Mo Alloy
1181096	Co-Cr-Mo Alloy
1181101	Co-Cr-Mo Alloy
1181102	Co-Cr-Mo Alloy
1181103	Co-Cr-Mo Alloy
1181104	Co-Cr-Mo Alloy
1181105	Co-Cr-Mo Alloy
1181106	Co-Cr-Mo Alloy
1181107	Co-Cr-Mo Alloy
1181111	Co-Cr-Mo Alloy

1181112	Co-Cr-Mo Alloy
1181113	Co-Cr-Mo Alloy
1181114	Co-Cr-Mo Alloy
1181115	Co-Cr-Mo Alloy
1181116	Co-Cr-Mo Alloy
1181117	Co-Cr-Mo Alloy
1181262	Co-Cr-Mo Alloy
1181263	Co-Cr-Mo Alloy
1181264	Co-Cr-Mo Alloy
1181265	Co-Cr-Mo Alloy
1181266	Co-Cr-Mo Alloy
1181267	Co-Cr-Mo Alloy
1181272	Co-Cr-Mo Alloy
1181273	Co-Cr-Mo Alloy
1181274	Co-Cr-Mo Alloy
1181275	Co-Cr-Mo Alloy
1181276	Co-Cr-Mo Alloy
1181277	Co-Cr-Mo Alloy
1182001	Co-Cr-Mo Alloy
1182002	Co-Cr-Mo Alloy
1182003	Co-Cr-Mo Alloy
1182004	Co-Cr-Mo Alloy
1182005	Co-Cr-Mo Alloy
1182006	Co-Cr-Mo Alloy
1182011	Co-Cr-Mo Alloy
1182012	Co-Cr-Mo Alloy
1182013	Co-Cr-Mo Alloy
1182014	Co-Cr-Mo Alloy
1182015	Co-Cr-Mo Alloy
1182016	Co-Cr-Mo Alloy
1183001	Co-Cr-Mo Alloy + TiNbN
1183002	Co-Cr-Mo Alloy + TiNbN
1183003	Co-Cr-Mo Alloy + TiNbN
1183004	Co-Cr-Mo Alloy + TiNbN
1183005	Co-Cr-Mo Alloy + TiNbN
1183006	Co-Cr-Mo Alloy + TiNbN
1183011	Co-Cr-Mo Alloy + TiNbN
1183012	Co-Cr-Mo Alloy + TiNbN
1183013	Co-Cr-Mo Alloy + TiNbN
1183014	Co-Cr-Mo Alloy + TiNbN
1183015	Co-Cr-Mo Alloy + TiNbN
1183016	Co-Cr-Mo Alloy + TiNbN
1188901	Co-Cr-Mo Alloy + TiNbN

1188902	Co-Cr-Mo Alloy + TiNbN
1188903	Co-Cr-Mo Alloy + TiNbN
1188904	Co-Cr-Mo Alloy + TiNbN
1188905	Co-Cr-Mo Alloy + TiNbN
1188906	Co-Cr-Mo Alloy + TiNbN
1188907	Co-Cr-Mo Alloy + TiNbN
1188911	Co-Cr-Mo Alloy + TiNbN
1188912	Co-Cr-Mo Alloy + TiNbN
1188913	Co-Cr-Mo Alloy + TiNbN
1188914	Co-Cr-Mo Alloy + TiNbN
1188915	Co-Cr-Mo Alloy + TiNbN
1188916	Co-Cr-Mo Alloy + TiNbN
1188917	Co-Cr-Mo Alloy + TiNbN
1188922	Co-Cr-Mo Alloy + TiNbN
1188923	Co-Cr-Mo Alloy + TiNbN
1188924	Co-Cr-Mo Alloy + TiNbN
1188925	Co-Cr-Mo Alloy + TiNbN
1188926	Co-Cr-Mo Alloy + TiNbN
1188927	Co-Cr-Mo Alloy + TiNbN
1188932	Co-Cr-Mo Alloy + TiNbN
1188933	Co-Cr-Mo Alloy + TiNbN
1188934	Co-Cr-Mo Alloy + TiNbN
1188935	Co-Cr-Mo Alloy + TiNbN
1188936	Co-Cr-Mo Alloy + TiNbN
1188937	Co-Cr-Mo Alloy + TiNbN
1188951	Co-Cr-Mo Alloy + TiNbN
1188952	Co-Cr-Mo Alloy + TiNbN
1188953	Co-Cr-Mo Alloy + TiNbN
1188954	Co-Cr-Mo Alloy + TiNbN
1188955	Co-Cr-Mo Alloy + TiNbN
1188956	Co-Cr-Mo Alloy + TiNbN
1188961	Co-Cr-Mo Alloy + TiNbN
1188962	Co-Cr-Mo Alloy + TiNbN
1188963	Co-Cr-Mo Alloy + TiNbN
1188964	Co-Cr-Mo Alloy + TiNbN
1188965	Co-Cr-Mo Alloy + TiNbN
1188966	Co-Cr-Mo Alloy + TiNbN
1188972	Co-Cr-Mo Alloy + TiNbN
1188973	Co-Cr-Mo Alloy + TiNbN
1188974	Co-Cr-Mo Alloy + TiNbN
1188975	Co-Cr-Mo Alloy + TiNbN
1188976	Co-Cr-Mo Alloy + TiNbN
1188982	Co-Cr-Mo Alloy + TiNbN

1188983	Co-Cr-Mo Alloy + TiNbN
1188984	Co-Cr-Mo Alloy + TiNbN
1188985	Co-Cr-Mo Alloy + TiNbN
1188986	Co-Cr-Mo Alloy + TiNbN
1189002	Co-Cr-Mo Alloy + TiNbN
1189003	Co-Cr-Mo Alloy + TiNbN
1189004	Co-Cr-Mo Alloy + TiNbN
1189005	Co-Cr-Mo Alloy + TiNbN
1189006	Co-Cr-Mo Alloy + TiNbN
1189007	Co-Cr-Mo Alloy + TiNbN
1189012	Co-Cr-Mo Alloy + TiNbN
1189013	Co-Cr-Mo Alloy + TiNbN
1189014	Co-Cr-Mo Alloy + TiNbN
1189015	Co-Cr-Mo Alloy + TiNbN
1189016	Co-Cr-Mo Alloy + TiNbN
1189017	Co-Cr-Mo Alloy + TiNbN
1210001	Ti6Al4V Alloy
1210002	Ti6Al4V Alloy
1210003	Ti6Al4V Alloy
1210011	Ti6Al4V Alloy
1210012	Ti6Al4V Alloy
1210013	Ti6Al4V Alloy
1210021	Ti6Al4V Alloy
1210022	Ti6Al4V Alloy
1210023	Ti6Al4V Alloy
1210101	Co-Cr-Mo Alloy
1210102	Co-Cr-Mo Alloy
1210103	Co-Cr-Mo Alloy
1210111	Co-Cr-Mo Alloy
1210112	Co-Cr-Mo Alloy
1210113	Co-Cr-Mo Alloy
1210121	Co-Cr-Mo Alloy
1210122	Co-Cr-Mo Alloy
1210123	Co-Cr-Mo Alloy
1210131	Co-Cr-Mo Alloy
1210132	Co-Cr-Mo Alloy
1210133	Co-Cr-Mo Alloy
1210151	Co-Cr-Mo Alloy + TiNbN
1210152	Co-Cr-Mo Alloy + TiNbN
1210153	Co-Cr-Mo Alloy + TiNbN
1210161	Co-Cr-Mo Alloy + TiNbN
1210162	Co-Cr-Mo Alloy + TiNbN
1210163	Co-Cr-Mo Alloy + TiNbN

1210171	Co-Cr-Mo Alloy + TiNbN
1210172	Co-Cr-Mo Alloy + TiNbN
1210173	Co-Cr-Mo Alloy + TiNbN
1210181	Co-Cr-Mo Alloy + TiNbN
1210182	Co-Cr-Mo Alloy + TiNbN
1210183	Co-Cr-Mo Alloy + TiNbN
1210210	Co-Cr-Mo Alloy + UHMWPE
1210220	Co-Cr-Mo Alloy + UHMWPE
1210230	Co-Cr-Mo Alloy + UHMWPE
1210240	Co-Cr-Mo Alloy + UHMWPE
1210250	Co-Cr-Mo Alloy + UHMWPE
1210260	Co-Cr-Mo Alloy + UHMWPE
1210301	Ti6Al4V Alloy
1210302	Ti6Al4V Alloy
1210303	Ti6Al4V Alloy
1210304	Ti6Al4V Alloy
1210305	Ti6Al4V Alloy
1210306	Ti6Al4V Alloy
1210307	Ti6Al4V Alloy
1210308	Ti6Al4V Alloy
1210309	Ti6Al4V Alloy
1210310	Ti6Al4V Alloy
1210311	Ti6Al4V Alloy
1210312	Ti6Al4V Alloy
1210401	Ti6Al4V Alloy
1210402	Ti6Al4V Alloy
1210403	Ti6Al4V Alloy
1210404	Ti6Al4V Alloy
1210411	Ti6Al4V Alloy
1210412	Ti6Al4V Alloy
1210413	Ti6Al4V Alloy
1210414	Ti6Al4V Alloy
1210421	Ti6Al4V Alloy
1210422	Ti6Al4V Alloy
1210423	Ti6Al4V Alloy
1210424	Ti6Al4V Alloy
1210501	Ti6Al4V Alloy
1210502	Ti6Al4V Alloy
1210503	Ti6Al4V Alloy
1210504	Ti6Al4V Alloy
1210511	Ti6Al4V Alloy
1210512	Ti6Al4V Alloy
1210513	Ti6Al4V Alloy

1210514	Ti6Al4V Alloy
1210521	Ti6Al4V Alloy
1210522	Ti6Al4V Alloy
1210523	Ti6Al4V Alloy
1210524	Ti6Al4V Alloy
1210601	Ti6Al4V Alloy
1210602	Ti6Al4V Alloy
1210603	Ti6Al4V Alloy
1210604	Ti6Al4V Alloy
1210701	Ti6Al4V Alloy
1210702	Ti6Al4V Alloy
1210703	Ti6Al4V Alloy
1210704	Ti6Al4V Alloy
1210711	Ti6Al4V Alloy
1210712	Ti6Al4V Alloy
1210713	Ti6Al4V Alloy
1210714	Ti6Al4V Alloy
1210721	Ti6Al4V Alloy
1210722	Ti6Al4V Alloy
1210723	Ti6Al4V Alloy
1210724	Ti6Al4V Alloy
1210801	Ti6Al4V Alloy
1210802	Ti6Al4V Alloy
1210803	Ti6Al4V Alloy
1210804	Ti6Al4V Alloy
1210811	Ti6Al4V Alloy
1210812	Ti6Al4V Alloy
1210813	Ti6Al4V Alloy
1210814	Ti6Al4V Alloy
1210821	Ti6Al4V Alloy
1210822	Ti6Al4V Alloy
1210823	Ti6Al4V Alloy
1210824	Ti6Al4V Alloy
1210901	Ti6Al4V Alloy
1210902	Ti6Al4V Alloy
1210903	Ti6Al4V Alloy
1210904	Ti6Al4V Alloy
1211001	Co-Cr-Mo Alloy
1211002	Co-Cr-Mo Alloy
1211003	Co-Cr-Mo Alloy
1211021	Co-Cr-Mo Alloy
1211022	Co-Cr-Mo Alloy
1211023	Co-Cr-Mo Alloy

1211081	Co-Cr-Mo Alloy + TiNbN
1211082	Co-Cr-Mo Alloy + TiNbN
1211083	Co-Cr-Mo Alloy + TiNbN
1211091	Co-Cr-Mo Alloy + TiNbN
1211092	Co-Cr-Mo Alloy + TiNbN
1211093	Co-Cr-Mo Alloy + TiNbN
1211101	Ti6Al4V Alloy
1211102	Ti6Al4V Alloy
1211103	Ti6Al4V Alloy
1211210	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211220	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211230	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211240	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211250	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211260	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1211301	Ti6Al4V Alloy
1211302	Ti6Al4V Alloy
1211303	Ti6Al4V Alloy
1211304	Ti6Al4V Alloy
1211305	Ti6Al4V Alloy
1211306	Ti6Al4V Alloy
1211307	Ti6Al4V Alloy
1211308	Ti6Al4V Alloy
1211309	Ti6Al4V Alloy
1211310	Ti6Al4V Alloy
1211311	Ti6Al4V Alloy
1211312	Ti6Al4V Alloy
1211401	Ti6Al4V Alloy
1211402	Ti6Al4V Alloy
1211403	Ti6Al4V Alloy
1211411	Ti6Al4V Alloy
1211412	Ti6Al4V Alloy
1211413	Ti6Al4V Alloy
1211421	Ti6Al4V Alloy
1211422	Ti6Al4V Alloy
1211423	Ti6Al4V Alloy
1211510	Co-Cr-Mo Alloy
1211520	Co-Cr-Mo Alloy
1211530	Co-Cr-Mo Alloy
1211540	Co-Cr-Mo Alloy
1211550	Co-Cr-Mo Alloy
1211560	Co-Cr-Mo Alloy
1211610	Co-Cr-Mo Alloy + TiNbN

1211620	Co-Cr-Mo Alloy + TiNbN
1211630	Co-Cr-Mo Alloy + TiNbN
1211640	Co-Cr-Mo Alloy + TiNbN
1211650	Co-Cr-Mo Alloy + TiNbN
1211660	Co-Cr-Mo Alloy + TiNbN
1212210	Co-Cr-Mo Alloy + UHMWPE
1212220	Co-Cr-Mo Alloy + UHMWPE
1212230	Co-Cr-Mo Alloy + UHMWPE
1212240	Co-Cr-Mo Alloy + UHMWPE
1212250	Co-Cr-Mo Alloy + UHMWPE
1212260	Co-Cr-Mo Alloy + UHMWPE
1212327	UHMWPE + Ti6Al4V Alloy
1212330	UHMWPE + Ti6Al4V Alloy
1212333	UHMWPE + Ti6Al4V Alloy
1212336	UHMWPE + Ti6Al4V Alloy
1212339	UHMWPE + Ti6Al4V Alloy
1212801	Ti6Al4V Alloy
1212802	Ti6Al4V Alloy
1212803	Ti6Al4V Alloy
1212804	Ti6Al4V Alloy
1212811	Ti6Al4V Alloy
1212812	Ti6Al4V Alloy
1212813	Ti6Al4V Alloy
1212814	Ti6Al4V Alloy
1212821	Ti6Al4V Alloy
1212822	Ti6Al4V Alloy
1212823	Ti6Al4V Alloy
1212824	Ti6Al4V Alloy
1212901	Ti6Al4V Alloy
1212902	Ti6Al4V Alloy
1212903	Ti6Al4V Alloy
1212904	Ti6Al4V Alloy
1212911	Ti6Al4V Alloy
1212912	Ti6Al4V Alloy
1212913	Ti6Al4V Alloy
1212914	Ti6Al4V Alloy
1212921	Ti6Al4V Alloy
1212922	Ti6Al4V Alloy
1212923	Ti6Al4V Alloy
1212924	Ti6Al4V Alloy
1213101	Ti6Al4V Alloy
1213102	Ti6Al4V Alloy
1213103	Ti6Al4V Alloy

1213111	Ti6Al4V Alloy
1213112	Ti6Al4V Alloy
1213113	Ti6Al4V Alloy
1213121	Ti6Al4V Alloy
1213122	Ti6Al4V Alloy
1213123	Ti6Al4V Alloy
1213210	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213220	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213230	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213240	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213250	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213260	Co-Cr-Mo Alloy + UHMWPE + TiNbN
1213301	Ti6Al4V Alloy
1213302	Ti6Al4V Alloy
1213303	Ti6Al4V Alloy
1213401	Ti6Al4V Alloy
1213402	Ti6Al4V Alloy
1213403	Ti6Al4V Alloy
1213411	Ti6Al4V Alloy
1213412	Ti6Al4V Alloy
1213413	Ti6Al4V Alloy
1213421	Ti6Al4V Alloy
1213422	Ti6Al4V Alloy
1213423	Ti6Al4V Alloy
1213529	UHMWPE
1213532	UHMWPE
1213535	UHMWPE
1213538	UHMWPE
1213540	UHMWPE
1213627	UHMWPE
1213630	UHMWPE
1213633	UHMWPE
1213636	UHMWPE
1213639	UHMWPE
1213701	Ti6Al4V Alloy
1213702	Ti6Al4V Alloy
1213703	Ti6Al4V Alloy
1213704	Ti6Al4V Alloy
1213705	Ti6Al4V Alloy
1213706	Ti6Al4V Alloy
1213707	Ti6Al4V Alloy
1213708	Ti6Al4V Alloy
1213709	Ti6Al4V Alloy

1213710	Ti6Al4V Alloy
1213711	Ti6Al4V Alloy
1213712	Ti6Al4V Alloy
1213901	Ti6Al4V Alloy
1213902	Ti6Al4V Alloy
1213903	Ti6Al4V Alloy
1213904	Ti6Al4V Alloy
1213905	Ti6Al4V Alloy
1213906	Ti6Al4V Alloy
1213907	Ti6Al4V Alloy
1213908	Ti6Al4V Alloy
1213909	Ti6Al4V Alloy
1213910	Ti6Al4V Alloy
1213911	Ti6Al4V Alloy
1213912	Ti6Al4V Alloy
1214001	Ti6Al4V Alloy
1214002	Ti6Al4V Alloy
1214003	Ti6Al4V Alloy
1214011	Ti6Al4V Alloy
1214012	Ti6Al4V Alloy
1214013	Ti6Al4V Alloy
1214021	Ti6Al4V Alloy
1214022	Ti6Al4V Alloy
1214023	Ti6Al4V Alloy
1214031	Ti6Al4V Alloy
1214032	Ti6Al4V Alloy
1214033	Ti6Al4V Alloy
1214051	Ti6Al4V Alloy
1214052	Ti6Al4V Alloy
1214053	Ti6Al4V Alloy
1214061	Ti6Al4V Alloy
1214062	Ti6Al4V Alloy
1214063	Ti6Al4V Alloy
1214071	Ti6Al4V Alloy
1214072	Ti6Al4V Alloy
1214073	Ti6Al4V Alloy
1214081	Ti6Al4V Alloy
1214082	Ti6Al4V Alloy
1214083	Ti6Al4V Alloy
1215001	Ti6Al4V Alloy + Hydroxyapatite
1215002	Ti6Al4V Alloy + Hydroxyapatite
1215003	Ti6Al4V Alloy + Hydroxyapatite
1215011	Ti6Al4V Alloy + Hydroxyapatite

1215012	Ti6Al4V Alloy + Hydroxyapatite
1215013	Ti6Al4V Alloy + Hydroxyapatite
1215021	Ti6Al4V Alloy + Hydroxyapatite
1215022	Ti6Al4V Alloy + Hydroxyapatite
1215023	Ti6Al4V Alloy + Hydroxyapatite
1215031	Ti6Al4V Alloy + Hydroxyapatite
1215032	Ti6Al4V Alloy + Hydroxyapatite
1215033	Ti6Al4V Alloy + Hydroxyapatite
1215051	Ti6Al4V Alloy + Hydroxyapatite
1215052	Ti6Al4V Alloy + Hydroxyapatite
1215053	Ti6Al4V Alloy + Hydroxyapatite
1215061	Ti6Al4V Alloy + Hydroxyapatite
1215062	Ti6Al4V Alloy + Hydroxyapatite
1215063	Ti6Al4V Alloy + Hydroxyapatite
1215071	Ti6Al4V Alloy + Hydroxyapatite
1215072	Ti6Al4V Alloy + Hydroxyapatite
1215073	Ti6Al4V Alloy + Hydroxyapatite
1215081	Ti6Al4V Alloy + Hydroxyapatite
1215082	Ti6Al4V Alloy + Hydroxyapatite
1215083	Ti6Al4V Alloy + Hydroxyapatite
1216001	Co-Cr-Mo Alloy
1216002	Co-Cr-Mo Alloy
1216003	Co-Cr-Mo Alloy
1216011	Co-Cr-Mo Alloy
1216012	Co-Cr-Mo Alloy
1216013	Co-Cr-Mo Alloy
1216021	Co-Cr-Mo Alloy
1216022	Co-Cr-Mo Alloy
1216023	Co-Cr-Mo Alloy
1216031	Co-Cr-Mo Alloy
1216032	Co-Cr-Mo Alloy
1216033	Co-Cr-Mo Alloy
1216051	Co-Cr-Mo Alloy
1216052	Co-Cr-Mo Alloy
1216053	Co-Cr-Mo Alloy
1216061	Co-Cr-Mo Alloy
1216062	Co-Cr-Mo Alloy
1216063	Co-Cr-Mo Alloy
1216071	Co-Cr-Mo Alloy
1216072	Co-Cr-Mo Alloy
1216073	Co-Cr-Mo Alloy
1216081	Co-Cr-Mo Alloy
1216082	Co-Cr-Mo Alloy

1216083	Co-Cr-Mo Alloy
1217001	Co-Cr-Mo Alloy + Hydroxyapatite
1217002	Co-Cr-Mo Alloy + Hydroxyapatite
1217003	Co-Cr-Mo Alloy + Hydroxyapatite
1217011	Co-Cr-Mo Alloy + Hydroxyapatite
1217012	Co-Cr-Mo Alloy + Hydroxyapatite
1217013	Co-Cr-Mo Alloy + Hydroxyapatite
1217021	Co-Cr-Mo Alloy + Hydroxyapatite
1217022	Co-Cr-Mo Alloy + Hydroxyapatite
1217023	Co-Cr-Mo Alloy + Hydroxyapatite
1217031	Co-Cr-Mo Alloy + Hydroxyapatite
1217032	Co-Cr-Mo Alloy + Hydroxyapatite
1217033	Co-Cr-Mo Alloy + Hydroxyapatite
1217051	Co-Cr-Mo Alloy + Hydroxyapatite
1217052	Co-Cr-Mo Alloy + Hydroxyapatite
1217053	Co-Cr-Mo Alloy + Hydroxyapatite
1217061	Co-Cr-Mo Alloy + Hydroxyapatite
1217062	Co-Cr-Mo Alloy + Hydroxyapatite
1217063	Co-Cr-Mo Alloy + Hydroxyapatite
1217071	Co-Cr-Mo Alloy + Hydroxyapatite
1217072	Co-Cr-Mo Alloy + Hydroxyapatite
1217073	Co-Cr-Mo Alloy + Hydroxyapatite
1217081	Co-Cr-Mo Alloy + Hydroxyapatite
1217082	Co-Cr-Mo Alloy + Hydroxyapatite
1217083	Co-Cr-Mo Alloy + Hydroxyapatite
1218000	Ti6Al4V Alloy
1218050	Ti6Al4V Alloy
1218055	Ti6Al4V Alloy
1218060	Ti6Al4V Alloy
1218065	Ti6Al4V Alloy
1218500	Ti6Al4V Alloy
1219001	Ti6Al4V Alloy
1219002	Ti6Al4V Alloy
1219003	Ti6Al4V Alloy
1219004	Ti6Al4V Alloy
1219005	Ti6Al4V Alloy
1219006	Ti6Al4V Alloy
1219007	Ti6Al4V Alloy
1219008	Ti6Al4V Alloy
1219009	Ti6Al4V Alloy
1219010	Ti6Al4V Alloy
1219011	Ti6Al4V Alloy
1219012	Ti6Al4V Alloy

1220001	Ti6Al4V Alloy
1220002	Ti6Al4V Alloy
1220003	Ti6Al4V Alloy
1220011	Ti6Al4V Alloy
1220012	Ti6Al4V Alloy
1220013	Ti6Al4V Alloy
1220031	Ti6Al4V Alloy
1220032	Ti6Al4V Alloy
1220033	Ti6Al4V Alloy
1220101	Ti6Al4V Alloy
1220102	Ti6Al4V Alloy
1220103	Ti6Al4V Alloy
1220131	Ti6Al4V Alloy
1220132	Ti6Al4V Alloy
1220133	Ti6Al4V Alloy
1220201	Ti6Al4V Alloy
1220202	Ti6Al4V Alloy
1220203	Ti6Al4V Alloy
1220204	Ti6Al4V Alloy
1220205	Ti6Al4V Alloy
1220206	Ti6Al4V Alloy
1220207	Ti6Al4V Alloy
1220208	Ti6Al4V Alloy
1220209	Ti6Al4V Alloy
1220210	Ti6Al4V Alloy
1220211	Ti6Al4V Alloy
1220212	Ti6Al4V Alloy
1220301	Ti6Al4V Alloy
1220302	Ti6Al4V Alloy
1220303	Ti6Al4V Alloy
1220304	Ti6Al4V Alloy
1220305	Ti6Al4V Alloy
1220306	Ti6Al4V Alloy
1220307	Ti6Al4V Alloy
1220308	Ti6Al4V Alloy
1220309	Ti6Al4V Alloy
1220310	Ti6Al4V Alloy
1220311	Ti6Al4V Alloy
1220312	Ti6Al4V Alloy
1220401	Ti6Al4V Alloy
1220402	Ti6Al4V Alloy
1220403	Ti6Al4V Alloy
1220404	Ti6Al4V Alloy

1220405	Ti6Al4V Alloy
1220406	Ti6Al4V Alloy
1220407	Ti6Al4V Alloy
1220408	Ti6Al4V Alloy
1220409	Ti6Al4V Alloy
1220410	Ti6Al4V Alloy
1220411	Ti6Al4V Alloy
1220412	Ti6Al4V Alloy
1220501	Ti6Al4V Alloy
1220502	Ti6Al4V Alloy
1220503	Ti6Al4V Alloy
1220504	Ti6Al4V Alloy
1220505	Ti6Al4V Alloy
1220506	Ti6Al4V Alloy
1220507	Ti6Al4V Alloy
1220508	Ti6Al4V Alloy
1220509	Ti6Al4V Alloy
1220510	Ti6Al4V Alloy
1220511	Ti6Al4V Alloy
1220512	Ti6Al4V Alloy
1220601	Ti6Al4V Alloy
1220602	Ti6Al4V Alloy
1220603	Ti6Al4V Alloy
1220604	Ti6Al4V Alloy
1220605	Ti6Al4V Alloy
1220606	Ti6Al4V Alloy
1220607	Ti6Al4V Alloy
1220608	Ti6Al4V Alloy
1220609	Ti6Al4V Alloy
1220610	Ti6Al4V Alloy
1220611	Ti6Al4V Alloy
1220612	Ti6Al4V Alloy
1220701	Ti6Al4V Alloy
1220702	Ti6Al4V Alloy
1220703	Ti6Al4V Alloy
1220704	Ti6Al4V Alloy
1220705	Ti6Al4V Alloy
1220706	Ti6Al4V Alloy
1220707	Ti6Al4V Alloy
1220708	Ti6Al4V Alloy
1220709	Ti6Al4V Alloy
1220710	Ti6Al4V Alloy
1220711	Ti6Al4V Alloy

1220712	Ti6Al4V Alloy
1220801	Ti6Al4V Alloy
1220802	Ti6Al4V Alloy
1220803	Ti6Al4V Alloy
1220804	Ti6Al4V Alloy
1220805	Ti6Al4V Alloy
1220806	Ti6Al4V Alloy
1220807	Ti6Al4V Alloy
1220808	Ti6Al4V Alloy
1220809	Ti6Al4V Alloy
1220810	Ti6Al4V Alloy
1220811	Ti6Al4V Alloy
1220812	Ti6Al4V Alloy
1220901	Ti6Al4V Alloy
1220902	Ti6Al4V Alloy
1220903	Ti6Al4V Alloy
1220904	Ti6Al4V Alloy
1220905	Ti6Al4V Alloy
1220906	Ti6Al4V Alloy
1220907	Ti6Al4V Alloy
1220908	Ti6Al4V Alloy
1220909	Ti6Al4V Alloy
1220910	Ti6Al4V Alloy
1220911	Ti6Al4V Alloy
1220912	Ti6Al4V Alloy
1221001	Ti6Al4V Alloy
1221002	Ti6Al4V Alloy
1221003	Ti6Al4V Alloy
1221004	Ti6Al4V Alloy
1221005	Ti6Al4V Alloy
1221006	Ti6Al4V Alloy
1221007	Ti6Al4V Alloy
1221008	Ti6Al4V Alloy
1221009	Ti6Al4V Alloy
1221010	Ti6Al4V Alloy
1221011	Ti6Al4V Alloy
1221012	Ti6Al4V Alloy
1221101	Ti6Al4V Alloy
1221102	Ti6Al4V Alloy
1221103	Ti6Al4V Alloy
1221104	Ti6Al4V Alloy
1221105	Ti6Al4V Alloy
1221106	Ti6Al4V Alloy

1221107	Ti6Al4V Alloy
1221108	Ti6Al4V Alloy
1221109	Ti6Al4V Alloy
1221110	Ti6Al4V Alloy
1221111	Ti6Al4V Alloy
1221112	Ti6Al4V Alloy
1221201	Ti6Al4V Alloy
1221202	Ti6Al4V Alloy
1221203	Ti6Al4V Alloy
1221204	Ti6Al4V Alloy
1221205	Ti6Al4V Alloy
1221206	Ti6Al4V Alloy
1221207	Ti6Al4V Alloy
1221208	Ti6Al4V Alloy
1221209	Ti6Al4V Alloy
1221210	Ti6Al4V Alloy
1221211	Ti6Al4V Alloy
1221212	Ti6Al4V Alloy
1221301	Ti6Al4V Alloy
1221302	Ti6Al4V Alloy
1221303	Ti6Al4V Alloy
1221304	Ti6Al4V Alloy
1221305	Ti6Al4V Alloy
1221306	Ti6Al4V Alloy
1221307	Ti6Al4V Alloy
1221308	Ti6Al4V Alloy
1221309	Ti6Al4V Alloy
1221310	Ti6Al4V Alloy
1221311	Ti6Al4V Alloy
1221312	Ti6Al4V Alloy
1221401	Ti6Al4V Alloy
1221402	Ti6Al4V Alloy
1221403	Ti6Al4V Alloy
1221404	Ti6Al4V Alloy
1221405	Ti6Al4V Alloy
1221406	Ti6Al4V Alloy
1221407	Ti6Al4V Alloy
1221408	Ti6Al4V Alloy
1221409	Ti6Al4V Alloy
1221410	Ti6Al4V Alloy
1221411	Ti6Al4V Alloy
1221412	Ti6Al4V Alloy
1221501	Ti6Al4V Alloy

1221502	Ti6Al4V Alloy
1221503	Ti6Al4V Alloy
1221504	Ti6Al4V Alloy
1221505	Ti6Al4V Alloy
1221506	Ti6Al4V Alloy
1221507	Ti6Al4V Alloy
1221508	Ti6Al4V Alloy
1221509	Ti6Al4V Alloy
1221510	Ti6Al4V Alloy
1221511	Ti6Al4V Alloy
1221512	Ti6Al4V Alloy
1221601	Ti6Al4V Alloy
1221602	Ti6Al4V Alloy
1221603	Ti6Al4V Alloy
1221604	Ti6Al4V Alloy
1221605	Ti6Al4V Alloy
1221606	Ti6Al4V Alloy
1221607	Ti6Al4V Alloy
1221608	Ti6Al4V Alloy
1221609	Ti6Al4V Alloy
1221610	Ti6Al4V Alloy
1221611	Ti6Al4V Alloy
1221612	Ti6Al4V Alloy
1221701	Ti6Al4V Alloy
1221702	Ti6Al4V Alloy
1221703	Ti6Al4V Alloy
1221704	Ti6Al4V Alloy
1221705	Ti6Al4V Alloy
1221706	Ti6Al4V Alloy
1221707	Ti6Al4V Alloy
1221708	Ti6Al4V Alloy
1221709	Ti6Al4V Alloy
1221710	Ti6Al4V Alloy
1221711	Ti6Al4V Alloy
1221712	Ti6Al4V Alloy
1221801	Ti6Al4V Alloy
1221802	Ti6Al4V Alloy
1221803	Ti6Al4V Alloy
1221804	Ti6Al4V Alloy
1221805	Ti6Al4V Alloy
1221806	Ti6Al4V Alloy
1221807	Ti6Al4V Alloy
1221808	Ti6Al4V Alloy

1221809	Ti6Al4V Alloy
1221810	Ti6Al4V Alloy
1221811	Ti6Al4V Alloy
1221812	Ti6Al4V Alloy
1221901	Ti6Al4V Alloy
1221902	Ti6Al4V Alloy
1221903	Ti6Al4V Alloy
1221904	Ti6Al4V Alloy
1221905	Ti6Al4V Alloy
1221906	Ti6Al4V Alloy
1221907	Ti6Al4V Alloy
1221908	Ti6Al4V Alloy
1221909	Ti6Al4V Alloy
1221910	Ti6Al4V Alloy
1221911	Ti6Al4V Alloy
1221912	Ti6Al4V Alloy
1222001	Ti6Al4V Alloy
1222002	Ti6Al4V Alloy
1222003	Ti6Al4V Alloy
1222004	Ti6Al4V Alloy
1222005	Ti6Al4V Alloy
1222006	Ti6Al4V Alloy
1222007	Ti6Al4V Alloy
1222008	Ti6Al4V Alloy
1222009	Ti6Al4V Alloy
1222010	Ti6Al4V Alloy
1222011	Ti6Al4V Alloy
1222012	Ti6Al4V Alloy
1222101	Ti6Al4V Alloy
1222102	Ti6Al4V Alloy
1222103	Ti6Al4V Alloy
1222104	Ti6Al4V Alloy
1222105	Ti6Al4V Alloy
1222106	Ti6Al4V Alloy
1222107	Ti6Al4V Alloy
1222108	Ti6Al4V Alloy
1222109	Ti6Al4V Alloy
1222110	Ti6Al4V Alloy
1222111	Ti6Al4V Alloy
1222112	Ti6Al4V Alloy
1226001	Co-Cr-Mo Alloy
1226002	Co-Cr-Mo Alloy
1226003	Co-Cr-Mo Alloy

1226021	Co-Cr-Mo Alloy
1226022	Co-Cr-Mo Alloy
1226023	Co-Cr-Mo Alloy
1227001	Co-Cr-Mo Alloy + Hydroxyapatite
1227002	Co-Cr-Mo Alloy + Hydroxyapatite
1227003	Co-Cr-Mo Alloy + Hydroxyapatite
1227021	Co-Cr-Mo Alloy + Hydroxyapatite
1227022	Co-Cr-Mo Alloy + Hydroxyapatite
1227023	Co-Cr-Mo Alloy + Hydroxyapatite
1228001	Co-Cr-Mo Alloy
1228002	Co-Cr-Mo Alloy
1228003	Co-Cr-Mo Alloy
1228021	Co-Cr-Mo Alloy
1228022	Co-Cr-Mo Alloy
1228023	Co-Cr-Mo Alloy
1229001	Co-Cr-Mo Alloy + Hydroxyapatite
1229002	Co-Cr-Mo Alloy + Hydroxyapatite
1229003	Co-Cr-Mo Alloy + Hydroxyapatite
1229021	Co-Cr-Mo Alloy + Hydroxyapatite
1229022	Co-Cr-Mo Alloy + Hydroxyapatite
1229023	Co-Cr-Mo Alloy + Hydroxyapatite
1301000	Ti6Al4V Alloy
1301001	Ti6Al4V Alloy
1301002	Ti6Al4V Alloy
1301003	Ti6Al4V Alloy
1301004	Ti6Al4V Alloy
1301100	Ti6Al4V Alloy
1301101	Ti6Al4V Alloy
1301102	Ti6Al4V Alloy
1301103	Ti6Al4V Alloy
1301104	Ti6Al4V Alloy
1301200	Ti6Al4V Alloy
1301201	Ti6Al4V Alloy
1301202	Ti6Al4V Alloy
1301203	Ti6Al4V Alloy
1301204	Ti6Al4V Alloy
1301300	Ti6Al4V Alloy
1301301	Ti6Al4V Alloy
1301302	Ti6Al4V Alloy
1301303	Ti6Al4V Alloy
1301304	Ti6Al4V Alloy
1301400	Ti6Al4V Alloy
1301401	Ti6Al4V Alloy

1301402	Ti6Al4V Alloy
1301403	Ti6Al4V Alloy
1301404	Ti6Al4V Alloy
1301500	Ti6Al4V Alloy
1301501	Ti6Al4V Alloy
1301502	Ti6Al4V Alloy
1301503	Ti6Al4V Alloy
1301504	Ti6Al4V Alloy
1301600	Ti6Al4V Alloy
1301601	Ti6Al4V Alloy
1301602	Ti6Al4V Alloy
1301603	Ti6Al4V Alloy
1301604	Ti6Al4V Alloy
1301800	Ti6Al4V Alloy
1301801	Ti6Al4V Alloy
1301802	Ti6Al4V Alloy
1301803	Ti6Al4V Alloy
1301804	Ti6Al4V Alloy
1302000	Ti6Al4V Alloy
1302001	Ti6Al4V Alloy
1302002	Ti6Al4V Alloy
1302003	Ti6Al4V Alloy
1302004	Ti6Al4V Alloy
1302200	Ti6Al4V Alloy
1302201	Ti6Al4V Alloy
1302202	Ti6Al4V Alloy
1302203	Ti6Al4V Alloy
1302204	Ti6Al4V Alloy
1302400	Ti6Al4V Alloy
1302401	Ti6Al4V Alloy
1302402	Ti6Al4V Alloy
1302403	Ti6Al4V Alloy
1302404	Ti6Al4V Alloy
1311000	Ti6Al4V Alloy
1311001	Ti6Al4V Alloy
1311002	Ti6Al4V Alloy
1311003	Ti6Al4V Alloy
1311004	Ti6Al4V Alloy
1311100	Ti6Al4V Alloy
1311101	Ti6Al4V Alloy
1311102	Ti6Al4V Alloy
1311103	Ti6Al4V Alloy
1311104	Ti6Al4V Alloy

1311200	Ti6Al4V Alloy
1311201	Ti6Al4V Alloy
1311202	Ti6Al4V Alloy
1311203	Ti6Al4V Alloy
1311204	Ti6Al4V Alloy
1311300	Ti6Al4V Alloy
1311301	Ti6Al4V Alloy
1311302	Ti6Al4V Alloy
1311303	Ti6Al4V Alloy
1311304	Ti6Al4V Alloy
1311400	Ti6Al4V Alloy
1311401	Ti6Al4V Alloy
1311402	Ti6Al4V Alloy
1311403	Ti6Al4V Alloy
1311404	Ti6Al4V Alloy
1311500	Ti6Al4V Alloy
1311501	Ti6Al4V Alloy
1311502	Ti6Al4V Alloy
1311503	Ti6Al4V Alloy
1311504	Ti6Al4V Alloy
1311600	Ti6Al4V Alloy
1311601	Ti6Al4V Alloy
1311602	Ti6Al4V Alloy
1311603	Ti6Al4V Alloy
1311604	Ti6Al4V Alloy
1311800	Ti6Al4V Alloy
1311801	Ti6Al4V Alloy
1311802	Ti6Al4V Alloy
1311803	Ti6Al4V Alloy
1311804	Ti6Al4V Alloy
1312000	Ti6Al4V Alloy
1312001	Ti6Al4V Alloy
1312002	Ti6Al4V Alloy
1312003	Ti6Al4V Alloy
1312004	Ti6Al4V Alloy
1312200	Ti6Al4V Alloy
1312201	Ti6Al4V Alloy
1312202	Ti6Al4V Alloy
1312203	Ti6Al4V Alloy
1312204	Ti6Al4V Alloy
1312400	Ti6Al4V Alloy
1312401	Ti6Al4V Alloy
1312402	Ti6Al4V Alloy

1312403	Ti6Al4V Alloy
1312404	Ti6Al4V Alloy
1381000	Ti6Al4V Alloy
1381001	Ti6Al4V Alloy
1381002	Ti6Al4V Alloy
1381003	Ti6Al4V Alloy
1381004	Ti6Al4V Alloy
1381200	Ti6Al4V Alloy
1381201	Ti6Al4V Alloy
1381202	Ti6Al4V Alloy
1381203	Ti6Al4V Alloy
1381204	Ti6Al4V Alloy
1381400	Ti6Al4V Alloy
1381401	Ti6Al4V Alloy
1381402	Ti6Al4V Alloy
1381403	Ti6Al4V Alloy
1381404	Ti6Al4V Alloy
1381600	Ti6Al4V Alloy
1381601	Ti6Al4V Alloy
1381602	Ti6Al4V Alloy
1381603	Ti6Al4V Alloy
1381604	Ti6Al4V Alloy
1381800	Ti6Al4V Alloy
1381801	Ti6Al4V Alloy
1381802	Ti6Al4V Alloy
1381803	Ti6Al4V Alloy
1381804	Ti6Al4V Alloy
1382000	Ti6Al4V Alloy
1382001	Ti6Al4V Alloy
1382002	Ti6Al4V Alloy
1382003	Ti6Al4V Alloy
1382004	Ti6Al4V Alloy
1400120	UHMWPE
1400123	UHMWPE
1400126	UHMWPE
1400327	UHMWPE
1400330	UHMWPE
1400333	UHMWPE
1400336	UHMWPE
1400339	UHMWPE
1401120	UHMWPE
1401123	UHMWPE
1401126	UHMWPE

1480027	UHMWPE
1480030	UHMWPE
1480033	UHMWPE
1480036	UHMWPE
1480039	UHMWPE
1480127	UHMWPE
1480130	UHMWPE
1480133	UHMWPE
1480136	UHMWPE
1480139	UHMWPE
1480220	UHMWPE
1480223	UHMWPE
1480226	UHMWPE
1611001	Co-Cr-Mo Alloy + Hydroxyapatite
1611002	Co-Cr-Mo Alloy + Hydroxyapatite
1611003	Co-Cr-Mo Alloy + Hydroxyapatite
1611004	Co-Cr-Mo Alloy + Hydroxyapatite
1611005	Co-Cr-Mo Alloy + Hydroxyapatite
1611006	Co-Cr-Mo Alloy + Hydroxyapatite
1611007	Co-Cr-Mo Alloy + Hydroxyapatite
1611008	Co-Cr-Mo Alloy + Hydroxyapatite
1611009	Co-Cr-Mo Alloy + Hydroxyapatite
1611010	Co-Cr-Mo Alloy + Hydroxyapatite
1611011	Co-Cr-Mo Alloy + Hydroxyapatite
1611021	Co-Cr-Mo Alloy + Hydroxyapatite
1611022	Co-Cr-Mo Alloy + Hydroxyapatite
1611023	Co-Cr-Mo Alloy + Hydroxyapatite
1611024	Co-Cr-Mo Alloy + Hydroxyapatite
1611025	Co-Cr-Mo Alloy + Hydroxyapatite
1611026	Co-Cr-Mo Alloy + Hydroxyapatite
1611027	Co-Cr-Mo Alloy + Hydroxyapatite
1611028	Co-Cr-Mo Alloy + Hydroxyapatite
1611029	Co-Cr-Mo Alloy + Hydroxyapatite
1611030	Co-Cr-Mo Alloy + Hydroxyapatite
1611031	Co-Cr-Mo Alloy + Hydroxyapatite
1612001	Co-Cr-Mo Alloy
1612002	Co-Cr-Mo Alloy
1612003	Co-Cr-Mo Alloy
1612004	Co-Cr-Mo Alloy
1612005	Co-Cr-Mo Alloy
1612006	Co-Cr-Mo Alloy
1612007	Co-Cr-Mo Alloy
1612008	Co-Cr-Mo Alloy

1612009	Co-Cr-Mo Alloy
1612010	Co-Cr-Mo Alloy
1612011	Co-Cr-Mo Alloy
1612021	Co-Cr-Mo Alloy
1612022	Co-Cr-Mo Alloy
1612023	Co-Cr-Mo Alloy
1612024	Co-Cr-Mo Alloy
1612025	Co-Cr-Mo Alloy
1612026	Co-Cr-Mo Alloy
1612027	Co-Cr-Mo Alloy
1612028	Co-Cr-Mo Alloy
1612029	Co-Cr-Mo Alloy
1612030	Co-Cr-Mo Alloy
1612031	Co-Cr-Mo Alloy
1613001	Co-Cr-Mo Alloy + Hydroxyapatite
1613002	Co-Cr-Mo Alloy + Hydroxyapatite
1613003	Co-Cr-Mo Alloy + Hydroxyapatite
1613004	Co-Cr-Mo Alloy + Hydroxyapatite
1613005	Co-Cr-Mo Alloy + Hydroxyapatite
1613006	Co-Cr-Mo Alloy + Hydroxyapatite
1613007	Co-Cr-Mo Alloy + Hydroxyapatite
1613008	Co-Cr-Mo Alloy + Hydroxyapatite
1613009	Co-Cr-Mo Alloy + Hydroxyapatite
1613010	Co-Cr-Mo Alloy + Hydroxyapatite
1613011	Co-Cr-Mo Alloy + Hydroxyapatite
1613021	Co-Cr-Mo Alloy + Hydroxyapatite
1613022	Co-Cr-Mo Alloy + Hydroxyapatite
1613023	Co-Cr-Mo Alloy + Hydroxyapatite
1613024	Co-Cr-Mo Alloy + Hydroxyapatite
1613025	Co-Cr-Mo Alloy + Hydroxyapatite
1613026	Co-Cr-Mo Alloy + Hydroxyapatite
1613027	Co-Cr-Mo Alloy + Hydroxyapatite
1613028	Co-Cr-Mo Alloy + Hydroxyapatite
1613029	Co-Cr-Mo Alloy + Hydroxyapatite
1613030	Co-Cr-Mo Alloy + Hydroxyapatite
1613031	Co-Cr-Mo Alloy + Hydroxyapatite
1614001	Co-Cr-Mo Alloy
1614002	Co-Cr-Mo Alloy
1614003	Co-Cr-Mo Alloy
1614004	Co-Cr-Mo Alloy
1614005	Co-Cr-Mo Alloy
1614006	Co-Cr-Mo Alloy
1614007	Co-Cr-Mo Alloy

1614008	Co-Cr-Mo Alloy
1614009	Co-Cr-Mo Alloy
1614010	Co-Cr-Mo Alloy
1614011	Co-Cr-Mo Alloy
1614021	Co-Cr-Mo Alloy
1614022	Co-Cr-Mo Alloy
1614023	Co-Cr-Mo Alloy
1614024	Co-Cr-Mo Alloy
1614025	Co-Cr-Mo Alloy
1614026	Co-Cr-Mo Alloy
1614027	Co-Cr-Mo Alloy
1614028	Co-Cr-Mo Alloy
1614029	Co-Cr-Mo Alloy
1614030	Co-Cr-Mo Alloy
1614031	Co-Cr-Mo Alloy
1615001	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615002	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615003	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615004	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615005	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615006	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615007	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615008	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615009	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615010	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615011	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615021	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615022	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615023	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615024	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615025	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615026	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615027	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615028	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615029	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615030	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1615031	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1616001	Ti6Al4V Alloy + TiNbN
1616002	Ti6Al4V Alloy + TiNbN
1616003	Ti6Al4V Alloy + TiNbN
1616004	Ti6Al4V Alloy + TiNbN
1616005	Ti6Al4V Alloy + TiNbN
1616006	Ti6Al4V Alloy + TiNbN

1616007	Ti6Al4V Alloy + TiNbN
1616008	Ti6Al4V Alloy + TiNbN
1616009	Ti6Al4V Alloy + TiNbN
1616010	Ti6Al4V Alloy + TiNbN
1616011	Ti6Al4V Alloy + TiNbN
1616021	Ti6Al4V Alloy + TiNbN
1616022	Ti6Al4V Alloy + TiNbN
1616023	Ti6Al4V Alloy + TiNbN
1616024	Ti6Al4V Alloy + TiNbN
1616025	Ti6Al4V Alloy + TiNbN
1616026	Ti6Al4V Alloy + TiNbN
1616027	Ti6Al4V Alloy + TiNbN
1616028	Ti6Al4V Alloy + TiNbN
1616029	Ti6Al4V Alloy + TiNbN
1616030	Ti6Al4V Alloy + TiNbN
1616031	Ti6Al4V Alloy + TiNbN
1618001	Co-Cr-Mo Alloy
1618002	Co-Cr-Mo Alloy
1618003	Co-Cr-Mo Alloy
1618004	Co-Cr-Mo Alloy
1618005	Co-Cr-Mo Alloy
1618006	Co-Cr-Mo Alloy
1618007	Co-Cr-Mo Alloy
1618008	Co-Cr-Mo Alloy
1618009	Co-Cr-Mo Alloy
1618010	Co-Cr-Mo Alloy
1618011	Co-Cr-Mo Alloy
1618021	Co-Cr-Mo Alloy
1618022	Co-Cr-Mo Alloy
1618023	Co-Cr-Mo Alloy
1618024	Co-Cr-Mo Alloy
1618025	Co-Cr-Mo Alloy
1618026	Co-Cr-Mo Alloy
1618027	Co-Cr-Mo Alloy
1618028	Co-Cr-Mo Alloy
1618029	Co-Cr-Mo Alloy
1618030	Co-Cr-Mo Alloy
1618031	Co-Cr-Mo Alloy
1618901	Co-Cr-Mo Alloy + TiNbN
1618902	Co-Cr-Mo Alloy + TiNbN
1618903	Co-Cr-Mo Alloy + TiNbN
1618904	Co-Cr-Mo Alloy + TiNbN
1618905	Co-Cr-Mo Alloy + TiNbN

1618906	Co-Cr-Mo Alloy + TiNbN
1618907	Co-Cr-Mo Alloy + TiNbN
1618908	Co-Cr-Mo Alloy + TiNbN
1618909	Co-Cr-Mo Alloy + TiNbN
1618910	Co-Cr-Mo Alloy + TiNbN
1618911	Co-Cr-Mo Alloy + TiNbN
1618921	Co-Cr-Mo Alloy + TiNbN
1618922	Co-Cr-Mo Alloy + TiNbN
1618923	Co-Cr-Mo Alloy + TiNbN
1618924	Co-Cr-Mo Alloy + TiNbN
1618925	Co-Cr-Mo Alloy + TiNbN
1618926	Co-Cr-Mo Alloy + TiNbN
1618927	Co-Cr-Mo Alloy + TiNbN
1618928	Co-Cr-Mo Alloy + TiNbN
1618929	Co-Cr-Mo Alloy + TiNbN
1618930	Co-Cr-Mo Alloy + TiNbN
1618931	Co-Cr-Mo Alloy + TiNbN
1621001	UHMWPE
1621002	UHMWPE
1621003	UHMWPE
1621004	UHMWPE
1621005	UHMWPE
1621006	UHMWPE
1621007	UHMWPE
1621008	UHMWPE
1621009	UHMWPE
1621010	UHMWPE
1621011	UHMWPE
1621021	UHMWPE
1621022	UHMWPE
1621023	UHMWPE
1621024	UHMWPE
1621025	UHMWPE
1621026	UHMWPE
1621027	UHMWPE
1621028	UHMWPE
1621029	UHMWPE
1621030	UHMWPE
1621031	UHMWPE
1621041	UHMWPE
1621042	UHMWPE
1621043	UHMWPE
1621044	UHMWPE

1621045	UHMWPE
1621046	UHMWPE
1621047	UHMWPE
1621048	UHMWPE
1621049	UHMWPE
1621050	UHMWPE
1621051	UHMWPE
1621101	UHMWPE
1621102	UHMWPE
1621103	UHMWPE
1621104	UHMWPE
1621105	UHMWPE
1621106	UHMWPE
1621107	UHMWPE
1621108	UHMWPE
1621109	UHMWPE
1621110	UHMWPE
1621111	UHMWPE
1621121	UHMWPE
1621122	UHMWPE
1621123	UHMWPE
1621124	UHMWPE
1621125	UHMWPE
1621126	UHMWPE
1621127	UHMWPE
1621128	UHMWPE
1621129	UHMWPE
1621130	UHMWPE
1621131	UHMWPE
1621141	UHMWPE
1621142	UHMWPE
1621143	UHMWPE
1621144	UHMWPE
1621145	UHMWPE
1621146	UHMWPE
1621147	UHMWPE
1621148	UHMWPE
1621149	UHMWPE
1621150	UHMWPE
1621151	UHMWPE
1621201	UHMWPE
1621202	UHMWPE
1621203	UHMWPE

1621204	UHMWPE
1621205	UHMWPE
1621206	UHMWPE
1621207	UHMWPE
1621208	UHMWPE
1621209	UHMWPE
1621210	UHMWPE
1621211	UHMWPE
1621221	UHMWPE
1621222	UHMWPE
1621223	UHMWPE
1621224	UHMWPE
1621225	UHMWPE
1621226	UHMWPE
1621227	UHMWPE
1621228	UHMWPE
1621229	UHMWPE
1621230	UHMWPE
1621231	UHMWPE
1621241	UHMWPE
1621242	UHMWPE
1621243	UHMWPE
1621244	UHMWPE
1621245	UHMWPE
1621246	UHMWPE
1621247	UHMWPE
1621248	UHMWPE
1621249	UHMWPE
1621250	UHMWPE
1621251	UHMWPE
1621301	UHMWPE
1621302	UHMWPE
1621303	UHMWPE
1621304	UHMWPE
1621305	UHMWPE
1621306	UHMWPE
1621307	UHMWPE
1621308	UHMWPE
1621309	UHMWPE
1621310	UHMWPE
1621311	UHMWPE
1621321	UHMWPE
1621322	UHMWPE

1621323	UHMWPE
1621324	UHMWPE
1621325	UHMWPE
1621326	UHMWPE
1621327	UHMWPE
1621328	UHMWPE
1621329	UHMWPE
1621330	UHMWPE
1621331	UHMWPE
1621341	UHMWPE
1621342	UHMWPE
1621343	UHMWPE
1621344	UHMWPE
1621345	UHMWPE
1621346	UHMWPE
1621347	UHMWPE
1621348	UHMWPE
1621349	UHMWPE
1621350	UHMWPE
1621351	UHMWPE
1621401	UHMWPE
1621402	UHMWPE
1621403	UHMWPE
1621404	UHMWPE
1621405	UHMWPE
1621406	UHMWPE
1621407	UHMWPE
1621408	UHMWPE
1621409	UHMWPE
1621410	UHMWPE
1621411	UHMWPE
1621421	UHMWPE
1621422	UHMWPE
1621423	UHMWPE
1621424	UHMWPE
1621425	UHMWPE
1621426	UHMWPE
1621427	UHMWPE
1621428	UHMWPE
1621429	UHMWPE
1621430	UHMWPE
1621431	UHMWPE
1621441	UHMWPE

1621442	UHMWPE
1621443	UHMWPE
1621444	UHMWPE
1621445	UHMWPE
1621446	UHMWPE
1621447	UHMWPE
1621448	UHMWPE
1621449	UHMWPE
1621450	UHMWPE
1621451	UHMWPE
1621501	UHMWPE
1621502	UHMWPE
1621503	UHMWPE
1621504	UHMWPE
1621505	UHMWPE
1621506	UHMWPE
1621507	UHMWPE
1621508	UHMWPE
1621509	UHMWPE
1621510	UHMWPE
1621511	UHMWPE
1621521	UHMWPE
1621522	UHMWPE
1621523	UHMWPE
1621524	UHMWPE
1621525	UHMWPE
1621526	UHMWPE
1621527	UHMWPE
1621528	UHMWPE
1621529	UHMWPE
1621530	UHMWPE
1621531	UHMWPE
1621541	UHMWPE
1621542	UHMWPE
1621543	UHMWPE
1621544	UHMWPE
1621545	UHMWPE
1621546	UHMWPE
1621547	UHMWPE
1621548	UHMWPE
1621549	UHMWPE
1621550	UHMWPE
1621551	UHMWPE

1621601	UHMWPE
1621602	UHMWPE
1621603	UHMWPE
1621604	UHMWPE
1621605	UHMWPE
1621606	UHMWPE
1621607	UHMWPE
1621608	UHMWPE
1621609	UHMWPE
1621610	UHMWPE
1621611	UHMWPE
1621621	UHMWPE
1621622	UHMWPE
1621623	UHMWPE
1621624	UHMWPE
1621625	UHMWPE
1621626	UHMWPE
1621627	UHMWPE
1621628	UHMWPE
1621629	UHMWPE
1621630	UHMWPE
1621631	UHMWPE
1621641	UHMWPE
1621642	UHMWPE
1621643	UHMWPE
1621644	UHMWPE
1621645	UHMWPE
1621646	UHMWPE
1621647	UHMWPE
1621648	UHMWPE
1621649	UHMWPE
1621650	UHMWPE
1621651	UHMWPE
1621701	UHMWPE
1621702	UHMWPE
1621703	UHMWPE
1621704	UHMWPE
1621705	UHMWPE
1621706	UHMWPE
1621707	UHMWPE
1621708	UHMWPE
1621709	UHMWPE
1621710	UHMWPE

1621711	UHMWPE
1621721	UHMWPE
1621722	UHMWPE
1621723	UHMWPE
1621724	UHMWPE
1621725	UHMWPE
1621726	UHMWPE
1621727	UHMWPE
1621728	UHMWPE
1621729	UHMWPE
1621730	UHMWPE
1621731	UHMWPE
1621741	UHMWPE
1621742	UHMWPE
1621743	UHMWPE
1621744	UHMWPE
1621745	UHMWPE
1621746	UHMWPE
1621747	UHMWPE
1621748	UHMWPE
1621749	UHMWPE
1621750	UHMWPE
1621751	UHMWPE
1622001	UHMWPE
1622002	UHMWPE
1622003	UHMWPE
1622004	UHMWPE
1622005	UHMWPE
1622006	UHMWPE
1622007	UHMWPE
1622008	UHMWPE
1622009	UHMWPE
1622010	UHMWPE
1622011	UHMWPE
1622021	UHMWPE
1622022	UHMWPE
1622023	UHMWPE
1622024	UHMWPE
1622025	UHMWPE
1622026	UHMWPE
1622027	UHMWPE
1622028	UHMWPE
1622029	UHMWPE

1622030	UHMWPE
1622031	UHMWPE
1622041	UHMWPE
1622042	UHMWPE
1622043	UHMWPE
1622044	UHMWPE
1622045	UHMWPE
1622046	UHMWPE
1622047	UHMWPE
1622048	UHMWPE
1622049	UHMWPE
1622050	UHMWPE
1622051	UHMWPE
1622401	UHMWPE
1622402	UHMWPE
1622403	UHMWPE
1622404	UHMWPE
1622405	UHMWPE
1622406	UHMWPE
1622407	UHMWPE
1622408	UHMWPE
1622409	UHMWPE
1622410	UHMWPE
1622411	UHMWPE
1622421	UHMWPE
1622422	UHMWPE
1622423	UHMWPE
1622424	UHMWPE
1622425	UHMWPE
1622426	UHMWPE
1622427	UHMWPE
1622428	UHMWPE
1622429	UHMWPE
1622430	UHMWPE
1622431	UHMWPE
1622441	UHMWPE
1622442	UHMWPE
1622443	UHMWPE
1622444	UHMWPE
1622445	UHMWPE
1622446	UHMWPE
1622447	UHMWPE
1622448	UHMWPE

1622449	UHMWPE
1622450	UHMWPE
1622451	UHMWPE
1631001	UHMWPE
1631002	UHMWPE
1631003	UHMWPE
1631004	UHMWPE
1631005	UHMWPE
1631006	UHMWPE
1631007	UHMWPE
1631008	UHMWPE
1631009	UHMWPE
1631010	UHMWPE
1631011	UHMWPE
1631021	UHMWPE
1631022	UHMWPE
1631023	UHMWPE
1631024	UHMWPE
1631025	UHMWPE
1631026	UHMWPE
1631027	UHMWPE
1631028	UHMWPE
1631029	UHMWPE
1631030	UHMWPE
1631031	UHMWPE
1631041	UHMWPE
1631042	UHMWPE
1631043	UHMWPE
1631044	UHMWPE
1631045	UHMWPE
1631046	UHMWPE
1631047	UHMWPE
1631048	UHMWPE
1631049	UHMWPE
1631050	UHMWPE
1631051	UHMWPE
1631101	UHMWPE
1631102	UHMWPE
1631103	UHMWPE
1631104	UHMWPE
1631105	UHMWPE
1631106	UHMWPE
1631107	UHMWPE

1631108	UHMWPE
1631109	UHMWPE
1631110	UHMWPE
1631111	UHMWPE
1631121	UHMWPE
1631122	UHMWPE
1631123	UHMWPE
1631124	UHMWPE
1631125	UHMWPE
1631126	UHMWPE
1631127	UHMWPE
1631128	UHMWPE
1631129	UHMWPE
1631130	UHMWPE
1631131	UHMWPE
1631141	UHMWPE
1631142	UHMWPE
1631143	UHMWPE
1631144	UHMWPE
1631145	UHMWPE
1631146	UHMWPE
1631147	UHMWPE
1631148	UHMWPE
1631149	UHMWPE
1631150	UHMWPE
1631151	UHMWPE
1631201	UHMWPE
1631202	UHMWPE
1631203	UHMWPE
1631204	UHMWPE
1631205	UHMWPE
1631206	UHMWPE
1631207	UHMWPE
1631208	UHMWPE
1631209	UHMWPE
1631210	UHMWPE
1631211	UHMWPE
1631221	UHMWPE
1631222	UHMWPE
1631223	UHMWPE
1631224	UHMWPE
1631225	UHMWPE
1631226	UHMWPE

1631227	UHMWPE
1631228	UHMWPE
1631229	UHMWPE
1631230	UHMWPE
1631231	UHMWPE
1631241	UHMWPE
1631242	UHMWPE
1631243	UHMWPE
1631244	UHMWPE
1631245	UHMWPE
1631246	UHMWPE
1631247	UHMWPE
1631248	UHMWPE
1631249	UHMWPE
1631250	UHMWPE
1631251	UHMWPE
1631301	UHMWPE
1631302	UHMWPE
1631303	UHMWPE
1631304	UHMWPE
1631305	UHMWPE
1631306	UHMWPE
1631307	UHMWPE
1631308	UHMWPE
1631309	UHMWPE
1631310	UHMWPE
1631311	UHMWPE
1631321	UHMWPE
1631322	UHMWPE
1631323	UHMWPE
1631324	UHMWPE
1631325	UHMWPE
1631326	UHMWPE
1631327	UHMWPE
1631328	UHMWPE
1631329	UHMWPE
1631330	UHMWPE
1631331	UHMWPE
1631341	UHMWPE
1631342	UHMWPE
1631343	UHMWPE
1631344	UHMWPE
1631345	UHMWPE

1631346	UHMWPE
1631347	UHMWPE
1631348	UHMWPE
1631349	UHMWPE
1631350	UHMWPE
1631351	UHMWPE
1631401	UHMWPE
1631402	UHMWPE
1631403	UHMWPE
1631404	UHMWPE
1631405	UHMWPE
1631406	UHMWPE
1631407	UHMWPE
1631408	UHMWPE
1631409	UHMWPE
1631410	UHMWPE
1631411	UHMWPE
1631421	UHMWPE
1631422	UHMWPE
1631423	UHMWPE
1631424	UHMWPE
1631425	UHMWPE
1631426	UHMWPE
1631427	UHMWPE
1631428	UHMWPE
1631429	UHMWPE
1631430	UHMWPE
1631431	UHMWPE
1631441	UHMWPE
1631442	UHMWPE
1631443	UHMWPE
1631444	UHMWPE
1631445	UHMWPE
1631446	UHMWPE
1631447	UHMWPE
1631448	UHMWPE
1631449	UHMWPE
1631450	UHMWPE
1631451	UHMWPE
1631501	UHMWPE
1631502	UHMWPE
1631503	UHMWPE
1631504	UHMWPE

1631505	UHMWPE
1631506	UHMWPE
1631507	UHMWPE
1631508	UHMWPE
1631509	UHMWPE
1631510	UHMWPE
1631511	UHMWPE
1631521	UHMWPE
1631522	UHMWPE
1631523	UHMWPE
1631524	UHMWPE
1631525	UHMWPE
1631526	UHMWPE
1631527	UHMWPE
1631528	UHMWPE
1631529	UHMWPE
1631530	UHMWPE
1631531	UHMWPE
1631541	UHMWPE
1631542	UHMWPE
1631543	UHMWPE
1631544	UHMWPE
1631545	UHMWPE
1631546	UHMWPE
1631547	UHMWPE
1631548	UHMWPE
1631549	UHMWPE
1631550	UHMWPE
1631551	UHMWPE
1631601	UHMWPE
1631602	UHMWPE
1631603	UHMWPE
1631604	UHMWPE
1631605	UHMWPE
1631606	UHMWPE
1631607	UHMWPE
1631608	UHMWPE
1631609	UHMWPE
1631610	UHMWPE
1631611	UHMWPE
1631621	UHMWPE
1631622	UHMWPE
1631623	UHMWPE

1631624	UHMWPE
1631625	UHMWPE
1631626	UHMWPE
1631627	UHMWPE
1631628	UHMWPE
1631629	UHMWPE
1631630	UHMWPE
1631631	UHMWPE
1631641	UHMWPE
1631642	UHMWPE
1631643	UHMWPE
1631644	UHMWPE
1631645	UHMWPE
1631646	UHMWPE
1631647	UHMWPE
1631648	UHMWPE
1631649	UHMWPE
1631650	UHMWPE
1631651	UHMWPE
1631701	UHMWPE
1631702	UHMWPE
1631703	UHMWPE
1631704	UHMWPE
1631705	UHMWPE
1631706	UHMWPE
1631707	UHMWPE
1631708	UHMWPE
1631709	UHMWPE
1631710	UHMWPE
1631711	UHMWPE
1631721	UHMWPE
1631722	UHMWPE
1631723	UHMWPE
1631724	UHMWPE
1631725	UHMWPE
1631726	UHMWPE
1631727	UHMWPE
1631728	UHMWPE
1631729	UHMWPE
1631730	UHMWPE
1631731	UHMWPE
1631741	UHMWPE
1631742	UHMWPE

1631743	UHMWPE
1631744	UHMWPE
1631745	UHMWPE
1631746	UHMWPE
1631747	UHMWPE
1631748	UHMWPE
1631749	UHMWPE
1631750	UHMWPE
1631751	UHMWPE
1632001	UHMWPE
1632002	UHMWPE
1632003	UHMWPE
1632004	UHMWPE
1632005	UHMWPE
1632006	UHMWPE
1632007	UHMWPE
1632008	UHMWPE
1632009	UHMWPE
1632010	UHMWPE
1632011	UHMWPE
1632021	UHMWPE
1632022	UHMWPE
1632023	UHMWPE
1632024	UHMWPE
1632025	UHMWPE
1632026	UHMWPE
1632027	UHMWPE
1632028	UHMWPE
1632029	UHMWPE
1632030	UHMWPE
1632031	UHMWPE
1632041	UHMWPE
1632042	UHMWPE
1632043	UHMWPE
1632044	UHMWPE
1632045	UHMWPE
1632046	UHMWPE
1632047	UHMWPE
1632048	UHMWPE
1632049	UHMWPE
1632050	UHMWPE
1632051	UHMWPE
1632401	UHMWPE

1632402	UHMWPE
1632403	UHMWPE
1632404	UHMWPE
1632405	UHMWPE
1632406	UHMWPE
1632407	UHMWPE
1632408	UHMWPE
1632409	UHMWPE
1632410	UHMWPE
1632411	UHMWPE
1632421	UHMWPE
1632422	UHMWPE
1632423	UHMWPE
1632424	UHMWPE
1632425	UHMWPE
1632426	UHMWPE
1632427	UHMWPE
1632428	UHMWPE
1632429	UHMWPE
1632430	UHMWPE
1632431	UHMWPE
1632441	UHMWPE
1632442	UHMWPE
1632443	UHMWPE
1632444	UHMWPE
1632445	UHMWPE
1632446	UHMWPE
1632447	UHMWPE
1632448	UHMWPE
1632449	UHMWPE
1632450	UHMWPE
1632451	UHMWPE
1648001	Co-Cr-Mo Alloy
1648002	Co-Cr-Mo Alloy
1648003	Co-Cr-Mo Alloy
1648004	Co-Cr-Mo Alloy
1648005	Co-Cr-Mo Alloy
1648006	Co-Cr-Mo Alloy
1648007	Co-Cr-Mo Alloy
1648027	Co-Cr-Mo Alloy
1648051	Co-Cr-Mo Alloy
1648052	Co-Cr-Mo Alloy
1648053	Co-Cr-Mo Alloy

1648054	Co-Cr-Mo Alloy
1648055	Co-Cr-Mo Alloy
1648056	Co-Cr-Mo Alloy
1648901	Co-Cr-Mo Alloy + Hydroxyapatite
1648902	Co-Cr-Mo Alloy + Hydroxyapatite
1648903	Co-Cr-Mo Alloy + Hydroxyapatite
1648904	Co-Cr-Mo Alloy + Hydroxyapatite
1648905	Co-Cr-Mo Alloy + Hydroxyapatite
1648906	Co-Cr-Mo Alloy + Hydroxyapatite
1648907	Co-Cr-Mo Alloy + Hydroxyapatite
1648927	Co-Cr-Mo Alloy + Hydroxyapatite
1648951	Co-Cr-Mo Alloy + Hydroxyapatite
1648952	Co-Cr-Mo Alloy + Hydroxyapatite
1648953	Co-Cr-Mo Alloy + Hydroxyapatite
1648954	Co-Cr-Mo Alloy + Hydroxyapatite
1648955	Co-Cr-Mo Alloy + Hydroxyapatite
1648956	Co-Cr-Mo Alloy + Hydroxyapatite
1651000	Ti6Al4V Alloy + Hydroxyapatite
1651001	Ti6Al4V Alloy + Hydroxyapatite
1651002	Ti6Al4V Alloy + Hydroxyapatite
1651003	Ti6Al4V Alloy + Hydroxyapatite
1651004	Ti6Al4V Alloy + Hydroxyapatite
1651005	Ti6Al4V Alloy + Hydroxyapatite
1651006	Ti6Al4V Alloy + Hydroxyapatite
1651007	Ti6Al4V Alloy + Hydroxyapatite
1651051	Ti6Al4V Alloy + Hydroxyapatite
1651052	Ti6Al4V Alloy + Hydroxyapatite
1651053	Ti6Al4V Alloy + Hydroxyapatite
1651054	Ti6Al4V Alloy + Hydroxyapatite
1651055	Ti6Al4V Alloy + Hydroxyapatite
1651056	Ti6Al4V Alloy + Hydroxyapatite
1652000	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652001	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652002	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652003	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652004	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652005	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652006	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652007	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652051	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652052	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652053	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652054	Ti6Al4V Alloy + TiNbN + Hydroxyapatite

1652055	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1652056	Ti6Al4V Alloy + TiNbN + Hydroxyapatite
1660001	Co-Cr-Mo Alloy
1660002	Co-Cr-Mo Alloy
1660003	Co-Cr-Mo Alloy
1660004	Co-Cr-Mo Alloy
1660005	Co-Cr-Mo Alloy
1660006	Co-Cr-Mo Alloy
1660007	Co-Cr-Mo Alloy
1660012	Co-Cr-Mo Alloy
1660023	Co-Cr-Mo Alloy
1660034	Co-Cr-Mo Alloy
1660045	Co-Cr-Mo Alloy
1660056	Co-Cr-Mo Alloy
1660067	Co-Cr-Mo Alloy
1661001	Co-Cr-Mo Alloy + Hydroxyapatite
1661002	Co-Cr-Mo Alloy + Hydroxyapatite
1661003	Co-Cr-Mo Alloy + Hydroxyapatite
1661004	Co-Cr-Mo Alloy + Hydroxyapatite
1661005	Co-Cr-Mo Alloy + Hydroxyapatite
1661006	Co-Cr-Mo Alloy + Hydroxyapatite
1661007	Co-Cr-Mo Alloy + Hydroxyapatite
1661012	Co-Cr-Mo Alloy + Hydroxyapatite
1661023	Co-Cr-Mo Alloy + Hydroxyapatite
1661034	Co-Cr-Mo Alloy + Hydroxyapatite
1661045	Co-Cr-Mo Alloy + Hydroxyapatite
1661056	Co-Cr-Mo Alloy + Hydroxyapatite
1661067	Co-Cr-Mo Alloy + Hydroxyapatite
1662001	Co-Cr-Mo Alloy
1662002	Co-Cr-Mo Alloy
1662003	Co-Cr-Mo Alloy
1662004	Co-Cr-Mo Alloy
1662005	Co-Cr-Mo Alloy
1662006	Co-Cr-Mo Alloy
1662007	Co-Cr-Mo Alloy
1662012	Co-Cr-Mo Alloy
1662023	Co-Cr-Mo Alloy
1662034	Co-Cr-Mo Alloy
1662045	Co-Cr-Mo Alloy
1662056	Co-Cr-Mo Alloy
1662067	Co-Cr-Mo Alloy
1663001	Co-Cr-Mo Alloy + Hydroxyapatite
1663002	Co-Cr-Mo Alloy + Hydroxyapatite

1663003	Co-Cr-Mo Alloy + Hydroxyapatite
1663004	Co-Cr-Mo Alloy + Hydroxyapatite
1663005	Co-Cr-Mo Alloy + Hydroxyapatite
1663006	Co-Cr-Mo Alloy + Hydroxyapatite
1663007	Co-Cr-Mo Alloy + Hydroxyapatite
1663012	Co-Cr-Mo Alloy + Hydroxyapatite
1663023	Co-Cr-Mo Alloy + Hydroxyapatite
1663034	Co-Cr-Mo Alloy + Hydroxyapatite
1663045	Co-Cr-Mo Alloy + Hydroxyapatite
1663056	Co-Cr-Mo Alloy + Hydroxyapatite
1663067	Co-Cr-Mo Alloy + Hydroxyapatite
1668001	Co-Cr-Mo Alloy
1668002	Co-Cr-Mo Alloy
1668003	Co-Cr-Mo Alloy
1668004	Co-Cr-Mo Alloy
1668005	Co-Cr-Mo Alloy
1668006	Co-Cr-Mo Alloy
1668007	Co-Cr-Mo Alloy
1668012	Co-Cr-Mo Alloy
1668023	Co-Cr-Mo Alloy
1668034	Co-Cr-Mo Alloy
1668045	Co-Cr-Mo Alloy
1668056	Co-Cr-Mo Alloy
1668067	Co-Cr-Mo Alloy
1668901	Co-Cr-Mo Alloy + TiNbN
1668902	Co-Cr-Mo Alloy + TiNbN
1668903	Co-Cr-Mo Alloy + TiNbN
1668904	Co-Cr-Mo Alloy + TiNbN
1668905	Co-Cr-Mo Alloy + TiNbN
1668906	Co-Cr-Mo Alloy + TiNbN
1668907	Co-Cr-Mo Alloy + TiNbN
1668912	Co-Cr-Mo Alloy + TiNbN
1668923	Co-Cr-Mo Alloy + TiNbN
1668934	Co-Cr-Mo Alloy + TiNbN
1668945	Co-Cr-Mo Alloy + TiNbN
1668956	Co-Cr-Mo Alloy + TiNbN
1668967	Co-Cr-Mo Alloy + TiNbN
1668990	Co-Cr-Mo Alloy + TiNbN
1669090	Co-Cr-Mo Alloy
1700002	Ti6Al4V Alloy
1700004	Ti6Al4V Alloy
1700005	Ti6Al4V Alloy
1700006	Ti6Al4V Alloy

1700007	Ti6Al4V Alloy
1700010	Ti6Al4V Alloy
1700011	Ti6Al4V Alloy
1700012	Ti6Al4V Alloy
1700013	Ti6Al4V Alloy
1700020	Ti6Al4V Alloy
1700021	Ti6Al4V Alloy
1700022	Ti6Al4V Alloy
1701000	Co-Cr-Mo Alloy + Ti
1701001	Co-Cr-Mo Alloy + Ti
1701002	Co-Cr-Mo Alloy + Ti
1701003	Co-Cr-Mo Alloy + Ti
1701004	Co-Cr-Mo Alloy + Ti
1701005	Co-Cr-Mo Alloy + Ti
1701006	Co-Cr-Mo Alloy + Ti
1701007	Co-Cr-Mo Alloy + Ti
1701051	Co-Cr-Mo Alloy + Ti
1701052	Co-Cr-Mo Alloy + Ti
1701053	Co-Cr-Mo Alloy + Ti
1701054	Co-Cr-Mo Alloy + Ti
1701055	Co-Cr-Mo Alloy + Ti
1701056	Co-Cr-Mo Alloy + Ti
1702000	Co-Cr-Mo Alloy
1702001	Co-Cr-Mo Alloy
1702002	Co-Cr-Mo Alloy
1702003	Co-Cr-Mo Alloy
1702004	Co-Cr-Mo Alloy
1702005	Co-Cr-Mo Alloy
1702006	Co-Cr-Mo Alloy
1702007	Co-Cr-Mo Alloy
1702051	Co-Cr-Mo Alloy
1702052	Co-Cr-Mo Alloy
1702053	Co-Cr-Mo Alloy
1702054	Co-Cr-Mo Alloy
1702055	Co-Cr-Mo Alloy
1702056	Co-Cr-Mo Alloy
1703000	Co-Cr-Mo Alloy
1703001	Co-Cr-Mo Alloy
1703002	Co-Cr-Mo Alloy
1703003	Co-Cr-Mo Alloy
1703004	Co-Cr-Mo Alloy
1703005	Co-Cr-Mo Alloy
1703006	Co-Cr-Mo Alloy

1703007	Co-Cr-Mo Alloy
1703051	Co-Cr-Mo Alloy
1703052	Co-Cr-Mo Alloy
1703053	Co-Cr-Mo Alloy
1703054	Co-Cr-Mo Alloy
1703055	Co-Cr-Mo Alloy
1703056	Co-Cr-Mo Alloy
1704000	Co-Cr-Mo Alloy + Hydroxyapatite
1704001	Co-Cr-Mo Alloy + Hydroxyapatite
1704002	Co-Cr-Mo Alloy + Hydroxyapatite
1704003	Co-Cr-Mo Alloy + Hydroxyapatite
1704004	Co-Cr-Mo Alloy + Hydroxyapatite
1704005	Co-Cr-Mo Alloy + Hydroxyapatite
1704006	Co-Cr-Mo Alloy + Hydroxyapatite
1704007	Co-Cr-Mo Alloy + Hydroxyapatite
1704051	Co-Cr-Mo Alloy + Hydroxyapatite
1704052	Co-Cr-Mo Alloy + Hydroxyapatite
1704053	Co-Cr-Mo Alloy + Hydroxyapatite
1704054	Co-Cr-Mo Alloy + Hydroxyapatite
1704055	Co-Cr-Mo Alloy + Hydroxyapatite
1704056	Co-Cr-Mo Alloy + Hydroxyapatite
1705000	Co-Cr-Mo Alloy + Hydroxyapatite
1705001	Co-Cr-Mo Alloy + Hydroxyapatite
1705002	Co-Cr-Mo Alloy + Hydroxyapatite
1705003	Co-Cr-Mo Alloy + Hydroxyapatite
1705004	Co-Cr-Mo Alloy + Hydroxyapatite
1705005	Co-Cr-Mo Alloy + Hydroxyapatite
1705006	Co-Cr-Mo Alloy + Hydroxyapatite
1705007	Co-Cr-Mo Alloy + Hydroxyapatite
1705051	Co-Cr-Mo Alloy + Hydroxyapatite
1705052	Co-Cr-Mo Alloy + Hydroxyapatite
1705053	Co-Cr-Mo Alloy + Hydroxyapatite
1705054	Co-Cr-Mo Alloy + Hydroxyapatite
1705055	Co-Cr-Mo Alloy + Hydroxyapatite
1705056	Co-Cr-Mo Alloy + Hydroxyapatite
1710000	Ti6Al4V Alloy
1710020	Ti6Al4V Alloy
1710025	Ti6Al4V Alloy
1710030	Ti6Al4V Alloy
1710035	Ti6Al4V Alloy
1710040	Ti6Al4V Alloy
1710045	Ti6Al4V Alloy
1710050	Ti6Al4V Alloy

1710101	Ti6Al4V Alloy
1710102	Ti6Al4V Alloy
1710103	Ti6Al4V Alloy
1711001	Co-Cr-Mo Alloy + Ti
1711002	Co-Cr-Mo Alloy + Ti
1711003	Co-Cr-Mo Alloy + Ti
1711004	Co-Cr-Mo Alloy + Ti
1711005	Co-Cr-Mo Alloy + Ti
1711006	Co-Cr-Mo Alloy + Ti
1711011	Co-Cr-Mo Alloy + Ti
1711012	Co-Cr-Mo Alloy + Ti
1711013	Co-Cr-Mo Alloy + Ti
1711014	Co-Cr-Mo Alloy + Ti
1711015	Co-Cr-Mo Alloy + Ti
1711016	Co-Cr-Mo Alloy + Ti
1712001	Co-Cr-Mo Alloy
1712002	Co-Cr-Mo Alloy
1712003	Co-Cr-Mo Alloy
1712004	Co-Cr-Mo Alloy
1712005	Co-Cr-Mo Alloy
1712006	Co-Cr-Mo Alloy
1712011	Co-Cr-Mo Alloy
1712012	Co-Cr-Mo Alloy
1712013	Co-Cr-Mo Alloy
1712014	Co-Cr-Mo Alloy
1712015	Co-Cr-Mo Alloy
1712016	Co-Cr-Mo Alloy
1713001	Co-Cr-Mo Alloy
1713002	Co-Cr-Mo Alloy
1713003	Co-Cr-Mo Alloy
1713004	Co-Cr-Mo Alloy
1713005	Co-Cr-Mo Alloy
1713006	Co-Cr-Mo Alloy
1713011	Co-Cr-Mo Alloy
1713012	Co-Cr-Mo Alloy
1713013	Co-Cr-Mo Alloy
1713014	Co-Cr-Mo Alloy
1713015	Co-Cr-Mo Alloy
1713016	Co-Cr-Mo Alloy
1714001	Co-Cr-Mo Alloy + Hydroxyapatite
1714002	Co-Cr-Mo Alloy + Hydroxyapatite
1714003	Co-Cr-Mo Alloy + Hydroxyapatite
1714004	Co-Cr-Mo Alloy + Hydroxyapatite

1714005	Co-Cr-Mo Alloy + Hydroxyapatite
1714006	Co-Cr-Mo Alloy + Hydroxyapatite
1714011	Co-Cr-Mo Alloy + Hydroxyapatite
1714012	Co-Cr-Mo Alloy + Hydroxyapatite
1714013	Co-Cr-Mo Alloy + Hydroxyapatite
1714014	Co-Cr-Mo Alloy + Hydroxyapatite
1714015	Co-Cr-Mo Alloy + Hydroxyapatite
1714016	Co-Cr-Mo Alloy + Hydroxyapatite
1715001	Co-Cr-Mo Alloy + Hydroxyapatite
1715002	Co-Cr-Mo Alloy + Hydroxyapatite
1715003	Co-Cr-Mo Alloy + Hydroxyapatite
1715004	Co-Cr-Mo Alloy + Hydroxyapatite
1715005	Co-Cr-Mo Alloy + Hydroxyapatite
1715006	Co-Cr-Mo Alloy + Hydroxyapatite
1715011	Co-Cr-Mo Alloy + Hydroxyapatite
1715012	1715011
1715013	1715012
1715014	1715013
1715015	1715014
1715016	1715015
1716001	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716002	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716003	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716004	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716005	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716006	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716011	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716012	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716013	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716014	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716015	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1716016	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1721001	Co-Cr-Mo Alloy + Ti
1721002	Co-Cr-Mo Alloy + Ti
1721003	Co-Cr-Mo Alloy + Ti
1721004	Co-Cr-Mo Alloy + Ti
1721005	Co-Cr-Mo Alloy + Ti
1721006	Co-Cr-Mo Alloy + Ti
1721011	Co-Cr-Mo Alloy + Ti
1721012	Co-Cr-Mo Alloy + Ti
1721013	Co-Cr-Mo Alloy + Ti
1721014	Co-Cr-Mo Alloy + Ti
1721015	Co-Cr-Mo Alloy + Ti

1721016	Co-Cr-Mo Alloy + Ti
1722001	Co-Cr-Mo Alloy
1722002	Co-Cr-Mo Alloy
1722003	Co-Cr-Mo Alloy
1722004	Co-Cr-Mo Alloy
1722005	Co-Cr-Mo Alloy
1722006	Co-Cr-Mo Alloy
1722011	Co-Cr-Mo Alloy
1722012	Co-Cr-Mo Alloy
1722013	Co-Cr-Mo Alloy
1722014	Co-Cr-Mo Alloy
1722015	Co-Cr-Mo Alloy
1722016	Co-Cr-Mo Alloy
1723001	Co-Cr-Mo Alloy
1723002	Co-Cr-Mo Alloy
1723003	Co-Cr-Mo Alloy
1723004	Co-Cr-Mo Alloy
1723005	Co-Cr-Mo Alloy
1723006	Co-Cr-Mo Alloy
1723011	Co-Cr-Mo Alloy
1723012	Co-Cr-Mo Alloy
1723013	Co-Cr-Mo Alloy
1723014	Co-Cr-Mo Alloy
1723015	Co-Cr-Mo Alloy
1723016	Co-Cr-Mo Alloy
1724001	Co-Cr-Mo Alloy + Hydroxyapatite
1724002	Co-Cr-Mo Alloy + Hydroxyapatite
1724003	Co-Cr-Mo Alloy + Hydroxyapatite
1724004	Co-Cr-Mo Alloy + Hydroxyapatite
1724005	Co-Cr-Mo Alloy + Hydroxyapatite
1724006	Co-Cr-Mo Alloy + Hydroxyapatite
1724011	Co-Cr-Mo Alloy + Hydroxyapatite
1724012	Co-Cr-Mo Alloy + Hydroxyapatite
1724013	Co-Cr-Mo Alloy + Hydroxyapatite
1724014	Co-Cr-Mo Alloy + Hydroxyapatite
1724015	Co-Cr-Mo Alloy + Hydroxyapatite
1724016	Co-Cr-Mo Alloy + Hydroxyapatite
1725001	Co-Cr-Mo Alloy + Hydroxyapatite
1725002	Co-Cr-Mo Alloy + Hydroxyapatite
1725003	Co-Cr-Mo Alloy + Hydroxyapatite
1725004	Co-Cr-Mo Alloy + Hydroxyapatite
1725005	Co-Cr-Mo Alloy + Hydroxyapatite
1725006	Co-Cr-Mo Alloy + Hydroxyapatite

1725011	Co-Cr-Mo Alloy + Hydroxyapatite
1725012	Co-Cr-Mo Alloy + Hydroxyapatite
1725013	Co-Cr-Mo Alloy + Hydroxyapatite
1725014	Co-Cr-Mo Alloy + Hydroxyapatite
1725015	Co-Cr-Mo Alloy + Hydroxyapatite
1725016	Co-Cr-Mo Alloy + Hydroxyapatite
1726001	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726002	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726003	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726004	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726005	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726006	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726011	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726012	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726013	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726014	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726015	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1726016	Co-Cr-Mo Alloy + Ti + Hydroxyapatite
1731001	Ti6Al4V Alloy
1731002	Ti6Al4V Alloy
1731003	Ti6Al4V Alloy
1731004	Ti6Al4V Alloy
1731005	Ti6Al4V Alloy
1731006	Ti6Al4V Alloy
1731011	Ti6Al4V Alloy
1731012	Ti6Al4V Alloy
1731013	Ti6Al4V Alloy
1731014	Ti6Al4V Alloy
1731015	Ti6Al4V Alloy
1731016	Ti6Al4V Alloy
1732001	Ti6Al4V Alloy
1732002	Ti6Al4V Alloy
1732003	Ti6Al4V Alloy
1732004	Ti6Al4V Alloy
1732005	Ti6Al4V Alloy
1732006	Ti6Al4V Alloy
1732011	Ti6Al4V Alloy
1732012	Ti6Al4V Alloy
1732013	Ti6Al4V Alloy
1732014	Ti6Al4V Alloy
1732015	Ti6Al4V Alloy
1732016	Ti6Al4V Alloy
1733001	Ti6Al4V Alloy

1733002	Ti6Al4V Alloy
1733003	Ti6Al4V Alloy
1733004	Ti6Al4V Alloy
1733005	Ti6Al4V Alloy
1733006	Ti6Al4V Alloy
1733011	Ti6Al4V Alloy
1733012	Ti6Al4V Alloy
1733013	Ti6Al4V Alloy
1733014	Ti6Al4V Alloy
1733015	Ti6Al4V Alloy
1733016	Ti6Al4V Alloy
1734001	Ti6Al4V Alloy + Hydroxyapatite
1734002	Ti6Al4V Alloy + Hydroxyapatite
1734003	Ti6Al4V Alloy + Hydroxyapatite
1734004	Ti6Al4V Alloy + Hydroxyapatite
1734005	Ti6Al4V Alloy + Hydroxyapatite
1734006	Ti6Al4V Alloy + Hydroxyapatite
1734011	Ti6Al4V Alloy + Hydroxyapatite
1734012	Ti6Al4V Alloy + Hydroxyapatite
1734013	Ti6Al4V Alloy + Hydroxyapatite
1734014	Ti6Al4V Alloy + Hydroxyapatite
1734015	Ti6Al4V Alloy + Hydroxyapatite
1734016	Ti6Al4V Alloy + Hydroxyapatite
1735001	Ti6Al4V Alloy + Hydroxyapatite
1735002	Ti6Al4V Alloy + Hydroxyapatite
1735003	Ti6Al4V Alloy + Hydroxyapatite
1735004	Ti6Al4V Alloy + Hydroxyapatite
1735005	Ti6Al4V Alloy + Hydroxyapatite
1735006	Ti6Al4V Alloy + Hydroxyapatite
1735011	Ti6Al4V Alloy + Hydroxyapatite
1735012	Ti6Al4V Alloy + Hydroxyapatite
1735013	Ti6Al4V Alloy + Hydroxyapatite
1735014	Ti6Al4V Alloy + Hydroxyapatite
1735015	Ti6Al4V Alloy + Hydroxyapatite
1735016	Ti6Al4V Alloy + Hydroxyapatite
1736001	Ti6Al4V Alloy + Hydroxyapatite
1736002	Ti6Al4V Alloy + Hydroxyapatite
1736003	Ti6Al4V Alloy + Hydroxyapatite
1736004	Ti6Al4V Alloy + Hydroxyapatite
1736005	Ti6Al4V Alloy + Hydroxyapatite
1736006	Ti6Al4V Alloy + Hydroxyapatite
1736011	Ti6Al4V Alloy + Hydroxyapatite
1736012	Ti6Al4V Alloy + Hydroxyapatite

1736013	Ti6Al4V Alloy + Hydroxyapatite
1736014	Ti6Al4V Alloy + Hydroxyapatite
1736015	Ti6Al4V Alloy + Hydroxyapatite
1736016	Ti6Al4V Alloy + Hydroxyapatite
1740901	UHMWPE
1740902	UHMWPE
1740903	UHMWPE
1740904	UHMWPE
1740905	UHMWPE
1740906	UHMWPE
1740911	UHMWPE
1740912	UHMWPE
1740913	UHMWPE
1740914	UHMWPE
1740915	UHMWPE
1740916	UHMWPE
1741001	UHMWPE
1741002	UHMWPE
1741003	UHMWPE
1741004	UHMWPE
1741005	UHMWPE
1741006	UHMWPE
1741011	UHMWPE
1741012	UHMWPE
1741013	UHMWPE
1741014	UHMWPE
1741015	UHMWPE
1741016	UHMWPE
1741101	UHMWPE
1741102	UHMWPE
1741103	UHMWPE
1741104	UHMWPE
1741105	UHMWPE
1741106	UHMWPE
1741111	UHMWPE
1741112	UHMWPE
1741113	UHMWPE
1741114	UHMWPE
1741115	UHMWPE
1741116	UHMWPE
1741201	UHMWPE
1741202	UHMWPE
1741203	UHMWPE

1741204	UHMWPE
1741205	UHMWPE
1741206	UHMWPE
1741211	UHMWPE
1741212	UHMWPE
1741213	UHMWPE
1741214	UHMWPE
1741215	UHMWPE
1741216	UHMWPE
1741301	UHMWPE
1741302	UHMWPE
1741303	UHMWPE
1741304	UHMWPE
1741305	UHMWPE
1741306	UHMWPE
1741311	UHMWPE
1741312	UHMWPE
1741313	UHMWPE
1741314	UHMWPE
1741315	UHMWPE
1741316	UHMWPE
1741401	UHMWPE
1741402	UHMWPE
1741403	UHMWPE
1741404	UHMWPE
1741405	UHMWPE
1741406	UHMWPE
1741411	UHMWPE
1741412	UHMWPE
1741413	UHMWPE
1741414	UHMWPE
1741415	UHMWPE
1741416	UHMWPE
1741601	UHMWPE
1741602	UHMWPE
1741603	UHMWPE
1741604	UHMWPE
1741605	UHMWPE
1741606	UHMWPE
1741611	UHMWPE
1741612	UHMWPE
1741613	UHMWPE
1741614	UHMWPE

1741615	UHMWPE
1741616	UHMWPE
1750001	Co-Cr-Mo Alloy + Hydroxyapatite
1750002	Co-Cr-Mo Alloy + Hydroxyapatite
1750003	Co-Cr-Mo Alloy + Hydroxyapatite
1750004	Co-Cr-Mo Alloy + Hydroxyapatite
1750005	Co-Cr-Mo Alloy + Hydroxyapatite
1750006	Co-Cr-Mo Alloy + Hydroxyapatite
1750007	Co-Cr-Mo Alloy + Hydroxyapatite
1750027	Co-Cr-Mo Alloy + Hydroxyapatite
1750051	Co-Cr-Mo Alloy + Hydroxyapatite
1750052	Co-Cr-Mo Alloy + Hydroxyapatite
1750055	Co-Cr-Mo Alloy + Hydroxyapatite
1750056	Co-Cr-Mo Alloy + Hydroxyapatite
1750083	Co-Cr-Mo Alloy + Hydroxyapatite
1750094	Co-Cr-Mo Alloy + Hydroxyapatite
1751001	Ti6Al4V Alloy
1751002	Ti6Al4V Alloy
1751003	Ti6Al4V Alloy
1751004	Ti6Al4V Alloy
1751005	Ti6Al4V Alloy
1751006	Ti6Al4V Alloy
1751007	Ti6Al4V Alloy
1751027	Ti6Al4V Alloy
1751051	Ti6Al4V Alloy
1751052	Ti6Al4V Alloy
1751055	Ti6Al4V Alloy
1751056	Ti6Al4V Alloy
1751083	Ti6Al4V Alloy
1751094	Ti6Al4V Alloy
1752001	Ti6Al4V Alloy + Hydroxyapatite
1752002	Ti6Al4V Alloy + Hydroxyapatite
1752003	Ti6Al4V Alloy + Hydroxyapatite
1752004	Ti6Al4V Alloy + Hydroxyapatite
1752005	Ti6Al4V Alloy + Hydroxyapatite
1752006	Ti6Al4V Alloy + Hydroxyapatite
1752007	Ti6Al4V Alloy + Hydroxyapatite
1752027	Ti6Al4V Alloy + Hydroxyapatite
1752051	Ti6Al4V Alloy + Hydroxyapatite
1752052	Ti6Al4V Alloy + Hydroxyapatite
1752055	Ti6Al4V Alloy + Hydroxyapatite
1752056	Ti6Al4V Alloy + Hydroxyapatite
1752083	Ti6Al4V Alloy + Hydroxyapatite

1752094	Ti6Al4V Alloy + Hydroxyapatite
1753001	Ti6Al4V Alloy
1753002	Ti6Al4V Alloy
1753003	Ti6Al4V Alloy
1753004	Ti6Al4V Alloy
1753005	Ti6Al4V Alloy
1753006	Ti6Al4V Alloy
1753007	Ti6Al4V Alloy
1753027	Ti6Al4V Alloy
1753051	Ti6Al4V Alloy
1753052	Ti6Al4V Alloy
1753055	Ti6Al4V Alloy
1753056	Ti6Al4V Alloy
1753083	Ti6Al4V Alloy
1753094	Ti6Al4V Alloy
1758001	Ti6Al4V Alloy
1758002	Ti6Al4V Alloy
1758003	Ti6Al4V Alloy
1758004	Ti6Al4V Alloy
1758005	Ti6Al4V Alloy
1758006	Ti6Al4V Alloy
1758007	Ti6Al4V Alloy
1758027	Ti6Al4V Alloy
1758051	Ti6Al4V Alloy
1758052	Ti6Al4V Alloy
1758055	Ti6Al4V Alloy
1758056	Ti6Al4V Alloy
1758083	Ti6Al4V Alloy
1758094	Ti6Al4V Alloy
1759000	Ti6Al4V Alloy
1760001	Ti6Al4V Alloy
1760002	Ti6Al4V Alloy
1760003	Ti6Al4V Alloy
1760004	Ti6Al4V Alloy
1760005	Ti6Al4V Alloy
1760006	Ti6Al4V Alloy
1760007	Ti6Al4V Alloy
1760012	Ti6Al4V Alloy
1760023	Ti6Al4V Alloy
1760034	Ti6Al4V Alloy
1760045	Ti6Al4V Alloy
1760056	Ti6Al4V Alloy
1760067	Ti6Al4V Alloy

1761001	Ti6Al4V Alloy + Hydroxyapatite
1761002	Ti6Al4V Alloy + Hydroxyapatite
1761003	Ti6Al4V Alloy + Hydroxyapatite
1761004	Ti6Al4V Alloy + Hydroxyapatite
1761005	Ti6Al4V Alloy + Hydroxyapatite
1761006	Ti6Al4V Alloy + Hydroxyapatite
1761007	Ti6Al4V Alloy + Hydroxyapatite
1761012	Ti6Al4V Alloy + Hydroxyapatite
1761023	Ti6Al4V Alloy + Hydroxyapatite
1761034	Ti6Al4V Alloy + Hydroxyapatite
1761045	Ti6Al4V Alloy + Hydroxyapatite
1761056	Ti6Al4V Alloy + Hydroxyapatite
1761067	Ti6Al4V Alloy + Hydroxyapatite
1762001	Ti6Al4V Alloy
1762002	Ti6Al4V Alloy
1762003	Ti6Al4V Alloy
1762004	Ti6Al4V Alloy
1762005	Ti6Al4V Alloy
1762006	Ti6Al4V Alloy
1762007	Ti6Al4V Alloy
1762012	Ti6Al4V Alloy
1762023	Ti6Al4V Alloy
1762034	Ti6Al4V Alloy
1762045	Ti6Al4V Alloy
1762056	Ti6Al4V Alloy
1762067	Ti6Al4V Alloy
1763001	Ti6Al4V Alloy + Hydroxyapatite
1763002	Ti6Al4V Alloy + Hydroxyapatite
1763003	Ti6Al4V Alloy + Hydroxyapatite
1763004	Ti6Al4V Alloy + Hydroxyapatite
1763005	Ti6Al4V Alloy + Hydroxyapatite
1763006	Ti6Al4V Alloy + Hydroxyapatite
1763007	Ti6Al4V Alloy + Hydroxyapatite
1763012	Ti6Al4V Alloy + Hydroxyapatite
1763023	Ti6Al4V Alloy + Hydroxyapatite
1763034	Ti6Al4V Alloy + Hydroxyapatite
1763045	Ti6Al4V Alloy + Hydroxyapatite
1763056	Ti6Al4V Alloy + Hydroxyapatite
1763067	Ti6Al4V Alloy + Hydroxyapatite
1768001	Ti6Al4V Alloy
1768002	Ti6Al4V Alloy
1768003	Ti6Al4V Alloy
1768004	Ti6Al4V Alloy

1768005	Ti6Al4V Alloy
1768006	Ti6Al4V Alloy
1768007	Ti6Al4V Alloy
1768012	Ti6Al4V Alloy
1768023	Ti6Al4V Alloy
1768034	Ti6Al4V Alloy
1768045	Ti6Al4V Alloy
1768056	Ti6Al4V Alloy
1768067	Ti6Al4V Alloy
1771001	Ti6Al4V Alloy
1771002	Ti6Al4V Alloy
1771003	Ti6Al4V Alloy
1771004	Ti6Al4V Alloy
1771005	Ti6Al4V Alloy
1771006	Ti6Al4V Alloy
1771007	Ti6Al4V Alloy
1771011	Ti6Al4V Alloy
1771012	Ti6Al4V Alloy
1771013	Ti6Al4V Alloy
1771014	Ti6Al4V Alloy
1771015	Ti6Al4V Alloy
1771016	Ti6Al4V Alloy
1771017	Ti6Al4V Alloy
1771051	Ti6Al4V Alloy
1771061	Ti6Al4V Alloy
1778001	Ti6Al4V Alloy
1778002	Ti6Al4V Alloy
1778003	Ti6Al4V Alloy
1778004	Ti6Al4V Alloy
1778005	Ti6Al4V Alloy
1778006	Ti6Al4V Alloy
1778007	Ti6Al4V Alloy
1778011	Ti6Al4V Alloy
1778012	Ti6Al4V Alloy
1778013	Ti6Al4V Alloy
1778014	Ti6Al4V Alloy
1778015	Ti6Al4V Alloy
1778016	Ti6Al4V Alloy
1778017	Ti6Al4V Alloy
1778051	Ti6Al4V Alloy
1778061	Ti6Al4V Alloy
1779000	Ti6Al4V Alloy
1779005	Ti6Al4V Alloy

1779010	Ti6Al4V Alloy
1779015	Ti6Al4V Alloy
1779020	Ti6Al4V Alloy
1781000	Co-Cr-Mo Alloy
1781001	Co-Cr-Mo Alloy
1781002	Co-Cr-Mo Alloy
1781003	Co-Cr-Mo Alloy
1781004	Co-Cr-Mo Alloy
1781005	Co-Cr-Mo Alloy
1781006	Co-Cr-Mo Alloy
1781007	Co-Cr-Mo Alloy
1781051	Co-Cr-Mo Alloy
1781052	Co-Cr-Mo Alloy
1781053	Co-Cr-Mo Alloy
1781054	Co-Cr-Mo Alloy
1781055	Co-Cr-Mo Alloy
1781056	Co-Cr-Mo Alloy
1782001	Co-Cr-Mo Alloy
1782002	Co-Cr-Mo Alloy
1782003	Co-Cr-Mo Alloy
1782004	Co-Cr-Mo Alloy
1782005	Co-Cr-Mo Alloy
1782006	Co-Cr-Mo Alloy
1782011	Co-Cr-Mo Alloy
1782012	Co-Cr-Mo Alloy
1782013	Co-Cr-Mo Alloy
1782014	Co-Cr-Mo Alloy
1782015	Co-Cr-Mo Alloy
1782016	Co-Cr-Mo Alloy
1783001	Co-Cr-Mo Alloy
1783002	Co-Cr-Mo Alloy
1783003	Co-Cr-Mo Alloy
1783004	Co-Cr-Mo Alloy
1783005	Co-Cr-Mo Alloy
1783006	Co-Cr-Mo Alloy
1783011	Co-Cr-Mo Alloy
1783012	Co-Cr-Mo Alloy
1783013	Co-Cr-Mo Alloy
1783014	Co-Cr-Mo Alloy
1783015	Co-Cr-Mo Alloy
1783016	Co-Cr-Mo Alloy
1784001	Ti6Al4V Alloy
1784002	Ti6Al4V Alloy

1784003	Ti6Al4V Alloy
1784004	Ti6Al4V Alloy
1784005	Ti6Al4V Alloy
1784006	Ti6Al4V Alloy
1784011	Ti6Al4V Alloy
1784012	Ti6Al4V Alloy
1784013	Ti6Al4V Alloy
1784014	Ti6Al4V Alloy
1784015	Ti6Al4V Alloy
1784016	Ti6Al4V Alloy
1788900	Co-Cr-Mo Alloy + TiNbN
1788901	Co-Cr-Mo Alloy + TiNbN
1788902	Co-Cr-Mo Alloy + TiNbN
1788903	Co-Cr-Mo Alloy + TiNbN
1788904	Co-Cr-Mo Alloy + TiNbN
1788905	Co-Cr-Mo Alloy + TiNbN
1788906	Co-Cr-Mo Alloy + TiNbN
1788907	Co-Cr-Mo Alloy + TiNbN
1788951	Co-Cr-Mo Alloy + TiNbN
1788952	Co-Cr-Mo Alloy + TiNbN
1788953	Co-Cr-Mo Alloy + TiNbN
1788954	Co-Cr-Mo Alloy + TiNbN
1788955	Co-Cr-Mo Alloy + TiNbN
1788956	Co-Cr-Mo Alloy + TiNbN
1790002	Ti6Al4V Alloy
1790004	Ti6Al4V Alloy
1790006	Ti6Al4V Alloy
1790011	Ti6Al4V Alloy
1790012	Ti6Al4V Alloy
1790013	Ti6Al4V Alloy
1790014	Ti6Al4V Alloy
1790015	Ti6Al4V Alloy
1790016	Ti6Al4V Alloy
1790017	Ti6Al4V Alloy
1790018	Ti6Al4V Alloy
1790019	Ti6Al4V Alloy
1790020	Ti6Al4V Alloy
1790021	Ti6Al4V Alloy
1790111	Ti6Al4V Alloy
1790112	Ti6Al4V Alloy
1790113	Ti6Al4V Alloy
1790114	Ti6Al4V Alloy
1790115	Ti6Al4V Alloy

1790116	Ti6Al4V Alloy
1790117	Ti6Al4V Alloy
1790118	Ti6Al4V Alloy
1790119	Ti6Al4V Alloy
1790120	Ti6Al4V Alloy
1790121	Ti6Al4V Alloy
1790211	Ti6Al4V Alloy
1790212	Ti6Al4V Alloy
1790213	Ti6Al4V Alloy
1790214	Ti6Al4V Alloy
1790215	Ti6Al4V Alloy
1790216	Ti6Al4V Alloy
1790217	Ti6Al4V Alloy
1790218	Ti6Al4V Alloy
1790219	Ti6Al4V Alloy
1790220	Ti6Al4V Alloy
1790221	Ti6Al4V Alloy
1790311	Ti6Al4V Alloy
1790312	Ti6Al4V Alloy
1790313	Ti6Al4V Alloy
1790314	Ti6Al4V Alloy
1790315	Ti6Al4V Alloy
1790316	Ti6Al4V Alloy
1790317	Ti6Al4V Alloy
1790318	Ti6Al4V Alloy
1790319	Ti6Al4V Alloy
1790320	Ti6Al4V Alloy
1790321	Ti6Al4V Alloy
1790411	Ti6Al4V Alloy
1790412	Ti6Al4V Alloy
1790413	Ti6Al4V Alloy
1790414	Ti6Al4V Alloy
1790415	Ti6Al4V Alloy
1790416	Ti6Al4V Alloy
1790417	Ti6Al4V Alloy
1790418	Ti6Al4V Alloy
1790419	Ti6Al4V Alloy
1790420	Ti6Al4V Alloy
1790421	Ti6Al4V Alloy
1790511	Ti6Al4V Alloy
1790512	Ti6Al4V Alloy
1790513	Ti6Al4V Alloy
1790514	Ti6Al4V Alloy

1790515	Ti6Al4V Alloy
1790516	Ti6Al4V Alloy
1790517	Ti6Al4V Alloy
1790518	Ti6Al4V Alloy
1790519	Ti6Al4V Alloy
1790520	Ti6Al4V Alloy
1790521	Ti6Al4V Alloy
1791011	Ti6Al4V Alloy
1791012	Ti6Al4V Alloy
1791013	Ti6Al4V Alloy
1791014	Ti6Al4V Alloy
1791015	Ti6Al4V Alloy
1791016	Ti6Al4V Alloy
1791017	Ti6Al4V Alloy
1791018	Ti6Al4V Alloy
1791019	Ti6Al4V Alloy
1791020	Ti6Al4V Alloy
1791021	Ti6Al4V Alloy
1791111	Ti6Al4V Alloy
1791112	Ti6Al4V Alloy
1791113	Ti6Al4V Alloy
1791114	Ti6Al4V Alloy
1791115	Ti6Al4V Alloy
1791116	Ti6Al4V Alloy
1791117	Ti6Al4V Alloy
1791118	Ti6Al4V Alloy
1791119	Ti6Al4V Alloy
1791120	Ti6Al4V Alloy
1791121	Ti6Al4V Alloy
1791211	Ti6Al4V Alloy
1791212	Ti6Al4V Alloy
1791213	Ti6Al4V Alloy
1791214	Ti6Al4V Alloy
1791215	Ti6Al4V Alloy
1791216	Ti6Al4V Alloy
1791217	Ti6Al4V Alloy
1791218	Ti6Al4V Alloy
1791219	Ti6Al4V Alloy
1791220	Ti6Al4V Alloy
1791221	Ti6Al4V Alloy
1791311	Ti6Al4V Alloy
1791312	Ti6Al4V Alloy
1791313	Ti6Al4V Alloy

1791314	Ti6Al4V Alloy
1791315	Ti6Al4V Alloy
1791316	Ti6Al4V Alloy
1791317	Ti6Al4V Alloy
1791318	Ti6Al4V Alloy
1791319	Ti6Al4V Alloy
1791320	Ti6Al4V Alloy
1791321	Ti6Al4V Alloy
1791411	Ti6Al4V Alloy
1791412	Ti6Al4V Alloy
1791413	Ti6Al4V Alloy
1791414	Ti6Al4V Alloy
1791415	Ti6Al4V Alloy
1791416	Ti6Al4V Alloy
1791417	Ti6Al4V Alloy
1791418	Ti6Al4V Alloy
1791419	Ti6Al4V Alloy
1791420	Ti6Al4V Alloy
1791421	Ti6Al4V Alloy
1791511	Ti6Al4V Alloy
1791512	Ti6Al4V Alloy
1791513	Ti6Al4V Alloy
1791514	Ti6Al4V Alloy
1791515	Ti6Al4V Alloy
1791516	Ti6Al4V Alloy
1791517	Ti6Al4V Alloy
1791518	Ti6Al4V Alloy
1791519	Ti6Al4V Alloy
1791520	Ti6Al4V Alloy
1791521	Ti6Al4V Alloy
1791615	Ti6Al4V Alloy
1791617	Ti6Al4V Alloy
1791619	Ti6Al4V Alloy
1791621	Ti6Al4V Alloy
1791623	Ti6Al4V Alloy
1791715	Ti6Al4V Alloy
1791717	Ti6Al4V Alloy
1791719	Ti6Al4V Alloy
1791721	Ti6Al4V Alloy
1791723	Ti6Al4V Alloy
1791725	Ti6Al4V Alloy
1791817	Ti6Al4V Alloy
1791819	Ti6Al4V Alloy

1791821	Ti6Al4V Alloy
1791823	Ti6Al4V Alloy
1791825	Ti6Al4V Alloy
1794010	Ti6Al4V Alloy
1794011	Ti6Al4V Alloy
1794012	Ti6Al4V Alloy
1794013	Ti6Al4V Alloy
1794014	Ti6Al4V Alloy
1794015	Ti6Al4V Alloy
1794016	Ti6Al4V Alloy
1794017	Ti6Al4V Alloy
1794018	Ti6Al4V Alloy
1794019	Ti6Al4V Alloy
1794020	Ti6Al4V Alloy
1794022	Ti6Al4V Alloy
1794024	Ti6Al4V Alloy
1794110	Ti6Al4V Alloy
1794111	Ti6Al4V Alloy
1794112	Ti6Al4V Alloy
1794113	Ti6Al4V Alloy
1794114	Ti6Al4V Alloy
1794115	Ti6Al4V Alloy
1794116	Ti6Al4V Alloy
1794117	Ti6Al4V Alloy
1794118	Ti6Al4V Alloy
1794119	Ti6Al4V Alloy
1794120	Ti6Al4V Alloy
1794122	Ti6Al4V Alloy
1794124	Ti6Al4V Alloy
1794210	Ti6Al4V Alloy
1794211	Ti6Al4V Alloy
1794212	Ti6Al4V Alloy
1794213	Ti6Al4V Alloy
1794214	Ti6Al4V Alloy
1794215	Ti6Al4V Alloy
1794216	Ti6Al4V Alloy
1794217	Ti6Al4V Alloy
1794218	Ti6Al4V Alloy
1794219	Ti6Al4V Alloy
1794220	Ti6Al4V Alloy
1794222	Ti6Al4V Alloy
1794224	Ti6Al4V Alloy
1794410	Ti6Al4V Alloy

1794411	Ti6Al4V Alloy
1794412	Ti6Al4V Alloy
1794413	Ti6Al4V Alloy
1794414	Ti6Al4V Alloy
1794415	Ti6Al4V Alloy
1794416	Ti6Al4V Alloy
1794417	Ti6Al4V Alloy
1794418	Ti6Al4V Alloy
1794419	Ti6Al4V Alloy
1794420	Ti6Al4V Alloy
1794422	Ti6Al4V Alloy
1794424	Ti6Al4V Alloy
1794510	Ti6Al4V Alloy
1794511	Ti6Al4V Alloy
1794512	Ti6Al4V Alloy
1794513	Ti6Al4V Alloy
1794514	Ti6Al4V Alloy
1794515	Ti6Al4V Alloy
1794516	Ti6Al4V Alloy
1794517	Ti6Al4V Alloy
1794518	Ti6Al4V Alloy
1794519	Ti6Al4V Alloy
1794520	Ti6Al4V Alloy
1794522	Ti6Al4V Alloy
1794524	Ti6Al4V Alloy
1795010	Ti6Al4V Alloy
1795011	Ti6Al4V Alloy
1795012	Ti6Al4V Alloy
1795013	Ti6Al4V Alloy
1795014	Ti6Al4V Alloy
1795015	Ti6Al4V Alloy
1795016	Ti6Al4V Alloy
1795017	Ti6Al4V Alloy
1795018	Ti6Al4V Alloy
1795019	Ti6Al4V Alloy
1795020	Ti6Al4V Alloy
1795022	Ti6Al4V Alloy
1795024	Ti6Al4V Alloy
1795110	Ti6Al4V Alloy
1795111	Ti6Al4V Alloy
1795112	Ti6Al4V Alloy
1795113	Ti6Al4V Alloy
1795114	Ti6Al4V Alloy

1795115	Ti6Al4V Alloy
1795116	Ti6Al4V Alloy
1795117	Ti6Al4V Alloy
1795118	Ti6Al4V Alloy
1795119	Ti6Al4V Alloy
1795120	Ti6Al4V Alloy
1795122	Ti6Al4V Alloy
1795124	Ti6Al4V Alloy
1795210	Ti6Al4V Alloy
1795211	Ti6Al4V Alloy
1795212	Ti6Al4V Alloy
1795213	Ti6Al4V Alloy
1795214	Ti6Al4V Alloy
1795215	Ti6Al4V Alloy
1795216	Ti6Al4V Alloy
1795217	Ti6Al4V Alloy
1795218	Ti6Al4V Alloy
1795219	Ti6Al4V Alloy
1795220	Ti6Al4V Alloy
1795222	Ti6Al4V Alloy
1795224	Ti6Al4V Alloy
1795410	Ti6Al4V Alloy
1795411	Ti6Al4V Alloy
1795412	Ti6Al4V Alloy
1795413	Ti6Al4V Alloy
1795414	Ti6Al4V Alloy
1795415	Ti6Al4V Alloy
1795416	Ti6Al4V Alloy
1795417	Ti6Al4V Alloy
1795418	Ti6Al4V Alloy
1795419	Ti6Al4V Alloy
1795420	Ti6Al4V Alloy
1795422	Ti6Al4V Alloy
1795424	Ti6Al4V Alloy
1795510	Ti6Al4V Alloy
1795511	Ti6Al4V Alloy
1795512	Ti6Al4V Alloy
1795513	Ti6Al4V Alloy
1795514	Ti6Al4V Alloy
1795515	Ti6Al4V Alloy
1795516	Ti6Al4V Alloy
1795517	Ti6Al4V Alloy
1795518	Ti6Al4V Alloy

1795519	Ti6Al4V Alloy
1795520	Ti6Al4V Alloy
1795522	Ti6Al4V Alloy
1795524	Ti6Al4V Alloy
1801001	UHMWPE
1801002	UHMWPE
1801003	UHMWPE
1801004	UHMWPE
1801005	UHMWPE
1801006	UHMWPE
1801007	UHMWPE
1801011	UHMWPE
1801012	UHMWPE
1801013	UHMWPE
1801014	UHMWPE
1801015	UHMWPE
1801016	UHMWPE
1801017	UHMWPE
1801031	UHMWPE
1801032	UHMWPE
1801033	UHMWPE
1801034	UHMWPE
1801035	UHMWPE
1801036	UHMWPE
1801037	UHMWPE
1801041	UHMWPE
1801042	UHMWPE
1801043	UHMWPE
1801044	UHMWPE
1801045	UHMWPE
1801046	UHMWPE
1801047	UHMWPE
1801071	UHMWPE
1801072	UHMWPE
1801073	UHMWPE
1801074	UHMWPE
1801075	UHMWPE
1801076	UHMWPE
1801077	UHMWPE
1801081	UHMWPE
1801082	UHMWPE
1801083	UHMWPE
1801084	UHMWPE

1801085	UHMWPE
1801086	UHMWPE
1801087	UHMWPE
1801101	UHMWPE
1801102	UHMWPE
1801103	UHMWPE
1801104	UHMWPE
1801105	UHMWPE
1801106	UHMWPE
1801107	UHMWPE
1801111	UHMWPE
1801112	UHMWPE
1801113	UHMWPE
1801114	UHMWPE
1801115	UHMWPE
1801116	UHMWPE
1801117	UHMWPE
1801131	UHMWPE
1801132	UHMWPE
1801133	UHMWPE
1801134	UHMWPE
1801135	UHMWPE
1801136	UHMWPE
1801137	UHMWPE
1801141	UHMWPE
1801142	UHMWPE
1801143	UHMWPE
1801144	UHMWPE
1801145	UHMWPE
1801146	UHMWPE
1801147	UHMWPE
1801171	UHMWPE
1801172	UHMWPE
1801173	UHMWPE
1801174	UHMWPE
1801175	UHMWPE
1801176	UHMWPE
1801177	UHMWPE
1801181	UHMWPE
1801182	UHMWPE
1801183	UHMWPE
1801184	UHMWPE
1801185	UHMWPE

1801186	UHMWPE
1801187	UHMWPE
1801201	UHMWPE
1801202	UHMWPE
1801203	UHMWPE
1801204	UHMWPE
1801205	UHMWPE
1801206	UHMWPE
1801207	UHMWPE
1801211	UHMWPE
1801212	UHMWPE
1801213	UHMWPE
1801214	UHMWPE
1801215	UHMWPE
1801216	UHMWPE
1801217	UHMWPE
1801231	UHMWPE
1801232	UHMWPE
1801233	UHMWPE
1801234	UHMWPE
1801235	UHMWPE
1801236	UHMWPE
1801237	UHMWPE
1801241	UHMWPE
1801242	UHMWPE
1801243	UHMWPE
1801244	UHMWPE
1801245	UHMWPE
1801246	UHMWPE
1801247	UHMWPE
1801271	UHMWPE
1801272	UHMWPE
1801273	UHMWPE
1801274	UHMWPE
1801275	UHMWPE
1801276	UHMWPE
1801277	UHMWPE
1801281	UHMWPE
1801282	UHMWPE
1801283	UHMWPE
1801284	UHMWPE
1801285	UHMWPE
1801286	UHMWPE

1801287	UHMWPE
1801301	UHMWPE
1801302	UHMWPE
1801303	UHMWPE
1801304	UHMWPE
1801305	UHMWPE
1801306	UHMWPE
1801307	UHMWPE
1801311	UHMWPE
1801312	UHMWPE
1801313	UHMWPE
1801314	UHMWPE
1801315	UHMWPE
1801316	UHMWPE
1801317	UHMWPE
1801331	UHMWPE
1801332	UHMWPE
1801333	UHMWPE
1801334	UHMWPE
1801335	UHMWPE
1801336	UHMWPE
1801337	UHMWPE
1801341	UHMWPE
1801342	UHMWPE
1801343	UHMWPE
1801344	UHMWPE
1801345	UHMWPE
1801346	UHMWPE
1801347	UHMWPE
1801371	UHMWPE
1801372	UHMWPE
1801373	UHMWPE
1801374	UHMWPE
1801375	UHMWPE
1801376	UHMWPE
1801377	UHMWPE
1801381	UHMWPE
1801382	UHMWPE
1801383	UHMWPE
1801384	UHMWPE
1801385	UHMWPE
1801386	UHMWPE
1801387	UHMWPE

1801401	UHMWPE
1801402	UHMWPE
1801403	UHMWPE
1801404	UHMWPE
1801405	UHMWPE
1801406	UHMWPE
1801407	UHMWPE
1801411	UHMWPE
1801412	UHMWPE
1801413	UHMWPE
1801414	UHMWPE
1801415	UHMWPE
1801416	UHMWPE
1801417	UHMWPE
1801431	UHMWPE
1801432	UHMWPE
1801433	UHMWPE
1801434	UHMWPE
1801435	UHMWPE
1801436	UHMWPE
1801437	UHMWPE
1801441	UHMWPE
1801442	UHMWPE
1801443	UHMWPE
1801444	UHMWPE
1801445	UHMWPE
1801446	UHMWPE
1801447	UHMWPE
1801471	UHMWPE
1801472	UHMWPE
1801473	UHMWPE
1801474	UHMWPE
1801475	UHMWPE
1801476	UHMWPE
1801477	UHMWPE
1801481	UHMWPE
1801482	UHMWPE
1801483	UHMWPE
1801484	UHMWPE
1801485	UHMWPE
1801486	UHMWPE
1801487	UHMWPE
1801501	UHMWPE

1801502	UHMWPE
1801503	UHMWPE
1801504	UHMWPE
1801505	UHMWPE
1801506	UHMWPE
1801507	UHMWPE
1801511	UHMWPE
1801512	UHMWPE
1801513	UHMWPE
1801514	UHMWPE
1801515	UHMWPE
1801516	UHMWPE
1801517	UHMWPE
1801531	UHMWPE
1801532	UHMWPE
1801533	UHMWPE
1801534	UHMWPE
1801535	UHMWPE
1801536	UHMWPE
1801537	UHMWPE
1801541	UHMWPE
1801542	UHMWPE
1801543	UHMWPE
1801544	UHMWPE
1801545	UHMWPE
1801546	UHMWPE
1801547	UHMWPE
1801571	UHMWPE
1801572	UHMWPE
1801573	UHMWPE
1801574	UHMWPE
1801575	UHMWPE
1801576	UHMWPE
1801577	UHMWPE
1801581	UHMWPE
1801582	UHMWPE
1801583	UHMWPE
1801584	UHMWPE
1801585	UHMWPE
1801586	UHMWPE
1801587	UHMWPE
1801601	UHMWPE
1801602	UHMWPE

1801603	UHMWPE
1801604	UHMWPE
1801605	UHMWPE
1801606	UHMWPE
1801607	UHMWPE
1801611	UHMWPE
1801612	UHMWPE
1801613	UHMWPE
1801614	UHMWPE
1801615	UHMWPE
1801616	UHMWPE
1801617	UHMWPE
1801631	UHMWPE
1801632	UHMWPE
1801633	UHMWPE
1801634	UHMWPE
1801635	UHMWPE
1801636	UHMWPE
1801637	UHMWPE
1801641	UHMWPE
1801642	UHMWPE
1801643	UHMWPE
1801644	UHMWPE
1801645	UHMWPE
1801646	UHMWPE
1801647	UHMWPE
1801671	UHMWPE
1801672	UHMWPE
1801673	UHMWPE
1801674	UHMWPE
1801675	UHMWPE
1801676	UHMWPE
1801677	UHMWPE
1801681	UHMWPE
1801682	UHMWPE
1801683	UHMWPE
1801684	UHMWPE
1801685	UHMWPE
1801686	UHMWPE
1801687	UHMWPE
1801701	UHMWPE
1801702	UHMWPE
1801703	UHMWPE

1801704	UHMWPE
1801705	UHMWPE
1801706	UHMWPE
1801707	UHMWPE
1801711	UHMWPE
1801712	UHMWPE
1801713	UHMWPE
1801714	UHMWPE
1801715	UHMWPE
1801716	UHMWPE
1801717	UHMWPE
1801731	UHMWPE
1801732	UHMWPE
1801733	UHMWPE
1801734	UHMWPE
1801735	UHMWPE
1801736	UHMWPE
1801737	UHMWPE
1801741	UHMWPE
1801742	UHMWPE
1801743	UHMWPE
1801744	UHMWPE
1801745	UHMWPE
1801746	UHMWPE
1801747	UHMWPE
1801771	UHMWPE
1801772	UHMWPE
1801773	UHMWPE
1801774	UHMWPE
1801775	UHMWPE
1801776	UHMWPE
1801777	UHMWPE
1801781	UHMWPE
1801782	UHMWPE
1801783	UHMWPE
1801784	UHMWPE
1801785	UHMWPE
1801786	UHMWPE
1801787	UHMWPE
1802001	UHMWPE
1802002	UHMWPE
1802003	UHMWPE
1802004	UHMWPE

1802005	UHMWPE
1802006	UHMWPE
1802007	UHMWPE
1802011	UHMWPE
1802012	UHMWPE
1802013	UHMWPE
1802014	UHMWPE
1802015	UHMWPE
1802016	UHMWPE
1802017	UHMWPE
1802031	UHMWPE
1802032	UHMWPE
1802033	UHMWPE
1802034	UHMWPE
1802035	UHMWPE
1802036	UHMWPE
1802037	UHMWPE
1802041	UHMWPE
1802042	UHMWPE
1802043	UHMWPE
1802044	UHMWPE
1802045	UHMWPE
1802046	UHMWPE
1802047	UHMWPE
1802071	UHMWPE
1802072	UHMWPE
1802073	UHMWPE
1802074	UHMWPE
1802075	UHMWPE
1802076	UHMWPE
1802077	UHMWPE
1802081	UHMWPE
1802082	UHMWPE
1802083	UHMWPE
1802084	UHMWPE
1802085	UHMWPE
1802086	UHMWPE
1802087	UHMWPE
1802401	UHMWPE
1802402	UHMWPE
1802403	UHMWPE
1802404	UHMWPE
1802405	UHMWPE

1802406	UHMWPE
1802407	UHMWPE
1802411	UHMWPE
1802412	UHMWPE
1802413	UHMWPE
1802414	UHMWPE
1802415	UHMWPE
1802416	UHMWPE
1802417	UHMWPE
1802431	UHMWPE
1802432	UHMWPE
1802433	UHMWPE
1802434	UHMWPE
1802435	UHMWPE
1802436	UHMWPE
1802437	UHMWPE
1802441	UHMWPE
1802442	UHMWPE
1802443	UHMWPE
1802444	UHMWPE
1802445	UHMWPE
1802446	UHMWPE
1802447	UHMWPE
1802471	UHMWPE
1802472	UHMWPE
1802473	UHMWPE
1802474	UHMWPE
1802475	UHMWPE
1802476	UHMWPE
1802477	UHMWPE
1802481	UHMWPE
1802482	UHMWPE
1802483	UHMWPE
1802484	UHMWPE
1802485	UHMWPE
1802486	UHMWPE
1802487	UHMWPE
1811001	UHMWPE
1811002	UHMWPE
1811003	UHMWPE
1811004	UHMWPE
1811005	UHMWPE
1811006	UHMWPE

1811007	UHMWPE
1811011	UHMWPE
1811012	UHMWPE
1811013	UHMWPE
1811014	UHMWPE
1811015	UHMWPE
1811016	UHMWPE
1811017	UHMWPE
1811201	UHMWPE
1811202	UHMWPE
1811203	UHMWPE
1811204	UHMWPE
1811205	UHMWPE
1811206	UHMWPE
1811207	UHMWPE
1811211	UHMWPE
1811212	UHMWPE
1811213	UHMWPE
1811214	UHMWPE
1811215	UHMWPE
1811216	UHMWPE
1811217	UHMWPE
1811401	UHMWPE
1811402	UHMWPE
1811403	UHMWPE
1811404	UHMWPE
1811405	UHMWPE
1811406	UHMWPE
1811407	UHMWPE
1811411	UHMWPE
1811412	UHMWPE
1811413	UHMWPE
1811414	UHMWPE
1811415	UHMWPE
1811416	UHMWPE
1811417	UHMWPE
1811701	UHMWPE
1811702	UHMWPE
1811703	UHMWPE
1811704	UHMWPE
1811705	UHMWPE
1811706	UHMWPE
1811707	UHMWPE

1811711	UHMWPE
1811712	UHMWPE
1811713	UHMWPE
1811714	UHMWPE
1811715	UHMWPE
1811716	UHMWPE
1811717	UHMWPE
1812001	UHMWPE
1812002	UHMWPE
1812003	UHMWPE
1812004	UHMWPE
1812005	UHMWPE
1812006	UHMWPE
1812007	UHMWPE
1812011	UHMWPE
1812012	UHMWPE
1812013	UHMWPE
1812014	UHMWPE
1812015	UHMWPE
1812016	UHMWPE
1812017	UHMWPE
1812401	UHMWPE
1812402	UHMWPE
1812403	UHMWPE
1812404	UHMWPE
1812405	UHMWPE
1812406	UHMWPE
1812407	UHMWPE
1812411	UHMWPE
1812412	UHMWPE
1812413	UHMWPE
1812414	UHMWPE
1812415	UHMWPE
1812416	UHMWPE
1812417	UHMWPE
1828001	Co-Cr-Mo Alloy
1828002	Co-Cr-Mo Alloy
1828003	Co-Cr-Mo Alloy
1828004	Co-Cr-Mo Alloy
1828005	Co-Cr-Mo Alloy
1828006	Co-Cr-Mo Alloy
1828007	Co-Cr-Mo Alloy
1828012	Co-Cr-Mo Alloy

1828023	Co-Cr-Mo Alloy
1828034	Co-Cr-Mo Alloy
1828045	Co-Cr-Mo Alloy
1828056	Co-Cr-Mo Alloy
1828067	Co-Cr-Mo Alloy
1828401	Co-Cr-Mo Alloy + TiNbN
1828402	Co-Cr-Mo Alloy + TiNbN
1828403	Co-Cr-Mo Alloy + TiNbN
1828404	Co-Cr-Mo Alloy + TiNbN
1828405	Co-Cr-Mo Alloy + TiNbN
1828406	Co-Cr-Mo Alloy + TiNbN
1828407	Co-Cr-Mo Alloy + TiNbN
1828412	Co-Cr-Mo Alloy + TiNbN
1828423	Co-Cr-Mo Alloy + TiNbN
1828434	Co-Cr-Mo Alloy + TiNbN
1828445	Co-Cr-Mo Alloy + TiNbN
1828456	Co-Cr-Mo Alloy + TiNbN
1828467	Co-Cr-Mo Alloy + TiNbN
1828501	Co-Cr-Mo Alloy
1828502	Co-Cr-Mo Alloy
1828503	Co-Cr-Mo Alloy
1828504	Co-Cr-Mo Alloy
1828505	Co-Cr-Mo Alloy
1828506	Co-Cr-Mo Alloy
1828507	Co-Cr-Mo Alloy
1828512	Co-Cr-Mo Alloy
1828523	Co-Cr-Mo Alloy
1828534	Co-Cr-Mo Alloy
1828545	Co-Cr-Mo Alloy
1828556	Co-Cr-Mo Alloy
1828567	Co-Cr-Mo Alloy
1828901	Co-Cr-Mo Alloy + TiNbN
1828902	Co-Cr-Mo Alloy + TiNbN
1828903	Co-Cr-Mo Alloy + TiNbN
1828904	Co-Cr-Mo Alloy + TiNbN
1828905	Co-Cr-Mo Alloy + TiNbN
1828906	Co-Cr-Mo Alloy + TiNbN
1828907	Co-Cr-Mo Alloy + TiNbN
1828912	Co-Cr-Mo Alloy + TiNbN
1828923	Co-Cr-Mo Alloy + TiNbN
1828934	Co-Cr-Mo Alloy + TiNbN
1828945	Co-Cr-Mo Alloy + TiNbN
1828956	Co-Cr-Mo Alloy + TiNbN

1828967	Co-Cr-Mo Alloy + TiNbN
1831001	UHMWPE
1831002	UHMWPE
1831003	UHMWPE
1831004	UHMWPE
1831005	UHMWPE
1831006	UHMWPE
1831007	UHMWPE
1831031	UHMWPE
1831032	UHMWPE
1831033	UHMWPE
1831034	UHMWPE
1831035	UHMWPE
1831036	UHMWPE
1831037	UHMWPE
1831041	UHMWPE
1831042	UHMWPE
1831043	UHMWPE
1831044	UHMWPE
1831045	UHMWPE
1831046	UHMWPE
1831047	UHMWPE
1831101	UHMWPE
1831102	UHMWPE
1831103	UHMWPE
1831104	UHMWPE
1831105	UHMWPE
1831106	UHMWPE
1831107	UHMWPE
1831131	UHMWPE
1831132	UHMWPE
1831133	UHMWPE
1831134	UHMWPE
1831135	UHMWPE
1831136	UHMWPE
1831137	UHMWPE
1831141	UHMWPE
1831142	UHMWPE
1831143	UHMWPE
1831144	UHMWPE
1831145	UHMWPE
1831146	UHMWPE
1831147	UHMWPE

1831201	UHMWPE
1831202	UHMWPE
1831203	UHMWPE
1831204	UHMWPE
1831205	UHMWPE
1831206	UHMWPE
1831207	UHMWPE
1831231	UHMWPE
1831232	UHMWPE
1831233	UHMWPE
1831234	UHMWPE
1831235	UHMWPE
1831236	UHMWPE
1831237	UHMWPE
1831241	UHMWPE
1831242	UHMWPE
1831243	UHMWPE
1831244	UHMWPE
1831245	UHMWPE
1831246	UHMWPE
1831247	UHMWPE
1831301	UHMWPE
1831302	UHMWPE
1831303	UHMWPE
1831304	UHMWPE
1831305	UHMWPE
1831306	UHMWPE
1831307	UHMWPE
1831331	UHMWPE
1831332	UHMWPE
1831333	UHMWPE
1831334	UHMWPE
1831335	UHMWPE
1831336	UHMWPE
1831337	UHMWPE
1831341	UHMWPE
1831342	UHMWPE
1831343	UHMWPE
1831344	UHMWPE
1831345	UHMWPE
1831346	UHMWPE
1831347	UHMWPE
1831401	UHMWPE

1831402	UHMWPE
1831403	UHMWPE
1831404	UHMWPE
1831405	UHMWPE
1831406	UHMWPE
1831407	UHMWPE
1831431	UHMWPE
1831432	UHMWPE
1831433	UHMWPE
1831434	UHMWPE
1831435	UHMWPE
1831436	UHMWPE
1831437	UHMWPE
1831441	UHMWPE
1831442	UHMWPE
1831443	UHMWPE
1831444	UHMWPE
1831445	UHMWPE
1831446	UHMWPE
1831447	UHMWPE
1831501	UHMWPE
1831502	UHMWPE
1831503	UHMWPE
1831504	UHMWPE
1831505	UHMWPE
1831506	UHMWPE
1831507	UHMWPE
1831531	UHMWPE
1831532	UHMWPE
1831533	UHMWPE
1831534	UHMWPE
1831535	UHMWPE
1831536	UHMWPE
1831537	UHMWPE
1831541	UHMWPE
1831542	UHMWPE
1831543	UHMWPE
1831544	UHMWPE
1831545	UHMWPE
1831546	UHMWPE
1831547	UHMWPE
1831601	UHMWPE
1831602	UHMWPE

1831603	UHMWPE
1831604	UHMWPE
1831605	UHMWPE
1831606	UHMWPE
1831607	UHMWPE
1831631	UHMWPE
1831632	UHMWPE
1831633	UHMWPE
1831634	UHMWPE
1831635	UHMWPE
1831636	UHMWPE
1831637	UHMWPE
1831641	UHMWPE
1831642	UHMWPE
1831643	UHMWPE
1831644	UHMWPE
1831645	UHMWPE
1831646	UHMWPE
1831647	UHMWPE
1831701	UHMWPE
1831702	UHMWPE
1831703	UHMWPE
1831704	UHMWPE
1831705	UHMWPE
1831706	UHMWPE
1831707	UHMWPE
1831731	UHMWPE
1831732	UHMWPE
1831733	UHMWPE
1831734	UHMWPE
1831735	UHMWPE
1831736	UHMWPE
1831737	UHMWPE
1831741	UHMWPE
1831742	UHMWPE
1831743	UHMWPE
1831744	UHMWPE
1831745	UHMWPE
1831746	UHMWPE
1831747	UHMWPE
1832001	UHMWPE
1832002	UHMWPE
1832003	UHMWPE

1832004	UHMWPE
1832005	UHMWPE
1832006	UHMWPE
1832007	UHMWPE
1832031	UHMWPE
1832032	UHMWPE
1832033	UHMWPE
1832034	UHMWPE
1832035	UHMWPE
1832036	UHMWPE
1832037	UHMWPE
1832041	UHMWPE
1832042	UHMWPE
1832043	UHMWPE
1832044	UHMWPE
1832045	UHMWPE
1832046	UHMWPE
1832047	UHMWPE
1832401	UHMWPE
1832402	UHMWPE
1832403	UHMWPE
1832404	UHMWPE
1832405	UHMWPE
1832406	UHMWPE
1832407	UHMWPE
1832431	UHMWPE
1832432	UHMWPE
1832433	UHMWPE
1832434	UHMWPE
1832435	UHMWPE
1832436	UHMWPE
1832437	UHMWPE
1832441	UHMWPE
1832442	UHMWPE
1832443	UHMWPE
1832444	UHMWPE
1832445	UHMWPE
1832446	UHMWPE
1832447	UHMWPE
1840901	UHMWPE
1840902	UHMWPE
1840903	UHMWPE
1840904	UHMWPE

1840905	UHMWPE
1840906	UHMWPE
1840911	UHMWPE
1840912	UHMWPE
1840913	UHMWPE
1840914	UHMWPE
1840915	UHMWPE
1840916	UHMWPE
1841001	UHMWPE
1841002	UHMWPE
1841003	UHMWPE
1841004	UHMWPE
1841005	UHMWPE
1841006	UHMWPE
1841011	UHMWPE
1841012	UHMWPE
1841013	UHMWPE
1841014	UHMWPE
1841015	UHMWPE
1841016	UHMWPE
1841101	UHMWPE
1841102	UHMWPE
1841103	UHMWPE
1841104	UHMWPE
1841105	UHMWPE
1841106	UHMWPE
1841111	UHMWPE
1841112	UHMWPE
1841113	UHMWPE
1841114	UHMWPE
1841115	UHMWPE
1841116	UHMWPE
1841201	UHMWPE
1841202	UHMWPE
1841203	UHMWPE
1841204	UHMWPE
1841205	UHMWPE
1841206	UHMWPE
1841211	UHMWPE
1841212	UHMWPE
1841213	UHMWPE
1841214	UHMWPE
1841215	UHMWPE

1841216	UHMWPE
1841301	UHMWPE
1841302	UHMWPE
1841303	UHMWPE
1841304	UHMWPE
1841305	UHMWPE
1841306	UHMWPE
1841311	UHMWPE
1841312	UHMWPE
1841313	UHMWPE
1841314	UHMWPE
1841315	UHMWPE
1841316	UHMWPE
1841401	UHMWPE
1841402	UHMWPE
1841403	UHMWPE
1841404	UHMWPE
1841405	UHMWPE
1841406	UHMWPE
1841411	UHMWPE
1841412	UHMWPE
1841413	UHMWPE
1841414	UHMWPE
1841415	UHMWPE
1841416	UHMWPE
1841601	UHMWPE
1841602	UHMWPE
1841603	UHMWPE
1841604	UHMWPE
1841605	UHMWPE
1841606	UHMWPE
1841611	UHMWPE
1841612	UHMWPE
1841613	UHMWPE
1841614	UHMWPE
1841615	UHMWPE
1841616	UHMWPE
1850901	UHMWPE
1850902	UHMWPE
1850903	UHMWPE
1850904	UHMWPE
1850905	UHMWPE
1850906	UHMWPE

1850911	UHMWPE
1850912	UHMWPE
1850913	UHMWPE
1850914	UHMWPE
1850915	UHMWPE
1850916	UHMWPE
1851001	UHMWPE
1851002	UHMWPE
1851003	UHMWPE
1851004	UHMWPE
1851005	UHMWPE
1851006	UHMWPE
1851011	UHMWPE
1851012	UHMWPE
1851013	UHMWPE
1851014	UHMWPE
1851015	UHMWPE
1851016	UHMWPE
1851101	UHMWPE
1851102	UHMWPE
1851103	UHMWPE
1851104	UHMWPE
1851105	UHMWPE
1851106	UHMWPE
1851111	UHMWPE
1851112	UHMWPE
1851113	UHMWPE
1851114	UHMWPE
1851115	UHMWPE
1851116	UHMWPE
1851201	UHMWPE
1851202	UHMWPE
1851203	UHMWPE
1851204	UHMWPE
1851205	UHMWPE
1851206	UHMWPE
1851211	UHMWPE
1851212	UHMWPE
1851213	UHMWPE
1851214	UHMWPE
1851215	UHMWPE
1851216	UHMWPE
1851301	UHMWPE

1851302	UHMWPE
1851303	UHMWPE
1851304	UHMWPE
1851305	UHMWPE
1851306	UHMWPE
1851311	UHMWPE
1851312	UHMWPE
1851313	UHMWPE
1851314	UHMWPE
1851315	UHMWPE
1851316	UHMWPE
1851401	UHMWPE
1851402	UHMWPE
1851403	UHMWPE
1851404	UHMWPE
1851405	UHMWPE
1851406	UHMWPE
1851411	UHMWPE
1851412	UHMWPE
1851413	UHMWPE
1851414	UHMWPE
1851415	UHMWPE
1851416	UHMWPE
1851601	UHMWPE
1851602	UHMWPE
1851603	UHMWPE
1851604	UHMWPE
1851605	UHMWPE
1851606	UHMWPE
1851611	UHMWPE
1851612	UHMWPE
1851613	UHMWPE
1851614	UHMWPE
1851615	UHMWPE
1851616	UHMWPE
1861001	UHMWPE
1861002	UHMWPE
1861003	UHMWPE
1861004	UHMWPE
1861005	UHMWPE
1861006	UHMWPE
1861007	UHMWPE
1861031	UHMWPE

1861032	UHMWPE
1861033	UHMWPE
1861034	UHMWPE
1861035	UHMWPE
1861036	UHMWPE
1861037	UHMWPE
1861041	UHMWPE
1861042	UHMWPE
1861043	UHMWPE
1861044	UHMWPE
1861045	UHMWPE
1861046	UHMWPE
1861047	UHMWPE
1861101	UHMWPE
1861102	UHMWPE
1861103	UHMWPE
1861104	UHMWPE
1861105	UHMWPE
1861106	UHMWPE
1861107	UHMWPE
1861131	UHMWPE
1861132	UHMWPE
1861133	UHMWPE
1861134	UHMWPE
1861135	UHMWPE
1861136	UHMWPE
1861137	UHMWPE
1861141	UHMWPE
1861142	UHMWPE
1861143	UHMWPE
1861144	UHMWPE
1861145	UHMWPE
1861146	UHMWPE
1861147	UHMWPE
1861201	UHMWPE
1861202	UHMWPE
1861203	UHMWPE
1861204	UHMWPE
1861205	UHMWPE
1861206	UHMWPE
1861207	UHMWPE
1861231	UHMWPE
1861232	UHMWPE

1861233	UHMWPE
1861234	UHMWPE
1861235	UHMWPE
1861236	UHMWPE
1861237	UHMWPE
1861241	UHMWPE
1861242	UHMWPE
1861243	UHMWPE
1861244	UHMWPE
1861245	UHMWPE
1861246	UHMWPE
1861247	UHMWPE
1861301	UHMWPE
1861302	UHMWPE
1861303	UHMWPE
1861304	UHMWPE
1861305	UHMWPE
1861306	UHMWPE
1861307	UHMWPE
1861331	UHMWPE
1861332	UHMWPE
1861333	UHMWPE
1861334	UHMWPE
1861335	UHMWPE
1861336	UHMWPE
1861337	UHMWPE
1861341	UHMWPE
1861342	UHMWPE
1861343	UHMWPE
1861344	UHMWPE
1861345	UHMWPE
1861346	UHMWPE
1861347	UHMWPE
1861401	UHMWPE
1861402	UHMWPE
1861403	UHMWPE
1861404	UHMWPE
1861405	UHMWPE
1861406	UHMWPE
1861407	UHMWPE
1861431	UHMWPE
1861432	UHMWPE
1861433	UHMWPE

1861434	UHMWPE
1861435	UHMWPE
1861436	UHMWPE
1861437	UHMWPE
1861441	UHMWPE
1861442	UHMWPE
1861443	UHMWPE
1861444	UHMWPE
1861445	UHMWPE
1861446	UHMWPE
1861447	UHMWPE
1861501	UHMWPE
1861502	UHMWPE
1861503	UHMWPE
1861504	UHMWPE
1861505	UHMWPE
1861506	UHMWPE
1861507	UHMWPE
1861531	UHMWPE
1861532	UHMWPE
1861533	UHMWPE
1861534	UHMWPE
1861535	UHMWPE
1861536	UHMWPE
1861537	UHMWPE
1861541	UHMWPE
1861542	UHMWPE
1861543	UHMWPE
1861544	UHMWPE
1861545	UHMWPE
1861546	UHMWPE
1861547	UHMWPE
1861601	UHMWPE
1861602	UHMWPE
1861603	UHMWPE
1861604	UHMWPE
1861605	UHMWPE
1861606	UHMWPE
1861607	UHMWPE
1861631	UHMWPE
1861632	UHMWPE
1861633	UHMWPE
1861634	UHMWPE

1861635	UHMWPE
1861636	UHMWPE
1861637	UHMWPE
1861641	UHMWPE
1861642	UHMWPE
1861643	UHMWPE
1861644	UHMWPE
1861645	UHMWPE
1861646	UHMWPE
1861647	UHMWPE
1861701	UHMWPE
1861702	UHMWPE
1861703	UHMWPE
1861704	UHMWPE
1861705	UHMWPE
1861706	UHMWPE
1861707	UHMWPE
1861731	UHMWPE
1861732	UHMWPE
1861733	UHMWPE
1861734	UHMWPE
1861735	UHMWPE
1861736	UHMWPE
1861737	UHMWPE
1861741	UHMWPE
1861742	UHMWPE
1861743	UHMWPE
1861744	UHMWPE
1861745	UHMWPE
1861746	UHMWPE
1861747	UHMWPE
1862001	UHMWPE
1862002	UHMWPE
1862003	UHMWPE
1862004	UHMWPE
1862005	UHMWPE
1862006	UHMWPE
1862007	UHMWPE
1862031	UHMWPE
1862032	UHMWPE
1862033	UHMWPE
1862034	UHMWPE
1862035	UHMWPE

1862036	UHMWPE
1862037	UHMWPE
1862041	UHMWPE
1862042	UHMWPE
1862043	UHMWPE
1862044	UHMWPE
1862045	UHMWPE
1862046	UHMWPE
1862047	UHMWPE
1862401	UHMWPE
1862402	UHMWPE
1862403	UHMWPE
1862404	UHMWPE
1862405	UHMWPE
1862406	UHMWPE
1862407	UHMWPE
1862431	UHMWPE
1862432	UHMWPE
1862433	UHMWPE
1862434	UHMWPE
1862435	UHMWPE
1862436	UHMWPE
1862437	UHMWPE
1862441	UHMWPE
1862442	UHMWPE
1862443	UHMWPE
1862444	UHMWPE
1862445	UHMWPE
1862446	UHMWPE
1862447	UHMWPE
1871001	UHMWPE
1871002	UHMWPE
1871003	UHMWPE
1871004	UHMWPE
1871005	UHMWPE
1871006	UHMWPE
1871021	UHMWPE
1871022	UHMWPE
1871023	UHMWPE
1871024	UHMWPE
1871025	UHMWPE
1871026	UHMWPE
1871027	UHMWPE

1871031	UHMWPE
1871032	UHMWPE
1871033	UHMWPE
1871034	UHMWPE
1871035	UHMWPE
1871036	UHMWPE
1871041	UHMWPE
1871042	UHMWPE
1871043	UHMWPE
1871044	UHMWPE
1871045	UHMWPE
1871046	UHMWPE
1871051	UHMWPE
1871052	UHMWPE
1871053	UHMWPE
1871054	UHMWPE
1871055	UHMWPE
1871056	UHMWPE
1871061	UHMWPE
1871062	UHMWPE
1871063	UHMWPE
1871064	UHMWPE
1871065	UHMWPE
1871066	UHMWPE
1871071	UHMWPE
1871072	UHMWPE
1871073	UHMWPE
1871074	UHMWPE
1871075	UHMWPE
1871076	UHMWPE
1871081	UHMWPE
1871082	UHMWPE
1871083	UHMWPE
1871084	UHMWPE
1871085	UHMWPE
1871086	UHMWPE
1871087	UHMWPE
1871091	UHMWPE
1871092	UHMWPE
1871093	UHMWPE
1871094	UHMWPE
1871095	UHMWPE
1871096	UHMWPE

1871097	UHMWPE
1871121	UHMWPE
1871122	UHMWPE
1871123	UHMWPE
1871124	UHMWPE
1871125	UHMWPE
1871126	UHMWPE
1871127	UHMWPE
1871181	UHMWPE
1871182	UHMWPE
1871183	UHMWPE
1871184	UHMWPE
1871185	UHMWPE
1871186	UHMWPE
1871187	UHMWPE
1871191	UHMWPE
1871192	UHMWPE
1871193	UHMWPE
1871194	UHMWPE
1871195	UHMWPE
1871196	UHMWPE
1871197	UHMWPE
1871201	UHMWPE
1871202	UHMWPE
1871203	UHMWPE
1871204	UHMWPE
1871205	UHMWPE
1871206	UHMWPE
1871221	UHMWPE
1871222	UHMWPE
1871223	UHMWPE
1871224	UHMWPE
1871225	UHMWPE
1871226	UHMWPE
1871227	UHMWPE
1871231	UHMWPE
1871232	UHMWPE
1871233	UHMWPE
1871234	UHMWPE
1871235	UHMWPE
1871236	UHMWPE
1871241	UHMWPE
1871242	UHMWPE

1871243	UHMWPE
1871244	UHMWPE
1871245	UHMWPE
1871246	UHMWPE
1871251	UHMWPE
1871252	UHMWPE
1871253	UHMWPE
1871254	UHMWPE
1871255	UHMWPE
1871256	UHMWPE
1871261	UHMWPE
1871262	UHMWPE
1871263	UHMWPE
1871264	UHMWPE
1871265	UHMWPE
1871266	UHMWPE
1871271	UHMWPE
1871272	UHMWPE
1871273	UHMWPE
1871274	UHMWPE
1871275	UHMWPE
1871276	UHMWPE
1871281	UHMWPE
1871282	UHMWPE
1871283	UHMWPE
1871284	UHMWPE
1871285	UHMWPE
1871286	UHMWPE
1871287	UHMWPE
1871291	UHMWPE
1871292	UHMWPE
1871293	UHMWPE
1871294	UHMWPE
1871295	UHMWPE
1871296	UHMWPE
1871297	UHMWPE
1871321	UHMWPE
1871322	UHMWPE
1871323	UHMWPE
1871324	UHMWPE
1871325	UHMWPE
1871326	UHMWPE
1871327	UHMWPE

1871381	UHMWPE
1871382	UHMWPE
1871383	UHMWPE
1871384	UHMWPE
1871385	UHMWPE
1871386	UHMWPE
1871387	UHMWPE
1871391	UHMWPE
1871392	UHMWPE
1871393	UHMWPE
1871394	UHMWPE
1871395	UHMWPE
1871396	UHMWPE
1871397	UHMWPE
1871401	UHMWPE
1871402	UHMWPE
1871403	UHMWPE
1871404	UHMWPE
1871405	UHMWPE
1871406	UHMWPE
1871421	UHMWPE
1871422	UHMWPE
1871423	UHMWPE
1871424	UHMWPE
1871425	UHMWPE
1871426	UHMWPE
1871427	UHMWPE
1871431	UHMWPE
1871432	UHMWPE
1871433	UHMWPE
1871434	UHMWPE
1871435	UHMWPE
1871436	UHMWPE
1871441	UHMWPE
1871442	UHMWPE
1871443	UHMWPE
1871444	UHMWPE
1871445	UHMWPE
1871446	UHMWPE
1871451	UHMWPE
1871452	UHMWPE
1871453	UHMWPE
1871454	UHMWPE

1871455	UHMWPE
1871456	UHMWPE
1871461	UHMWPE
1871462	UHMWPE
1871463	UHMWPE
1871464	UHMWPE
1871465	UHMWPE
1871466	UHMWPE
1871471	UHMWPE
1871472	UHMWPE
1871473	UHMWPE
1871474	UHMWPE
1871475	UHMWPE
1871476	UHMWPE
1871481	UHMWPE
1871482	UHMWPE
1871483	UHMWPE
1871484	UHMWPE
1871485	UHMWPE
1871486	UHMWPE
1871487	UHMWPE
1871491	UHMWPE
1871492	UHMWPE
1871493	UHMWPE
1871494	UHMWPE
1871495	UHMWPE
1871496	UHMWPE
1871497	UHMWPE
1871521	UHMWPE
1871522	UHMWPE
1871523	UHMWPE
1871524	UHMWPE
1871525	UHMWPE
1871526	UHMWPE
1871527	UHMWPE
1871581	UHMWPE
1871582	UHMWPE
1871583	UHMWPE
1871584	UHMWPE
1871585	UHMWPE
1871586	UHMWPE
1871587	UHMWPE
1871591	UHMWPE

1871592	UHMWPE
1871593	UHMWPE
1871594	UHMWPE
1871595	UHMWPE
1871596	UHMWPE
1871597	UHMWPE
1871601	UHMWPE
1871602	UHMWPE
1871603	UHMWPE
1871604	UHMWPE
1871605	UHMWPE
1871606	UHMWPE
1871621	UHMWPE
1871622	UHMWPE
1871623	UHMWPE
1871624	UHMWPE
1871625	UHMWPE
1871626	UHMWPE
1871627	UHMWPE
1871631	UHMWPE
1871632	UHMWPE
1871633	UHMWPE
1871634	UHMWPE
1871635	UHMWPE
1871636	UHMWPE
1871641	UHMWPE
1871642	UHMWPE
1871643	UHMWPE
1871644	UHMWPE
1871645	UHMWPE
1871646	UHMWPE
1871651	UHMWPE
1871652	UHMWPE
1871653	UHMWPE
1871654	UHMWPE
1871655	UHMWPE
1871656	UHMWPE
1871661	UHMWPE
1871662	UHMWPE
1871663	UHMWPE
1871664	UHMWPE
1871665	UHMWPE
1871666	UHMWPE

1871671	UHMWPE
1871672	UHMWPE
1871673	UHMWPE
1871674	UHMWPE
1871675	UHMWPE
1871676	UHMWPE
1871681	UHMWPE
1871682	UHMWPE
1871683	UHMWPE
1871684	UHMWPE
1871685	UHMWPE
1871686	UHMWPE
1871687	UHMWPE
1871691	UHMWPE
1871692	UHMWPE
1871693	UHMWPE
1871694	UHMWPE
1871695	UHMWPE
1871696	UHMWPE
1871697	UHMWPE
1871701	UHMWPE
1871702	UHMWPE
1871703	UHMWPE
1871704	UHMWPE
1871705	UHMWPE
1871706	UHMWPE
1871721	UHMWPE
1871722	UHMWPE
1871723	UHMWPE
1871724	UHMWPE
1871725	UHMWPE
1871726	UHMWPE
1871727	UHMWPE
1871731	UHMWPE
1871732	UHMWPE
1871733	UHMWPE
1871734	UHMWPE
1871735	UHMWPE
1871736	UHMWPE
1871741	UHMWPE
1871742	UHMWPE
1871743	UHMWPE
1871744	UHMWPE

1871745	UHMWPE
1871746	UHMWPE
1871751	UHMWPE
1871752	UHMWPE
1871753	UHMWPE
1871754	UHMWPE
1871755	UHMWPE
1871756	UHMWPE
1871761	UHMWPE
1871762	UHMWPE
1871763	UHMWPE
1871764	UHMWPE
1871765	UHMWPE
1871766	UHMWPE
1871771	UHMWPE
1871772	UHMWPE
1871773	UHMWPE
1871774	UHMWPE
1871775	UHMWPE
1871776	UHMWPE
1871781	UHMWPE
1871782	UHMWPE
1871783	UHMWPE
1871784	UHMWPE
1871785	UHMWPE
1871786	UHMWPE
1871787	UHMWPE
1871791	UHMWPE
1871792	UHMWPE
1871793	UHMWPE
1871794	UHMWPE
1871795	UHMWPE
1871796	UHMWPE
1871797	UHMWPE
1871801	UHMWPE
1871802	UHMWPE
1871803	UHMWPE
1871804	UHMWPE
1871805	UHMWPE
1871806	UHMWPE
1871831	UHMWPE
1871832	UHMWPE
1871833	UHMWPE

1871834	UHMWPE
1871835	UHMWPE
1871836	UHMWPE
1871841	UHMWPE
1871842	UHMWPE
1871843	UHMWPE
1871844	UHMWPE
1871845	UHMWPE
1871846	UHMWPE
1871851	UHMWPE
1871852	UHMWPE
1871853	UHMWPE
1871854	UHMWPE
1871855	UHMWPE
1871856	UHMWPE
1871861	UHMWPE
1871862	UHMWPE
1871863	UHMWPE
1871864	UHMWPE
1871865	UHMWPE
1871866	UHMWPE
1871871	UHMWPE
1871872	UHMWPE
1871873	UHMWPE
1871874	UHMWPE
1871875	UHMWPE
1871876	UHMWPE
1872001	UHMWPE
1872002	UHMWPE
1872003	UHMWPE
1872004	UHMWPE
1872005	UHMWPE
1872006	UHMWPE
1872021	UHMWPE
1872022	UHMWPE
1872023	UHMWPE
1872024	UHMWPE
1872025	UHMWPE
1872026	UHMWPE
1872027	UHMWPE
1872031	UHMWPE
1872032	UHMWPE
1872033	UHMWPE

1872034	UHMWPE
1872035	UHMWPE
1872036	UHMWPE
1872041	UHMWPE
1872042	UHMWPE
1872043	UHMWPE
1872044	UHMWPE
1872045	UHMWPE
1872046	UHMWPE
1872051	UHMWPE
1872052	UHMWPE
1872053	UHMWPE
1872054	UHMWPE
1872055	UHMWPE
1872056	UHMWPE
1872061	UHMWPE
1872062	UHMWPE
1872063	UHMWPE
1872064	UHMWPE
1872065	UHMWPE
1872066	UHMWPE
1872071	UHMWPE
1872072	UHMWPE
1872073	UHMWPE
1872074	UHMWPE
1872075	UHMWPE
1872076	UHMWPE
1872081	UHMWPE
1872082	UHMWPE
1872083	UHMWPE
1872084	UHMWPE
1872085	UHMWPE
1872086	UHMWPE
1872087	UHMWPE
1872091	UHMWPE
1872092	UHMWPE
1872093	UHMWPE
1872094	UHMWPE
1872095	UHMWPE
1872096	UHMWPE
1872097	UHMWPE
1872201	UHMWPE
1872202	UHMWPE

1872203	UHMWPE
1872204	UHMWPE
1872205	UHMWPE
1872206	UHMWPE
1872231	UHMWPE
1872232	UHMWPE
1872233	UHMWPE
1872234	UHMWPE
1872235	UHMWPE
1872236	UHMWPE
1872241	UHMWPE
1872242	UHMWPE
1872243	UHMWPE
1872244	UHMWPE
1872245	UHMWPE
1872246	UHMWPE
1872251	UHMWPE
1872252	UHMWPE
1872253	UHMWPE
1872254	UHMWPE
1872255	UHMWPE
1872256	UHMWPE
1872261	UHMWPE
1872262	UHMWPE
1872263	UHMWPE
1872264	UHMWPE
1872265	UHMWPE
1872266	UHMWPE
1872271	UHMWPE
1872272	UHMWPE
1872273	UHMWPE
1872274	UHMWPE
1872275	UHMWPE
1872276	UHMWPE
1872401	UHMWPE
1872402	UHMWPE
1872403	UHMWPE
1872404	UHMWPE
1872405	UHMWPE
1872406	UHMWPE
1872421	UHMWPE
1872422	UHMWPE
1872423	UHMWPE

1872424	UHMWPE
1872425	UHMWPE
1872426	UHMWPE
1872427	UHMWPE
1872431	UHMWPE
1872432	UHMWPE
1872433	UHMWPE
1872434	UHMWPE
1872435	UHMWPE
1872436	UHMWPE
1872441	UHMWPE
1872442	UHMWPE
1872443	UHMWPE
1872444	UHMWPE
1872445	UHMWPE
1872446	UHMWPE
1872451	UHMWPE
1872452	UHMWPE
1872453	UHMWPE
1872454	UHMWPE
1872455	UHMWPE
1872456	UHMWPE
1872461	UHMWPE
1872462	UHMWPE
1872463	UHMWPE
1872464	UHMWPE
1872465	UHMWPE
1872466	UHMWPE
1872471	UHMWPE
1872472	UHMWPE
1872473	UHMWPE
1872474	UHMWPE
1872475	UHMWPE
1872476	UHMWPE
1872481	UHMWPE
1872482	UHMWPE
1872483	UHMWPE
1872484	UHMWPE
1872485	UHMWPE
1872486	UHMWPE
1872487	UHMWPE
1872491	UHMWPE
1872492	UHMWPE

1872493	UHMWPE
1872494	UHMWPE
1872495	UHMWPE
1872496	UHMWPE
1872497	UHMWPE
1872601	UHMWPE
1872602	UHMWPE
1872603	UHMWPE
1872604	UHMWPE
1872605	UHMWPE
1872606	UHMWPE
1872631	UHMWPE
1872632	UHMWPE
1872633	UHMWPE
1872634	UHMWPE
1872635	UHMWPE
1872636	UHMWPE
1872641	UHMWPE
1872642	UHMWPE
1872643	UHMWPE
1872644	UHMWPE
1872645	UHMWPE
1872646	UHMWPE
1872651	UHMWPE
1872652	UHMWPE
1872653	UHMWPE
1872654	UHMWPE
1872655	UHMWPE
1872656	UHMWPE
1872661	UHMWPE
1872662	UHMWPE
1872663	UHMWPE
1872664	UHMWPE
1872665	UHMWPE
1872666	UHMWPE
1872671	UHMWPE
1872672	UHMWPE
1872673	UHMWPE
1872674	UHMWPE
1872675	UHMWPE
1872676	UHMWPE
1872801	UHMWPE
1872802	UHMWPE

1872803	UHMWPE
1872804	UHMWPE
1872805	UHMWPE
1872806	UHMWPE
1872831	UHMWPE
1872832	UHMWPE
1872833	UHMWPE
1872834	UHMWPE
1872835	UHMWPE
1872836	UHMWPE
1872841	UHMWPE
1872842	UHMWPE
1872843	UHMWPE
1872844	UHMWPE
1872845	UHMWPE
1872846	UHMWPE
1872851	UHMWPE
1872852	UHMWPE
1872853	UHMWPE
1872854	UHMWPE
1872855	UHMWPE
1872856	UHMWPE
1872861	UHMWPE
1872862	UHMWPE
1872863	UHMWPE
1872864	UHMWPE
1872865	UHMWPE
1872866	UHMWPE
1872871	UHMWPE
1872872	UHMWPE
1872873	UHMWPE
1872874	UHMWPE
1872875	UHMWPE
1872876	UHMWPE
1873001	UHMWPE
1873002	UHMWPE
1873003	UHMWPE
1873004	UHMWPE
1873005	UHMWPE
1873006	UHMWPE
1873031	UHMWPE
1873032	UHMWPE
1873033	UHMWPE

1873034	UHMWPE
1873035	UHMWPE
1873036	UHMWPE
1873041	UHMWPE
1873042	UHMWPE
1873043	UHMWPE
1873044	UHMWPE
1873045	UHMWPE
1873046	UHMWPE
1873051	UHMWPE
1873052	UHMWPE
1873053	UHMWPE
1873054	UHMWPE
1873055	UHMWPE
1873056	UHMWPE
1873061	UHMWPE
1873062	UHMWPE
1873063	UHMWPE
1873064	UHMWPE
1873065	UHMWPE
1873066	UHMWPE
1873071	UHMWPE
1873072	UHMWPE
1873073	UHMWPE
1873074	UHMWPE
1873075	UHMWPE
1873076	UHMWPE
1881001	UHMWPE
1881002	UHMWPE
1881003	UHMWPE
1881004	UHMWPE
1881005	UHMWPE
1881006	UHMWPE
1881031	UHMWPE
1881032	UHMWPE
1881033	UHMWPE
1881034	UHMWPE
1881035	UHMWPE
1881036	UHMWPE
1881041	UHMWPE
1881042	UHMWPE
1881043	UHMWPE
1881044	UHMWPE

1881045	UHMWPE
1881046	UHMWPE
1881051	UHMWPE
1881052	UHMWPE
1881053	UHMWPE
1881054	UHMWPE
1881055	UHMWPE
1881056	UHMWPE
1881201	UHMWPE
1881202	UHMWPE
1881203	UHMWPE
1881204	UHMWPE
1881205	UHMWPE
1881206	UHMWPE
1881231	UHMWPE
1881232	UHMWPE
1881233	UHMWPE
1881234	UHMWPE
1881235	UHMWPE
1881236	UHMWPE
1881241	UHMWPE
1881242	UHMWPE
1881243	UHMWPE
1881244	UHMWPE
1881245	UHMWPE
1881246	UHMWPE
1881251	UHMWPE
1881252	UHMWPE
1881253	UHMWPE
1881254	UHMWPE
1881255	UHMWPE
1881256	UHMWPE
1881401	UHMWPE
1881402	UHMWPE
1881403	UHMWPE
1881404	UHMWPE
1881405	UHMWPE
1881406	UHMWPE
1881431	UHMWPE
1881432	UHMWPE
1881433	UHMWPE
1881434	UHMWPE
1881435	UHMWPE

1881436	UHMWPE
1881441	UHMWPE
1881442	UHMWPE
1881443	UHMWPE
1881444	UHMWPE
1881445	UHMWPE
1881446	UHMWPE
1881451	UHMWPE
1881452	UHMWPE
1881453	UHMWPE
1881454	UHMWPE
1881455	UHMWPE
1881456	UHMWPE
1881601	UHMWPE
1881602	UHMWPE
1881603	UHMWPE
1881604	UHMWPE
1881605	UHMWPE
1881606	UHMWPE
1881631	UHMWPE
1881632	UHMWPE
1881633	UHMWPE
1881634	UHMWPE
1881635	UHMWPE
1881636	UHMWPE
1881641	UHMWPE
1881642	UHMWPE
1881643	UHMWPE
1881644	UHMWPE
1881645	UHMWPE
1881646	UHMWPE
1881701	UHMWPE
1881702	UHMWPE
1881703	UHMWPE
1881704	UHMWPE
1881705	UHMWPE
1881706	UHMWPE
1881731	UHMWPE
1881732	UHMWPE
1881733	UHMWPE
1881734	UHMWPE
1881735	UHMWPE
1881736	UHMWPE

1881741	UHMWPE
1881742	UHMWPE
1881743	UHMWPE
1881744	UHMWPE
1881745	UHMWPE
1881746	UHMWPE
1881751	UHMWPE
1881752	UHMWPE
1881753	UHMWPE
1881754	UHMWPE
1881755	UHMWPE
1881756	UHMWPE
1881801	UHMWPE
1881802	UHMWPE
1881803	UHMWPE
1881804	UHMWPE
1881805	UHMWPE
1881806	UHMWPE
1881831	UHMWPE
1881832	UHMWPE
1881833	UHMWPE
1881834	UHMWPE
1881835	UHMWPE
1881836	UHMWPE
1881841	UHMWPE
1881842	UHMWPE
1881843	UHMWPE
1881844	UHMWPE
1881845	UHMWPE
1881846	UHMWPE
1882001	UHMWPE
1882002	UHMWPE
1882003	UHMWPE
1882004	UHMWPE
1882005	UHMWPE
1882006	UHMWPE
1882031	UHMWPE
1882032	UHMWPE
1882033	UHMWPE
1882034	UHMWPE
1882035	UHMWPE
1882036	UHMWPE
1882041	UHMWPE

1882042	UHMWPE
1882043	UHMWPE
1882044	UHMWPE
1882045	UHMWPE
1882046	UHMWPE
1882051	UHMWPE
1882052	UHMWPE
1882053	UHMWPE
1882054	UHMWPE
1882055	UHMWPE
1882056	UHMWPE
1882201	UHMWPE
1882202	UHMWPE
1882203	UHMWPE
1882204	UHMWPE
1882205	UHMWPE
1882206	UHMWPE
1882231	UHMWPE
1882232	UHMWPE
1882233	UHMWPE
1882234	UHMWPE
1882235	UHMWPE
1882236	UHMWPE
1882241	UHMWPE
1882242	UHMWPE
1882243	UHMWPE
1882244	UHMWPE
1882245	UHMWPE
1882246	UHMWPE
1882401	UHMWPE
1882402	UHMWPE
1882403	UHMWPE
1882404	UHMWPE
1882405	UHMWPE
1882406	UHMWPE
1882431	UHMWPE
1882432	UHMWPE
1882433	UHMWPE
1882434	UHMWPE
1882435	UHMWPE
1882436	UHMWPE
1882441	UHMWPE
1882442	UHMWPE

1882443	UHMWPE
1882444	UHMWPE
1882445	UHMWPE
1882446	UHMWPE
1882451	UHMWPE
1882452	UHMWPE
1882453	UHMWPE
1882454	UHMWPE
1882455	UHMWPE
1882456	UHMWPE
1882601	UHMWPE
1882602	UHMWPE
1882603	UHMWPE
1882604	UHMWPE
1882605	UHMWPE
1882606	UHMWPE
1882631	UHMWPE
1882632	UHMWPE
1882633	UHMWPE
1882634	UHMWPE
1882635	UHMWPE
1882636	UHMWPE
1882641	UHMWPE
1882642	UHMWPE
1882643	UHMWPE
1882644	UHMWPE
1882645	UHMWPE
1882646	UHMWPE
1882801	UHMWPE
1882802	UHMWPE
1882803	UHMWPE
1882804	UHMWPE
1882805	UHMWPE
1882806	UHMWPE
1882831	UHMWPE
1882832	UHMWPE
1882833	UHMWPE
1882834	UHMWPE
1882835	UHMWPE
1882836	UHMWPE
1882841	UHMWPE
1882842	UHMWPE
1882843	UHMWPE

1882844	UHMWPE
1882845	UHMWPE
1882846	UHMWPE
1883001	UHMWPE
1883002	UHMWPE
1883003	UHMWPE
1883004	UHMWPE
1883005	UHMWPE
1883006	UHMWPE
1883031	UHMWPE
1883032	UHMWPE
1883033	UHMWPE
1883034	UHMWPE
1883035	UHMWPE
1883036	UHMWPE
1883041	UHMWPE
1883042	UHMWPE
1883043	UHMWPE
1883044	UHMWPE
1883045	UHMWPE
1883046	UHMWPE
1883100	Ti6Al4V Alloy
1884002	Ti6Al4V Alloy
1884004	Ti6Al4V Alloy
1884006	Ti6Al4V Alloy
1885000	UHMWPE
1885010	Ti6Al4V Alloy
1885012	Ti6Al4V Alloy
1885014	Ti6Al4V Alloy
1885017	Ti6Al4V Alloy
1885020	Ti6Al4V Alloy
1885024	Ti6Al4V Alloy
1885050	Ti6Al4V Alloy
1885060	Ti6Al4V Alloy
1890001	Ti6Al4V Alloy
1890002	Ti6Al4V Alloy
1890003	Ti6Al4V Alloy
1890004	Ti6Al4V Alloy
1890005	Ti6Al4V Alloy
1890006	Ti6Al4V Alloy
1890007	Ti6Al4V Alloy
1890012	Ti6Al4V Alloy
1890023	Ti6Al4V Alloy

1890034	Ti6Al4V Alloy
1890045	Ti6Al4V Alloy
1890056	Ti6Al4V Alloy
1890067	Ti6Al4V Alloy
1891001	Ti6Al4V Alloy + Hydroxyapatite
1891002	Ti6Al4V Alloy + Hydroxyapatite
1891003	Ti6Al4V Alloy + Hydroxyapatite
1891004	Ti6Al4V Alloy + Hydroxyapatite
1891005	Ti6Al4V Alloy + Hydroxyapatite
1891006	Ti6Al4V Alloy + Hydroxyapatite
1891007	Ti6Al4V Alloy + Hydroxyapatite
1891012	Ti6Al4V Alloy + Hydroxyapatite
1891023	Ti6Al4V Alloy + Hydroxyapatite
1891034	Ti6Al4V Alloy + Hydroxyapatite
1891045	Ti6Al4V Alloy + Hydroxyapatite
1891056	Ti6Al4V Alloy + Hydroxyapatite
1891067	Ti6Al4V Alloy + Hydroxyapatite
1892001	Ti6Al4V Alloy
1892002	Ti6Al4V Alloy
1892003	Ti6Al4V Alloy
1892004	Ti6Al4V Alloy
1892005	Ti6Al4V Alloy
1892006	Ti6Al4V Alloy
1892007	Ti6Al4V Alloy
1892012	Ti6Al4V Alloy
1892023	Ti6Al4V Alloy
1892034	Ti6Al4V Alloy
1892045	Ti6Al4V Alloy
1892056	Ti6Al4V Alloy
1892067	Ti6Al4V Alloy
1893001	Ti6Al4V Alloy + Hydroxyapatite
1893002	Ti6Al4V Alloy + Hydroxyapatite
1893003	Ti6Al4V Alloy + Hydroxyapatite
1893004	Ti6Al4V Alloy + Hydroxyapatite
1893005	Ti6Al4V Alloy + Hydroxyapatite
1893006	Ti6Al4V Alloy + Hydroxyapatite
1893007	Ti6Al4V Alloy + Hydroxyapatite
1893012	Ti6Al4V Alloy + Hydroxyapatite
1893023	Ti6Al4V Alloy + Hydroxyapatite
1893034	Ti6Al4V Alloy + Hydroxyapatite
1893045	Ti6Al4V Alloy + Hydroxyapatite
1893056	Ti6Al4V Alloy + Hydroxyapatite
1893067	Ti6Al4V Alloy + Hydroxyapatite

1898001	Ti6Al4V Alloy
1898002	Ti6Al4V Alloy
1898003	Ti6Al4V Alloy
1898004	Ti6Al4V Alloy
1898005	Ti6Al4V Alloy
1898006	Ti6Al4V Alloy
1898007	Ti6Al4V Alloy
1898012	Ti6Al4V Alloy
1898023	Ti6Al4V Alloy
1898034	Ti6Al4V Alloy
1898045	Ti6Al4V Alloy
1898056	Ti6Al4V Alloy
1898067	Ti6Al4V Alloy
1901001	UHMWPE
1901002	UHMWPE
1901003	UHMWPE
1901004	UHMWPE
1901005	UHMWPE
1901006	UHMWPE
1901007	UHMWPE
1901031	UHMWPE
1901032	UHMWPE
1901033	UHMWPE
1901034	UHMWPE
1901035	UHMWPE
1901036	UHMWPE
1901037	UHMWPE
1901041	UHMWPE
1901042	UHMWPE
1901043	UHMWPE
1901044	UHMWPE
1901045	UHMWPE
1901046	UHMWPE
1901047	UHMWPE
1901053	UHMWPE
1901054	UHMWPE
1901055	UHMWPE
1901056	UHMWPE
1901057	UHMWPE
1901058	UHMWPE
1901101	UHMWPE
1901102	UHMWPE
1901103	UHMWPE

1901104	UHMWPE
1901105	UHMWPE
1901106	UHMWPE
1901107	UHMWPE
1901131	UHMWPE
1901132	UHMWPE
1901133	UHMWPE
1901134	UHMWPE
1901135	UHMWPE
1901136	UHMWPE
1901137	UHMWPE
1901141	UHMWPE
1901142	UHMWPE
1901143	UHMWPE
1901144	UHMWPE
1901145	UHMWPE
1901146	UHMWPE
1901147	UHMWPE
1901153	UHMWPE
1901154	UHMWPE
1901155	UHMWPE
1901156	UHMWPE
1901157	UHMWPE
1901158	UHMWPE
1901201	UHMWPE
1901202	UHMWPE
1901203	UHMWPE
1901204	UHMWPE
1901205	UHMWPE
1901206	UHMWPE
1901207	UHMWPE
1901231	UHMWPE
1901232	UHMWPE
1901233	UHMWPE
1901234	UHMWPE
1901235	UHMWPE
1901236	UHMWPE
1901237	UHMWPE
1901241	UHMWPE
1901242	UHMWPE
1901243	UHMWPE
1901244	UHMWPE
1901245	UHMWPE

1901246	UHMWPE
1901247	UHMWPE
1901253	UHMWPE
1901254	UHMWPE
1901255	UHMWPE
1901256	UHMWPE
1901257	UHMWPE
1901258	UHMWPE
1901301	UHMWPE
1901302	UHMWPE
1901303	UHMWPE
1901304	UHMWPE
1901305	UHMWPE
1901306	UHMWPE
1901307	UHMWPE
1901331	UHMWPE
1901332	UHMWPE
1901333	UHMWPE
1901334	UHMWPE
1901335	UHMWPE
1901336	UHMWPE
1901337	UHMWPE
1901341	UHMWPE
1901342	UHMWPE
1901343	UHMWPE
1901344	UHMWPE
1901345	UHMWPE
1901346	UHMWPE
1901347	UHMWPE
1901353	UHMWPE
1901354	UHMWPE
1901355	UHMWPE
1901356	UHMWPE
1901357	UHMWPE
1901358	UHMWPE
1901401	UHMWPE
1901402	UHMWPE
1901403	UHMWPE
1901404	UHMWPE
1901405	UHMWPE
1901406	UHMWPE
1901407	UHMWPE
1901431	UHMWPE

1901432	UHMWPE
1901433	UHMWPE
1901434	UHMWPE
1901435	UHMWPE
1901436	UHMWPE
1901437	UHMWPE
1901441	UHMWPE
1901442	UHMWPE
1901443	UHMWPE
1901444	UHMWPE
1901445	UHMWPE
1901446	UHMWPE
1901447	UHMWPE
1901453	UHMWPE
1901454	UHMWPE
1901455	UHMWPE
1901456	UHMWPE
1901457	UHMWPE
1901458	UHMWPE
1901501	UHMWPE
1901502	UHMWPE
1901503	UHMWPE
1901504	UHMWPE
1901505	UHMWPE
1901506	UHMWPE
1901507	UHMWPE
1901531	UHMWPE
1901532	UHMWPE
1901533	UHMWPE
1901534	UHMWPE
1901535	UHMWPE
1901536	UHMWPE
1901537	UHMWPE
1901541	UHMWPE
1901542	UHMWPE
1901543	UHMWPE
1901544	UHMWPE
1901545	UHMWPE
1901546	UHMWPE
1901547	UHMWPE
1901553	UHMWPE
1901554	UHMWPE
1901555	UHMWPE

1901556	UHMWPE
1901557	UHMWPE
1901558	UHMWPE
1901601	UHMWPE
1901602	UHMWPE
1901603	UHMWPE
1901604	UHMWPE
1901605	UHMWPE
1901606	UHMWPE
1901607	UHMWPE
1901631	UHMWPE
1901632	UHMWPE
1901633	UHMWPE
1901634	UHMWPE
1901635	UHMWPE
1901636	UHMWPE
1901637	UHMWPE
1901641	UHMWPE
1901642	UHMWPE
1901643	UHMWPE
1901644	UHMWPE
1901645	UHMWPE
1901646	UHMWPE
1901647	UHMWPE
1901653	UHMWPE
1901654	UHMWPE
1901655	UHMWPE
1901656	UHMWPE
1901657	UHMWPE
1901658	UHMWPE
1901701	UHMWPE
1901702	UHMWPE
1901703	UHMWPE
1901704	UHMWPE
1901705	UHMWPE
1901706	UHMWPE
1901707	UHMWPE
1901731	UHMWPE
1901732	UHMWPE
1901733	UHMWPE
1901734	UHMWPE
1901735	UHMWPE
1901736	UHMWPE

1901737	UHMWPE
1901741	UHMWPE
1901742	UHMWPE
1901743	UHMWPE
1901744	UHMWPE
1901745	UHMWPE
1901746	UHMWPE
1901747	UHMWPE
1901753	UHMWPE
1901754	UHMWPE
1901755	UHMWPE
1901756	UHMWPE
1901757	UHMWPE
1901758	UHMWPE
1902001	UHMWPE
1902002	UHMWPE
1902003	UHMWPE
1902004	UHMWPE
1902005	UHMWPE
1902006	UHMWPE
1902007	UHMWPE
1902031	UHMWPE
1902032	UHMWPE
1902033	UHMWPE
1902034	UHMWPE
1902035	UHMWPE
1902036	UHMWPE
1902037	UHMWPE
1902041	UHMWPE
1902042	UHMWPE
1902043	UHMWPE
1902044	UHMWPE
1902045	UHMWPE
1902046	UHMWPE
1902047	UHMWPE
1902053	UHMWPE
1902054	UHMWPE
1902055	UHMWPE
1902056	UHMWPE
1902057	UHMWPE
1902058	UHMWPE
1902401	UHMWPE
1902402	UHMWPE

1902403	UHMWPE
1902404	UHMWPE
1902405	UHMWPE
1902406	UHMWPE
1902407	UHMWPE
1902431	UHMWPE
1902432	UHMWPE
1902433	UHMWPE
1902434	UHMWPE
1902435	UHMWPE
1902436	UHMWPE
1902437	UHMWPE
1902441	UHMWPE
1902442	UHMWPE
1902443	UHMWPE
1902444	UHMWPE
1902445	UHMWPE
1902446	UHMWPE
1902447	UHMWPE
1902453	UHMWPE
1902454	UHMWPE
1902455	UHMWPE
1902456	UHMWPE
1902457	UHMWPE
1902458	UHMWPE
1911001	UHMWPE
1911002	UHMWPE
1911003	UHMWPE
1911004	UHMWPE
1911005	UHMWPE
1911006	UHMWPE
1911007	UHMWPE
1911031	UHMWPE
1911032	UHMWPE
1911033	UHMWPE
1911034	UHMWPE
1911035	UHMWPE
1911036	UHMWPE
1911037	UHMWPE
1911041	UHMWPE
1911042	UHMWPE
1911043	UHMWPE
1911044	UHMWPE

1911045	UHMWPE
1911046	UHMWPE
1911047	UHMWPE
1911053	UHMWPE
1911054	UHMWPE
1911055	UHMWPE
1911056	UHMWPE
1911057	UHMWPE
1911058	UHMWPE
1911101	UHMWPE
1911102	UHMWPE
1911103	UHMWPE
1911104	UHMWPE
1911105	UHMWPE
1911106	UHMWPE
1911107	UHMWPE
1911131	UHMWPE
1911132	UHMWPE
1911133	UHMWPE
1911134	UHMWPE
1911135	UHMWPE
1911136	UHMWPE
1911137	UHMWPE
1911141	UHMWPE
1911142	UHMWPE
1911143	UHMWPE
1911144	UHMWPE
1911145	UHMWPE
1911146	UHMWPE
1911147	UHMWPE
1911153	UHMWPE
1911154	UHMWPE
1911155	UHMWPE
1911156	UHMWPE
1911157	UHMWPE
1911158	UHMWPE
1911201	UHMWPE
1911202	UHMWPE
1911203	UHMWPE
1911204	UHMWPE
1911205	UHMWPE
1911206	UHMWPE
1911207	UHMWPE

1911231	UHMWPE
1911232	UHMWPE
1911233	UHMWPE
1911234	UHMWPE
1911235	UHMWPE
1911236	UHMWPE
1911237	UHMWPE
1911241	UHMWPE
1911242	UHMWPE
1911243	UHMWPE
1911244	UHMWPE
1911245	UHMWPE
1911246	UHMWPE
1911247	UHMWPE
1911253	UHMWPE
1911254	UHMWPE
1911255	UHMWPE
1911256	UHMWPE
1911257	UHMWPE
1911258	UHMWPE
1911301	UHMWPE
1911302	UHMWPE
1911303	UHMWPE
1911304	UHMWPE
1911305	UHMWPE
1911306	UHMWPE
1911307	UHMWPE
1911331	UHMWPE
1911332	UHMWPE
1911333	UHMWPE
1911334	UHMWPE
1911335	UHMWPE
1911336	UHMWPE
1911337	UHMWPE
1911341	UHMWPE
1911342	UHMWPE
1911343	UHMWPE
1911344	UHMWPE
1911345	UHMWPE
1911346	UHMWPE
1911347	UHMWPE
1911353	UHMWPE
1911354	UHMWPE

1911355	UHMWPE
1911356	UHMWPE
1911357	UHMWPE
1911358	UHMWPE
1911401	UHMWPE
1911402	UHMWPE
1911403	UHMWPE
1911404	UHMWPE
1911405	UHMWPE
1911406	UHMWPE
1911407	UHMWPE
1911431	UHMWPE
1911432	UHMWPE
1911433	UHMWPE
1911434	UHMWPE
1911435	UHMWPE
1911436	UHMWPE
1911437	UHMWPE
1911441	UHMWPE
1911442	UHMWPE
1911443	UHMWPE
1911444	UHMWPE
1911445	UHMWPE
1911446	UHMWPE
1911447	UHMWPE
1911453	UHMWPE
1911454	UHMWPE
1911455	UHMWPE
1911456	UHMWPE
1911457	UHMWPE
1911458	UHMWPE
1911501	UHMWPE
1911502	UHMWPE
1911503	UHMWPE
1911504	UHMWPE
1911505	UHMWPE
1911506	UHMWPE
1911507	UHMWPE
1911531	UHMWPE
1911532	UHMWPE
1911533	UHMWPE
1911534	UHMWPE
1911535	UHMWPE

1911536	UHMWPE
1911537	UHMWPE
1911541	UHMWPE
1911542	UHMWPE
1911543	UHMWPE
1911544	UHMWPE
1911545	UHMWPE
1911546	UHMWPE
1911547	UHMWPE
1911553	UHMWPE
1911554	UHMWPE
1911555	UHMWPE
1911556	UHMWPE
1911557	UHMWPE
1911558	UHMWPE
1911601	UHMWPE
1911602	UHMWPE
1911603	UHMWPE
1911604	UHMWPE
1911605	UHMWPE
1911606	UHMWPE
1911607	UHMWPE
1911631	UHMWPE
1911632	UHMWPE
1911633	UHMWPE
1911634	UHMWPE
1911635	UHMWPE
1911636	UHMWPE
1911637	UHMWPE
1911641	UHMWPE
1911642	UHMWPE
1911643	UHMWPE
1911644	UHMWPE
1911645	UHMWPE
1911646	UHMWPE
1911647	UHMWPE
1911653	UHMWPE
1911654	UHMWPE
1911655	UHMWPE
1911656	UHMWPE
1911657	UHMWPE
1911658	UHMWPE
1911701	UHMWPE

1911702	UHMWPE
1911703	UHMWPE
1911704	UHMWPE
1911705	UHMWPE
1911706	UHMWPE
1911707	UHMWPE
1911731	UHMWPE
1911732	UHMWPE
1911733	UHMWPE
1911734	UHMWPE
1911735	UHMWPE
1911736	UHMWPE
1911737	UHMWPE
1911741	UHMWPE
1911742	UHMWPE
1911743	UHMWPE
1911744	UHMWPE
1911745	UHMWPE
1911746	UHMWPE
1911747	UHMWPE
1911753	UHMWPE
1911754	UHMWPE
1911755	UHMWPE
1911756	UHMWPE
1911757	UHMWPE
1911758	UHMWPE
1912001	UHMWPE
1912002	UHMWPE
1912003	UHMWPE
1912004	UHMWPE
1912005	UHMWPE
1912006	UHMWPE
1912007	UHMWPE
1912031	UHMWPE
1912032	UHMWPE
1912033	UHMWPE
1912034	UHMWPE
1912035	UHMWPE
1912036	UHMWPE
1912037	UHMWPE
1912041	UHMWPE
1912042	UHMWPE
1912043	UHMWPE

1912044	UHMWPE
1912045	UHMWPE
1912046	UHMWPE
1912047	UHMWPE
1912053	UHMWPE
1912054	UHMWPE
1912055	UHMWPE
1912056	UHMWPE
1912057	UHMWPE
1912058	UHMWPE
1912401	UHMWPE
1912402	UHMWPE
1912403	UHMWPE
1912404	UHMWPE
1912405	UHMWPE
1912406	UHMWPE
1912407	UHMWPE
1912431	UHMWPE
1912432	UHMWPE
1912433	UHMWPE
1912434	UHMWPE
1912435	UHMWPE
1912436	UHMWPE
1912437	UHMWPE
1912441	UHMWPE
1912442	UHMWPE
1912443	UHMWPE
1912444	UHMWPE
1912445	UHMWPE
1912446	UHMWPE
1912447	UHMWPE
1912453	UHMWPE
1912454	UHMWPE
1912455	UHMWPE
1912456	UHMWPE
1912457	UHMWPE
1912458	UHMWPE
1921001	UHMWPE
1921002	UHMWPE
1921003	UHMWPE
1921004	UHMWPE
1921005	UHMWPE
1921006	UHMWPE

1921007	UHMWPE
1921031	UHMWPE
1921032	UHMWPE
1921033	UHMWPE
1921034	UHMWPE
1921035	UHMWPE
1921036	UHMWPE
1921037	UHMWPE
1921041	UHMWPE
1921042	UHMWPE
1921043	UHMWPE
1921044	UHMWPE
1921045	UHMWPE
1921046	UHMWPE
1921047	UHMWPE
1921053	UHMWPE
1921054	UHMWPE
1921055	UHMWPE
1921056	UHMWPE
1921057	UHMWPE
1921058	UHMWPE
1921101	UHMWPE
1921102	UHMWPE
1921103	UHMWPE
1921104	UHMWPE
1921105	UHMWPE
1921106	UHMWPE
1921107	UHMWPE
1921131	UHMWPE
1921132	UHMWPE
1921133	UHMWPE
1921134	UHMWPE
1921135	UHMWPE
1921136	UHMWPE
1921137	UHMWPE
1921141	UHMWPE
1921142	UHMWPE
1921143	UHMWPE
1921144	UHMWPE
1921145	UHMWPE
1921146	UHMWPE
1921147	UHMWPE
1921153	UHMWPE

1921154	UHMWPE
1921155	UHMWPE
1921156	UHMWPE
1921157	UHMWPE
1921158	UHMWPE
1921201	UHMWPE
1921202	UHMWPE
1921203	UHMWPE
1921204	UHMWPE
1921205	UHMWPE
1921206	UHMWPE
1921207	UHMWPE
1921231	UHMWPE
1921232	UHMWPE
1921233	UHMWPE
1921234	UHMWPE
1921235	UHMWPE
1921236	UHMWPE
1921237	UHMWPE
1921241	UHMWPE
1921242	UHMWPE
1921243	UHMWPE
1921244	UHMWPE
1921245	UHMWPE
1921246	UHMWPE
1921247	UHMWPE
1921253	UHMWPE
1921254	UHMWPE
1921255	UHMWPE
1921256	UHMWPE
1921257	UHMWPE
1921258	UHMWPE
1921301	UHMWPE
1921302	UHMWPE
1921303	UHMWPE
1921304	UHMWPE
1921305	UHMWPE
1921306	UHMWPE
1921307	UHMWPE
1921331	UHMWPE
1921332	UHMWPE
1921333	UHMWPE
1921334	UHMWPE

1921335	UHMWPE
1921336	UHMWPE
1921337	UHMWPE
1921341	UHMWPE
1921342	UHMWPE
1921343	UHMWPE
1921344	UHMWPE
1921345	UHMWPE
1921346	UHMWPE
1921347	UHMWPE
1921353	UHMWPE
1921354	UHMWPE
1921355	UHMWPE
1921356	UHMWPE
1921357	UHMWPE
1921358	UHMWPE
1921401	UHMWPE
1921402	UHMWPE
1921403	UHMWPE
1921404	UHMWPE
1921405	UHMWPE
1921406	UHMWPE
1921407	UHMWPE
1921431	UHMWPE
1921432	UHMWPE
1921433	UHMWPE
1921434	UHMWPE
1921435	UHMWPE
1921436	UHMWPE
1921437	UHMWPE
1921441	UHMWPE
1921442	UHMWPE
1921443	UHMWPE
1921444	UHMWPE
1921445	UHMWPE
1921446	UHMWPE
1921447	UHMWPE
1921453	UHMWPE
1921454	UHMWPE
1921455	UHMWPE
1921456	UHMWPE
1921457	UHMWPE
1921458	UHMWPE

1921501	UHMWPE
1921502	UHMWPE
1921503	UHMWPE
1921504	UHMWPE
1921505	UHMWPE
1921506	UHMWPE
1921507	UHMWPE
1921531	UHMWPE
1921532	UHMWPE
1921533	UHMWPE
1921534	UHMWPE
1921535	UHMWPE
1921536	UHMWPE
1921537	UHMWPE
1921541	UHMWPE
1921542	UHMWPE
1921543	UHMWPE
1921544	UHMWPE
1921545	UHMWPE
1921546	UHMWPE
1921547	UHMWPE
1921553	UHMWPE
1921554	UHMWPE
1921555	UHMWPE
1921556	UHMWPE
1921557	UHMWPE
1921558	UHMWPE
1921601	UHMWPE
1921602	UHMWPE
1921603	UHMWPE
1921604	UHMWPE
1921605	UHMWPE
1921606	UHMWPE
1921607	UHMWPE
1921631	UHMWPE
1921632	UHMWPE
1921633	UHMWPE
1921634	UHMWPE
1921635	UHMWPE
1921636	UHMWPE
1921637	UHMWPE
1921641	UHMWPE
1921642	UHMWPE

1921643	UHMWPE
1921644	UHMWPE
1921645	UHMWPE
1921646	UHMWPE
1921647	UHMWPE
1921653	UHMWPE
1921654	UHMWPE
1921655	UHMWPE
1921656	UHMWPE
1921657	UHMWPE
1921658	UHMWPE
1921701	UHMWPE
1921702	UHMWPE
1921703	UHMWPE
1921704	UHMWPE
1921705	UHMWPE
1921706	UHMWPE
1921707	UHMWPE
1921731	UHMWPE
1921732	UHMWPE
1921733	UHMWPE
1921734	UHMWPE
1921735	UHMWPE
1921736	UHMWPE
1921737	UHMWPE
1921741	UHMWPE
1921742	UHMWPE
1921743	UHMWPE
1921744	UHMWPE
1921745	UHMWPE
1921746	UHMWPE
1921747	UHMWPE
1921753	UHMWPE
1921754	UHMWPE
1921755	UHMWPE
1921756	UHMWPE
1921757	UHMWPE
1921758	UHMWPE
1922001	UHMWPE
1922002	UHMWPE
1922003	UHMWPE
1922004	UHMWPE
1922005	UHMWPE

1922006	UHMWPE
1922007	UHMWPE
1922031	UHMWPE
1922032	UHMWPE
1922033	UHMWPE
1922034	UHMWPE
1922035	UHMWPE
1922036	UHMWPE
1922037	UHMWPE
1922041	UHMWPE
1922042	UHMWPE
1922043	UHMWPE
1922044	UHMWPE
1922045	UHMWPE
1922046	UHMWPE
1922047	UHMWPE
1922053	UHMWPE
1922054	UHMWPE
1922055	UHMWPE
1922056	UHMWPE
1922057	UHMWPE
1922058	UHMWPE
1922401	UHMWPE
1922402	UHMWPE
1922403	UHMWPE
1922404	UHMWPE
1922405	UHMWPE
1922406	UHMWPE
1922407	UHMWPE
1922431	UHMWPE
1922432	UHMWPE
1922433	UHMWPE
1922434	UHMWPE
1922435	UHMWPE
1922436	UHMWPE
1922437	UHMWPE
1922441	UHMWPE
1922442	UHMWPE
1922443	UHMWPE
1922444	UHMWPE
1922445	UHMWPE
1922446	UHMWPE
1922447	UHMWPE

1922453	UHMWPE
1922454	UHMWPE
1922455	UHMWPE
1922456	UHMWPE
1922457	UHMWPE
1922458	UHMWPE
1931001	UHMWPE
1931002	UHMWPE
1931003	UHMWPE
1931004	UHMWPE
1931005	UHMWPE
1931006	UHMWPE
1931007	UHMWPE
1931011	UHMWPE
1931012	UHMWPE
1931013	UHMWPE
1931014	UHMWPE
1931015	UHMWPE
1931016	UHMWPE
1931017	UHMWPE
1931031	UHMWPE
1931032	UHMWPE
1931033	UHMWPE
1931034	UHMWPE
1931035	UHMWPE
1931036	UHMWPE
1931037	UHMWPE
1931041	UHMWPE
1931042	UHMWPE
1931043	UHMWPE
1931044	UHMWPE
1931045	UHMWPE
1931046	UHMWPE
1931047	UHMWPE
1931053	UHMWPE
1931054	UHMWPE
1931055	UHMWPE
1931056	UHMWPE
1931057	UHMWPE
1931058	UHMWPE
1931073	UHMWPE
1931074	UHMWPE
1931101	UHMWPE

1931102	UHMWPE
1931103	UHMWPE
1931104	UHMWPE
1931105	UHMWPE
1931106	UHMWPE
1931107	UHMWPE
1931131	UHMWPE
1931132	UHMWPE
1931133	UHMWPE
1931134	UHMWPE
1931135	UHMWPE
1931136	UHMWPE
1931137	UHMWPE
1931141	UHMWPE
1931142	UHMWPE
1931143	UHMWPE
1931144	UHMWPE
1931145	UHMWPE
1931146	UHMWPE
1931147	UHMWPE
1931153	UHMWPE
1931154	UHMWPE
1931155	UHMWPE
1931156	UHMWPE
1931157	UHMWPE
1931158	UHMWPE
1931201	UHMWPE
1931202	UHMWPE
1931203	UHMWPE
1931204	UHMWPE
1931205	UHMWPE
1931206	UHMWPE
1931207	UHMWPE
1931211	UHMWPE
1931212	UHMWPE
1931213	UHMWPE
1931214	UHMWPE
1931215	UHMWPE
1931216	UHMWPE
1931217	UHMWPE
1931231	UHMWPE
1931232	UHMWPE
1931233	UHMWPE

1931234	UHMWPE
1931235	UHMWPE
1931236	UHMWPE
1931237	UHMWPE
1931241	UHMWPE
1931242	UHMWPE
1931243	UHMWPE
1931244	UHMWPE
1931245	UHMWPE
1931246	UHMWPE
1931247	UHMWPE
1931253	UHMWPE
1931254	UHMWPE
1931255	UHMWPE
1931256	UHMWPE
1931257	UHMWPE
1931258	UHMWPE
1931273	UHMWPE
1931274	UHMWPE
1931301	UHMWPE
1931302	UHMWPE
1931303	UHMWPE
1931304	UHMWPE
1931305	UHMWPE
1931306	UHMWPE
1931307	UHMWPE
1931331	UHMWPE
1931332	UHMWPE
1931333	UHMWPE
1931334	UHMWPE
1931335	UHMWPE
1931336	UHMWPE
1931337	UHMWPE
1931341	UHMWPE
1931342	UHMWPE
1931343	UHMWPE
1931344	UHMWPE
1931345	UHMWPE
1931346	UHMWPE
1931347	UHMWPE
1931353	UHMWPE
1931354	UHMWPE
1931355	UHMWPE

1931356	UHMWPE
1931357	UHMWPE
1931358	UHMWPE
1931401	UHMWPE
1931402	UHMWPE
1931403	UHMWPE
1931404	UHMWPE
1931405	UHMWPE
1931406	UHMWPE
1931407	UHMWPE
1931411	UHMWPE
1931412	UHMWPE
1931413	UHMWPE
1931414	UHMWPE
1931415	UHMWPE
1931416	UHMWPE
1931417	UHMWPE
1931431	UHMWPE
1931432	UHMWPE
1931433	UHMWPE
1931434	UHMWPE
1931435	UHMWPE
1931436	UHMWPE
1931437	UHMWPE
1931441	UHMWPE
1931442	UHMWPE
1931443	UHMWPE
1931444	UHMWPE
1931445	UHMWPE
1931446	UHMWPE
1931447	UHMWPE
1931453	UHMWPE
1931454	UHMWPE
1931455	UHMWPE
1931456	UHMWPE
1931457	UHMWPE
1931458	UHMWPE
1931473	UHMWPE
1931474	UHMWPE
1931501	UHMWPE
1931502	UHMWPE
1931503	UHMWPE
1931504	UHMWPE

1931505	UHMWPE
1931506	UHMWPE
1931507	UHMWPE
1931531	UHMWPE
1931532	UHMWPE
1931533	UHMWPE
1931534	UHMWPE
1931535	UHMWPE
1931536	UHMWPE
1931537	UHMWPE
1931541	UHMWPE
1931542	UHMWPE
1931543	UHMWPE
1931544	UHMWPE
1931545	UHMWPE
1931546	UHMWPE
1931547	UHMWPE
1931553	UHMWPE
1931554	UHMWPE
1931555	UHMWPE
1931556	UHMWPE
1931557	UHMWPE
1931558	UHMWPE
1931601	UHMWPE
1931602	UHMWPE
1931603	UHMWPE
1931604	UHMWPE
1931605	UHMWPE
1931606	UHMWPE
1931607	UHMWPE
1931631	UHMWPE
1931632	UHMWPE
1931633	UHMWPE
1931634	UHMWPE
1931635	UHMWPE
1931636	UHMWPE
1931637	UHMWPE
1931641	UHMWPE
1931642	UHMWPE
1931643	UHMWPE
1931644	UHMWPE
1931645	UHMWPE
1931646	UHMWPE

1931647	UHMWPE
1931653	UHMWPE
1931654	UHMWPE
1931655	UHMWPE
1931656	UHMWPE
1931657	UHMWPE
1931658	UHMWPE
1931701	UHMWPE
1931702	UHMWPE
1931703	UHMWPE
1931704	UHMWPE
1931705	UHMWPE
1931706	UHMWPE
1931707	UHMWPE
1931711	UHMWPE
1931712	UHMWPE
1931713	UHMWPE
1931714	UHMWPE
1931715	UHMWPE
1931716	UHMWPE
1931717	UHMWPE
1931731	UHMWPE
1931732	UHMWPE
1931733	UHMWPE
1931734	UHMWPE
1931735	UHMWPE
1931736	UHMWPE
1931737	UHMWPE
1931741	UHMWPE
1931742	UHMWPE
1931743	UHMWPE
1931744	UHMWPE
1931745	UHMWPE
1931746	UHMWPE
1931747	UHMWPE
1931753	UHMWPE
1931754	UHMWPE
1931755	UHMWPE
1931756	UHMWPE
1931757	UHMWPE
1931758	UHMWPE
1931773	UHMWPE
1931774	UHMWPE

1932001	UHMWPE
1932002	UHMWPE
1932003	UHMWPE
1932004	UHMWPE
1932005	UHMWPE
1932006	UHMWPE
1932007	UHMWPE
1932011	UHMWPE
1932012	UHMWPE
1932013	UHMWPE
1932014	UHMWPE
1932015	UHMWPE
1932016	UHMWPE
1932017	UHMWPE
1932031	UHMWPE
1932032	UHMWPE
1932033	UHMWPE
1932034	UHMWPE
1932035	UHMWPE
1932036	UHMWPE
1932037	UHMWPE
1932041	UHMWPE
1932042	UHMWPE
1932043	UHMWPE
1932044	UHMWPE
1932045	UHMWPE
1932046	UHMWPE
1932047	UHMWPE
1932053	UHMWPE
1932054	UHMWPE
1932055	UHMWPE
1932056	UHMWPE
1932057	UHMWPE
1932058	UHMWPE
1932073	UHMWPE
1932074	UHMWPE
1932401	UHMWPE
1932402	UHMWPE
1932403	UHMWPE
1932404	UHMWPE
1932405	UHMWPE
1932406	UHMWPE
1932407	UHMWPE

1932411	UHMWPE
1932412	UHMWPE
1932413	UHMWPE
1932414	UHMWPE
1932415	UHMWPE
1932416	UHMWPE
1932417	UHMWPE
1932431	UHMWPE
1932432	UHMWPE
1932433	UHMWPE
1932434	UHMWPE
1932435	UHMWPE
1932436	UHMWPE
1932437	UHMWPE
1932441	UHMWPE
1932442	UHMWPE
1932443	UHMWPE
1932444	UHMWPE
1932445	UHMWPE
1932446	UHMWPE
1932447	UHMWPE
1932453	UHMWPE
1932454	UHMWPE
1932455	UHMWPE
1932456	UHMWPE
1932457	UHMWPE
1932458	UHMWPE
1932473	UHMWPE
1932474	UHMWPE
1941001	UHMWPE
1941002	UHMWPE
1941003	UHMWPE
1941004	UHMWPE
1941005	UHMWPE
1941006	UHMWPE
1941007	UHMWPE
1941031	UHMWPE
1941032	UHMWPE
1941033	UHMWPE
1941034	UHMWPE
1941035	UHMWPE
1941036	UHMWPE
1941037	UHMWPE

1941041	UHMWPE
1941042	UHMWPE
1941043	UHMWPE
1941044	UHMWPE
1941045	UHMWPE
1941046	UHMWPE
1941047	UHMWPE
1941053	UHMWPE
1941054	UHMWPE
1941055	UHMWPE
1941056	UHMWPE
1941057	UHMWPE
1941058	UHMWPE
1941101	UHMWPE
1941102	UHMWPE
1941103	UHMWPE
1941104	UHMWPE
1941105	UHMWPE
1941106	UHMWPE
1941107	UHMWPE
1941131	UHMWPE
1941132	UHMWPE
1941133	UHMWPE
1941134	UHMWPE
1941135	UHMWPE
1941136	UHMWPE
1941137	UHMWPE
1941141	UHMWPE
1941142	UHMWPE
1941143	UHMWPE
1941144	UHMWPE
1941145	UHMWPE
1941146	UHMWPE
1941147	UHMWPE
1941153	UHMWPE
1941154	UHMWPE
1941155	UHMWPE
1941156	UHMWPE
1941157	UHMWPE
1941158	UHMWPE
1941201	UHMWPE
1941202	UHMWPE
1941203	UHMWPE

1941204	UHMWPE
1941205	UHMWPE
1941206	UHMWPE
1941207	UHMWPE
1941231	UHMWPE
1941232	UHMWPE
1941233	UHMWPE
1941234	UHMWPE
1941235	UHMWPE
1941236	UHMWPE
1941237	UHMWPE
1941241	UHMWPE
1941242	UHMWPE
1941243	UHMWPE
1941244	UHMWPE
1941245	UHMWPE
1941246	UHMWPE
1941247	UHMWPE
1941253	UHMWPE
1941254	UHMWPE
1941255	UHMWPE
1941256	UHMWPE
1941257	UHMWPE
1941258	UHMWPE
1941301	UHMWPE
1941302	UHMWPE
1941303	UHMWPE
1941304	UHMWPE
1941305	UHMWPE
1941306	UHMWPE
1941307	UHMWPE
1941331	UHMWPE
1941332	UHMWPE
1941333	UHMWPE
1941334	UHMWPE
1941335	UHMWPE
1941336	UHMWPE
1941337	UHMWPE
1941341	UHMWPE
1941342	UHMWPE
1941343	UHMWPE
1941344	UHMWPE
1941345	UHMWPE

1941346	UHMWPE
1941347	UHMWPE
1941353	UHMWPE
1941354	UHMWPE
1941355	UHMWPE
1941356	UHMWPE
1941357	UHMWPE
1941358	UHMWPE
1941401	UHMWPE
1941402	UHMWPE
1941403	UHMWPE
1941404	UHMWPE
1941405	UHMWPE
1941406	UHMWPE
1941407	UHMWPE
1941431	UHMWPE
1941432	UHMWPE
1941433	UHMWPE
1941434	UHMWPE
1941435	UHMWPE
1941436	UHMWPE
1941437	UHMWPE
1941441	UHMWPE
1941442	UHMWPE
1941443	UHMWPE
1941444	UHMWPE
1941445	UHMWPE
1941446	UHMWPE
1941447	UHMWPE
1941453	UHMWPE
1941454	UHMWPE
1941455	UHMWPE
1941456	UHMWPE
1941457	UHMWPE
1941458	UHMWPE
1941501	UHMWPE
1941502	UHMWPE
1941503	UHMWPE
1941504	UHMWPE
1941505	UHMWPE
1941506	UHMWPE
1941507	UHMWPE
1941531	UHMWPE

1941532	UHMWPE
1941533	UHMWPE
1941534	UHMWPE
1941535	UHMWPE
1941536	UHMWPE
1941537	UHMWPE
1941541	UHMWPE
1941542	UHMWPE
1941543	UHMWPE
1941544	UHMWPE
1941545	UHMWPE
1941546	UHMWPE
1941547	UHMWPE
1941553	UHMWPE
1941554	UHMWPE
1941555	UHMWPE
1941556	UHMWPE
1941557	UHMWPE
1941558	UHMWPE
1941601	UHMWPE
1941602	UHMWPE
1941603	UHMWPE
1941604	UHMWPE
1941605	UHMWPE
1941606	UHMWPE
1941607	UHMWPE
1941631	UHMWPE
1941632	UHMWPE
1941633	UHMWPE
1941634	UHMWPE
1941635	UHMWPE
1941636	UHMWPE
1941637	UHMWPE
1941641	UHMWPE
1941642	UHMWPE
1941643	UHMWPE
1941644	UHMWPE
1941645	UHMWPE
1941646	UHMWPE
1941647	UHMWPE
1941653	UHMWPE
1941654	UHMWPE
1941655	UHMWPE

1941656	UHMWPE
1941657	UHMWPE
1941658	UHMWPE
1941701	UHMWPE
1941702	UHMWPE
1941703	UHMWPE
1941704	UHMWPE
1941705	UHMWPE
1941706	UHMWPE
1941707	UHMWPE
1941731	UHMWPE
1941732	UHMWPE
1941733	UHMWPE
1941734	UHMWPE
1941735	UHMWPE
1941736	UHMWPE
1941737	UHMWPE
1941741	UHMWPE
1941742	UHMWPE
1941743	UHMWPE
1941744	UHMWPE
1941745	UHMWPE
1941746	UHMWPE
1941747	UHMWPE
1941753	UHMWPE
1941754	UHMWPE
1941755	UHMWPE
1941756	UHMWPE
1941757	UHMWPE
1941758	UHMWPE
1942001	UHMWPE
1942002	UHMWPE
1942003	UHMWPE
1942004	UHMWPE
1942005	UHMWPE
1942006	UHMWPE
1942007	UHMWPE
1942031	UHMWPE
1942032	UHMWPE
1942033	UHMWPE
1942034	UHMWPE
1942035	UHMWPE
1942036	UHMWPE

1942037	UHMWPE
1942041	UHMWPE
1942042	UHMWPE
1942043	UHMWPE
1942044	UHMWPE
1942045	UHMWPE
1942046	UHMWPE
1942047	UHMWPE
1942053	UHMWPE
1942054	UHMWPE
1942055	UHMWPE
1942056	UHMWPE
1942057	UHMWPE
1942058	UHMWPE
1942401	UHMWPE
1942402	UHMWPE
1942403	UHMWPE
1942404	UHMWPE
1942405	UHMWPE
1942406	UHMWPE
1942407	UHMWPE
1942431	UHMWPE
1942432	UHMWPE
1942433	UHMWPE
1942434	UHMWPE
1942435	UHMWPE
1942436	UHMWPE
1942437	UHMWPE
1942441	UHMWPE
1942442	UHMWPE
1942443	UHMWPE
1942444	UHMWPE
1942445	UHMWPE
1942446	UHMWPE
1942447	UHMWPE
1942453	UHMWPE
1942454	UHMWPE
1942455	UHMWPE
1942456	UHMWPE
1942457	UHMWPE
1942458	UHMWPE
1951001	UHMWPE
1951002	UHMWPE

1951003	UHMWPE
1951004	UHMWPE
1951005	UHMWPE
1951006	UHMWPE
1951007	UHMWPE
1951011	UHMWPE
1951012	UHMWPE
1951013	UHMWPE
1951014	UHMWPE
1951015	UHMWPE
1951016	UHMWPE
1951017	UHMWPE
1951031	UHMWPE
1951032	UHMWPE
1951033	UHMWPE
1951034	UHMWPE
1951035	UHMWPE
1951036	UHMWPE
1951037	UHMWPE
1951041	UHMWPE
1951042	UHMWPE
1951043	UHMWPE
1951044	UHMWPE
1951045	UHMWPE
1951046	UHMWPE
1951047	UHMWPE
1951053	UHMWPE
1951054	UHMWPE
1951055	UHMWPE
1951056	UHMWPE
1951057	UHMWPE
1951058	UHMWPE
1951063	UHMWPE
1951064	UHMWPE
1951065	UHMWPE
1951066	UHMWPE
1951067	UHMWPE
1951068	UHMWPE
1951071	UHMWPE
1951072	UHMWPE
1951073	UHMWPE
1951074	UHMWPE
1951075	UHMWPE

1951076	UHMWPE
1951077	UHMWPE
1951081	UHMWPE
1951082	UHMWPE
1951083	UHMWPE
1951084	UHMWPE
1951085	UHMWPE
1951086	UHMWPE
1951087	UHMWPE
1951101	UHMWPE
1951102	UHMWPE
1951103	UHMWPE
1951104	UHMWPE
1951105	UHMWPE
1951106	UHMWPE
1951107	UHMWPE
1951111	UHMWPE
1951112	UHMWPE
1951113	UHMWPE
1951114	UHMWPE
1951115	UHMWPE
1951116	UHMWPE
1951117	UHMWPE
1951131	UHMWPE
1951132	UHMWPE
1951133	UHMWPE
1951134	UHMWPE
1951135	UHMWPE
1951136	UHMWPE
1951137	UHMWPE
1951141	UHMWPE
1951142	UHMWPE
1951143	UHMWPE
1951144	UHMWPE
1951145	UHMWPE
1951146	UHMWPE
1951147	UHMWPE
1951153	UHMWPE
1951154	UHMWPE
1951155	UHMWPE
1951156	UHMWPE
1951157	UHMWPE
1951158	UHMWPE

1951163	UHMWPE
1951164	UHMWPE
1951165	UHMWPE
1951166	UHMWPE
1951167	UHMWPE
1951168	UHMWPE
1951171	UHMWPE
1951172	UHMWPE
1951173	UHMWPE
1951174	UHMWPE
1951175	UHMWPE
1951176	UHMWPE
1951177	UHMWPE
1951181	UHMWPE
1951182	UHMWPE
1951183	UHMWPE
1951184	UHMWPE
1951185	UHMWPE
1951186	UHMWPE
1951187	UHMWPE
1951201	UHMWPE
1951202	UHMWPE
1951203	UHMWPE
1951204	UHMWPE
1951205	UHMWPE
1951206	UHMWPE
1951207	UHMWPE
1951211	UHMWPE
1951212	UHMWPE
1951213	UHMWPE
1951214	UHMWPE
1951215	UHMWPE
1951216	UHMWPE
1951217	UHMWPE
1951231	UHMWPE
1951232	UHMWPE
1951233	UHMWPE
1951234	UHMWPE
1951235	UHMWPE
1951236	UHMWPE
1951237	UHMWPE
1951241	UHMWPE
1951242	UHMWPE

1951243	UHMWPE
1951244	UHMWPE
1951245	UHMWPE
1951246	UHMWPE
1951247	UHMWPE
1951253	UHMWPE
1951254	UHMWPE
1951255	UHMWPE
1951256	UHMWPE
1951257	UHMWPE
1951258	UHMWPE
1951263	UHMWPE
1951264	UHMWPE
1951265	UHMWPE
1951266	UHMWPE
1951267	UHMWPE
1951268	UHMWPE
1951271	UHMWPE
1951272	UHMWPE
1951273	UHMWPE
1951274	UHMWPE
1951275	UHMWPE
1951276	UHMWPE
1951277	UHMWPE
1951281	UHMWPE
1951282	UHMWPE
1951283	UHMWPE
1951284	UHMWPE
1951285	UHMWPE
1951286	UHMWPE
1951287	UHMWPE
1951301	UHMWPE
1951302	UHMWPE
1951303	UHMWPE
1951304	UHMWPE
1951305	UHMWPE
1951306	UHMWPE
1951307	UHMWPE
1951311	UHMWPE
1951312	UHMWPE
1951313	UHMWPE
1951314	UHMWPE
1951315	UHMWPE

1951316	UHMWPE
1951317	UHMWPE
1951331	UHMWPE
1951332	UHMWPE
1951333	UHMWPE
1951334	UHMWPE
1951335	UHMWPE
1951336	UHMWPE
1951337	UHMWPE
1951341	UHMWPE
1951342	UHMWPE
1951343	UHMWPE
1951344	UHMWPE
1951345	UHMWPE
1951346	UHMWPE
1951347	UHMWPE
1951353	UHMWPE
1951354	UHMWPE
1951355	UHMWPE
1951356	UHMWPE
1951357	UHMWPE
1951358	UHMWPE
1951363	UHMWPE
1951364	UHMWPE
1951365	UHMWPE
1951366	UHMWPE
1951367	UHMWPE
1951368	UHMWPE
1951371	UHMWPE
1951372	UHMWPE
1951373	UHMWPE
1951374	UHMWPE
1951375	UHMWPE
1951376	UHMWPE
1951377	UHMWPE
1951381	UHMWPE
1951382	UHMWPE
1951383	UHMWPE
1951384	UHMWPE
1951385	UHMWPE
1951386	UHMWPE
1951387	UHMWPE
1951401	UHMWPE

1951402	UHMWPE
1951403	UHMWPE
1951404	UHMWPE
1951405	UHMWPE
1951406	UHMWPE
1951407	UHMWPE
1951411	UHMWPE
1951412	UHMWPE
1951413	UHMWPE
1951414	UHMWPE
1951415	UHMWPE
1951416	UHMWPE
1951417	UHMWPE
1951431	UHMWPE
1951432	UHMWPE
1951433	UHMWPE
1951434	UHMWPE
1951435	UHMWPE
1951436	UHMWPE
1951437	UHMWPE
1951441	UHMWPE
1951442	UHMWPE
1951443	UHMWPE
1951444	UHMWPE
1951445	UHMWPE
1951446	UHMWPE
1951447	UHMWPE
1951453	UHMWPE
1951454	UHMWPE
1951455	UHMWPE
1951456	UHMWPE
1951457	UHMWPE
1951458	UHMWPE
1951463	UHMWPE
1951464	UHMWPE
1951465	UHMWPE
1951466	UHMWPE
1951467	UHMWPE
1951468	UHMWPE
1951471	UHMWPE
1951472	UHMWPE
1951473	UHMWPE
1951474	UHMWPE

1951475	UHMWPE
1951476	UHMWPE
1951477	UHMWPE
1951481	UHMWPE
1951482	UHMWPE
1951483	UHMWPE
1951484	UHMWPE
1951485	UHMWPE
1951486	UHMWPE
1951487	UHMWPE
1951501	UHMWPE
1951502	UHMWPE
1951503	UHMWPE
1951504	UHMWPE
1951505	UHMWPE
1951506	UHMWPE
1951507	UHMWPE
1951511	UHMWPE
1951512	UHMWPE
1951513	UHMWPE
1951514	UHMWPE
1951515	UHMWPE
1951516	UHMWPE
1951517	UHMWPE
1951531	UHMWPE
1951532	UHMWPE
1951533	UHMWPE
1951534	UHMWPE
1951535	UHMWPE
1951536	UHMWPE
1951537	UHMWPE
1951541	UHMWPE
1951542	UHMWPE
1951543	UHMWPE
1951544	UHMWPE
1951545	UHMWPE
1951546	UHMWPE
1951547	UHMWPE
1951553	UHMWPE
1951554	UHMWPE
1951555	UHMWPE
1951556	UHMWPE
1951557	UHMWPE

1951558	UHMWPE
1951563	UHMWPE
1951564	UHMWPE
1951565	UHMWPE
1951566	UHMWPE
1951567	UHMWPE
1951568	UHMWPE
1951571	UHMWPE
1951572	UHMWPE
1951573	UHMWPE
1951574	UHMWPE
1951575	UHMWPE
1951576	UHMWPE
1951577	UHMWPE
1951581	UHMWPE
1951582	UHMWPE
1951583	UHMWPE
1951584	UHMWPE
1951585	UHMWPE
1951586	UHMWPE
1951587	UHMWPE
1951601	UHMWPE
1951602	UHMWPE
1951603	UHMWPE
1951604	UHMWPE
1951605	UHMWPE
1951606	UHMWPE
1951607	UHMWPE
1951611	UHMWPE
1951612	UHMWPE
1951613	UHMWPE
1951614	UHMWPE
1951615	UHMWPE
1951616	UHMWPE
1951617	UHMWPE
1951631	UHMWPE
1951632	UHMWPE
1951633	UHMWPE
1951634	UHMWPE
1951635	UHMWPE
1951636	UHMWPE
1951637	UHMWPE
1951641	UHMWPE

1951642	UHMWPE
1951643	UHMWPE
1951644	UHMWPE
1951645	UHMWPE
1951646	UHMWPE
1951647	UHMWPE
1951653	UHMWPE
1951654	UHMWPE
1951655	UHMWPE
1951656	UHMWPE
1951657	UHMWPE
1951658	UHMWPE
1951663	UHMWPE
1951664	UHMWPE
1951665	UHMWPE
1951666	UHMWPE
1951667	UHMWPE
1951668	UHMWPE
1951671	UHMWPE
1951672	UHMWPE
1951673	UHMWPE
1951674	UHMWPE
1951675	UHMWPE
1951676	UHMWPE
1951677	UHMWPE
1951681	UHMWPE
1951682	UHMWPE
1951683	UHMWPE
1951684	UHMWPE
1951685	UHMWPE
1951686	UHMWPE
1951687	UHMWPE
1951701	UHMWPE
1951702	UHMWPE
1951703	UHMWPE
1951704	UHMWPE
1951705	UHMWPE
1951706	UHMWPE
1951707	UHMWPE
1951711	UHMWPE
1951712	UHMWPE
1951713	UHMWPE
1951714	UHMWPE

1951715	UHMWPE
1951716	UHMWPE
1951717	UHMWPE
1951731	UHMWPE
1951732	UHMWPE
1951733	UHMWPE
1951734	UHMWPE
1951735	UHMWPE
1951736	UHMWPE
1951737	UHMWPE
1951741	UHMWPE
1951742	UHMWPE
1951743	UHMWPE
1951744	UHMWPE
1951745	UHMWPE
1951746	UHMWPE
1951747	UHMWPE
1951753	UHMWPE
1951754	UHMWPE
1951755	UHMWPE
1951756	UHMWPE
1951757	UHMWPE
1951758	UHMWPE
1951763	UHMWPE
1951764	UHMWPE
1951765	UHMWPE
1951766	UHMWPE
1951767	UHMWPE
1951768	UHMWPE
1951771	UHMWPE
1951772	UHMWPE
1951773	UHMWPE
1951774	UHMWPE
1951775	UHMWPE
1951776	UHMWPE
1951777	UHMWPE
1951781	UHMWPE
1951782	UHMWPE
1951783	UHMWPE
1951784	UHMWPE
1951785	UHMWPE
1951786	UHMWPE
1951787	UHMWPE

1952001	UHMWPE
1952002	UHMWPE
1952003	UHMWPE
1952004	UHMWPE
1952005	UHMWPE
1952006	UHMWPE
1952007	UHMWPE
1952011	UHMWPE
1952012	UHMWPE
1952013	UHMWPE
1952014	UHMWPE
1952015	UHMWPE
1952016	UHMWPE
1952017	UHMWPE
1952031	UHMWPE
1952032	UHMWPE
1952033	UHMWPE
1952034	UHMWPE
1952035	UHMWPE
1952036	UHMWPE
1952037	UHMWPE
1952041	UHMWPE
1952042	UHMWPE
1952043	UHMWPE
1952044	UHMWPE
1952045	UHMWPE
1952046	UHMWPE
1952047	UHMWPE
1952053	UHMWPE
1952054	UHMWPE
1952055	UHMWPE
1952056	UHMWPE
1952057	UHMWPE
1952058	UHMWPE
1952063	UHMWPE
1952064	UHMWPE
1952065	UHMWPE
1952066	UHMWPE
1952067	UHMWPE
1952068	UHMWPE
1952071	UHMWPE
1952072	UHMWPE
1952073	UHMWPE

1952074	UHMWPE
1952075	UHMWPE
1952076	UHMWPE
1952077	UHMWPE
1952081	UHMWPE
1952082	UHMWPE
1952083	UHMWPE
1952084	UHMWPE
1952085	UHMWPE
1952086	UHMWPE
1952087	UHMWPE
1952401	UHMWPE
1952402	UHMWPE
1952403	UHMWPE
1952404	UHMWPE
1952405	UHMWPE
1952406	UHMWPE
1952407	UHMWPE
1952411	UHMWPE
1952412	UHMWPE
1952413	UHMWPE
1952414	UHMWPE
1952415	UHMWPE
1952416	UHMWPE
1952417	UHMWPE
1952431	UHMWPE
1952432	UHMWPE
1952433	UHMWPE
1952434	UHMWPE
1952435	UHMWPE
1952436	UHMWPE
1952437	UHMWPE
1952441	UHMWPE
1952442	UHMWPE
1952443	UHMWPE
1952444	UHMWPE
1952445	UHMWPE
1952446	UHMWPE
1952447	UHMWPE
1952453	UHMWPE
1952454	UHMWPE
1952455	UHMWPE
1952456	UHMWPE

1952457	UHMWPE
1952458	UHMWPE
1952463	UHMWPE
1952464	UHMWPE
1952465	UHMWPE
1952466	UHMWPE
1952467	UHMWPE
1952468	UHMWPE
1952471	UHMWPE
1952472	UHMWPE
1952473	UHMWPE
1952474	UHMWPE
1952475	UHMWPE
1952476	UHMWPE
1952477	UHMWPE
1952481	UHMWPE
1952482	UHMWPE
1952483	UHMWPE
1952484	UHMWPE
1952485	UHMWPE
1952486	UHMWPE
1952487	UHMWPE
1961001	UHMWPE
1961002	UHMWPE
1961003	UHMWPE
1961004	UHMWPE
1961005	UHMWPE
1961006	UHMWPE
1961007	UHMWPE
1961053	UHMWPE
1961054	UHMWPE
1961201	UHMWPE
1961202	UHMWPE
1961203	UHMWPE
1961204	UHMWPE
1961205	UHMWPE
1961206	UHMWPE
1961207	UHMWPE
1961253	UHMWPE
1961254	UHMWPE
1961401	UHMWPE
1961402	UHMWPE
1961403	UHMWPE

1961404	UHMWPE
1961405	UHMWPE
1961406	UHMWPE
1961407	UHMWPE
1961453	UHMWPE
1961454	UHMWPE
1961701	UHMWPE
1961702	UHMWPE
1961703	UHMWPE
1961704	UHMWPE
1961705	UHMWPE
1961706	UHMWPE
1961707	UHMWPE
1961753	UHMWPE
1961754	UHMWPE
1962001	UHMWPE
1962002	UHMWPE
1962003	UHMWPE
1962004	UHMWPE
1962005	UHMWPE
1962006	UHMWPE
1962007	UHMWPE
1962053	UHMWPE
1962054	UHMWPE
1962401	UHMWPE
1962402	UHMWPE
1962403	UHMWPE
1962404	UHMWPE
1962405	UHMWPE
1962406	UHMWPE
1962407	UHMWPE
1962453	UHMWPE
1962454	UHMWPE
1971001	UHMWPE
1971002	UHMWPE
1971003	UHMWPE
1971004	UHMWPE
1971005	UHMWPE
1971006	UHMWPE
1971007	UHMWPE
1971011	UHMWPE
1971012	UHMWPE
1971013	UHMWPE

1971014	UHMWPE
1971015	UHMWPE
1971016	UHMWPE
1971017	UHMWPE
1971053	UHMWPE
1971054	UHMWPE
1971073	UHMWPE
1971074	UHMWPE
1971201	UHMWPE
1971202	UHMWPE
1971203	UHMWPE
1971204	UHMWPE
1971205	UHMWPE
1971206	UHMWPE
1971207	UHMWPE
1971211	UHMWPE
1971212	UHMWPE
1971213	UHMWPE
1971214	UHMWPE
1971215	UHMWPE
1971216	UHMWPE
1971217	UHMWPE
1971253	UHMWPE
1971254	UHMWPE
1971273	UHMWPE
1971274	UHMWPE
1971401	UHMWPE
1971402	UHMWPE
1971403	UHMWPE
1971404	UHMWPE
1971405	UHMWPE
1971406	UHMWPE
1971407	UHMWPE
1971411	UHMWPE
1971412	UHMWPE
1971413	UHMWPE
1971414	UHMWPE
1971415	UHMWPE
1971416	UHMWPE
1971417	UHMWPE
1971453	UHMWPE
1971454	UHMWPE
1971473	UHMWPE

1971474	UHMWPE
1971701	UHMWPE
1971702	UHMWPE
1971703	UHMWPE
1971704	UHMWPE
1971705	UHMWPE
1971706	UHMWPE
1971707	UHMWPE
1971711	UHMWPE
1971712	UHMWPE
1971713	UHMWPE
1971714	UHMWPE
1971715	UHMWPE
1971716	UHMWPE
1971717	UHMWPE
1971753	UHMWPE
1971754	UHMWPE
1971773	UHMWPE
1971774	UHMWPE
1972001	UHMWPE
1972002	UHMWPE
1972003	UHMWPE
1972004	UHMWPE
1972005	UHMWPE
1972006	UHMWPE
1972007	UHMWPE
1972011	UHMWPE
1972012	UHMWPE
1972013	UHMWPE
1972014	UHMWPE
1972015	UHMWPE
1972016	UHMWPE
1972017	UHMWPE
1972053	UHMWPE
1972054	UHMWPE
1972073	UHMWPE
1972074	UHMWPE
1972401	UHMWPE
1972402	UHMWPE
1972403	UHMWPE
1972404	UHMWPE
1972405	UHMWPE
1972406	UHMWPE

1972407	UHMWPE
1972411	UHMWPE
1972412	UHMWPE
1972413	UHMWPE
1972414	UHMWPE
1972415	UHMWPE
1972416	UHMWPE
1972417	UHMWPE
1972453	UHMWPE
1972454	UHMWPE
1972473	UHMWPE
1972474	UHMWPE
1981001	UHMWPE
1981002	UHMWPE
1981003	UHMWPE
1981004	UHMWPE
1981005	UHMWPE
1981006	UHMWPE
1981007	UHMWPE
1981011	UHMWPE
1981012	UHMWPE
1981013	UHMWPE
1981014	UHMWPE
1981015	UHMWPE
1981016	UHMWPE
1981017	UHMWPE
1981031	UHMWPE
1981032	UHMWPE
1981033	UHMWPE
1981034	UHMWPE
1981035	UHMWPE
1981036	UHMWPE
1981037	UHMWPE
1981041	UHMWPE
1981042	UHMWPE
1981043	UHMWPE
1981044	UHMWPE
1981045	UHMWPE
1981046	UHMWPE
1981047	UHMWPE
1981053	UHMWPE
1981054	UHMWPE
1981063	UHMWPE

1981064	UHMWPE
1981071	UHMWPE
1981072	UHMWPE
1981073	UHMWPE
1981074	UHMWPE
1981075	UHMWPE
1981076	UHMWPE
1981077	UHMWPE
1981081	UHMWPE
1981082	UHMWPE
1981083	UHMWPE
1981084	UHMWPE
1981085	UHMWPE
1981086	UHMWPE
1981087	UHMWPE
1981101	UHMWPE
1981102	UHMWPE
1981103	UHMWPE
1981104	UHMWPE
1981105	UHMWPE
1981106	UHMWPE
1981107	UHMWPE
1981111	UHMWPE
1981112	UHMWPE
1981113	UHMWPE
1981114	UHMWPE
1981115	UHMWPE
1981116	UHMWPE
1981117	UHMWPE
1981131	UHMWPE
1981132	UHMWPE
1981133	UHMWPE
1981134	UHMWPE
1981135	UHMWPE
1981136	UHMWPE
1981137	UHMWPE
1981141	UHMWPE
1981142	UHMWPE
1981143	UHMWPE
1981144	UHMWPE
1981145	UHMWPE
1981146	UHMWPE
1981147	UHMWPE

1981153	UHMWPE
1981154	UHMWPE
1981163	UHMWPE
1981164	UHMWPE
1981171	UHMWPE
1981172	UHMWPE
1981173	UHMWPE
1981174	UHMWPE
1981175	UHMWPE
1981176	UHMWPE
1981177	UHMWPE
1981181	UHMWPE
1981182	UHMWPE
1981183	UHMWPE
1981184	UHMWPE
1981185	UHMWPE
1981186	UHMWPE
1981187	UHMWPE
1981201	UHMWPE
1981202	UHMWPE
1981203	UHMWPE
1981204	UHMWPE
1981205	UHMWPE
1981206	UHMWPE
1981207	UHMWPE
1981211	UHMWPE
1981212	UHMWPE
1981213	UHMWPE
1981214	UHMWPE
1981215	UHMWPE
1981216	UHMWPE
1981217	UHMWPE
1981231	UHMWPE
1981232	UHMWPE
1981233	UHMWPE
1981234	UHMWPE
1981235	UHMWPE
1981236	UHMWPE
1981237	UHMWPE
1981241	UHMWPE
1981242	UHMWPE
1981243	UHMWPE
1981244	UHMWPE

1981245	UHMWPE
1981246	UHMWPE
1981247	UHMWPE
1981253	UHMWPE
1981254	UHMWPE
1981263	UHMWPE
1981264	UHMWPE
1981271	UHMWPE
1981272	UHMWPE
1981273	UHMWPE
1981274	UHMWPE
1981275	UHMWPE
1981276	UHMWPE
1981277	UHMWPE
1981281	UHMWPE
1981282	UHMWPE
1981283	UHMWPE
1981284	UHMWPE
1981285	UHMWPE
1981286	UHMWPE
1981287	UHMWPE
1981301	UHMWPE
1981302	UHMWPE
1981303	UHMWPE
1981304	UHMWPE
1981305	UHMWPE
1981306	UHMWPE
1981307	UHMWPE
1981311	UHMWPE
1981312	UHMWPE
1981313	UHMWPE
1981314	UHMWPE
1981315	UHMWPE
1981316	UHMWPE
1981317	UHMWPE
1981331	UHMWPE
1981332	UHMWPE
1981333	UHMWPE
1981334	UHMWPE
1981335	UHMWPE
1981336	UHMWPE
1981337	UHMWPE
1981341	UHMWPE

1981342	UHMWPE
1981343	UHMWPE
1981344	UHMWPE
1981345	UHMWPE
1981346	UHMWPE
1981347	UHMWPE
1981353	UHMWPE
1981354	UHMWPE
1981363	UHMWPE
1981364	UHMWPE
1981371	UHMWPE
1981372	UHMWPE
1981373	UHMWPE
1981374	UHMWPE
1981375	UHMWPE
1981376	UHMWPE
1981377	UHMWPE
1981381	UHMWPE
1981382	UHMWPE
1981383	UHMWPE
1981384	UHMWPE
1981385	UHMWPE
1981386	UHMWPE
1981387	UHMWPE
1981401	UHMWPE
1981402	UHMWPE
1981403	UHMWPE
1981404	UHMWPE
1981405	UHMWPE
1981406	UHMWPE
1981407	UHMWPE
1981411	UHMWPE
1981412	UHMWPE
1981413	UHMWPE
1981414	UHMWPE
1981415	UHMWPE
1981416	UHMWPE
1981417	UHMWPE
1981431	UHMWPE
1981432	UHMWPE
1981433	UHMWPE
1981434	UHMWPE
1981435	UHMWPE

1981436	UHMWPE
1981437	UHMWPE
1981441	UHMWPE
1981442	UHMWPE
1981443	UHMWPE
1981444	UHMWPE
1981445	UHMWPE
1981446	UHMWPE
1981447	UHMWPE
1981453	UHMWPE
1981454	UHMWPE
1981463	UHMWPE
1981464	UHMWPE
1981471	UHMWPE
1981472	UHMWPE
1981473	UHMWPE
1981474	UHMWPE
1981475	UHMWPE
1981476	UHMWPE
1981477	UHMWPE
1981481	UHMWPE
1981482	UHMWPE
1981483	UHMWPE
1981484	UHMWPE
1981485	UHMWPE
1981486	UHMWPE
1981487	UHMWPE
1981501	UHMWPE
1981502	UHMWPE
1981503	UHMWPE
1981504	UHMWPE
1981505	UHMWPE
1981506	UHMWPE
1981507	UHMWPE
1981511	UHMWPE
1981512	UHMWPE
1981513	UHMWPE
1981514	UHMWPE
1981515	UHMWPE
1981516	UHMWPE
1981517	UHMWPE
1981531	UHMWPE
1981532	UHMWPE

1981533	UHMWPE
1981534	UHMWPE
1981535	UHMWPE
1981536	UHMWPE
1981537	UHMWPE
1981541	UHMWPE
1981542	UHMWPE
1981543	UHMWPE
1981544	UHMWPE
1981545	UHMWPE
1981546	UHMWPE
1981547	UHMWPE
1981553	UHMWPE
1981554	UHMWPE
1981563	UHMWPE
1981564	UHMWPE
1981571	UHMWPE
1981572	UHMWPE
1981573	UHMWPE
1981574	UHMWPE
1981575	UHMWPE
1981576	UHMWPE
1981577	UHMWPE
1981581	UHMWPE
1981582	UHMWPE
1981583	UHMWPE
1981584	UHMWPE
1981585	UHMWPE
1981586	UHMWPE
1981587	UHMWPE
1981601	UHMWPE
1981602	UHMWPE
1981603	UHMWPE
1981604	UHMWPE
1981605	UHMWPE
1981606	UHMWPE
1981607	UHMWPE
1981611	UHMWPE
1981612	UHMWPE
1981613	UHMWPE
1981614	UHMWPE
1981615	UHMWPE
1981616	UHMWPE

1981617	UHMWPE
1981631	UHMWPE
1981632	UHMWPE
1981633	UHMWPE
1981634	UHMWPE
1981635	UHMWPE
1981636	UHMWPE
1981637	UHMWPE
1981641	UHMWPE
1981642	UHMWPE
1981643	UHMWPE
1981644	UHMWPE
1981645	UHMWPE
1981646	UHMWPE
1981647	UHMWPE
1981653	UHMWPE
1981654	UHMWPE
1981663	UHMWPE
1981664	UHMWPE
1981671	UHMWPE
1981672	UHMWPE
1981673	UHMWPE
1981674	UHMWPE
1981675	UHMWPE
1981676	UHMWPE
1981677	UHMWPE
1981681	UHMWPE
1981682	UHMWPE
1981683	UHMWPE
1981684	UHMWPE
1981685	UHMWPE
1981686	UHMWPE
1981687	UHMWPE
1981701	UHMWPE
1981702	UHMWPE
1981703	UHMWPE
1981704	UHMWPE
1981705	UHMWPE
1981706	UHMWPE
1981707	UHMWPE
1981711	UHMWPE
1981712	UHMWPE
1981713	UHMWPE

1981714	UHMWPE
1981715	UHMWPE
1981716	UHMWPE
1981717	UHMWPE
1981731	UHMWPE
1981732	UHMWPE
1981733	UHMWPE
1981734	UHMWPE
1981735	UHMWPE
1981736	UHMWPE
1981737	UHMWPE
1981741	UHMWPE
1981742	UHMWPE
1981743	UHMWPE
1981744	UHMWPE
1981745	UHMWPE
1981746	UHMWPE
1981747	UHMWPE
1981753	UHMWPE
1981754	UHMWPE
1981763	UHMWPE
1981764	UHMWPE
1981771	UHMWPE
1981772	UHMWPE
1981773	UHMWPE
1981774	UHMWPE
1981775	UHMWPE
1981776	UHMWPE
1981777	UHMWPE
1981781	UHMWPE
1981782	UHMWPE
1981783	UHMWPE
1981784	UHMWPE
1981785	UHMWPE
1981786	UHMWPE
1981787	UHMWPE
1982001	UHMWPE
1982002	UHMWPE
1982003	UHMWPE
1982004	UHMWPE
1982005	UHMWPE
1982006	UHMWPE
1982007	UHMWPE

1982011	UHMWPE
1982012	UHMWPE
1982013	UHMWPE
1982014	UHMWPE
1982015	UHMWPE
1982016	UHMWPE
1982017	UHMWPE
1982031	UHMWPE
1982032	UHMWPE
1982033	UHMWPE
1982034	UHMWPE
1982035	UHMWPE
1982036	UHMWPE
1982037	UHMWPE
1982041	UHMWPE
1982042	UHMWPE
1982043	UHMWPE
1982044	UHMWPE
1982045	UHMWPE
1982046	UHMWPE
1982047	UHMWPE
1982053	UHMWPE
1982054	UHMWPE
1982063	UHMWPE
1982064	UHMWPE
1982071	UHMWPE
1982072	UHMWPE
1982073	UHMWPE
1982074	UHMWPE
1982075	UHMWPE
1982076	UHMWPE
1982077	UHMWPE
1982081	UHMWPE
1982082	UHMWPE
1982083	UHMWPE
1982084	UHMWPE
1982085	UHMWPE
1982086	UHMWPE
1982087	UHMWPE
1982401	UHMWPE
1982402	UHMWPE
1982403	UHMWPE
1982404	UHMWPE

1982405	UHMWPE
1982406	UHMWPE
1982407	UHMWPE
1982411	UHMWPE
1982412	UHMWPE
1982413	UHMWPE
1982414	UHMWPE
1982415	UHMWPE
1982416	UHMWPE
1982417	UHMWPE
1982431	UHMWPE
1982432	UHMWPE
1982433	UHMWPE
1982434	UHMWPE
1982435	UHMWPE
1982436	UHMWPE
1982437	UHMWPE
1982441	UHMWPE
1982442	UHMWPE
1982443	UHMWPE
1982444	UHMWPE
1982445	UHMWPE
1982446	UHMWPE
1982447	UHMWPE
1982453	UHMWPE
1982454	UHMWPE
1982463	UHMWPE
1982464	UHMWPE
1982471	UHMWPE
1982472	UHMWPE
1982473	UHMWPE
1982474	UHMWPE
1982475	UHMWPE
1982476	UHMWPE
1982477	UHMWPE
1982481	UHMWPE
1982482	UHMWPE
1982483	UHMWPE
1982484	UHMWPE
1982485	UHMWPE
1982486	UHMWPE
1982487	UHMWPE
2200321	Ti6Al4V Alloy

2200324	Ti6Al4V Alloy
2200327	Ti6Al4V Alloy
2200330	Ti6Al4V Alloy
2200333	Ti6Al4V Alloy
2200336	Ti6Al4V Alloy
2200338	Ti6Al4V Alloy
2200340	Ti6Al4V Alloy
2200421	Ti6Al4V Alloy
2200424	Ti6Al4V Alloy
2200427	Ti6Al4V Alloy
2200430	Ti6Al4V Alloy
2200433	Ti6Al4V Alloy
2200436	Ti6Al4V Alloy
2200438	Ti6Al4V Alloy
2200440	Ti6Al4V Alloy
2200521	Ti6Al4V Alloy
2200524	Ti6Al4V Alloy
2200527	Ti6Al4V Alloy
2200530	Ti6Al4V Alloy
2200533	Ti6Al4V Alloy
2200536	Ti6Al4V Alloy
2200538	Ti6Al4V Alloy
2200540	Ti6Al4V Alloy
2200721	Ti6Al4V Alloy
2200724	Ti6Al4V Alloy
2200727	Ti6Al4V Alloy
2200730	Ti6Al4V Alloy
2200733	Ti6Al4V Alloy
2200736	Ti6Al4V Alloy
2200738	Ti6Al4V Alloy
2200740	Ti6Al4V Alloy
2200921	Ti6Al4V Alloy
2200924	Ti6Al4V Alloy
2200927	Ti6Al4V Alloy
2200930	Ti6Al4V Alloy
2200933	Ti6Al4V Alloy
2200936	Ti6Al4V Alloy
2200938	Ti6Al4V Alloy
2200940	Ti6Al4V Alloy
2210321	INOX STEEL AISI 316L
2210324	INOX STEEL AISI 316L
2210327	INOX STEEL AISI 316L
2210330	INOX STEEL AISI 316L

2210333	INOX STEEL AISI 316L
2210336	INOX STEEL AISI 316L
2210338	INOX STEEL AISI 316L
2210340	INOX STEEL AISI 316L
2210421	INOX STEEL AISI 316L
2210424	INOX STEEL AISI 316L
2210427	INOX STEEL AISI 316L
2210430	INOX STEEL AISI 316L
2210433	INOX STEEL AISI 316L
2210436	INOX STEEL AISI 316L
2210438	INOX STEEL AISI 316L
2210440	INOX STEEL AISI 316L
2210521	INOX STEEL AISI 316L
2210524	INOX STEEL AISI 316L
2210527	INOX STEEL AISI 316L
2210530	INOX STEEL AISI 316L
2210533	INOX STEEL AISI 316L
2210536	INOX STEEL AISI 316L
2210538	INOX STEEL AISI 316L
2210540	INOX STEEL AISI 316L
2210721	INOX STEEL AISI 316L
2210724	INOX STEEL AISI 316L
2210727	INOX STEEL AISI 316L
2210730	INOX STEEL AISI 316L
2210733	INOX STEEL AISI 316L
2210736	INOX STEEL AISI 316L
2210738	INOX STEEL AISI 316L
2210740	INOX STEEL AISI 316L
2210921	INOX STEEL AISI 316L
2210924	INOX STEEL AISI 316L
2210927	INOX STEEL AISI 316L
2210930	INOX STEEL AISI 316L
2210933	INOX STEEL AISI 316L
2210936	INOX STEEL AISI 316L
2210938	INOX STEEL AISI 316L
2210940	INOX STEEL AISI 316L
2300825	Ti6Al4V Alloy
2300830	Ti6Al4V Alloy
2300835	Ti6Al4V Alloy
2300840	Ti6Al4V Alloy
2300845	Ti6Al4V Alloy
2300850	Ti6Al4V Alloy
2306525	Ti6Al4V Alloy

2306530	Ti6Al4V Alloy
2306535	Ti6Al4V Alloy
2306540	Ti6Al4V Alloy
2310825	INOX STEEL AISI 316L
2310830	INOX STEEL AISI 316L
2310835	INOX STEEL AISI 316L
2310840	INOX STEEL AISI 316L
2310845	INOX STEEL AISI 316L
2310850	INOX STEEL AISI 316L
2316525	INOX STEEL AISI 316L
2316530	INOX STEEL AISI 316L
2316535	INOX STEEL AISI 316L
2316540	INOX STEEL AISI 316L
2500010	Co-Cr-Mo Alloy + Hydroxyapatite
2500020	Co-Cr-Mo Alloy + Hydroxyapatite
2500030	Co-Cr-Mo Alloy + Hydroxyapatite
2500110	Ti6Al4V Alloy + TiNbN
2500120	Ti6Al4V Alloy + TiNbN
2500130	Ti6Al4V Alloy + TiNbN
2501105	UHMWPE
2501106	UHMWPE
2501107	UHMWPE
2501108	UHMWPE
2501125	UHMWPE
2501126	UHMWPE
2501127	UHMWPE
2501128	UHMWPE
2501135	UHMWPE
2501136	UHMWPE
2501137	UHMWPE
2501138	UHMWPE
2501205	UHMWPE
2501206	UHMWPE
2501207	UHMWPE
2501208	UHMWPE
2501225	UHMWPE
2501226	UHMWPE
2501227	UHMWPE
2501228	UHMWPE
2501235	UHMWPE
2501236	UHMWPE
2501237	UHMWPE
2501238	UHMWPE

2501305	UHMWPE
2501306	UHMWPE
2501307	UHMWPE
2501308	UHMWPE
2501325	UHMWPE
2501326	UHMWPE
2501327	UHMWPE
2501328	UHMWPE
2501335	UHMWPE
2501336	UHMWPE
2501337	UHMWPE
2501338	UHMWPE
2502010	Co-Cr-Mo Alloy + Hydroxyapatite
2502020	Co-Cr-Mo Alloy + Hydroxyapatite
2502030	Co-Cr-Mo Alloy + Hydroxyapatite
2502110	Ti6Al4V Alloy + TiNbN
2502120	Ti6Al4V Alloy + TiNbN
2502130	Ti6Al4V Alloy + TiNbN
2800010	Ti6Al4V Alloy
2800110	Ti6Al4V Alloy + TiNbN
2800120	Ti6Al4V Alloy + TiNbN
2800210	UHMWPE
2800220	UHMWPE
2800230	UHMWPE
2800310	UHMWPE
2800320	UHMWPE
2800330	UHMWPE
2800420	UHMWPE
2800430	UHMWPE
2800460	UHMWPE
2800470	UHMWPE
2800506	Ti6Al4V Alloy
2800510	Ti6Al4V Alloy
2800610	Ti6Al4V Alloy + TiNbN
2800620	Ti6Al4V Alloy + TiNbN
2800630	Ti6Al4V Alloy + TiNbN
2800710	Ti6Al4V Alloy + TiNbN
2800720	Ti6Al4V Alloy + TiNbN
2800730	Ti6Al4V Alloy + TiNbN
2800810	Ti6Al4V Alloy
2800820	Ti6Al4V Alloy
2800830	Ti6Al4V Alloy
2800915	Ti6Al4V Alloy

2800920	Ti6Al4V Alloy
2800925	Ti6Al4V Alloy
2800930	Ti6Al4V Alloy
2800935	Ti6Al4V Alloy
2800940	Ti6Al4V Alloy
2800950	Ti6Al4V Alloy
2801006	Ti6Al4V Alloy
2801008	Ti6Al4V Alloy
2801106	Ti6Al4V Alloy
2801108	Ti6Al4V Alloy
2801110	Ti6Al4V Alloy
2801112	Ti6Al4V Alloy
2801114	Ti6Al4V Alloy
2801116	Ti6Al4V Alloy
2801210	Ti6Al4V Alloy
2801310	Ti6Al4V Alloy
2801320	Ti6Al4V Alloy
2801330	Ti6Al4V Alloy
2801406	Ti6Al4V Alloy
2801408	Ti6Al4V Alloy
2801410	Ti6Al4V Alloy
2801412	Ti6Al4V Alloy
2801414	Ti6Al4V Alloy
2801416	Ti6Al4V Alloy
2801506	Ti6Al4V Alloy
2801508	Ti6Al4V Alloy
2801510	Ti6Al4V Alloy
2801512	Ti6Al4V Alloy
2801514	Ti6Al4V Alloy
2801516	Ti6Al4V Alloy
2801600	Ti6Al4V Alloy
2801650	Ti6Al4V Alloy
2802035	Ti6Al4V Alloy
2802045	Ti6Al4V Alloy
2802055	Ti6Al4V Alloy
2802065	Ti6Al4V Alloy
2802075	Ti6Al4V Alloy
2802085	Ti6Al4V Alloy
2802095	Ti6Al4V Alloy
2802105	Ti6Al4V Alloy
2881406	Ti6Al4V Alloy
2881408	Ti6Al4V Alloy
2881410	Ti6Al4V Alloy

2881412	Ti6Al4V Alloy
2881414	Ti6Al4V Alloy
2881416	Ti6Al4V Alloy
3100000	Ti6Al4V Alloy
3100001	Ti6Al4V Alloy
3100002	Ti6Al4V Alloy
3100003	Ti6Al4V Alloy
3100004	Ti6Al4V Alloy
3100005	Ti6Al4V Alloy
3100006	Ti6Al4V Alloy
3100007	Ti6Al4V Alloy
3100010	Ti6Al4V Alloy
3100011	Ti6Al4V Alloy
3100012	Ti6Al4V Alloy
3100013	Ti6Al4V Alloy
3100014	Ti6Al4V Alloy
3100015	Ti6Al4V Alloy
3100016	Ti6Al4V Alloy
3100017	Ti6Al4V Alloy
3100105	Ti6Al4V Alloy
3100115	Ti6Al4V Alloy
3100125	Ti6Al4V Alloy
3100135	Ti6Al4V Alloy
3100145	Ti6Al4V Alloy
3100155	Ti6Al4V Alloy
3100165	Ti6Al4V Alloy
3100205	Ti6Al4V Alloy
3100215	Ti6Al4V Alloy
3100225	Ti6Al4V Alloy
3100235	Ti6Al4V Alloy
3100245	Ti6Al4V Alloy
3100255	Ti6Al4V Alloy
3100265	Ti6Al4V Alloy
3180001	Ti6Al4V Alloy
3180002	Ti6Al4V Alloy
3180003	Ti6Al4V Alloy
3180004	Ti6Al4V Alloy
3180005	Ti6Al4V Alloy

3180006	Ti6Al4V Alloy
3180007	Ti6Al4V Alloy
3180011	Ti6Al4V Alloy
3180012	Ti6Al4V Alloy
3180013	Ti6Al4V Alloy
3180014	Ti6Al4V Alloy
3180015	Ti6Al4V Alloy
3180016	Ti6Al4V Alloy
3180017	Ti6Al4V Alloy
3200101	Ti6Al4V Alloy + TiNbN
3200102	Ti6Al4V Alloy + TiNbN
3200103	Ti6Al4V Alloy + TiNbN
3200104	Ti6Al4V Alloy + TiNbN
3200105	Ti6Al4V Alloy + TiNbN
3200106	Ti6Al4V Alloy + TiNbN
3300001	UHMWPE + Ti6Al4V Alloy + TiNbN
3300002	UHMWPE + Ti6Al4V Alloy + TiNbN
3300003	UHMWPE + Ti6Al4V Alloy + TiNbN
3300004	UHMWPE + Ti6Al4V Alloy + TiNbN
3300005	UHMWPE + Ti6Al4V Alloy + TiNbN
3300006	UHMWPE + Ti6Al4V Alloy + TiNbN
3300101	UHMWPE + Ti6Al4V Alloy + TiNbN
3300102	UHMWPE + Ti6Al4V Alloy + TiNbN
3300103	UHMWPE + Ti6Al4V Alloy + TiNbN
3300104	UHMWPE + Ti6Al4V Alloy + TiNbN
3300105	UHMWPE + Ti6Al4V Alloy + TiNbN
3300106	UHMWPE + Ti6Al4V Alloy + TiNbN
3300201	UHMWPE + Ti6Al4V Alloy + TiNbN
3300202	UHMWPE + Ti6Al4V Alloy + TiNbN
3300203	UHMWPE + Ti6Al4V Alloy + TiNbN
3300204	UHMWPE + Ti6Al4V Alloy + TiNbN
3300205	UHMWPE + Ti6Al4V Alloy + TiNbN
3300206	UHMWPE + Ti6Al4V Alloy + TiNbN
3300301	UHMWPE + Ti6Al4V Alloy + TiNbN
3300302	UHMWPE + Ti6Al4V Alloy + TiNbN
3300303	UHMWPE + Ti6Al4V Alloy + TiNbN
3300304	UHMWPE + Ti6Al4V Alloy + TiNbN
3300305	UHMWPE + Ti6Al4V Alloy + TiNbN
3300306	UHMWPE + Ti6Al4V Alloy + TiNbN
3301001	Co-Cr-Mo Alloy
3301002	Co-Cr-Mo Alloy
3301003	Co-Cr-Mo Alloy
3301004	Co-Cr-Mo Alloy

3301005	Co-Cr-Mo Alloy
3301006	Co-Cr-Mo Alloy
3301101	Co-Cr-Mo Alloy
3301102	Co-Cr-Mo Alloy
3301103	Co-Cr-Mo Alloy
3301104	Co-Cr-Mo Alloy
3301105	Co-Cr-Mo Alloy
3301106	Co-Cr-Mo Alloy
3301151	Ti6Al4V Alloy + TiNbN
3301152	Ti6Al4V Alloy + TiNbN
3301153	Ti6Al4V Alloy + TiNbN
3301154	Ti6Al4V Alloy + TiNbN
3301155	Ti6Al4V Alloy + TiNbN
3301156	Ti6Al4V Alloy + TiNbN
3301201	Co-Cr-Mo Alloy
3301202	Co-Cr-Mo Alloy
3301203	Co-Cr-Mo Alloy
3301204	Co-Cr-Mo Alloy
3301205	Co-Cr-Mo Alloy
3301206	Co-Cr-Mo Alloy
3301251	Ti6Al4V Alloy + TiNbN
3301252	Ti6Al4V Alloy + TiNbN
3301253	Ti6Al4V Alloy + TiNbN
3301254	Ti6Al4V Alloy + TiNbN
3301255	Ti6Al4V Alloy + TiNbN
3301256	Ti6Al4V Alloy + TiNbN
3301301	Co-Cr-Mo Alloy
3301302	Co-Cr-Mo Alloy
3301303	Co-Cr-Mo Alloy
3301304	Co-Cr-Mo Alloy
3301305	Co-Cr-Mo Alloy
3301306	Co-Cr-Mo Alloy
3301351	Ti6Al4V Alloy + TiNbN
3301352	Ti6Al4V Alloy + TiNbN
3301353	Ti6Al4V Alloy + TiNbN
3301354	Ti6Al4V Alloy + TiNbN
3301355	Ti6Al4V Alloy + TiNbN
3301356	Ti6Al4V Alloy + TiNbN
3301451	Ti6Al4V Alloy + TiNbN
3301452	Ti6Al4V Alloy + TiNbN
3301453	Ti6Al4V Alloy + TiNbN
3301454	Ti6Al4V Alloy + TiNbN
3301455	Ti6Al4V Alloy + TiNbN

3301456	Ti6Al4V Alloy + TiNbN
3301551	Ti6Al4V Alloy + TiNbN
3301552	Ti6Al4V Alloy + TiNbN
3301553	Ti6Al4V Alloy + TiNbN
3301554	Ti6Al4V Alloy + TiNbN
3301555	Ti6Al4V Alloy + TiNbN
3301556	Ti6Al4V Alloy + TiNbN
3301651	Ti6Al4V Alloy + TiNbN
3301652	Ti6Al4V Alloy + TiNbN
3301653	Ti6Al4V Alloy + TiNbN
3301654	Ti6Al4V Alloy + TiNbN
3301655	Ti6Al4V Alloy + TiNbN
3301656	Ti6Al4V Alloy + TiNbN
3400000	Ti6Al4V Alloy + TiNbN
3400001	Ti6Al4V Alloy + TiNbN
3400002	Ti6Al4V Alloy + TiNbN
3400003	Ti6Al4V Alloy + TiNbN
3400004	Ti6Al4V Alloy + TiNbN
3400009	Ti6Al4V Alloy + TiNbN
3400010	Ti6Al4V Alloy + TiNbN
3400011	Ti6Al4V Alloy + TiNbN
3400012	Ti6Al4V Alloy + TiNbN
3400013	Ti6Al4V Alloy + TiNbN
3400014	Ti6Al4V Alloy + TiNbN
3400019	Ti6Al4V Alloy + TiNbN
3400100	Ti6Al4V Alloy + TiNbN
3400101	Ti6Al4V Alloy + TiNbN
3400102	Ti6Al4V Alloy + TiNbN
3400103	Ti6Al4V Alloy + TiNbN
3400104	Ti6Al4V Alloy + TiNbN
3400109	Ti6Al4V Alloy + TiNbN
3400110	Ti6Al4V Alloy + TiNbN
3400111	Ti6Al4V Alloy + TiNbN
3400112	Ti6Al4V Alloy + TiNbN
3400113	Ti6Al4V Alloy + TiNbN
3400114	Ti6Al4V Alloy + TiNbN
3400119	Ti6Al4V Alloy + TiNbN
3400202	Ti6Al4V Alloy + TiNbN
3400203	Ti6Al4V Alloy + TiNbN
3400204	Ti6Al4V Alloy + TiNbN